









Chemistry 344: Spectroscopy Problem Set 2

Name (print): _____

(Not for credit)


TA Name (print): _____


- I. For each of the following molecules and fragments below, predict the multiplicity of each of the signals that you would expect to see in an $^1\text{H-NMR}$ spectrum. Include an approximation of the expected coupling value in Hz. See the example below.



- II. For each of the $^1\text{H-NMR}$ spectra below, determine the structure of the molecule responsible for it and assign the protons in the molecule to its corresponding $^1\text{H-NMR}$ signal.

300 MHz $^1\text{H NMR}$
In CDCl_3







The signals at δ 6.60 and δ 5.10 ppm are shown below with an axis Hz to assist in determining the coupling constant.




- III. The three $^1\text{H-NMR}$ spectra below from 7.0 – 7.7 ppm correspond to three chloronitrobenzenes. Assign each one to its molecule using the coupling patterns.



A.




B.



C.

IV. While $^1\text{H-NMR}$ is a powerful tool in structure determination, it can also be used to determine the relative ratio of two molecules in a mixture. This is particularly useful for determining product ratios in chemical reactions where more than one product is generated. The following $^1\text{H-NMR}$ spectrum is of a crude mixture of two isomers.

- A. Given that the signal at δ 8.10 ppm (36.61 H) is due to two H-atoms in the major product and the signals at δ 7.85 (5.67 H) is due to a single H-atom in the minor product, assign the remaining signals to either the major or minor isomer.



- B. Based upon your assignments, what is the product ratio of major to minor isomer?

- V. Determine the structures of the following compounds based on the $^1\text{H-NMR}$ and other spectroscopy data provided. At first, it may help to sketch the spectrum. (300 MHz in CDCl_3)
- $\text{C}_4\text{H}_8\text{O}_2$** : δ 1.15 (triplet, 3H), δ 2.33 (quartet, 2H), and δ 3.67 (singlet, 3H)
 - $\text{C}_9\text{H}_{10}\text{O}$** : δ 1.22 (triplet, 3H), δ 2.98 (quartet, 2H), δ 7.43 (multiplet, 2H), δ 7.53 (triplet of triplets, 1H), and δ 7.94 (multiplet, 2H)
 - $\text{C}_{10}\text{H}_{14}$** : δ 1.3 (singlet, 9H), δ 7.17 (triplet of triplets, 1H), δ 7.29 (triplet, 2H), and δ 7.38 (doublet, 2H)
 - $\text{C}_{10}\text{H}_{12}\text{O}_2$** : δ 2.01 (singlet, 3H), δ 2.92 (triplet, 2H), δ 4.27 (triplet, 2H), and δ 7.21 (multiplet, 5H)
IR spectrum contains a strong absorption at 1740 cm^{-1} MS contains a strong signal at $m/z = 43$
 - $\text{C}_8\text{H}_7\text{N}$** : δ 3.70 (singlet, 2H), and δ 7.3 (multiplet, 5H) IR spectrum contains a strong absorption at $\sim 2250\text{ cm}^{-1}$
 - $\text{C}_4\text{H}_6\text{Cl}_2\text{O}_2$** : δ 1.36 (triplet, 3H), δ 4.34 (quartet, 2H), and δ 5.96 (singlet, 1H)
 - $\text{C}_7\text{H}_{14}\text{O}$** : δ 0.91 (triplet, 6H), δ 1.60 (sextet, 4H), and δ 2.37 (triplet, 4H)
 - $\text{C}_5\text{H}_{10}\text{O}_2$** : δ 1.23 (doublet, 6H), δ 2.02 (singlet, 3H), and δ 4.99 (septet, 1H)
 - $\text{C}_6\text{H}_{12}\text{O}_2$** : δ 1.44 (singlet, 9H), and δ 1.96 (singlet, 3H)
 - $\text{C}_9\text{H}_{10}\text{O}_2$** : δ 2.08 (singlet, 3H), δ 5.09 (singlet, 2H), and δ 7.34 (multiplet, 5H)
 - $\text{C}_{10}\text{H}_{12}\text{O}_2$** : δ 1.23 (triplet, 3H), δ 3.60 (singlet, 2H), δ 4.13 (quartet, 2H), and δ 7.28 (multiplet, 5H)
 - $\text{C}_{10}\text{H}_{12}\text{O}_2$** : δ 1.52 (doublet, 3H), δ 2.05 (singlet, 3H), δ 5.87 (quartet, 1H), and δ 7.30 (multiplet, 5H) IR spectrum contains a strong absorption at 1742 cm^{-1}
 - $\text{C}_4\text{H}_8\text{O}_2$** : δ 1.39 (doublet, 3H), δ 2.22 (singlet, 3H), δ 3.91 (broad singlet, 1H), and δ 4.27 (quartet, 1H) IR spectrum contains a strong absorption at 3451 cm^{-1}
 - $\text{C}_{10}\text{H}_{12}\text{O}$** : δ 1.01 (triplet, 3H), δ 2.45 (quartet, 2H), δ 3.67 (singlet, 2H), and δ 7.24 (multiplet, 5H)
 - $\text{C}_2\text{H}_4\text{Br}_2$** : δ 2.46 (doublet, 3H), and δ 5.84 (quartet, 1H)
 - $\text{C}_3\text{H}_6\text{Br}_2$** : δ 2.35 (quintet, 2H), and δ 3.56 (triplet, 4H)
 - $\text{C}_8\text{H}_9\text{Br}$** : δ 2.02 (doublet, 3H), δ 5.18 (quartet, 1H), δ 7.29 (asymmetrical doublet, 2H), and δ 7.41 (multiplet, 3H)
 - $\text{C}_{14}\text{H}_{14}$** : δ 2.90 (singlet, 4H), δ 7.17 (multiplet, 6H), and δ 7.26 (multiplet, 4H)
 - $\text{C}_{11}\text{H}_{17}\text{N}$** : δ 1.04 (triplet, 6H), δ 2.51 (quartet, 4H), δ 3.56 (singlet, 2H), and δ 7.33 (multiplet, 5H)
 - $\text{C}_3\text{H}_5\text{ClO}_2$** : δ 1.74 (doublet, 3H), δ 4.45 (quartet, 1H), and δ 12.2 (singlet, 1H)
 - $\text{C}_3\text{H}_5\text{ClO}_2$** : δ 2.87 (triplet, 2H), δ 3.76 (triplet, 2H), and δ 11.8 (singlet, 1H)

- v) $\text{C}_{10}\text{H}_{14}$: δ 1.22 (doublet, 6H), δ 2.30 (singlet, 3H), δ 2.86 (septet, 1H), and δ 7.0 (symmetrical multiplet, 4H)
- w) $\text{C}_7\text{H}_{12}\text{O}_4$: δ 1.29 (triplet, 6H), δ 3.36 (singlet, 2H), and δ 4.22 (quartet, 4H)
- x) $\text{C}_9\text{H}_8\text{O}$: δ 3.54 (singlet, 4H), and δ 7.26 (symmetrical multiplet, 4H)
- y) C_9H_{10} : C_9H_{10} : δ 2.04 (quintet, 2H), δ 2.90 (triplet, 4H), δ 7.10 (multiplet, 2H), and δ 7.20 (multiplet, 2H)
- z) $\text{C}_3\text{H}_4\text{O}_2$: δ 3.55 (triplet, 2H), and 4.28 (triplet, 2H)
- aa) $\text{C}_5\text{H}_9\text{BrO}_2$: δ 1.29 (triplet, 3H), δ 2.92 (triplet, 2H), δ 3.58 (triplet, 2H), and δ 4.19 (quartet, 2H)
- bb) $\text{C}_5\text{H}_9\text{BrO}_2$: δ 1.17 (triplet, 3H), δ 1.63 (doublet, 3H), δ 3.79 (quartet, 1H), and δ 3.91 (quartet, 2H)
- cc) $\text{C}_6\text{H}_{13}\text{NO}_2$: δ 1.29 (triplet, 3H), δ 2.36 (singlet, 6H), δ 3.16 (singlet, 2H), and δ 4.20 (quartet, 2H)
IR spectrum contains a strong absorption at 1749 cm^{-1}


VI. For each of the spectra on the subsequent pages, determine the structure of the molecule and assign the ^1H -NMR and ^{13}C -NMR signals where possible. These exercises are authentic spectra taken with the NMR instrument used in the organic laboratory courses. Careful analysis of these spectra is great preparation for your use of NMR spectroscopy this semester.

CHEM 344 Unknown W

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_21Jan2013

Pulse Sequence: s2pu1

Unknown W
 $C_3H_5O_2Br$ 1H -NMR



Unknown W

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_21Jan2013


Pulse Sequence: s2pu1

Solvent: cdcl3
Ambient temperature
File: 5502
Mercury-300BB "orglabnmr"

Relax. delay 1.000 sec
Pulse 45.0 degrees
Acq. time 1.000 sec
Width 18867.9 Hz
128 repetitions
OBSERVE C13, 75.4428085 MHz
DECOUPLE H1, 300.0323276 MHz
Power 38 dB
continuously on
WALTZ-16 modulated
DATA PROCESSING
Line broadening 1.0 Hz
FT size 65536
Total time 8 min, 49 sec

Unknown W
 $C_3H_5O_2Br$ ^{13}C -NMR

See Course Website or Loudon
page 622-629 for more
information on ^{13}C -NMR.




CHEM 344 Unknown K

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_13Jun2012

Pulse Sequence: s2pu1

Unknown K
 $C_{12}H_{14}O_4$ 1H -NMR




CHEM 344 Unknown K

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_13Jun2012

Pulse Sequence: s2pu1

Unknown K
 $C_{12}H_{14}O_4$ 1H -NMR




CHEM 344 Unknown I

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_13Jun2008

Unknown I
 $C_{11}H_{12}O_3$ 1H -NMR

Pulse Sequence: s2pu1




CHEM 344 Unknown I

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_13Jun2008

Pulse Sequence: s2pu1

Unknown I
 $C_{11}H_{12}O_3$ 1H -NMR





CHEM 344 Unknown T

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_20Sep2011

Pulse Sequence: s2pu1

Unknown T
 $C_{10}H_{10}O_4$ 1H -NMR





CHEM 344 Unknown P

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_21Jan2013

Pulse Sequence: s2pu1

Unknown P
 $C_9H_9O_3Br$ 1H -NMR





Unknown P

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_21Jan2013

Pulse Sequence: s2pu1
Solvent: C6D6
Ambient temperature
File: 5102
Mercury-300BB "orglabnmr"

Relax. delay 1.000 sec
Pulse 45.0 degrees
Acq. time 1.000 sec
Width 18867.9 Hz
256 repetitions
OBSERVE C13, 75.4428153 MHz
DECOUPLE H1, 300.0323546 MHz
Power 38 dB
continuously on
WALTZ-16 modulated
DATA PROCESSING
FT size 65536
Total time 8 min, 49 sec

Unknown P
 $C_9H_9O_3Br$ ^{13}C -NMR



Unknown P

Archive directory: /export/home/orglabTA/vnmrsys/data
Sample directory: auto_21Jan2013

Pulse Sequence: s2pu1
Solvent: C6D6
Ambient temperature
File: 5102
Mercury-300BB "orglabnmr"

Relax. delay 1.000 sec
Pulse 45.0 degrees
Acq. time 1.000 sec
Width 18867.9 Hz
256 repetitions
OBSERVE C13, 75.4428153 MHz
DECOUPLE H1, 300.0323546 MHz
Power 38 dB
continuously on
WALTZ-16 modulated
DATA PROCESSING
FT size 65536
Total time 8 min, 49 sec

Unknown P
 $C_9H_9O_3Br$ ^{13}C -NMR

