
## **Building a Photometer to Measure Absorbance College For Kids**

In a general sense, a photometer is an instrument that measures the power of a source of light. The measurement of light power is extremely useful in chemistry, since it allows us to measure light that has been absorbed, scattered, or emitted by a chemical sample, which in turn provides information about the structure and the amount of the molecules in the sample. In fact, the measurement of light power provides the basis for <u>every</u> spectroscopic instrument. In this experiment you are going to build your own photometer with an LED light source and a CdS photoconductive detector.

## **Construction of the Photometer**

Some of the photometer is already constructed for you but you do have to complete the construction. To complete the construction, solder in place: the two resistors, a switch, and the battery leads. After your soldering is complete attach the four legs to support the photometer. Plug an LED and CdS detector in the appropriate sockets and attach a battery.



## **Measurement of Absorbance**

Follow these steps to measure the absorbance of the solution:

- Insert the appropriate colored LED.
- Place water in the sample cell and measure the voltage across the CdS
- Place the solution in the sample cell and measure the new voltage across the CdS.
- Calculate the absorbance using the Excel Spreadsheet.

Measure the absorbance of the three food coloring solutions with the red, green, and blue LEDs.