

17. Heteronuclear Correlation on AM Consoles


by cg fry: created 2/22/94 – revised 7/31/95

I. Discussion

Heteronuclear correlation maps interactions from an X nucleus (e.g., ^{13}C) to ^1H . The nominal experiment, referred to here as the HETCOR experiment, is implemented on the Bruker spectrometers, and uses X-nucleus detection with ^1H decoupling. The inverse experiment (i.e., HMQC), currently only possible on the Unity 500, uses ^1H detection with X-nucleus decoupling, and has much improved sensitivity over the standard HETCOR experiment. HETCOR will work ok down to 20 mg material in a 5 mm tube; more dilute samples should be run on the Unity with the inverse experiment.

The HETCOR experiment is relatively straightforward, requiring as a minimum running standard ^1H and ^{13}C spectra with optimized sweep widths as setup. It is important to have correct pulse width calibrations for both ^1H and ^{13}C transmitters, however, for correct optimized use of the the experiment. An understanding of filter bandwidths and folding is helpful to properly set the sweep widths.

II. Critical Parameters



$$\mathbf{D3} \quad - \quad = \frac{1}{2J_{CH}} \sim 3.3 \text{ ms}$$

change if have unusual J_{CH} couplings

$$\mathbf{D4} \quad - \quad = \frac{1}{4J_{CH}}$$

P1 – ^1H 90° pulse width; see calibration sheet next to spectrometer, or recalibrate before experiment

P3 – ^{13}C 90° pulse width; see calibration sheet

S1 – 0H for full ^1H power

S2 – same decoupling power used for standard ^{13}C spectrum (typically 20H)

$$\mathbf{IN} \quad - \quad \text{increment for 2D; make sure} \quad \mathbf{IN} = \frac{1}{4 \cdot \mathbf{SWI}} = \frac{1}{2 \cdot \left(^1\text{H sweep width} \right)} = \mathbf{DW} \text{ (} ^1\text{H dwell)}$$

III. Acquiring 2D HETCOR Spectra

1. **Use job 1:** Obtain a routine 1H spectrum. Save this spectrum using **WR**, as *nameH.001* where *name* is a seven letter or less name used for all HETCOR data for this sample. Do not use alphanumeric in the suffix for this file, to stay compatible with the DISNMR plotting routine CP2P. Note **SR** after setting the reference (using, e.g., EP-G command).

Note SR(1H) and save 1H spectrum as *nameH.001*

2. Use **EP** and CNTL-O command to set **SW** and **O1** for region of interest. Take a full (e.g., SI=16 or 32K) 1H spectrum with new settings, being careful that foldovers do not interfere with important 1H resonances.

Note O1(1H), and DW(1H)

3. **Use job 2:** Obtain a routine ^{13}C spectrum using CPD (i.e. Waltz-16) decoupling with **O2** set to the **O1** value used in the 1H spectrum. Set the reference, and write down the value for **SR**.

Note SR(^{13}C)

4. Now obtain DEPT-45 to remove quats. Use **EP** and CNTL-O command to set **SW** and **O1** to include only protonated ^{13}C plus ~10%. Take full (e.g., SI=32 or 64K) ^{13}C spectrum with new settings. Make sure you have sufficient S/N to observe all carbons of interest; S/N~3 is minimal, S/N~8 should be ok (i.e. decrease/increase NS until S/N~8). Save a good S/N ^{13}C Dept-45 spectrum as *nameC.001*

Note DW(^{13}C) and O1(^{13}C) and save ^{13}C Dept-45 spectrum as *nameC.001*

5. **In job 3:** enter **PJ XHCORR.1DJ**

RJ2D XHCORR.2DJ

6. **Return to job 2**, and transfer carbon parameters to job 3: **TR <ret> 3 <ret>**

8. Set **SI=2K** (larger or smaller depending on ^{13}C SW; 4k for SW > 150 ppm, 1k < 70 ppm) and retake the ^{13}C Dept-45 spectrum. Make sure you have sufficient S/N to observe all carbons of interest; S/N~3 is minimal, S/N~8 should be OK (i.e. decrease/increase NS until S/N~8).

7. Enter the following parameters

NE=128 (increase if want better 1H resolution, and have sufficient spectrometer time)

ND0=2

IN=[DW(1H); see step 2]; enter IN as number followed by **U**; e.g.,

IN <ret> 120U for DW(1H)=120

8. Enter

ST2D

and make sure that:

TD1=NE=128 (or if enough time, multiple of 128)

SI1=TD1×2

TD2=SI2=TD=SI=2K

SW2=SW

Check that $SR=SR2=SR(^{13}C)$ and $SR1=SR(^1H)$

4. Enter: **XFB**

and be patient. 2k×256w will take ~8 min on DataStation.

5. Enter: **EP2D** (DataStation/AM's) or **AP2D** (AC's)

Set contour levels appropriately, and **ESC-X**

6. Enter: **CP2P**

```
FILENAME IN F1: nameH.001  
FILENAME IN F2: nameC.001  
NO. OF PENS: 1 (or 2 if have red pen)  
NO. OF CONTOUR LEVELS: 4 (1-7 levels)  
OUTLINE BOX: Y  
GRID: N or Y
```

7. Enter: **NP** to feed in new paper to plotter or printer