


### 13. $^1\text{H}$ Spin-Lattice Relaxation, $T_1$

by cg fry: created 12/11/93 – revised 4/1/96

#### I. Discussion

Spin-lattice relaxation is nominally an exponential process, and becomes important for any quantitative, exchange, or cross-relaxation type of experiment such as NOE. The relaxation from zero magnetization and from an inversion pulse are both shown in the figure below.



T1-RECV.CDR

Often a quick determination of the  $T_1$  is sufficient; for most liquids, use the rapid, inversion-recovery method in Section II below. This method involves finding  $\tau_n$  for the  $^1\text{H}$  of interest. For more accurate determinations, use the procedure described in Section IV.

#### II. Rapid Determination of $T_1$ by Inversion Recovery Null Method

(good for any sample—e.g. any liquid—but is semi-quantitative at best)

- Obtain a reasonable quality  $^1\text{H}$  spectrum. Locate protons of interest for visual observation.
- $180^\circ$  and  $90^\circ$  pulse lengths must be reasonably accurate. If your sample is highly polar (e.g., ionic water solution), you must calibrate beforehand. Otherwise, if tuned correctly, prior calibrations can be used.
- Decrease SI and TD (keep SI=TD) to 8k if possible.
- Set AI=1 (answer Y to REINITIALIZATION?), NS=1.
- Reacquire fid with ZG, EF, and phase manually in EP.
- **AS T1NULL.AU**

**P2=16**  $\mu\text{s}$  ( $180^\circ$  pulse length; see calibration sheet on wall)

**VD VDLIST.001**

**1 0.01** s (to start with)

**2 EN** (or )

**RD = 0**

**PW = 8**  $\mu\text{s}$  ( $90^\circ$  pulse width)

**DE** (set by spectrometer)

**NS = 1**

**DS = 0**

- **NE = 1**
- **AU** <ret> (or equivalently AU T1NULL.AU<ret>)
- **VD=0.01** should give large negative peaks. Lengthen **VD** until proton peaks of interest are nulled (slightly shorter **VD** should give negative peaks; slightly longer should give positive peaks). For the nulled **VD<sub>n</sub>**,  $T_1 \cong \text{VD}_n / \ln 2 = \text{VD}_n * 1.44$ .
- Find **VD<sub>n</sub>** for other protons of interest.

### III. Rapid, approximate measurement of T<sub>1</sub> by Repetition-Rate method

(only for samples where  $T_2 \ll T_1$ , e.g. solids [but not for nominal liquids; see previous section])

1. if present, **RJ H1REPRAT.005** or **.010** (for 5mm or 10mm probes, respectively), or:  
 set **NS = 1**, **DS = 0**, **RD = 0**, **SW = 20ppm** (10000 on AM500, 7200 on AM360), **SI = 16K**,  
**AQ = 0.1s**, and **PW = 90°** pulse length (see calibration sheet next to spectrometer; or recalibrate  
 90° - see section xxx)
2. set **O1** close to, or on peak of most interest (use **EP**, set cursor on peak and enter **O1** then **M**)
3. check **RG** (I like to reduce TD = 1k to see beginning of FID clearly, with RD set to compensate for shorter AQ), then wait at least 5T<sub>1</sub>;
4. set absolute intensity mode by typing **AI**, answer **Y** if queried REINITIALIZE?
5. **ZG**, then **EF**, go into **EP**, and **P**; after good phasing type **M** and <RET> out of **EP**
6. adjust vertical scaling so 2<(height of peak of interest)<4 divisions; note height and scaling
7. set **DS = 8**, and **RD** so that RD+AQ ~ T<sub>1</sub>/2
8. **ZG**, then **EFP**
9. adjust **RD** until height of peak of interest is half height noted in step 6 (i.e. same height with scaling doubled)
10.  $T_1 \cong (\text{AQ} + \text{RD}) / \ln 2$

### IV. Quantitative measurement of T<sub>1</sub> by Inversion Recovery Method

#### A. Comments

- A reasonable estimate of T<sub>1</sub> must be known to correctly setup: use the Inversion Recovery Null Method (Section II).
- 180° and 90° pulse lengths must be accurate. If your sample is highly polar (e.g., ionic water solution), you must calibrate beforehand. Otherwise, if tuned correctly, prior calibrations can be used.
- For precise measurements, use D1 = 10 T<sub>1</sub>; D1 = 5 T<sub>1</sub> will provide reasonable values if experimental times get long. The number of VD values can also be decreased to lessen experimental times  $\cong \text{NS} * \text{NE} * (2 * \text{D1}) * (\text{second \# in VC list})$

**B. Acquisition Set-up**

1. setup normally as for high-resolution  $^1H$  spectra; make sure RG is set correctly (not too big!), and NS = 8
2. set-up two lists used by the automation routine: a VC (variable counter) list, and a VD (variable delay) list
  - i) enter the variable delay editor by typing: **VD<RET>**
    - use the default name VDLIST.001, or enter a unique name you can recall later
    - enter the delays; Table I giving a reasonable grouping
    - enter a number followed by S, M or U for seconds, milliseconds, or microseconds
    - end the list by entering **ESC**

**Table I. A reasonable grouping of VD delays for Inversion-Recovery  $T_1$  sequence, where  $T_1^r$  is the  $T_1$  estimate (usually from repetition-rate experiment).**

$0.01T_1^r$	$0.5T_1^r$
$10T_1^r$	$10T_1^r$
$2T_1^r$	$5T_1^r$
$0.3T_1^r$	$0.1T_1^r$
$10T_1^r$	$0.01T_1^r$
$0.7T_1^r$	$1.4T_1^r$
$1.0T_1^r$	$10T_1^r$

- ii) enter the variable counter list editor by typing: **VC<RET>**
    - use the default name VCLIST.001, or enter a unique name you can recall later
    - the first line must contain the number of delays in the VD list = NE
    - the second line contains the number of complete cycles of NS and VD to do; for  $^1H T_1 = 1$
    - end the list by entering **ESC**
3. start automation setup for the inversion-recovery sequence: **AS INVREC.AU**
  - i) change values for the VC list, or enter **<CNTL>-L** to list
  - ii) **D1** =  $10T_1^r$
  - iii) **P2** =  $2*PW$ , where PW =  $90^\circ$  pulse length (see calibration sheet next to spectrometer, or recalibrate using procedure in section xxx)
  - iv) **RD** = 0, **DE** set by computer, **NS** = 8 (or multiple of), **DS** = 0
  - v) **NE** = number of delays in VD list
4. experiment time  $\cong$  NS\*NE\*(2\*D1)\*(second # in VC list)
5. start automation routine: **AU<RET><RET>**; enter unique filename with *no* suffix

### C. $T_1$ Analysis

1. Transfer data to data station using Lightnet (can use wildcard), or enter other session
2. check **NE** = size VD list
3. set **AI** = 1, Y if asks REINITIALIZE?
4. **RE** filename.002 (assuming .002 has longest VD), do **EF**, then **EP**, **P**, and when good phase, **M**
5. do all FT's by running automation routine: **AU FT.AU**  
enter filename of  $T_1$  set, then new name of FT'd spectra (e.g., T1SET and T1SETF)
6. if want plot, go through DPO (Y to both axes, N to parameter list, Y is want TI), CX=CY=10 to 20, MAXY = 10 to 20
7. go into EP, set cursor on top of peak of interest and enter **T**, then <RET> to accept T1PTS as filename  
go to next peak of interest and enter **T** (no <RET>); repeat for all peaks of interest
8. enter  $T_1$  routine: **T1<RET>**
  - i) enter **PD**, N to integrals list, N to VD list, Y to EP list
  - ii) when see display, enter **CT1** to calculate  $T_1$
  - iii) enter **PLTD** to plot (can change TI and DPO prior to plot)
  - iv) enter **NXTP** to go to next peak
  - v) enter **QUIT** to exit  $T_1$  routine