12. HOMODEC on AM/AC Spectrometers

by cg fry: created 12/20/94 - updated 3/22/96

I. Setting Up

- Obtain a normal ${}^{1}H$ NMR spectrum. Save file using **WR** command.
- Choose **all irradiation frequencies** that you want to irradiate, e.g., you have two such frequencies (*vide infra*).
- Frequencies can be selected by moving the cursor in EP mode to a desired position in the spectrum, followed by **O2** and **L**. On the first frequency, the spectrometer will offer FQLIST.001 as a default list name: you may hit **<ret>** to accept this name, or enter a new name of your own. All other frequencies, require just **O2** and **L**. Repeat this until you have chosen all the irradiation frequencies.
- <**ret**> to exit <u>EP</u> to Command Interpreter.
- AS HOMODEC.AU <ret>

D1 = 1 (sec)

O2 FQLIST.001

- 1 = (first irradition frequency) <**ret**>
- 2 = (second irradiation frequency) <**ret**>
- ...and so on, if you have more irradiation frequencies.

3 = END

(Note that this frequency list is automatically set if you have chosen the frequencies with O2 and L commands in <u>EP</u> subroutine.)

S3 = 20L (decoupler power; use 20L to 10L; the lower the number is, the higher the power)

 $\mathbf{RD} = \mathbf{2.0}$ (0 can be used, but better if 1 or 2 sec is used)

PW = **2.5** (whatever is used in job file)

DE = (dead time: set automatically by the spectrometer)

NS = 8 (number of scans)

DS = 2 (number of dummy scans; **0** works fine in most cases)

NE = 2 (number of experiments = number of irradiation frequencies)

II. Taking ¹H NMR Decoupling Spectrum

- **EXPT** does not work!! Calculate the experimental time as (RD+AQ)*(NS+DS)*NE*(~1.2). The factor 1.2 estimates the amount of disk access time spent. Use a larger number like 1.5 if RD+AQ<2 s, or 1.0 if RD+AQ>10 s. If the experiment is too long, adjust <u>NS</u>, <u>DS</u>, <u>AQ</u> or <u>RD</u> accordingly.
- ZE; AU <ret>

INPUT FREQ. LIST FILENAME #1: **FQLIST <ret>** INPUT FID FILENAME #2: **KL** (any filename with no extension) **<ret>**

• FIDs are now collected and stored automatically into FILENAME #2's with extensions 1 through the total number of frequencies (i.e., KL.1 and KL.2).

III. Processing FIDs

Work up the first irradiation spectrum as follows:

LB = 0; RE KL.1 (to read FILENAME #2.1); EF (do EM plus FT); EP; <ctrl> R; BP; <ctrl> B; phase the spectrum with C- and D-knobs; M; <ret>.

- Save the resulting irradiation spectrum, e.g., WR SKL.1.
- Process other irradiation spectra using <u>EFP</u> (<u>EM</u> + <u>FT</u> + <u>PK</u>) as follows:

RE KL.2 (to read FILENAME #2.2); EFP; WR SKL.2 (the second irradiation spectrum; any filename can be used).

Or set NE = # files, and use AU FT.AU to automatically process all the files.

- To plot these spectra, set DPO, CX, and CY. Then, PX. In WP console, USR FASTU to send the • stored spectrum to the PC data stations.
- PO (decoupler power off) before leaving the spectrometer. ٠

IV. Listing of Automation Routine

;

```
HOMODEC.AU
;
   HOMONUCLEAR DECOUPLING USING ONE FREQ. LIST AND ONE POWER SETTING.
   1 FL #1 / INPUT FREQ. LIST ; READ IN FREQ. LIST AND INITIALIZE POINTER
   2 ZE
   3 D1 HD O2 S3 ; TURN ON HOMODEC. WITH POWER S3 AND SET
                    ; O2 FREQ. FROM CURRENT FL LIST, INCREMENT
                    ; FREQ. POINTER.
   4 GO = 4
                   ; ACQUIRE DATA
   5 WR #2
                   ; INPUT FID
                   ; STORE FID
                ; INCREMENT FILE EXTENSION
; LOOP FOR NEXT EXPERIMENT
   б
      IF #2
   7
      IN = 2
                    ; LOOP FOR NEXT EXPERIMENT
   8
     D1 DO
      EXIT
                   ; EXIT WITH DEC. OFF
  PROGRAM REQUESTS A FILENAME FOR FREQ. LIST AND FID STORAGE
;
  D1 = 1 SEC TO SET FREQ.
;
  S3 = DESIRED DEC. POWER
;
  NE DEFINES NUMBER OF ITEMS IN FL LIST = NO. OF EXPERIMENTS
;
; USE RD AS DESIRED AND AT LEAST 2 DUMMY SCANS.
```