

5. AM/AC Spectrometer Software

by *cg fry*: created 12/10/94 – updated 20.Aug.2002

I. ADAKOS (Aspect Disk And Keyboard Operating System)

- ADAKOS is the base operating system of the ASPECT computer, the acquisition and analysis computer on the AM spectrometers. Some functions of ADAKOS include listing of directories, batch copying or deletion of files, renaming of files, initialization of new disks and restoration of deleted files. This program is started automatically whenever the computer is rebooted (<STOP><CLEAR><DISK>).

All ADAKOS commands have two or three letters, e.g., **DIR**, **COP**, and **DEL**.

MO (monitor) to exit to the ADAKOS subprogram from DISNMR.

* **DSP ON** ("*" is the AM-console prompt) to turn on the display.

DIR to list all files in the directory: **H** to halt and **C** to continue. <ctrl> **Q** to stop.

SDIR *.* for sorted directory (sorted by alphabet; can be quite slow).

- A filename has up to eight characters, followed by up to four characters. The wild card * represents only groups of four characters for *only* blocks 1-4, 5-8 or for the suffix. **DIR** *.* will only show filenames having four characters or less. **DIR** aa*.* is invalid, since * cannot represent the 3rd through 6th characters. ? can be used for any single character (only available in ADAKOS, unfortunately not in DISR).
- D1 is the system disk. Use F1 for the floppy disk. The disk designation must follow the filename unless it is the default drive (D1). e.g.,

DIR *.*.*=F1

to list all files on the floppy.

COP <input file>/<output file>=F1

to copy a file from disk onto a floppy.

COP <input file>/=D2

to copy with the same filename.

COP *.*.*=F1/=D1

to copy every file on a floppy to the hard disk.

DEL <filename>=D1

to delete a file from disk.

DISR91

to run DISNMR (**RUN DISNMR** on WP-console).

II. DISNMR

A. Control Keys and System Commands

<ctrl> **D** toggles the grid ON and OFF on the display;

<ctrl> **L** toggles between the lock display only, the data only, and both.

<ctrl> **H** halts following next acquisition; active only in current job

<ctrl> **E** halts immediately; does not save data; active only in current job

<ctrl> **K** halts all activity including acquisitions, plotting and printing; active in all jobs; use only as a *last resort* (try <ctrl>H first, the <ctrl>E)

<ctrl> **X** toggles between 1st and 2nd sessions

<ctrl> **Y** toggles display from interleaved to separate quadrature channels

- **OUT** <ret> for setting up output devices.
COMMAND INTERPRETER? **B** (both).
OUTPUT DEVICE? **P** (printer only).

B. Files in DISNMR

- In DISNMR, any file except the current one can be deleted in the same manner as above. **DIR** is also available, but the wild character ? can't be used. The * wildcard can only represent characters 1-4 or 5-8, or the suffix:

DI *.* will list all files having 4 characters or less in prefix
DI *.*.* will list all files
DI cgf1*.* will list all files starting with **cgf1**
DI cgf*.* *not allowed* since * must start in position 1 or 5 of prefix
DI cgf?*.* *not allowed* in DISNMR; ok in Adakos

- There are many types of files:

File Type Process	1D data files	2D data files	Job files	Shim files	Automation
Common suffix	**.*001	**.*SER	**.*1DJ **.*2DJ	**.*SHIM	**.*AU
Disk I/O	WR, RE	(.*AU writes), RE	WJ, RJ, PJ	WSH, RSH	EDIT, AU
Directory	DI:F (FIDs) DI:S (spectra)	DI *.*.*SER DI *.*.*SMX	DI:J	DI:Z	DI:A

- TI** will allow entry of up to 80 freeform characters, stored with data file. *Users should include a brief sample description, and notebook number and page in the TI area.*
- RE <filename>** to read a file ("?" following any command opens a helpfile). One can also save or read files from an external 8" or 3-1/2" floppy disk (F1) as follows:
WR <filename>=F1 to save a file to floppy
RE <filename>=F1 to read a file.
- COP** and **DEL** commands work as in Adakos.

C. Common Acquisition Commands (WJ, RJ)

AQ – acquisition time	NS – # scans
AS – automation setup	O1 – transmitter offset frequency
CO – continue acquisition	O2 – decoupler offset frequency
CPD – composite pulse decoupling (Waltz-16 decoupling)	PO – decoupler power off
DO – decoupler gated off	RD – relaxation delay
DE – pre-scan delay	RG – receiver gain (actually, receiver attenuation in dB)
DP – decoupler power (actually, decoupler attenuation in dB)	RGA – automatic receiver gain setting
DS – # dummy scans	SF – spectrometer frequency
GO – begin acquisition	SI – spectrum size
GS – GO setup mode	SW – sweep width
HD – pulsed homonuclear single frequency decoupling	TD – time domain size
II – initialize interface	ZE – zero data and initialize for acquisition
NE – # experiments	

D. Common Pulse Sequence Commands and Parameters (WJ, RJ)

- (**:**D) – () operates on decoupler
An – transmitter quadrature phase
Bn – decoupler quadrature phase
Dn – delay length for nth delay
 D0 is 2D t₂ delay incremented by
 IN (*not* the loop parameter)
FL – input frequency list
GO=j – perform GO, jump to line j when
 finished; continues NS times
- IF** – increment file number
IN=j – loop index jumps to line j;
 continues NE times
IN – time incr. in 2D set by SW1
O2 – read decoupler frequency,
 increment list pointer
Pn – pulse length for nth pulse
Sn – decoupler power setting

E. Common Data Processing Commands (WJ, PJ)


- AI** – absolute intensity
 = 0 (default) autoscales
 = 1 absolute intensity
AT – additive transfer of two spectra
BC – apply linear baseline correction
FT – Fourier transform
PK – apply phase correction last
 stored in **EP** mode
SR – spectrum reference frequency
- EM** – exponential multiply
 $E_i = \exp[-\pi t_i \cdot LB]$
 chose **LB** ~ linewidth
GM – gaussian multiply
 $G_i = \exp[-\pi \cdot t_i \cdot LB + \frac{\pi \cdot t_i^2 \cdot LB}{2AQ \cdot GB}]$
 chose **LB** ~ -linewidth
 chose **GB** ~ fraction AQ
 where S/N first = 0

F. Common EP Commands

- ^B** – re-enable knobs A and B
^F – display current F1,F2 region
^O – set SW,O1 to expanded region
^R – display whole spectrum
LINE FEED – baseline points selection mode
A – add phase corrections
B – phase on biggest peak
C – phase on cursor selected peak
D – dual display
D – decrease separation
I – increase separation
S – subtract spectra
E – toggle EP info. displays
F – enter frequency limits F1,F2
G – set spectrum reference (SR)
H – hard copy (print) cursor info.
I – enter integration mode
A – calibrate (normalize) integrals
M – toggle knobs C,D from cursor mode
 to slope and bias mode
S – section plot of displayed integral
X – full plot of displayed integral
Z – zero integral; set starting-end points
 of integrals
- J** – enter line fitting mode
K – enter interactive baseline correction
M – set minimum intensity for peak picking
P – phase on biggest peak displayed
O1 – printout absolute freq. of cursor
M – will store frequency as O1
O2 – printout absolute frequency of cursor
L – stores freq to specified freq list file
M – will store frequency as O2
R – use twice to set freq limits for expanded
 display (^B to reset A,B knobs)
S – section plot
T – enter cursor to T1PNTS file
U – update region for plotting
X – plot expanded region (PX)
Y – change CY,CX within EP
< – left shift one point (LS)
> – right shift one point (RS)
: – toggle
5 – peak peak and display
6 – move display down 1/8 screen
7 – move display up 1/8 screen
8 – display rms S/N
9 – toggle to dots/connected display

III. Plotting with DISNMR

A. Overview



Plotting on Aspect computers starts similarly to any computer system, with page and plotting regions. The plotting region depends on the particular plotting device; here we'll assume either an HP LaserJet printer, or an HP 7475 or 7550A plotter. Landscape plotting is default on both plotters. The spectral region plotted is set with the **EP U** (update) command. Plot parameters can be read in with the **PJ** command (these jobfiles may not be present; we use NUTS to plot nearly all spectra):

PJ A3PLOT.1DJ	11×17" landscape plotting
PJ A4PLOT.1DJ	8.5×11" landscape plotting
PJ A3STKPLT.1DJ	11×17" stack plotting with rotated (portrait) spectra
PJ A4STKPLT.1DJ	8.5×11" stack plotting with rotated (portrait) spectra

The parameters and commands that control plotting (see figure above):

Parameter/ Command	Description	Landscape 11×17" A3	Landscape 8.5×11" A4	Portrait 11×17" A3	Portrait 8.5×11" A4
X0	X offset for origin	0	0	0	0
Y0	Y offset for origin	0-4	0-2	0-4	0-4
CX	X-axis length of spectrum	35	24	24	16
CY	Height of reference (usually tallest) peak in spectrum	20	16	30	20
MAXY	Maximum Y deviation plotted	22	18	32	22
DPO rotate?	define plotter options	N	N	Y	Y

B. Standard Plots

1. **RE** spectrum in
2. **PJ** plot jobfile (see above)
3. check **DPO**, **CY**, and **MAXY**
4. **EP**, get spectrum on screen as wanting to plot and enter **U** to update plot region
5. use **X** and **S** inside **EP**, or exit **EP** and use **PX**
6. *remember to cap pens on pen plotters when finished*

C. Stack Plots

1. Start by **RE**'ing bottom spectrum
2. **PJ** plot jobfile (see above)
3. check **DPO** (plot x-axis now, and rotate if want portrait stack), **CY**, and **MAXY** (title height)
4. check **X0,Y0** for initial plot (use **OP** if using **STACK.AU**; manual plotting is recommended)
5. use **EP** and **U** to update plot region; exit **EP**
6. plot first, bottom spectrum with axis and title using **PX**
7. **RE** next spectrum
8. change **DPO** for no axis or title plotting
9. increment **X0,Y0**
10. plot using **PX**
11. goto 9 and continue for all spectra
12. *remember to cap pens on pen plotters when finished*

IV. Data Transfers**A. UW Network**

- Windows NT - Backbone file server network for PC's
- completely transparent to users, data is uplinked to a Snap server
- NMR-Link - Bruker file transfer network; bidirectional
- writes files in Bruker image format (widely read by NMR software)
 - uses PC AT's to automatically send FASTRAN data to WinNT server
- FTP - File transfer protocol; generic file transfer between many platforms
- NFS - Sun's Network File system (must use FTP in our facility to transfer)

(All FID files must be stored on AC spectrometer hard disk by using **WR <filename>** before the transfer.)

B. NMR-Link

- NMR-Link exports and imports Bruker binary (and with a non-default setup, ASCII) data. The two commands used in this facility are:

topc filename.ext to send files to the PC network

frompc filename.ext to receive files from the PC network

The above two commands are always issues while sitting at spectrometers or the datastation, with a PC-AT providing dedicated file servicing.

- ***NMR-Link at an AC Spectrometer***

When at a spectrometer, the commands are run from the 2nd session, toggle to/from with the <cntl>-X command.

- ***C. FTP to/from NMRSnap***

- Login to NMRSnap.chem.wisc.edu from any FTP-capable computer (see postings in lab for password).
- The only directory that has write privileges is \\nmrsnap\DataArchives\temp

V. Formatting Floppy Disks

- Insert a blank 3-1/2" floppy disk in the disk drive and turn the drive ON. **MO** to exit to the monitor from DISNMR. **DSP ON** will toggle display to the monitor. Then, do as follows:

***ADAKET** <ret>

SELECT FUNCTION? **F** (to format a floppy disk)

ERASE DATA? **Y** <ret> ...formatting begins and takes ~45 seconds...

M (to exit to the monitor)

DIR=F1 (to verify formatting)

LOAD ADAKOS

- Reboot the computer by **STOP-CLEAR-DISK**. If DISNMR is not loaded automatically, type **DISR91**, followed by <ret>.