


Calibration of 90° Pulsewidths

Measurement of pw90 can often be skipped on Bruker spectrometers equipped with automatic tune-and-match (ATM) probes. On manually tuned probes — e.g., the Inova600 and Avance360 — the procedure is highly recommended.

Calibration of ¹H 90° pulsewidths in TopSpin:¹

Measuring the length of the standard proton pulse to achieve a 90° nutation² is a primary calibration done for NMR experiments. This pulsewidth on the observe channel, **p1**, (or pw90³) is probe-dependent and is always matched to the power used, **plw1**. A series of spectra at increasing pulse length — which increasingly nutate the magnetization — is collected via Bruker's **popt** automation routine. The best measurement is made at the 360° null crossing, with $pw90 = pw360 / 4$.


- rpar PROTON all ↵ ns 1 ↵ ds 0 ↵ pulprog zg ↵**
- Do normal **lock**, **atma**, **getprosol**, etc.
- Changing PULPROG to **zg** uses **p1** directly, rather than **p1*0.33** as in **zg30**.
- efp ↵ apk ↵ .all ↵ abs n ↵**
- Set **PH_mod = pk** so spectra are phased during **popt**. Set at cmdnd line or from the PROCARS tab.
- Expand the spectrum about peaks in the midrange ppm's, (don't expand on solvent or impurity multiplets).
- Right-click on the spectrum, and choose:



Save Display Region To...

and use the top selection. *You must run this command prior to doing a **popt** experiment.*

- popt ↵** will open the table to the right. Enter PARAMETER **p1** and STARTVAL = **4**.
- Delete other values existing in the table.
- Set **INC = 4**, **NEXP = 16**. **ENDVAL** will automatically update.
- Click **START OPTIMIZE**. Give the appropriate responses [don't know what **APPEND** means; doesn't seem to matter].
- FIRST TIME ISSUE:** you may need to make the following changes that, once done, will be remembered by TopSpin (then again, you might prefer the original settings...)





To keep **popt** from constant changing away from the stacked spectrum presentation:

¹ A (relatively) new TopSpin command, **pulsecal**, works fine in samples that have a very strong singlet peak, but fails otherwise. It is therefore a very good alternative to the above for biological samples in 90:10% H₂O:D₂O. It is not recommended for samples in other solvents.

² *Nutation* in NMR refers to the rotation of magnetization about a rotating-frame axis, e.g. x' . For magnetization initially along the z axis (parallel to \mathbf{B}_0), then in a left-hand system, the magnetization nutates $z \rightarrow y' \rightarrow -z \rightarrow -y' \rightarrow z \rightarrow \dots$

³ I (cgf) often use the Varian symbol **pw90** for the 90° pulsewidth, as it is more directly descriptive; in vnmr/vnmrj the observe channel power is **tpwr**. **pw360** is not used in vnmr, but the extension for a 360° pulsewidth is obvious.

MANAGE → PREFERENCES → ACQUISITION uncheck
“Auto open acquisition window after “zg”

- (n) You should now be seeing a horizontal stack plot growing. Click  to get it scaled correctly, and  to get the stack mid-screen vertically.

Note the sinusoidal behavior:
90° at ~ 15 μs, 180° at 30 μs,
360° at 60 μs (on the 400;
differs elsewhere). The 360°
crossover is the best to measure.
90° is too flat; 180° too
dependent on repetition rate.

- (o) TopSpin automatically moves the spectra into the 999 processing folder. To get back to the acquisition folder, type:

re 1 1 ↵

where you'll see the last spectrum acquired, with the parameters used for that spectrum: in this example, the last spectrum has **p1=96 μs**.

- (p) Run **popt** again, narrowing in on the 360° region. A good choice for the above example is: STARTVAL = 53, INC = 0.4, NEXP = 16.
- (q) **p1 = pw360 /4**. To invoke the new pulselength throughout all calculated pulses (including shaped pulses), enter the following in each new experiment:

getprosol 1H <pw90> <power>

where **<power>** is stated in attenuation, not watts (e.g., -13.94 dBW, not 24.8 watts).

If **pw90 = 8.15 @ -13.94**, then after reading in a new experiment such as noesy1d, you would enter:

getprosol 1H 8.15 -13.94 ↵

This command will recalculate the shaped pulse based on the new pulselength.

