

# TopSpin Data Archiving to Castor

last update: 30 Apr 2014 (cgf)

- Automatic archiving occurs upon completion of **zg**. Usually the archived dataset (that on castor) will have raw data (fid or ser) only, and no spectrum exists on castor. IconNMR does process prior to archiving.
- Intermediate processing using: **tr** ↓ **efp** ↓ will create a processed spectrum of the data, but with only the number of scans from the **tr** command. This often poor spectrum (e.g., showing just NS=32 of the 640 actually completed) will be copied during auto-archiving. Process the fid or ser file, not the spectrum file (1r or 2rr)!
- Some acquisition types will not auto-archive at all: **multizg** and **multi\_zgvd** are two; others may exist.
- To use processing done on the *spectrometer* after **zg** completion, copy the data folder from *spect* to castor.
- On the spectrometer, open the file browser [in the Linux menu: **Applications** → **System Tools** → **File Browser**] and make 2 bookmarks:
  - spect* pointing to /home/username/persephone
  - castor* pointing to /nmrdata/username/public\_html/persephone

toggles to normal path display

click to show data on spectrometer

click to show data on castor

Always drag from *spect* to castor!

