

# Direct Observation of X-Nuclei within TopSpin

updated: 5 Sept 2014 (cgf)

## I. Sample Concentrations:


Direct observation of X-nuclei will succeed or fail primarily as a function of sample concentration. Other factors can be very important, especially for low- $\gamma$  nuclei. For such nuclei, coupling to an abundant spin —  $^1\text{H}$ ,  $^{31}\text{P}$ ,  $^{19}\text{F}$  being the most important — will often determine whether any observation can be made or not. Labeling will often be the only other choice for nuclei such as  $^{15}\text{N}$ . [Note the “Direct” component of this discussion. Many nuclei can often be straightforwardly observed via “indirect” detection, usually via  $^1\text{H}$ -detected 2D experiments.]

Minimum concentrations have not been determined for most nuclei. The following are guides: overnight runs for callisto, Prodigy, and BBFO+ (on Av500); 2 hr runs on Av400

$^{13}\text{C}$  nat. abundance: *callisto*  $\geq 300 \mu\text{M}$ , Prodigy  $\geq 1 \text{ mM}$ ; BBFO+  $\geq 3\text{mM}$  Av400  $\geq 15\text{mM}$

$^{19}\text{F}$ : Prodigy  $\geq 50\mu\text{M}$  (w no  $^1\text{H}$  coupling), BBFO+  $\geq 120 \mu\text{M}$ , Av400  $\geq 1\text{mM}$

## II. Summary of Commands:


<b>sx #</b>	[callisto: isx #]
<b>new</b> ↓ or <b>rpar</b> ↓	;choose parameter set [follow <b>rpar</b> ↓ with <b>getprosol</b> ↓]
<b>atma</b> ↓ or <b>atmm</b> ↓	;be patient, as can take up to 2 mins
<b>topshim gui</b> ↓	;add tunexy and/or convcomp
<b>ased</b> ↓	;optimize parameters, esp <b>d1</b> and <b>ns</b>
<b>rga</b> ↓	;needed only for $^1\text{H}$ , $^{19}\text{F}$ and highly-concentrated $^{31}\text{P}$
<b>zg</b> ↓	; <b>tr #</b> and <b>halt #</b> are useful
 or <b>o1p</b> ↓	;adjusts <b>sw</b> and <b>o1p</b> , and will change <b>aq</b> (always check)

→ redo **atma**↓ if **o1p** changes  $> 50 \text{ kHz}$  (check via **o1**)

→ start with **sw** large, and reduce for better resolution ( $\approx 1/\text{aq}$ )

## III. General Direct-X Experiment Guide:

1. Acquisition of a  $^1\text{H}$  spectrum is always recommended, prior to any X-direct experiment (including  $^{13}\text{C}$  experiments). The  $^1\text{H}$  spectrum will confirm shim quality, and can be used to check the referencing of the X experiment through the use of the **xref** au routine.
2. **new**↓ or **rpar**↓ and choose the proper parameter set. A few experiments can be found from the Show Recommended list.



- Most parameter lists are found with the “Show Recommended” unchecked, with a portion of list shown on the next page. This list can be reduced using the “enter any string, \*, ?” search box. For direct-X experiments, these are always listed as nucleus\_symbol-isotope\_number. Thus, enter **N15** to display direct-detection  $^{15}\text{N}$  parameter sets.


**Bruker-supplied parameter sets.** Direct-X parameter sets are circled: **green** circles indicate parameters sets OK for BBFO and Prodigy probes; **blue** circles show nuclei accessible only on the BBFO+ probe ( $^{35}\text{Cl}$  is too low in  $\gamma$  for the Prodigy probe, falling below  $^{15}\text{N}$ ); **red** circles show parameter sets that should never be used, and should be gotten from the UWChem folder.

File Name	Parameter 1	Parameter 2	Parameter 3	Parameter 4	Parameter 5
AL27ND	APSY_HNCA_32	APSY_HNCACB_32	APSY_HNCO_32	APSY_HNCOCA_42	
APSY_HNCOACB_32	APSY_HNCOCANH_62	ASSURE_13C	ASSURE_19F	ASSURE_1H	
ASSURE_31P	B_HNCA CBGP3D	B_HNCA CBIGP3D	B_HNCA COGP3D	B_HNCA COGP4D	
B_HNCA GP3D	B_HNCA IGP3D	B_HNCOACBGP3D	B_HNCOACBGP4D	B_HNCOAGP3D	
B_HNCOAGP4D	B_HNCOGP3D	B_HNCOIGP3D	B_HSQCETF3GPSI	B_TRHNCA CBGP3D	
B_TRHNCA CBGP3D	B_TRHNCA COGP3D	B_TRHNCA GP3D	B_TRHNCA IGP3D	B_TRHNCOACBGP3D	
B_TRHNCOAGP3D	B_TRHNCOGP3D	B_TRHNCOIGP3D	B_TROSYETF3GPSI	B_TROSYF3GPPH	
R11ZG	BESTPROFILE	C_CACO	C_CACO_IA	C_CACO_S3	
C_CAN_IASQ	C_CAN_MQ	C_CAN_MQ,2	C_CANCO_IA3D	C_CANCO_IA3D,2	
C_CANCO_IA3D	C_CBCACO_IA3D	C_CBCACO_S33D	C_CBCACON_IA3D	C_CBCANCO_IA3D	
C_CCO_IA3D	C_CCO_S33D	C_CCCON_IA3D	C_CCFLOPSY16	C_CCFLOPSY16_CT	
C_CCFLOPSY16_CTIA	C_CCFLOPSY16_IA	C_CCNYESY	C_CCNYESY_CT	C_CCNYESY2	
C_COCA	C_COCA_IA	C_COCA_MQ	C_COCA_MQ,2	C_CON_IASQ	
C_CON_MQ	C_CON_MQIA	C_CON_SQ	C_COSY	C_COSY_CT	
C_COSY2_CT	C_HACACO_3D	C_HCACO_IA3D	C_HCACO_S33D	C_HCAN_IA3D	
C_HCANCO_IA3D	C_HCANCO_IA3D	C_HCBCA_IA3D	C_HCBCACO_IA3D	C_HCBCACO_S33D	
C_HCBCAN_IA3D	C_HCCFLOPSY16_3D	C_HNCA_IA3D	C_HNCAO_IA3D	C_HNCAO_S33D	
C_HNCO_IA3D	C_HNCOCA_IA3D	C_HNCOCA2_IA3D	C13APT	C13CPD	
C13CPD32	C13CPDSN	C13DE45SN	C13DEPT135	C13DEPT135p	
C13DEPT45	C13DEPT90	C13GD	C13HUMP	C13IG	
C13MULT	C13MULT135	C13MULT90	C13MULTCOMP	C13OFF	
C13PPTI	C13RESOL	C13SENS	CBCACONHGP3D	CBCACONHGPWG3D	
CBCACONHGPWG3D,2	CBCACONHGPWG4D	CBCANHGP3D	CBCANHGPWG3D	CCACONHGP2H3D	
CCACONHGP3D	CCACONHGP3D,2	CCANHGP2H3D	CCANHGP3D	CCANHGP3D,2	
CCCONHGP2H3D	CCCONHGP3D	CD111ZG	CD113ZG	CL35ZG	
CL37ZG	CMCQ_PROTON	CMCQ_WET	CMCse_13C	CMCse_15NHMBcf2	
CMCse_15NHSQCf2	CMCse_1H	CMCse_COSY	CMCse_HMBC	CMCse_HSQC	
COSY45SW	COSY90SW	COSYCWGPPSQF	COSYCWPHPS	COSYDCPHWT	
COSYDFGPPH19	COSYDQFPHSW	COSYGPDPHPSW	COSYGFPIXSW	COSYGPMSFW	
COSYGPSW	COSYPHPR	DIPS2ESFBGPPH	DIPS2ESGPPH	DIPS2ETGPSI19	
DIPS2GPPH19	DIPS2ESGPPH	DIPSITRET3GPPH	F19	F19CPD	
FHSQCXCF3GPPH	FHSQCF3GPPH	GA71ZG	gradshim1d1h	gradshim1d1h_f	
gradshim1d2h	gradshim1d2h_f	gradshimdata	gradshimrcb3d	H2OSUPMLEV	
H2OSUPNOESY	H2OSUPROESY	HACAHB COSYGP3D	HACONHGPWG3D	HANHGPWG3D	
HBCBCGCDCEHGP	HBCBCGCDHDGP	HBHACBANHGPWG4D	HBHACONHGP3D	HBHACONHGPWG3D	
HBHACONHGPWG4D	HBHACONHGP3D	HBHACONHGPWG3D	HACOCB3D	HACOCB3D	

- Some parameter sets, such as F19, must be modified to match the specific hardware setup on the spectrometer. Or the parameter set may simply have more optimal parameters than Bruker provides (e.g., for  $^{11}\text{B}$  and  $^{10}\text{B}$ ). These parameter sets are kept in a folder specific to the Chemistry NMR Facility at `/home/topspin3.1/uwchem/par`.

**Some experiments must utilize parameter sets provided in the UWChem folder!** Any experiment involving  $^{19}\text{F}$  on Persephone falls in this class of experiments. It is therefore important to check the list of experiments available in that folder, and always use that set if it is available. Note the set of parameters shown in the screen capture below: it is up-to-date as of 23 June 2013.

**UWChem-supplied parameter sets from Persephone.** Note the three sets of parameters are all probe specific: BBFO, FCH, and Prodigy. Always use only those sets matching the probe currently installed (the on-line calendar has the schedule). Three  $^{19}\text{F}$ -only parameter sets are circled in blue.



3. Tune the probe as normal, using **atma**. If this fails for some reason, use **atmm**. *Be patient!* The routine accesses capacitor rods that make mechanical changes to the probe. In the Prodigy, the rods pass through vacuum and cryogenic areas; it can only go so fast. Impatience can lead to a stuck rod, and potentially to probe damage.


Note: **atma** may take up to 2 min to complete. Wait for the task to complete: this involves tuning X, and then  $^1\text{H}$ , for X{ $^1\text{H}$ } experiments. Typing in TopSpin during **atma** or **atmm** may freeze the software (close TopSpin, wait 5 s, re-open TopSpin will usually solve this).

4. Shim as normal, using **topshim gui**. Addition of **convcomp** is important on the cryoprobes for any volatile solvent (all common solvents other than  $\text{D}_2\text{O}$  and DMSO). Look at the Report panel, and wait until it finishes before switching windows or typing anything into TopSpin.

Rerun topshim if the “final Bo stdDev” > 1Hz. Try TUNE → Before on XYs as shown, as that may help. Filtering or centrifuging the sample may be important to obtaining good quality shims.

5. **ased** and adjust parameters as necessary; the most important are usually **d1** and **ns**.


Note: Quantitation is often nontrivial for X-nuclei. See section III for more information.



6. As always on Bruker spectrometers, important information is contained in the pulse sequence listing. Clicking on the **PulseProg** tab will provide a listing. Particularly useful is the bottom comments section of the listing, which provides parameters comments, such as minimum **ns** cycles for any sequence. **halt #** on a minimum **ns** cycle will assist in avoiding artifacts.



7. Note also that some time delays may be calculated, such as DELTA in the example shown below. These delays must not be set to negative values. TopSpin will warn about such occurrences, but only after a **zg**, and will then lock up the acquisition (yes, this definitely should be fixed by Bruker...). When this happens, exit TopSpin completely, restart it, then:
  - i) fix the parameter (in **zpgp30**, **d1** > 100m)
  - ii) **ii restart**.




8. **rga** only needs to be run on  $^1\text{H}$ ,  $^{19}\text{F}$ , or perhaps for  $^{31}\text{P}$  on very concentrated samples (or labeled samples). Won't hurt anything to run **rga**, however, which usually will end up at 203.
9. Background signals can be reduced using longer **de** values, if no other broad signals are observed in the sample. The tradeoff is more 1<sup>st</sup>-order phase correction required as **de** increases, which causes baseline rolling.  $^{19}\text{F}$  is the most common nuclei where such problems occur, but the larger **de** could also be used for  $^{29}\text{Si}$  and  $^{11}\text{B}$ . See staff for more information.
10. Acquire as normal with **zg**. Observe during the run with **tr** or **tr #** to transfer at scan #. **halt #** is recommended where # is a multiple of the minimum **ns** phase cycle (usually 8).
11. A significant issue in many X-direct experiments is **ns** and **d1** settings when observation is not quickly made. A standard sample at a known (high) concentration that can be easily observed is the best method for working through such issues.
11. Adjust the sweepwidth **sw** and spectrum center **o1p** as needed. Bruker digital filters will exclude observation of any signals outside the spectral window.

$$\text{sw} \leq 2000000$$

can be used, but the probe usually can observe signals only within a 200 kHz bandwidth (bw) about **o1p**. An example of a  $^{13}\text{C}$  tune is to the right, where an approx bw is indicated.

If **o1p** is changed by > 50 kHz, the probe should be retuned (e.g., with **atma**).



Reacquire overlapping spectral regions to cover the complete span of chemical shifts for any detection of a novel compound. See *cgf* for more information and assistance.

#### IV. Other comments:

1. Obtainable and digital resolution are often issues with X-spectra, especially when sweepwidths get large, or are changed to access different chemical shift regions. The following guidelines will help:

Always check **aq** after any change in **sw**.

$$\text{obtainable resolution} \sim 1/\mathbf{aq} \quad (1)$$

Eq (1) is not precisely correct, but the dependence is right. Long **aq** acquiring a lot of noise will not assist with resolution. But if one hopes to resolve peaks that are separated by 1Hz, then **aq**  $\geq$  1 sec is recommended.

$$\text{digital resolution} = \mathbf{sw} / \mathbf{si} \quad \text{where } \mathbf{si} \leq \mathbf{1-2 \times td} \quad (2)$$

For **si**: 1 (one full zero-fill) if  $\mathbf{td} = 2^n$ , but can be almost 2 if  $\mathbf{td} = 2^n + 1$ . In any event:

**si = 2 × td** is a good setting.


2. Matched filters are ok for X-spectra: **lb** = 1/**aq**. For **aq** > 1, usually start with **lb**=1, but reduce if helps with resolution.
3. Quantitation: Usually, acquire with default parameters (usually with a faster guess for **d1<sub>fast</sub>**) until a good s/n is obtained. Integrate(deconvolute) the peak(s) of interest.

Acquire another spectrum with **d1<sub>new</sub>**  $\geq$  **d1<sub>fast</sub>** × 2.

If the integral doesn't change, or an important ratio doesn't change, likely all is OK. But if it does change, repeat with yet a longer **d1**. Repeat until the integral/ratio doesn't change. This is a repetition rate method for estimating  $T_1$  values.

4.  $^{13}\text{C}$  DEPT, or INEPT, or APT. Take your choice, they're all OK.

When the J-value is known, INEPT has the best sensitivity. The `ineptrd` pulse sequence is preferred; good pulsewidth calibrations are important: see *cgf* for more info.

5. Indirect detection via  $^1\text{H}$  (or possibly  $^{31}\text{P}$ ) is often superior to the point of being the only option. This is true for  $^{15}\text{N}$  at natural abundance, and for most nuclei that resonate at frequencies below  $^{15}\text{N}$ .
6. Referencing in TopSpin is typically very good. The software automatically references the spectrum to the assumed position of the  $^2\text{H}$  lock signal (e.g.,  $\text{CDCl}_3$  at 7.2 ppm). A very strong recommendation, even so, is to *always* keep a  $^1\text{H}$  spectrum taken in the same conditions (same sample, same temp, etc) with the X spectrum. The automation routine **xref** in TopSpin, or the equivalent routine in MNova , can then be used to check the referencing of the spectrum. In TopSpin:

a) While in the X-spectrum experiment, type **xref**.

b) enter the filename+expno+procno of the referenced  $^1\text{H}$  spectrum

That's all there is to it. With the  $^1\text{H}$  spectrum, you can go back at any time (e.g., many years from now when you've moved elsewhere) and check a  $^{31}\text{P}$  or  $^{19}\text{F}$  (etc.) reference. This could avoid contention over published data, as one example, without having to remake the sample.