

¹¹B and ¹⁰B by Liquid NMR

updated: 2 Apr 2013 (cgf)

Boron is a good nucleus for direct detection in solution by NMR. Relaxation of the isotopes of boron varies considerably depending on the specific chemistry, and experimental optimization quite important. Typically, the relaxation is directly related to the observed linewidths, providing a fairly straightforward method for optimization of the repetition times.

¹¹B is most commonly observed due to its higher magnetogyric ratio (160.4 MHz on a 500) and natural abundance (80.4%). It is spin 3/2.

¹⁰B may be better for compounds where asymmetric bonding leads to broader linewidths, as noted in a recent publication.¹ Since ¹⁰B is spin-3, its relaxation in solid materials is so fast (linewidths then are very broad) as to make it undetectable. ¹⁰B spectra therefore have no background from solid probe components (typically boron nitride). This background can be a major issue in ¹¹B spectra.

¹⁰B is, however, at least 10× less sensitive than ¹¹B, due to its lower magnetogyric ratio (53.7 MHz on a 500) and natural abundance (19.6%). ¹⁰B spectra will therefore take ≥ 10² times longer to obtain than ¹¹B from the same sample, so sample concentrations must be relatively high for ¹⁰B spectra.

Experimental Setup:

Obtain spectra for similar materials, and measure the full width at half-height of the peaks, Δν. For the narrowest peak, estimate the spin-lattice relaxation, T₁, as follows:

$$T_1 \approx T_2 = \frac{1}{\pi \Delta\nu} \quad (1)$$

- (a) Set **aq** = **3×T₁** using the narrowest peak of interest in the spectrum of a similar compound.
- (b) Set **d1** = **30μs**
- (c) Set **ns** very large, and use **halt** when sufficient scans have been acquired.
- (d) If peaks get very narrow — less than 6 Hz (and note that LB adds to the linewidth directly) — switch the pulse sequence to **pulprog** = **zg30** and **aq** ~ **1×T₁**.
- (e) ¹H decoupling may be desirable for samples with narrow peaks: **pulprog** = **zgpg30**, **aq** ~ **1×T₁** and make sure to **getprosol** prior to **rga** or **zg**.

To remove background signals in ¹¹B spectra, acquire an identical spectrum using another tube containing the same amount of the same solvent only. Take the difference of this blank with the ¹¹B spectrum from the sample.

¹ Peter Kiraly, "Background-free solution boron NMR spectroscopy," Magn. Reson. Chem. 50(9) 2012 620-626.