

## Quick Notes for Use of Bruker Avance Spectrometers at NMRFAM:

[cgfry: 99.04.12]

1. Use the Console window to start up xwinnmr: **xwinnmr**
2. Use an Xshell window to look at pulse sequences:  
**cdpp**  
**vi pulsesequence.cgf**  
(it is a good idea to copy the sequence to your initials; that gives you privileges to comment and edit the sequence, since you will own it; use e.g., **cp pulsesequence.fa pulsesequence.cgf**)

One can also use editors within xwinnmr by typing: **edprog(?)**

The default editor that is used can be changed in **setres**.

3. Inside **setres** make sure to turn off *ZGsafe* (to 'no') as otherwise this will prevent **multizg** from continuing without user interaction (defeating the purpose of multizg).
4. To set the temperature, use:  
**te↵** (enter new temp in kelvin)  
**teset↵** (to force the change to the controller)
5. Lock using **lock**  
Select the solvent, and the spectrometer should lock up. Enter **lockdisp** to open the lock display window.
6. Once the lock settles (after a temp change), gradient shim by typing:  
**gradshim↵**

A gradient shimming window will open; make sure the user is your id (e.g., *cgfry*) and not gradshim.

Click on the **Start Gradient shimming** button, and wait until a graph comes up and the SGI finishes pushing the new shim values into the shim unit. At this point, the shims should be pretty good; you can play around with moving Z6 or Z5 significantly and re-gradient shimming to try to improve, or of course one could manually work on the shims.

**wsh** and **rsh** work similarly to the AC/AM's, for writing and saving shims. However, once again Bruker forces all files from all users into a common directory, so make certain to append the suffix *.cgf* (use your own initials). That way:

**rsh \*.cgf** will show just your shim files.

7. We have been using the following sequences:

zgpr.cgf	(poor quality presat—use $\leq 30\text{dB}$ —but ok for pw90 calibrations)
zg3919.cgf	(check quality and O1 setting for 3919 water suppression; p1 must be set = $90^\circ$ )
roesy3919fa.cgf	(seems to be ok; spinlock power 15dB ~ 4kHz)
tocsy3919fa.cgf	(seems a problem with phase cycle; will recommend another when find one that works better)

You can see all pulse sequences by going into the **ased** window and clicking on the down arrow next to the pulsesequence. You can also change it directly by typing **pulprog** and typing in the new name.

8. Parameter sets can be read in to match the sequences listed above by typing:

**rpar \*cgf** (lists just subset written by cgfry)  
**rpar** (lists all parameters sets saved by any user)

**wpar** will save your parameters; all files from all users go into a common directory, so append the suffix initials so you can isolate just your files with the first rpar command above

9. Old data sets can be read in by typing:

**dir** and clicking on the desired set.

In Unix, you can get to your dataset by typing: **cddata**

The path is /u/data/cgfry/nmr or from a remote computer for the 750  
 /dmx750/data/cgfry/nmr

Make new datasets, or new experiment numbers with:

**edc**

10. Check the acquisition parameters by using:

**ased**

ased bring up a panel with pulse acquisition parameters specific to the pulsesequence chosen. If you change the pulse program within ased, exit and re-enter ased before changing other parameters.

Then check other parameters for reasons unknown not included in ased using:

**eda**

o1, td1, nd0, and various other parameters have to be checked in eda or with the command line. Typically use **DQD**.

It is recommended that the pulsesequence be open and checked while going through the parameters for the first time.

11. Check the gradient sequence using **grdprog** at the command prompt. Hopefully the sequence will have the name documented in the comment lines.

Check **const21** etc. as required for the gradient sequence; follow comments in the pulseprograms.

12. Use **zg acqu** to start acquisition and switch to real-time fid observation to check the receiver gain level. Clipping occurs at the top and bottom of the unaltered acqu window ( $\pm 12,000$ ). For a 2D sequence, watch both the first and second rows, don't rely on just the first row as a check. For COSY-type sequences, change **d0** ~ 10ms so the magnetization is echoing back for an **rg** check (remember to set **d0** back after adjusting **rg**).
13. Backup data by using **tar** as usual, being careful to *not* include absolute paths in the input portion of the command:

**cddata**

**tar cvf /shetland\_data/dir/tarfile.tar datasetname**

From a remote computer, **ftp** or **telnet** can connect to **hereford** or **shetland**. Make sure to switch to binary transfer with **bin** before performing the **ftp get** command.

14. Use **edp** to process the data in **xwinnmr**. The comments in the pulse program will tell whether the data was acquired **TPPI**, **States-TPPI**, or **States**.

**xfb** will transform in both directions keeping all quadrants of the data.

**xfb n** will transform, but not keep the imaginary components; much faster, but you won't be able to phase without retransforming here.

Use the **+/-** buttons to turn on negative and both signs of 2d sets.

1d: use **PHASE CUR** and click middle button to define the toggle point, then click-hold the left button on **PH0** and phase at the toggle point. Click-hold the left button on **PH1** to 1st-order phase far away from the toggle point.

**pk** will apply the previous phase correction

2d: use **PHASE CUR ROW** and click middle button in 2d set; click on **1** to place in the 1 window; repeat for **2** and **3** windows. Then do **PH0** and **PH1**. Repeat complete process for **CUR COL** if **F1** needs phasing.

15. See the companion guide for information on how to translate the data into **vnmr**.