

VII. Deuterium NMR on Varian Spectrometers

[updated: 25 July 2010]

1. Working in fully protonated solvent usually works best. The primary problem will be shimming without a lock signal. Use of a deuterated solvent (10:90% ^2H : ^1H) can be used if long (>1hr) acquisition times are needed; this will enable magnetic field drift compensation using the ^2H lock. A ^1H spectrum will typically be acquired to provide referencing.
2. Hook up and view the LOCK-channel with the HP sweeper. Insert the proper cap rod for highest sensitivity ^2H detection [the lock channel can be used to simplify ^2H experiments—ask NMR staff for more info—if samples are at high concentration].
3. a) For bbext on the UNITY-500: detune the lock channel using the recessed screw—use **only** the non-magnetic screwdrivers in the lab!—close to the X-match rod (not at the outer rim, but toward the center of the probe)
b) For hcx and bbswgo on the INOVA-500, or other older probes: turn only the outer black knob/capacitor at the bottom of the probe.
4. Hook up and tune the X-channel on the HP scope to ^2H ; the resonance should not show asymmetry (or dual resonance) if the lock channel is sufficiently detuned.
5. For most samples containing protons, ^1H decoupling will be beneficial (for both NOE and decoupling). A ^1H spectrum will provide the best method of referencing the ^2H spectrum.
– Tune the ^1H channel on the probe.
6. a) In the ACQI panel, turn LOCK off and
b) set the LOCK POWER = 0
c) Detach the lock cable from the probe.

These steps should be performed, unless you need a lock on a ^2H solvent, for long experiments. If the lock is used, you will observe the lock transmitter signal, which can easily be mistaken for a real ^2H signal.

7. INOVA-500 only: remove the ^2H /lock-reject filter that is attached at the transmitter side of the X-rf cables from the probe: moving back from the probe along the X-cables is first the X-specific-filter, then another cable to two filters—one of which is the ^2H reject filter—attached to the preamp box next to the magnet leg.
8. Shim using one of the following methods:

Shim on the solute:

- a) Acquire a normal ^1H 1d and expand around a solute singlet peak (or *known* simple multiplet), using nt=1 ss=0 d1~ T_1 , pw=pw90/3, gain optimized. Do not shim on a solvent peak; these are broadened by relaxation dampening too much to be used for shimming. For per-deuterated samples, a ^2H singlet can be used for shimming.
- b) Expand so the singlet fills $\geq 30\%$ of the display window, and type **gf**.
- c) In acqi window, click on FID and shim on the spectrum.

Gradient Shimming is possible, and can be quicker once the interface is learned:

1st use requires a short setup, best checked by the NMR staff

[Unity: gzlvl=3200 d3=0,0.005 sw=80000 sqsinebell gzsize=5]

[Inova: read ^1H shimmap in for hcx probe]

- a) connect proton cable for ^1H detection
 - b) start PFG shimming software: **gmapsys**↵
 - c) Unity-500 only: menu button: **SETUP**
 menu button: **FIND GZLVL/GZWIN**
 - d) make a gradient map: menu button: **SHIMMAPS**
 Unity only: menu button: **MAKE SHIMMAP**
 Inova spec: menu button: **AUTOMAKE SHIMMAP**
 - e) let the shimmap finish (takes about 1min)
 - f) shim z-shims menu button: **AUTOSHIM**
 - g) always click on menu button: **QUIT**
 - h) may have to go into ACQI window and restart sample spinning
9. Acquire a ^1H spectrum. If ADC or receiver clipping occurs at **gain=0**, try **pw=1**. If clipping still occurs, set **tpwr=45**. You need only enough to reference the spectrum. Save it with **svf**, and archive it with the ^2H spectra.
 10. Move to a new experiment: e.g., if the ^1H spectrum is in exp1, do **jexp2**↵.
 11. Read in ^2H , SOLVENT parameters, and acquire the ^2H spectra. For quantitative spectra, set **dm='nny'** **d1=2** (likely long enough). For best sensitivity, set **dm='yyy'**.
 12. Reference the ^2H spectra using **xref**↵. Save as normal with **svf**.
 13. When finished, you *must*:
 - a) retune the lock channel on the probe
 - b) INOVA-500 only: re-insert the lock-reject filter on the X-channel
 - c) UNITY-500 only: reconnect the lock cabling with ^2H filter to the probe's lock port
 - d) insert and lock the CDCl_3 standard sample.