

## Acquiring Kinetics Data in VNMR

updated: 2011 Dec 1

### Improvements in $^1\text{H}$ Kinetics Experiments (new as of 22 Nov 2011)

Investigations carried out by Paul White and cgfry in fall 2011 demonstrate clearly the adverse effects too short a repetition delay ( $d1+at$ ) can have on kinetics data:

- Peak intensities and integrations are significantly affected, because of incomplete relaxation recovery. Measurement of  $T_1$  values may enable some correction of the results.
- In a standard experiment (e.g., `s2pul`) that does not include homospoil or gradient crushers, the peak intensities become increasingly dependent on chemical shift as the repetition delay is shortened. Somewhere between  $1\text{s} \leq d1+at \leq 3\text{s}$ , the data become unusable for quantitative studies; reasonably accurate post-acquisition corrections to the data cannot be made.

The following set of recommendations should be added to experimental setup protocols for  $^1\text{H}$  kinetics experiments; they minimize errors arising from small repetition delay times, and enable reasonable methods to be used to correct the data. They will always improve data quality, and should therefore *always* be used for quantitative work.

- Use gradient crushing as enabled by `sspul='y'` in the new sequence `s2sspul` for all  $^1\text{H}$  kinetics experiments. This change is effected on our Varian high-field spectrometers (but not yet on Hermes) by setting up the kinetics run as done previously (and described in the next section), and then typing into vnmr:

```
s2sspul↵
```

This macro changes the experiment (from `s2pul` to `s2sspul`) which adds gradient crushing before the `d1` period. All other parameters are left unchanged.




Figure 1. The `s2sspul` pulse sequence with  $G_z-90^\circ[\text{pss}]-G_z$  gradient crusher sequence preceding the `d1` period. The rest of the sequence is standard `s2pul` sequence.

- The gradient crushing period removes any `at` (acquisition time) dependence from the repetition delay: with `s2sspul`, the longitudinal magnetization utilized in the next experiment

is recovered only during **d1**. Recommendations for repetition delays in **s2sspul** experiments thus equate to recommendations for setting the **d1** value.

3. When **d1**  $\geq 3 \times T_1$ , use a 90° pulse: **pw=pw90**.

When **d1**  $< 3 \times T_1$ , use a 30° pulse. Estimations of systematic errors and corrections to the results can be made with knowledge/measurement of the proton  $T_1$  values. Note that  $T_1$  values can change during a reaction, so measurements of  $T_1$  both before and after the reaction may be required. **cgfry** can assist with how make such estimations and corrections.

## Experimental Setup Protocols

**Arrays in VNMR turn acquiring kinetics data into a fairly simple task (Section I). More complex situations can also be accommodated, involving relatively complex acquisition schemes with accurate time-stamps being generated with the spectra.**

### I. Typical (Simplest) Setup

- A. Set up a normal  $^1H$  (or  $^{13}C$ , or  $^{31}P$ , etc.) experiment, and be confident with your temperature calibrations. Acquire a spectrum with good sensitivity and good line shape.
- B. Make sure **d1+at**  $\geq 3 \times T_1$ ; this is not a strict requirement, since it is the *difference* in  $T_1$  values between the two (or more) interesting protons (or  $^{13}C$ , etc.) that matters; if the  $T_1$  values are identical, all values of **d1+at** will give quantitative results (albeit with different signal-to-noise). But some awareness of the  $T_1$  values is important to insure that the data is quantitative. Using  $3 \times T_1$  may introduce up to 5% systematic errors;  $5 \times T_1$  reduces these errors to  $\leq 2\%$ . **at** need be only long enough to provide *sufficient* resolution.
- C. To assist with fast kinetics, shim on a non-reacting sample having identical solution height—or keep the solution height  $\geq 3 \times$  coil length;  $3 \times 16\text{mm} = 48\text{mm}$  for the 500 MHz broadband and hcx probes—and good quality tubes (at least Wilmad 507, 528's being better). These conditions will give the best shims when the reacting sample is introduced, with careful work allowing immediate acquisition of reasonable quality spectra.
- D. For faster kinetics (half-life  $< 10$  min), keep **nt** as small as possible; for slower kinetics, keep **nt** equal to a multiple of 8. The time per experiment,  $t_i$ , will (typically) equal:

$$t_i = \text{nt} \times (\text{d1} + \text{at})$$

- E. If **nt=1** is used for fast kinetics, beware of the center glitch, especially on the UNITY-500. Change **tof** (using **movetof**) so no interesting peaks are at the center of the spectrum.
- F. Array **pad**, the *pre-acquisition delay* parameter, using something like:

**array('pad',50,30,0)**

which will give 50 values all with **pad**=30 secs. The timing in a proton experiment with

**nt=8, d1=0, at=2,**

and **pad** as set above will be as follows, where  $td_i$  = time at mid-point of each experiment:

1<sup>st</sup> pad:        30 s



**A. Non-linear pad arrays:**

A more efficient dataset can be obtained by setting **pad** such that more spectra are acquired early on, and fewer later in the run. In the example above, no data is acquired after 40 mins into the reaction. To acquire data out to multiple hours would require a very large array if **pad** values are kept constant at 30s; many more spectra will be acquired than would really be required. One solution would be to lengthen the **pad** values later in the array as follows:

**pad=0,1,1,1,1,1,2,2,2,3,3,4,5,6,7,8,10,12,14,16,18,20,23,26,29,32,35,40,45,50,55,60,65,  
70,80,90,100,120,140,180,220,300,400,600,800,1400,2000,3000,5000,10000**

which will acquire data out to ~ 7 h. These 50 spectra may much better determine the kinetics than the previous example. The **pad** array above, even so, was made up ad-hoc, not being well-thought out. To best optimize the kinetics run, some modeling of the kinetics is recommended, to estimate the optimum set of pad values.

**B. Arraying nt diagonally with a non-linear pad array:**

In the example above, issues are easy to identify at both extremes of the **pad** array:

- (i) For the early **pad** values, if **nt=8** is used, the time-resolution is 16 secs; reducing **nt**, with **nt=1**, will improve that to 2 secs (but, of course, only if the sensitivity allows).
- (ii) For the later **pad** values, keeping **nt=8** (or smaller) may limit the data quality. Much larger **nt** could be performed with no negative effects in time-resolution, since the reaction is proceeding more slowly.

*It therefore makes sense to increase both **nt** and **pad** throughout the kinetics run.* An example macro that sets **nt** and **pad**, and makes them track with each other (i.e., diagonalizes them) is given in `/vnmr/maclib/padnt.UW`. In this example, the kinetics are modeled by  $k = 0.003 \text{ s}^{-1}$  at the working temperature. **d1+at=3s** are estimated to be a sufficiently long repetition delay. 25 spectra will be acquired. If the sensitivity is allows **nt=1**, we can now accurately plot out the reaction time-course, and estimate the times at which equal amounts of reaction will occur. For 25 spectra, we want data after every 4% of reaction which occur at:

**Example 1:  $k = 0.003 \text{ s}^{-1}$  (**at=2, d1=1**).**

*optimal times = 0,14,28,43,58,74,91,110,129,149,170,193,218,245,274,305,341,  
380,424,476,536,611,707,842,1073*

If **d1+at=3**, and **nt=1** are OK, we can adjust the parameters to match as in **padnt.UW**:

**pad=0,8,8,9,9,10,11,10,7,8,9,11,13,15,11,7,12,15,20,28,24,27,48,63,135**  
**nt=1,2,2,2,2,2,2,4,4,4,4,4,4,8,8,8,8,8,16,16,16,32,32**  
**array=(pad,nt)'**

where the computed mid-spectrum times for each experiment are identical to the *optimal times*, except for the first spectrum with a mid-point at 1.5s. The last statement in the macro is required, which forces **pad** and **nt** values to increment together; this statement *diagonalizes* the square matrix.

Example 2:  $k = 0.0003 \text{ s}^{-1}$  (**at**=2, **d1**=3).

*optimal times* = 0,136,277,426,581,743,914,1095,1285,1487,1702,1932,2179,2446,  
2736,3054,3405,3798,4243,4757,5364,6108,7067,8419,10729

Here, the optimal times are a simple order-of-magnitude increase from Example 1, but the computed parameters are not. Data quality is improved by keeping **nt** a multiple of 8:

**pad**=0,56,62,68,75,83,91,100,71,42,55,70,87,107,130,157,192,152,125,194,288,264,  
319,391,1031

**nt**=8,16,16,16,16,16,16,32,32,32,32,32,32,32,32,64,64,64,64,128,128,256,256

**array**='(pad,nt)'

The *optimal times*, **pad**, and **nt** arrays are computed using an Excel spreadsheet; these two examples are given in /vnmr/maclib/**padnt.xls.zip**.

### C. Acquiring non-linear arrayed kinetics data sets with time-stamp data:

Once the padnt array is set, the data must still be acquired. Obtaining the data with accurate time-stamps is now very useful. To do this, first shim and set up the 1D parameters as discussed above (setting **pad** and **nt** by any method). Then enter:

**kintime**.**↓**

at the vnmr command prompt.

- i) **kintime** queries whether an ASCII file kintimes.txt from a previous run is backed up. kintimes.txt will be written to your current data location (**pwd** at the vnmr command line), and typically should be moved into the dataset folder after data acquisition is complete. For example, suppose I acquire a kinetics dataset and save it using **svf('kinexample')**. This will write the data to the folder

/b6hf/cic/fry/kinexample.fid

If I used **kintime** and **au** to acquire the dataset, a file /b6hf/cic/fry/kintimes.txt will have been written during the kinetics run. To properly save this data, it needs to be moved into the dataset folder using the following command in a Unix terminal window:

mv kintimes.txt kinexample.fid.**↓**

- ii) **kintime** sets the parameter: **wnt='kintimewr'**

Now you are ready to insert the reacting sample. Do so and enter:

**au**.**↓**

to start the acquisition.

**au** is a version of **go** that executes the parameter string **wnt** at the end of each arrayed experiment. **kintimewr** writes four values for each experiment: (i) the experiment #, (ii) the date, (iii) the time-stamp for the  $i^{\text{th}}$  experiment, and (iv) the length of the  $i^{\text{th}}$  experiment.

The time-stamp is the accurate *END-TIME* of the  $i^{\text{th}}$  experiment, not the middle. Taking this time and subtracting  $\frac{1}{2}$  the length of the  $i^{\text{th}}$  experiment —  $\mathbf{nt}_i \times (\mathbf{d1+at})/2$  — will give the correct mid-point time for each experiment.