

## II. X Acquisition (e.g., $^{31}\text{P}$ or $^{13}\text{C}$ )

(22-Jun-98)

### A. Discussion

- The spectrometer acquisition time of an X nucleus experiment should be estimated from similar experiments on the same equipment if possible using eqn (1) or (2) below.  $^{13}\text{C}$  and  $^{31}\text{P}$  experiment times on the 500's and 360 can be estimated from similar experiments on the AC-300 (Athena) using:

$$S/N \propto cB_o^{3/2}t^{1/2}\mathbf{x} \quad (1)$$


where  $c$  = concentration,  $B_o$  = magnetic field strength,  $t$  = experimental time, and  $\mathbf{x}$  = probe filling factor. To estimate the time under different conditions that would give identical  $S/N$ , then:

$$t_{\text{new}} = t_{\text{old}} \cdot \left[ \left( \frac{c_{\text{old}}}{c_{\text{new}}} \right)^2 \left( \frac{B_{o_{\text{old}}}}{B_{o_{\text{new}}}} \right)^3 \right] \quad (2)$$

A reasonable estimate (although direct comparison to experiments on similar compounds is preferred) on the UW Chem MRF equipment can be made starting with the observation that a 0.1 M solution will typically give publication quality spectra in 20 min on the AC-300 (Athena).

- Standard decoupled (NOE-enhanced), quantitative decoupled (no-NOE; Bruker's INVGATE.AU) and NOE-enhanced coupled spectra (Bruker's GATEDEC.AU) can be obtained through use of the S2PUL sequence.
- Nuclei with negative  $\gamma$  values, such as  $^{15}\text{N}$  and  $^{29}\text{Si}$ , are best acquired using polarization transfer sequences; DEPT in general is the preferred sequence over INEPT if more than one  $J_{\text{XH}}$  value is involved; quaternary moieties of negative  $\gamma$  nuclei are best obtained through long-range coupling unless the  $T_1(X)$  values are known to be reasonably short.

### 1d $^{13}\text{C}$ Acquisition (s2pul)



**B. Critical Parameters**

- d1** – relaxation delay; **d1** >  $T_1(^1H)$  to obtain optimum NOE
- dof** – usually =0 when **solvent** set correctly; should be within 5ppm of  $^1H$  coupled to X nucleus of interest
- dpwr** – decoupler power (larger number is higher power); *typically never* > 46
- dmf** – decoupling strength =  $[pw90(^1H \text{ at } dpwr)]^{-1}$
- dmm** – decoupler modulation mode; either **dmm='ccp'** **dseq='waltz16'** **dres=90** or **dmm='ccw'** [both are equivalent] is best for typical compounds
- dm** – decoupler on/off flag (see Table 6)


**C. 1d X{ $^1H$ } Acquisition**

- start by setting the probe parameter appropriately
- to setup parameters (method **b** is recommended)
  - a) either read in a data file (rt) or parameter file (rtp) for standard X{ $^1H$ } acquisition, or
  - ☞ b) click on MAIN MENU SETUP and select nucleus and solvent
  - c) alternatively, run s2pul macro (but normally will need to reset **spin='n'**, **gain=0**; check **temp** and **vttype**, and check **tpwr** (=58 typically) and **pw** (=5-8 $\mu$ s typically).
- use UWMACROS TUNE PROBE TUNE C13 [or **tunec** macro] and tune  $^{13}C$  channel; reattach  $^{13}C$  X cable to X Obs BNC
  - make sure ¼ wave cable is correct for  $^{13}C$
  - make sure low-pass (brown) filter is in-line
  - make sure correct probe cap is inserted (*none* in bbold probe)
  - make sure no external attenuator is in-line
- use UWMACROS TUNE PROBE TUNE H1 [or use **tuneh** macro (will have to reset gain after tunec if UWmacro is not use)] and tune decoupler channel; reattach  $^1H$  decoupler cable to decoupler BNC

**Table 6. Common Decoupler Parameter Settings**

Parameter Settings	Comments
<b>dm='yyy'</b> <b>dpwr=46</b> <b>su</b>	typical for $^{13}C$ acquisition; full on decoupler; keep <b>dpwr</b> <b>£48</b>
<b>dm='nnn'</b> <b>su</b>	turns decoupler off; always finish with this command + <b>su</b>
<b>dm='nny'</b> <b>su</b>	inverse gated (Bruker's INVGATE.AU) mode; gives quantitative (assuming <b>d1</b> >5 $T_1$ !) decoupled spectra
<b>dm='yyn'</b> <b>su</b>	gated spectra (Bruker's GATEDEC.AU); gives coupled spectra but with NOE buildup (DEPT is necessary for $^{29}Si$ , $^{15}N$ and other - $\gamma$ nuclei)
<b>dmm='p'</b>	normal setup for $^{13}C$ acquisition (8/1/97; <b>dmm='w'</b> works fine)
<b>dmm='ccp'</b>	normal setup for DEPT and INEPT most 2d $^{13}C$ experiments, does hard pulses followed by pulsed decoupling during acquisition
<b>dseq='waltz16'</b> <b>dres=90</b>	normal setup for $^{13}C$ (and all X) acquisitions; $^1H$ waltz-16 decoupling
<b>dseq='garp1'</b> <b>dres=1</b>	normal setup for all inverse experiments; X garp-1 decoupling
<b>dmf</b> = [1/pw90]	sets decoupler pulsewidth for composite pulses, where pw90 is the $^1H$ decoupler 90° pulse width at the <b>dpwr</b> setting used for decoupling
<b>dpwr</b> = [ typically ~42]	decoupler power in dB; typically want 90° pulsewidth = 100-150 $\mu$ s

## 1d <sup>13</sup>C Variations in VNMR



- use **dm='yyy'** **su** to turn on decoupler (or UWMACRO DECOUPLER ON; see Table 6 and figure below for more decoupler information)
- use **go** to start acquisition
  - **ga** will automatically apply a **wft** (weighted Fourier transform) following acquisition
  - use **wft** to manually transform
  - **dsc** will display the axis
  - **dsx** will apply **wft dsc(-3)**
- **gain = 30 to 40** should usually work for <sup>13</sup>C acquisition
  - if always getting ADC OVERLOAD beep, then receiver gain is too high
  - it is important the receiver gain be optimized for X nucleus experiments
  - see <sup>1</sup>H section above for how to set the gain accurately
- optimize **sw** by setting right (right mouse button) and left (left mouse button) cursors, and use **movesw**
- set reference by placing cursor close to peak, then type **nl rl** command, e.g., **rl(77p)**

### D. <sup>13</sup>C (X) pw90 and {<sup>1</sup>H} Decoupler Calibrations

- general rule of thumb for these calibrations is:
  - short runs use facility numbers
  - if probe problems are suspected, check pw90's of X and <sup>1</sup>H observe (not decouple)
    - if numbers are close to facility values, probe is likely OK
    - if pw90 is much less than facility value, you are doing something wrong (figure it out! :-)
    - if pw90 is much longer (>1μs) than facility value, find TA or facility staff
  - always perform calibrations for overnight or longer runs for PT-type experiments; for standard decoupled experiments, calibrations are rarely needed
- <sup>13</sup>C and other X nuclei pw90 calibrations require concentrated or labeled samples
  - for <sup>13</sup>C, use 50% benzene in acetone-d<sub>6</sub>; do not degas these samples
  - addition of GdCl<sub>2</sub> or Cr-acac can improve T<sub>1</sub>'s dramatically

- set **tpwr** to desired setting; see Table 2 (on lab wall for most current!!) for calibrations for probes in UW Chem MRF facility
- make sure you are on resonance (set cursor on multiplet, then **movetof**)
- set **d1=20s** or longer depending on sample  $T_1$
- use  $360^\circ$  pulse for final checks always!!
- once the observe pw90 is obtained, assuming in exp1 then **jexp2**
  - move  $^1\text{H}$  cable to observe port, and find resonance for benzene doublet (take one scan, move cursor to middle of doublet and use **movetof**); write this value of tof down as  $^1\text{H}$  dof
- assuming  $^{13}\text{C}$  still in exp1 and  $^1\text{H}$  in exp2, do **jexp1 mp(3) jexp3**
  - put  $^1\text{H}$  on resonance by setting **dof** =  $^1\text{H}$  dof from above
  - setup decoupler calibration experiment with UW macro **pwxdcc90**
  - set **pwxlvl** = **60** and find **pwxc** where antiphase doublet nulls; this is decoupler hard 90
  - set **pplvl** = **pwxlvl** and **pp** = **pwxc** for **hetcor** and **dept** experiments
  - set **pwxlvl** appropriately for decoupling (probe dependent; see Table 2 on wall) and find **pwxc** where antiphase doublet nulls; set **dpwr** = **pwxlvl** and **dmf** =  $1/(\text{pwxc} * 1\text{e-}6)$

### ***E. 1d $^1\text{H}$ Data Workup and Plotting***

- typically will want **lb** = **2** or **3** for  $^{13}\text{C}$  experiments
- see 1d  $^1\text{H}$  section E for plotting description