

## VII. COSY – 2d Homonuclear Correlation

(21-Jul-01)

### A. Discussion


- Absolute value mode COSY often provides sufficient homonuclear correlation data in a very efficient, robust experiment.
- $J_{HH} \geq 4$  Hz can typically be observed, but this is dependent on **sw** and **ni**

$$\text{digital resolution in } F2 = \frac{sw}{np} \quad \text{digital resolution in } F1 = \frac{sw1}{2ni}$$

assuming zero-filling only in F1 (not common in F2). For very large ni, very small couplings can be observed—often well below the linewidth of a high-resolution 1D experiment.

- A common rule-of-thumb is that  $J_{HH} \geq \text{dres1}/3$  will be observed in a COSY spectrum; this estimate is definitely a guideline only.  $\text{dres1} \equiv \text{digital resolution in } F1$  defined above.
- The most efficient method of observing small  $J_{HH}$  is with long-range COSY, which involves this sequence and **tau** > 0 (typically **tau** = 50-200 ms).
- **cosy** loses sensitivity from the strong resolution enhancing sinebell (and sinebell-squared) apodization functions used to overcome the absolute value processing. Even so, minimum phase cycling (**nt**=4) usually gives sufficient sensitivity. GCOSY is even faster (**nt**=1, **nt**=2 is better).
- DQF-COSY (next section) removes all singlets, and gives much cleaner diagonals. Use this sequence after a COSY if crosspeaks close to the diagonal or on a strong singlet region are needed.
- **gcosy** is preferred if a PFG (e.g., hcx or bbswg) probe is installed. **nt**=1 can be used with this sequence, reducing total acquisition time by 4.

### 2d Absolute Value COSY (relay)



### B. Critical Parameters

- p1** = 90° for maximum sensitivity; 45° (typical) to reduce intermultiplet crosspeaks
- ni** = usually set satisfactorily by **cosy** macro;  $sw1/2ni$  gives usable resolution; must have  $\leq 12$  Hz, but want  $\sim 6$  Hz if time allows for best results
- nt** = 4 is minimum phase cycle
- sw** = set using **movesw** macro (preceded by boxing selected area of spectrum)
- sw1** = must equal **sw** for **foldt** macro (commonly used; recommended for fast-COSY)

**tau** = set 50-150 ms typical for long -range COSY  
**d1** = best set  $\geq 2 \times T_1$ , but often run with **d1**  $\sim T_1$

### C. COSY 2d Acquisition

- for COSY, setup as in  $^1H$  1d and optimize **sw** with the **movesw** macro, then retake the data to ensure **sw** is set correctly
  - move the parameters: assuming you ran the  $^1H$  in exp1 and have nothing important in exp2, the command **mp(2) jexp2** entered from exp1 will work
  - from exp2 then use SETUP SEQUENCES COSY or enter the macro **cosy** or **gcosy**
  - check that **ni** ( $sw/2ni \leq 12\text{Hz/pt}$ ,  $\sim 6\text{Hz/pt}$  if time allows) **d1** and **np** are correct (check **time**), then enter **au** to run cosy

### D. Calibrations

- can use facility pulse width calibrations; COSY is very forgiving with pulse widths
- even so, performing a **pw90** calibration is always recommended for all 2d experiments

### E. 2d Data Workup and Plotting

- COSY data can be processed with the **do2d** macro, or use the **wft2d** command (see also **man('cosy')** or **man('gcosy')**
  - **do2d** applies symmetrization with the **foldt** macro
- general display and plot commands

**dconi** ; displays 2d's with color map (faster than contours)

plot projections before issuing the following command using menus (use the PLOT button)

**pcon** ; plots contours

**pconpos** ; plots phase sensitive contours with positive peaks having 10 contours and neg peaks having 1 contour

**pconneg** ; opposite of pconpos

**plot2dhr** ; plots with high-resolution traces (must have 1d in another experiment)