

IV. INEPT – Insensitive Nuclei Enhanced by Polarization Transfer

(22-Jun-98)


A. Discussion of INEPT

See the discussion section in the previous section on DEPT.

- INEPT should be used if one J value is involved and it's size is. Otherwise use of DEPT is recommended.
- The previous statement is not intended to dissuade students from experimenting with the different version of polarization transfer. (After all, the fearless facility manager has been known to give not-such-good advice on occasion. Please do let me know your observation/preferences if you try both INEPT and DEPT on same/similar compounds.)

Refocused Decoupled INEPT

(ineptrd.c [written by cgfry] from Bruker's INEPTRD.AU)



B. Critical Parameters

- d1** – relaxation delay; typically = 1-2s
- j** = 140Hz; *change if you want to observe X with $J_{\text{XH}} > 180$ Hz or < 110 Hz*
- INEPT is sensitive to accurate J coupling being entered; use DEPT unless you are fairly certain you know J reasonably accurately
- pw, tpwr** – observe X 90° pulse width **pw** at power level **tpwr**
- pp, pplvl** – high power ^1H 90° pulse width **pp** at power level **pplvl**
- mult** – multiplier for **d3** delay, where **mult**= **1** all positive, **2** -CH2- nulled, **3** -CH2- inverted from -CH3 and -CH<

C. INEPT Acquisition

- the INEPTRD macro (run by the menu SETUP SEQUENCES INEPT) is *not* the standard Varian sequence, but in my opinion is superior; it incorporates all aspects of Bruker's INEPTRD and INEPTP sequences for decoupled and coupled acquisition, respectively
- for short runs, use facility calibrations for **pw, tpwr, pp, pplvl, dpwr** (**£46**) and **dmf**

- FILE SETUP SEQUENCES INEPT will setup Inept correctly, including for non- ^{13}C acquisitions
 - **mult** is set by the number of coupled protons and **j** ; **d3** = **mult/4j**
 - the **d2** interpulse delay is set according to **j**
- for overnights or longer runs, recalibrate (at least) observe and (best) decoupler pulse widths
- set **mult** as needed; see below

D. Calibration

- see ^{13}C section for nominal ^{13}C (X) and ^1H decoupler calibrations
- The delay **d2** and **mult*d3** are calculated by the DEPT macro as follows:

$$\mathbf{d2} = \frac{1}{4J} \quad \text{echo period involving two } \mathbf{d2} \text{ delays creates } ^1\text{H} \text{ antiphase state}$$

$$\mathbf{d3} = \frac{\mathbf{mult}}{4J} = \frac{1}{2pJ} \left[\sin^{-1} \left(\frac{1}{\sqrt{n}} \right) (\text{rad}) \right] = \frac{1}{360^\circ J} \left[\sin^{-1} \left(\frac{1}{\sqrt{n}} \right) (\text{deg}) \right]$$

where $n = \# \text{ I nuclei coupled with } J_{XH}$

example: suppose have trimethylsilyl, and want to detect ^{29}Si , and $J_{\text{Si-H}} = 2 \text{ Hz}$,

then $\mathbf{d2} = 1/(4 \times 2 \text{ Hz}) = 125 \text{ msec}$

$$\text{and } \mathbf{d3} = \frac{1}{360^\circ \times 2\text{Hz}} \sin^{-1} \left(\frac{1}{\sqrt{9}} \right) = \frac{19.47^\circ}{360^\circ \times 2\text{Hz}} = 27 \text{ msec}$$

$$\text{or } \mathbf{mult} = \frac{2}{180^\circ} \sin^{-1} \left(\frac{1}{\sqrt{9}} \right) = \frac{2 \times 19.47^\circ}{180^\circ} = 0.216 \quad (\text{same as in DEPT setup})$$

These values are calculated internal in the /vnmr/psglib/ineptrd.c pulse sequence code.

E. Data Workup and Plotting

- workup is similar (identical) to ^{13}C 1d and DEPT