

### III. Performing $^1H T_1$ Estimates

[updated: 13 July 2010]

1. Acquire a normal  $^1H$  1d spectrum, as described in Section I (sbs\_1H1d).
2. Calibrate and set **pw90**, as described in Section II (sbs\_1Hpw90).
3. The macro **T1est** will automatically setup parameters as follows.

a) **p1=pw90\*2 d1=10 dps**

manual checks:

set **nt=1 ss=0 vp=80 d2=0.01**

**ga** should show inverted magnetization everywhere.

**d2=10 ga** should show positive magnetization everywhere, except perhaps for solvent (or other very slow relaxing protons) peaks.

Adjust **d2** to give crossovers for important multiplets (see 4 below).

b) Using an array for **d2** is OK, but presents the problem of what to set **d1** to

→ **d2=0.1,0.25,0.4,0.6,0.9,1.3,2.0,3.5 d1=10 go dsa** is usually a good setup for most compounds, but beware that **d1** may be too small, especially for samples sealed in non-O<sub>2</sub> containing atmospheres.

Acquire with: **go** then **dsa**

4. Note the crossovers for each important multiplet, and note these points as **d2<sub>null</sub>** times. Especially note the fastest and slowest relaxing protons of interest.
5.  $T_1 \approx 1.4 \times d2_{null}$  for each proton. If  $T_1 > 1/3 \times d1 = 3$  sec, then **d1** was too short during the experiment. Repeat with a longer **d1** value.
6. See the  $T_1$  section in VUG for more information, especially how to obtain quantitative  $T_1$  values (rather than the estimates provided by **d2<sub>null</sub>**). A brief outline on how to analyze and plot quantitative data in VNMR is:
  - a) set **Th** (button on the **ds** menu) ;all peaks above line will be found by **dpf**
  - b) **dpf** ;display peak frequencies
  - c) **fp** ;writes list of peak amplitudes
  - d) **t1** ;performs analysis and outputs T1 analysis
  - e) **expl** ;displays a plot of the T1 fits

7. To return to normal acquisitions, set:

**T1off**

or equivalently: **p1=0 d2=0 vp=20 dps**

→ useful related commands: **da ff dc dg dssh dssl pl('all')**