

## Inverse 2D Heteronuclear Correlation Experiments

### Step-by-step setup for inverse 2D heterocorrelation experiments:

1. Acquire a normal 1D proton spectrum in exp1.
  - Turn the spinner off, and tune both X (always first) and 1H channels of the probe.
  - Shim normally, including X and Y shims as required to achieve reasonable line shape.
  - Calibrate **pw=pw90** and check the setting of **gain**.
  - Set **d1**  $\geq T_1$ .
  - Reacquire the 1D  $^1H$  spectrum, and optimize **sw** with ~10% of the spectrum as good baseline on both edges. Use **movesw ga** to effect the change. Do not exclude any solute peaks except –OH or –NH, etc.

2. **jexp2** MAIN MENU SETUP NUC,SOLV [C13 solvent]

– **nt=1 ss=0 go** <wait for acquisition> **dsx**

The above will show a one-scan spectrum; if the solvent contains carbons, they will usually show, e.g.,  $CDCl_3$  at 77ppm, which you can then reference directly. **xref** is preferred for referencing (follow it with **rl**).

– **movesw ga** ;sets up the X-nucleus 1D **sw** and **tof** which will become **sw1** and **dof** for the 2D hsqc/hmbc

→ 160p to 0p is usually ok for hsqc

→ 200p to 0p usually ok for hmbc unless keto carbons are further downfield

→ **cr=160p delta=160p** will set first range exactly (although eyeball is ok)

This important step optimizes the sweepwidth in the indirect dimension, which optimizes the resolution for fixed experiment time; there should be ~10% baseline on either end of the spectrum, but no more if it can be avoided (err on the side of too large an **sw**, however, if you are uncertain about the range of chemical shifts).

3. **jexp3 mf(1,3) dsx** ;the  $^1H$  spectrum is used for setting up hsqc-hmbc exp types

Type in:	<b>HSQCAD</b>	(or <b>HSQC</b> )	for best and <b>mult=2</b> spectra, or
	<b>gHSQCAD</b>	(or <b>gHSQC</b> )	for <b>mult=0 nt=1</b> spectra, or
	<b>gHMBCAD</b>	(or <b>gHMBC</b> )	for n-bond (long-range) data

*See the discussion at the end of this document for more about these variants.*

Answer **2** (for this example) when asked for the location of the  $^{13}C$  experiment

→ does **sw** (exp2) → **sw1** (exp3)

**tof** (exp2) → **dof** (exp3)

and transfers referencing parameters

**4. Important checks to perform every experiment:**

- a) cables     $^1H$  on observe    (change needed on U500 only)  
                $^{13}C$  on decouple    (change needed on U500 only)  
                $^{13}C$  filter in-line (the low-pass filter will not work!)
- b) check **nt=2xj**    ;**nt=2** is minimum for **HSQC**, **nt=1** is ok for **gHSQC** if mult=0)  
                           ;**nt=8** is the recommended minimum phase cycle for **hmbc**  
 check **ni**            ;**dres1=sw1/(2xni)**    need enough to resolve  $^{13}C$  chem shifts  
 check **d1**  $\geq 1 \times T_1$     for  $^1H$  of interest  
 check **mult**  $\rightarrow$  = 0 all peaks are positive  
                           = 2  $-CH_2-$  are inverted wrt  $-CH<$  and methyls
- c) Check that the 1<sup>st</sup> row/spectrum; wait until FID gets to 2 in the STATUS panel, then:

**dsx dc va**

If you can see proton signals, then **nt** is OK; in fact, maybe 4x larger than necessary.

**Other notes:**

- check experimental time **time\_1** ; **nt**, **ni** and **d1** are the primary factors here
- **ff** may be needed on 1<sup>st</sup> processing (**prun**)
- **dc** will keep 1D traces from going below bottom of screen
- “**waveform timing error**” often displays upon **go**; this warning is OK, that the waveform cannot be produced exactly (but it’s close enough)

**5. Processing:**

- **prun dqcon plot2dhr** are the simplest/best methods to process, display and plot  
       **wft2da**    ; to transform data manually, see the **dqcosy** in VUG for details  
       **dconi**    ;to display or DISPLAY ...  
       **dqcon**    ;gives a nicer contour plot, but best done only on expansions (can be slow)  
       **pcon page**    ;simplest command for plotting (**plot2dhr** can also be used)  
       **plot2dhr**    ;queried macro for plotting 2Ds with high-resolution 1Ds as traces
- 2x [type **setLP1(2\*ni)**] to 4x [**setLP1(4\*ni)**] linear prediction in the indirect, F1 dimension should normally be applied during final processing of hsqc/hmbc type data; will improve sensitivity and resolution of the data. [use **prun(2)** or **prun(4)** ]

## HETEROCORRELATION 2D EXPERIMENT VARIANTS

**HSQC**, **HSQCAD**, **gHSQC**, and **gHSQCAD**, as well as **HMQC** and **gHMQC**, are all similar experiments, providing chemical-shift correlations between  $^1H$  and  $^{13}C$  having one-bond  $J$ -coupling of ~110 to 170 Hz (approx. value set by **j1xh**). **gHMBC** and **gHMBCAD** provide chemical-shift correlations using multiple-bond  $J$ -coupling (set by **jnxh**, typically between 2 to 25 Hz); one-bond correlations are filtered out (at size **j1xh**).

The versions starting with **g** use pulsed-field gradients (PFG) for *coherence selection*: the gradients reduce artifacts, but also reduce sensitivity by a factor of 2. Empirically, we find the higher sensitivity overcomes artifact noise in **HSQC** and **HSQCAD**, whereas artifact reduction is crucial to **gHMBC** and **gHMBCAD**; these four are the recommended experiments.

The versions ending in **AD** are relatively new experiments that improve data quality by compensating for mismatches in  $J$ -coupling (e.g., **j1xh** is unlikely to match all the 1-bond  $J_{CH}$  couplings in a sample), and by improving the X-nucleus [ $^{13}C$ ] 180° pulses. Initial experiments with **HSQCAD** have demonstrated clear improvements over **HSQC** and **gHSQC**, especially for *edited* experiments (a DEPT-135 analog; when **mult**=2).

**HSQCAD** with **nt=2** (**mult**=0 or 2) is the recommended 1-bond experiment.

**gHSQC** (and **gHSQCAD**) has a significant advantage when working with **nt=1**, halving the time from an **nt=2** HSQC experiment. The sample concentration must of course be sufficiently high to allow good signal-to-noise with a 1-scan experiment (20mM is enough). For *edited* hsqc (with **mult=2** rather than the default **mult=0**), **nt=2** must be used, and **HSQCAD** is strongly recommended. *Edited* hsqc is analogous to a DEPT-135, in presenting  $-CH_2-$  inverted with respect to  $-CH<$  and  $-CH_3$ , with some loss in sensitivity compared to the **mult=0** experiment (which is analogous to DEPT-45).

**HSQC** (and **HSQCAD**) gives narrower lines than **HMQC** by removing  $^1H-^1H$   $J$ -couplings during the F1 ( $^{13}C$ ) evolution.

**HSQC** is more sensitive to pulse-width errors than **HMQC**, since it contains more 180° pulses, and also gives a bit more phase distortion. **HMQC** is therefore preferred on older spectrometers (**HSQC** is preferred in our facility, even on the UNITY).

**HSQC** is better for high-MW (or similarly for low temperatures and/or viscous solvents) compounds since the internal times in the experiment are shorter than **HMQC**.

All six 1-bond experiments are processed with **gaussian** or **sqcosine** apodization.

All six 1-bond experiments can be processed with 2× (or 4×) linear prediction in the indirect, F1 dimension using **setLP1(2\*ni)** [see also the PROCESS tab]. Linear prediction can provide significant sensitivity and resolution enhancements, but may also generate(enhance) artifacts.

**gHMBCAD** is recommended for observing  $^1\text{H}$ - $^{13}\text{C}$  long-range J-couplings.

Since delays that are based on the long-range couplings (set by **jnxh**) are large, the sequence length is kept as short as possible. More importantly, the inherent variation in long-range  $^1\text{H}$ - $^{13}\text{C}$  couplings in any compound prevents refocusing of antiphase magnetization using a single delay period.  $^{13}\text{C}$  decoupling is therefore not possible, and all peaks will be  $^{13}\text{C}$ -coupled in the  $^1\text{H}$  (direct, F2) dimension. It is not uncommon that a **gHMBCAD** must be run twice, or even three times, to cover the range of possible couplings in a compound, e.g., with **jnxh** = 4, 9, and 16 (*not* as an array, but in three separate experiments). The elimination of artifacts by use of gradients (PFG) in **hmbc** is crucial to the success of the experiment; because of this, there is no non-gradient version of **gHMBCAD**.

**gHMBC** is absolute value in both dimensions, and processed with **sinebell** apodization, 2× [**setLP1(2\*ni)**] or 4× [**setLP1(4\*ni)**] linear prediction in F1 (**proc1='lp'**; see the PROCESS tab) is recommended. [**prun** does **sinebell** by default; use **prun(2)** or **prun(4)** to perform 2× or 4× linear prediction]

**gHMBCAD** is a newer version of the experiment that improves data quality by improving the X-nucleus [ $^{13}\text{C}$ ] refocusing pulse, and being phase-sensitive in F1 (remaining absolute value in F2).

2×F1 (to 4×F1) linear prediction is recommended, and is setup by typing:

**setLP1(2\*ni)**

Apodization reset (*after* setting up lp) by typing:

**sq sinebell('f2') sq cosine('f1')** [use **prun** or **prun(2)** ]