VNMR Pulse Sequences

Programming and Hardware Aspects Pub. No. 01-999014-00, Rev. A0398



VNMR Pulse Sequences Programming and Hardware Aspects Pub. No. 01-999014-00, Rev. A0398

By Rolf Kyburz

rolf@nmr.varian.ch Varian International AG Zug/Switzerland

Technical editor: James Welch

Copyright ©1996, 1998 by Varian, Inc. 3120 Hansen Way, Palo Alto, California 94304 http://www.varianinc.com All rights reserved. Printed in the United States.

The information in this document has been carefully checked and is believed to be entirely reliable. However, no responsibility is assumed for inaccuracies. Statements in this document are not intended to create any warranty, expressed or implied. Specifications and performance characteristics of the software described in this manual may be changed at any time without notice. Varian reserves the right to make changes in any products herein to improve reliability, function, or design. Varian does not assume any liability arising out of the application or use of any product or circuit described herein; neither does it convey any license under its patent rights nor the rights of others. Inclusion in this document does not imply that any particular feature is standard on the instrument.

The rights to use the information provided in this material is retained by Varian. All material is provided solely for your use under the VNMR software license and may not be transferred.

UNITY *INOVA*, UNITY *plus*, UNITY, *GEMINI 2000*, *GLIDE*, VXR, XL, VNMR, VnmrS, VnmrX, VnmrI, VnmrV, VnmrSGI, MAGICAL, AutoLock, AutoShim, AutoPhase, limNET, Ultra•nmr, Indirect•nmr, Auto•nmr, Triple•nmr, MagicAngle•nmr, Proton•nmr, Bioproton•nmr, ASM, and SMS are registered trademarks or trademarks of Varian, Inc. Sun, SunOS, Suninstall, SPARC, SPARCstation, Sun-3, Sun-4, SunCD, SunView, and NFS are registered trademarks of Sun Microsystems, Inc. and SPARC International. Oxford is a registered trademark of Oxford Instruments LTD. Ethernet is a registered trademarks or trademarks of their respective holders.

Overview of Contents

| Disclaimer | | 13 |
|-------------|-----------------------------------|-----|
| Foreword | | 13 |
| Acknowledg | ments | 14 |
| Conventions | in This Manual | 14 |
| Chapter 1. | Overview | 15 |
| Chapter 2. | Sequence Generation: seqgen | 17 |
| Chapter 3. | Object Library Generation: psggen | 29 |
| Chapter 4. | Гime Events | 37 |
| Chapter 5. | Submit to Acquisition: go | 55 |
| Chapter 6. | Acquisition Process | 59 |
| Chapter 7. | Digital Components | 65 |
| Chapter 8. | Acquisition CPU and Acode | 69 |
| Chapter 9. | Pulse Programmers | 85 |
| Chapter 10. | Phase Calculations | 95 |
| Chapter 11. | Phase Tables | 115 |
| Chapter 12. | AP Bus Traffic | 137 |
| Chapter 13. | Acquisition CPU Communication | 145 |
| Chapter 14. | Repeating Events | 147 |
| Chapter 15. | Decisions | 155 |
| Chapter 16. | Waveform Generators | 163 |
| Chapter 17. | Pulsed Field Gradients | 199 |
| Chapter 18. | Acquiring Data | 205 |
| Chapter 19. | Multidimensional Experiments | 215 |
| Chapter 20. | Solid-State NMR Experiments | 227 |
| Chapter 21. | (Micro)Imaging Experiments | 239 |
| Chapter 22. | Role of Macros and Parameters | 241 |
| Chapter 23. | Putting It All Together | 255 |
| Chapter 24. | Syntax Guidelines | 259 |
| Chapter 25. | Debugging a Pulse Sequence | 265 |
| Index | | 269 |

Table of Contents

| Disclaimer | 13 |
|--|----|
| Foreword | 13 |
| Acknowledgments | 14 |
| Conventions in This Manual | 14 |
| Chapter 1. Overview | 15 |
| 1.1 Pulse Sequence Execution | 15 |
| 1.2 What to Expect in This Manual | 16 |
| Chapter 2. Sequence Generation: seqgen | 17 |
| 2.1 Modifying the File for dps | 17 |
| 2.2 Running the make Program | 18 |
| 2.3 Calling the C Preprocessor | 19 |
| 2.4 Checking Syntax with lint | 22 |
| 2.5 Compiling and Linking | 24 |
| Chapter 3. Object Library Generation: psggen | 29 |
| 3.1 How Are Object Libraries Generated? | |
| 3.2 Adding Changes to the Object Libraries | |
| 3.3 Adding a New Precompiled Object | |
| Chapter 4. Time Events | 37 |
| 4.1 How Do Delays Work? | |
| Simple Delays | |
| Delays With Homospoil Pulse | 38 |
| Other Delays | 39 |
| 4.2 How Do Pulses Work? | 40 |
| Pulses on the Observe Channel | 40 |
| Simple Pulses on Other RF Channels | 48 |
| Simultaneous Pulses on Different RF Channels | 49 |
| Composite Pulses | |
| Considerations for the Delays Following the Last Pulse | 50 |
| 4.3 Other State-Related Pulse Sequence Statements | |
| Direct Gating | |
| Implicit Gating | 53 |
| 4.4 Basic Purpose of a Pulse Sequence | 53 |
| Chapter 5. Submit to Acquisition: go | 55 |
| 5.1 The Tasks for go | 55 |
| 5.2 Tasks for the Pulse Sequence Executable | 57 |

Table of Contents

| 5.3 Using go('acqi') | 58 |
|---|------------|
| Chapter 6. Acquisition Process | . 59 |
| 6.1 Starting the Acquisition Operating System | 59 |
| 6.2 Queuing and Starting the Acquisition | |
| 6.3 Downloading the FID | |
| 6.4 Controlling Acqproc | |
| Chapter 7. Digital Components | . 65 |
| 7.1 Main Boards | 66 |
| 7.2 Bus Structures | 67 |
| Chapter 8. Acquisition CPU and Acode | . 69 |
| 8.1 CPU Address Space | 69 |
| 8.2 Looking at Acode | |
| Methods of Interpreting the Contents of Acode Files | |
| 8.3 Structure of Acode Files | |
| Acode File Header LC Data Structure | |
| The AUTOD Data Structure | |
| The Instruction Section | |
| 8.4 Acode Interpretation | 82 |
| FIFO Flow | |
| Acode Size Limitations, Acode Buffering | 83 |
| Chapter 9. Pulse Programmers | . 85 |
| 9.1 Layout of the Pulse Programmer | 85 |
| 9.2 Fast Bits | 88 |
| 9.3 Timers and Timer Words | 89 |
| 9.4 Problems with Timer Word Errors | 92 |
| 9.5 Timer Words and Fast Bits in the Acode | 93 |
| Chapter 10. Phase Calculations | . 95 |
| 10.1 How Do Phase Calculations Work? | 95 |
| The Tools | |
| Phase Calculations in the Acode | |
| 10.2 Case 1: Decoding Phase Calculations | |
| 6 | 101 101 |
| | 101 |
| | 106 |
| | 106 |
| 10.5 Steady-State Phase Cycling | 109 |

| 10.6 C Constructs and Phase Calculations | 110 |
|---|-----|
| 10.7 Why Phase Calculations? | 111 |
| 10.8 Real-Time Random Numbers | 112 |
| Chapter 11. Phase Tables | 115 |
| 11.1 Basic Syntax | 115 |
| Shorthand Notation | 116 |
| Advanced Features | 117 |
| How Does a Table Work? | 119 |
| 11.2 Inline Phase Tables | |
| 11.3 Table Math | |
| 11.4 Phase Tables in the Acode | 122 |
| 11.5 Tables vs. Real-Time Calculations | |
| Point-to-Point Comparison | 123 |
| Comparison by Examples | 124 |
| 11.6 Combining the Best of the Two Worlds | |
| 11.7 Using Tables as Source for Random Numbers | 135 |
| Chapter 12. AP Bus Traffic | 137 |
| 12.1 What Is the AP Bus | 137 |
| 12.2 What Devices are Driven by the AP Bus? | 139 |
| 12.3 AP Bus Words in the Acode | |
| 12.4 Timing Considerations | 141 |
| Chapter 13. Acquisition CPU Communication | 145 |
| 13.1 Regular Pulse Sequence Communication | 145 |
| 13.2 Diagnostics and Error Output | |
| Chapter 14. Repeating Events | 147 |
| 14.1 C Loops | |
| 14.2 Real-Time Loops | 148 |
| 14.2 Real Time Loops 14.3 Hardware Loops | 150 |
| Chapter 15. Decisions | 155 |
| 15.1 Decisions and Branchings in C | 155 |
| Decisions Set by the status Statement | 159 |
| Checking Flag Parameters | 160 |
| 15.2 Real-Time Decisions | 160 |
| Programming Real-Time Decisions | 160 |
| Generating the Flag Variable | 162 |
| Chapter 16. Waveform Generators | 163 |
| 16.1 How Does a Waveform Generator Fit Into the System? | 163 |

Table of Contents

| 1 | 16.2 | How Does a Waveform Generator Work? | 165 |
|---|--|---|--|
| | | Sequence of Events in a Waveform Generator | 167 |
| | | How Are Patterns Stored in a Waveform Generator? | 168 |
| | | Waveform Generator Instruction Words | 169 |
| | | Waveform Generator Data File | 170 |
| | | Executing Waveform Generator Patterns | 172 |
| 1 | 16.3 | Using Waveform Generators for Shaped Pulses | 172 |
| | | Programming Shaped Pulses: An Example | 173 |
| 1 | 16.4 | Using Waveform Generators for Programmed Modulation | 177 |
| | | Programming Pattern Decoupling and Spinlock Experiments | 178 |
| | | How Does Pattern Modulation Work Internally? | 183 |
| 1 | 16.5 | What If a Waveform Generator Is Not Available | 189 |
| | | Programmed Decoupling | 189 |
| | | Shaped Pulses | 190 |
| 1 | 16.6 | Using a Waveform Generator for Shaping Gradient Pulses | 195 |
| | | | |
| Chap | ter : | 17. Pulsed Field Gradients | 199 |
| 1 | 17.1 | Pulse Sequence Statements for PFG Gradient Control | 199 |
| 1 | 17.2 | Shaping Pulsed Field Gradients | 202 |
| 1 | 17.3 | PFG Experiments Using Homospoil Pulses | 204 |
| | | | |
| | | | |
| Chap | ter : | 18. Acquiring Data | 205 |
| | | 18. Acquiring Data Implicit Acquisition | |
| 1 | 18.1 | | 205 |
| 1 | 18.1 18.2 | Implicit Acquisition | 205 208 |
| 1 | 18.1 18.2 18.3 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences | 205 208 209 |
| 1 | 18.1 18.2 18.3 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting | 205 208 209 |
| 1 | 18.1 18.2 18.3 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences | 205 208 209 210 |
| 1 | 18.1 18.2 18.3 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals | 205 208 209 210 210 |
| 1 | 18.1 18.2 18.3 18.4 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts | 205 208 209 210 210 212 |
| 1 | 18.1 18.2 18.3 18.4 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting | 205 208 209 210 210 212 213 |
| 1 | 18.1 18.2 18.3 18.4 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting | 205 208 209 210 210 212 213 213 214 |
| 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays | 205 208 209 210 212 213 214 215 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays | 205 208 209 210 212 213 214 215 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays 9. Multidimensional Experiments Indirect Time Domain Incrementation | 205 208 209 210 212 213 214 215 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays 19. Multidimensional Experiments Indirect Time Domain Incrementation nD Quadrature Detection | 205 208 209 210 212 213 214 215 215 217 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays 19. Multidimensional Experiments Indirect Time Domain Incrementation nD Quadrature Detection Absolute Value nD Experiments | 205 208 209 210 212 213 214 215 215 217 217 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays Indirect Time Domain Incrementation nD Quadrature Detection Absolute Value nD Experiments Phase-Sensitive nD Experiments: States/Haberkorn/Ruben | 205 208 209 210 212 213 214 215 215 217 217 217 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays 19. Multidimensional Experiments Indirect Time Domain Incrementation nD Quadrature Detection Absolute Value nD Experiments Phase-Sensitive nD Experiments: States/Haberkorn/Ruben Axial Peak Displacement (FAD) Phase-Sensitive nD Experiments: TPPI Phase-Sensitive nD Experiments: TPPI | 205 208 209 210 212 213 214 215 215 217 217 217 219 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays 19. Multidimensional Experiments Indirect Time Domain Incrementation nD Quadrature Detection Absolute Value nD Experiments Phase-Sensitive nD Experiments: States/Haberkorn/Ruben Axial Peak Displacement (FAD) Phase-Sensitive nD Experiments: TPPI Phase-Sensitive nD Experiments: Arrayed TPPI Folding in Indirect Dimensions | 205 208 209 210 212 213 214 215 217 217 217 217 219 219 |
| 1 1 1 1 1 1 1 1 1 1 1 1 | 18.1 18.2 18.3 18.4 18.5 ter : | Implicit Acquisition Explicit Acquisition Multi-FID Sequences Receiver Phase Shifting Detection of NMR signals Quadrature Receiver Phase Shifts Small Angle Receiver Phase Shifting Housekeeping Delays 19. Multidimensional Experiments Indirect Time Domain Incrementation nD Quadrature Detection Absolute Value nD Experiments Phase-Sensitive nD Experiments: States/Haberkorn/Ruben Axial Peak Displacement (FAD) Phase-Sensitive nD Experiments: TPPI Phase-Sensitive nD Experiments: TPPI | 205 208 209 210 212 213 214 215 215 217 217 217 217 219 220 |

| Chapter 20. Solid-State NMR Experiments | 227 |
|--|------------|
| 20.1 Cross-Polarization MAS Experiments | 227 |
| AP Bus Events in CP/MAS Experiments | 227 |
| Using a Waveform Generator in CP/MAS Experiments | 228 |
| 20.2 Sideband Suppression in MAS Experiments | 230 |
| 20.3 Rotor Synchronization | 233 |
| Measuring the Rotor Period Duration | 233 |
| Waiting for Triggers | |
| Rotor-Synchronized Experiments | 235 |
| 20.4 Multipulse Experiments | 235 |
| 20.5 Other Line-Narrowing Techniques | 236 |
| Chapter 21. (Micro)Imaging Experiments | 239 |
| Chapter 22. Role of Macros and Parameters | 241 |
| 22.1 Creating New Parameters in VNMR | 242 |
| 22.2 Using New Parameters in C | 246 |
| Numeric Parameters | 246 |
| String Parameters | 246 |
| 22.3 Adding New Parameters to the Display | 247 |
| 22.4 Doing It All by Macro | 248 |
| Macros for 1D Pulse Sequences | 248 |
| Macros for 2D Pulse Sequences | 251 |
| 22.5 Switching Between Similar Sequences | 252 |
| Chapter 23. Putting It All Together | 255 |
| 23.1 Starting a New Sequence | 255 |
| Programming by Modifying an Existing Pulse Sequence | 255 |
| Programming by the Top-Down Approach | 255 |
| 23.2 Testing a Sequence and Related Files | 256 |
| 23.3 Submitting a Pulse Sequence to the User Library | 257 |
| Chapter 24. Syntax Guidelines | 259 |
| 24.1 General C Syntax | |
| Comments | 259 |
| Indentation | 259 |
| Variables | 260 |
| 24.2 Outdated PSG Utilities | 260 |
| Device Addresses | 260 |
| Functions with Device Addresses | 261 |
| | - |
| Replacing power and pwrf Statements | 262 |
| C Constructs for Phase-Sensitive nD NMR | |
| | 262 263 |

| Using dps | 263 |
|--|-----|
| Chapter 25. Debugging a Pulse Sequence | 265 |
| 25.1 Debugging the Parameters | 266 |
| 25.2 Debugging the Software | 266 |
| 25.3 Debugging the Hardware | 267 |
| Index | 269 |

List of Figures

| Figure 1. Stages of pulse sequence execution | 16 |
|--|----|
| Figure 2. Compiling from the make utility | 19 |
| Figure 3. Diagram of the cpp program and include files | 20 |
| Figure 4. The lint check | 24 |
| Figure 5. Compiling and linking to form an executable program | 25 |
| Figure 6. Events associated with the psggen shell script | 30 |
| Figure 7. Timing diagram with transmitter, receiver, and transmitter phase | 42 |
| Figure 8. Switching mode for receiving an NMR signal | 43 |
| Figure 9. Maximum power directed into the probe | 43 |
| Figure 10. Timing diagram in "real life" | 44 |
| Figure 11. Phase shifting on UNITY plus systems | 46 |
| Figure 12. Phase shifting on systems prior to UNITYplus | 46 |
| Figure 13. Filter delay | 51 |
| Figure 14. Block diagram for a pulse programmer | 53 |
| Figure 15. Processes surrounding the go command | 55 |
| Figure 16. Acquisition control by Acqproc | 60 |
| Figure 17. Digital components of a spectrometer | 65 |
| Figure 18. Structure of the pulse programmer | 85 |
| Figure 19. UNITY plus waveform generator circuitry | 64 |
| Figure 20. UNITY programmable pulse modulator circuitry | 65 |
| Figure 21. Waveform generator board 10 | 66 |
| Figure 22. Detection of NMR signals | 10 |

List of Tables

| Table 1. Pulse programmer characteristics 87 |
|--|
| Table 2. Fast bit assignments, output boards and acquisition control boards |
| Table 3. Fast-bit assignments on pulse sequence control boards 90 |
| Table 4. Single- and double-precision timer word characteristics 92 |
| Table 5. Single- and double-precision timer word characteristics, output boards 92 |
| Table 6. Real-time math operators 96 |
| Table 7. Predefined AP bus delay constants. 142 |
| Table 8. PSG hardware flag and configuration variables 156 |
| Table 9. Waveform generator instruction words 169 |
| Table 10. Waveform generator gate control for pulse shapes 174 |
| Table 11. Comparison of waveform generator pattern words 196 |
| Table 12. Dealing with real and imaginary signal components 213 |
| Table 13. VNMR acquisition parameters used for nD experiments 215 |
| Table 14. Variables used in nD pulse sequences 215 |
| Table 15. VNMR parameter types and propertie 242 |
| Table 16. Predefined, indirect parameter limits 244 |
| Table 17. RF channel naming convention 261 |
| Table 18. Equivalent PSG functions with and without device address 261 |

Disclaimer

The information in this manual is intended to assist users in topics beyond the normal NMR spectrometer system hardware and software support provided by Varian. Some of this information is provided on an as-is basis, and Varian support and service personnel may be unable to answer questions related to information given in this manual.

Foreword

The main purpose of this manual is to complement the contents of the manual *VNMR User Programming* by providing comprehensive information on the *entire* subject of pulse sequence programming and execution as well as spectrometer control, trying to open up "black boxes" or "white spots" that might exist in a user's vision of a UNITY*plus*, UNITY, or VXR-S spectrometer. No attempt was made to describe all pulse sequence functions in detail—many of them will not even be mentioned here, but are covered by the manual *VNMR User Programming*.

Not all of the sections are equally important for the pulse sequence programmer. Most certainly it is possible to write correct pulse sequences without having access to the material presented here. It is the author's firm belief, however, that somebody who has more background knowledge on pulse sequence and spectrometer internals (as presented in this manual) will be more efficient and will less likely get stuck with the question of "How do I implement this on my spectrometer?".

The chapters on seqgen and psggen (chapters 2 and 3) are mainly written for people that want to learn about the programming internals of a pulse sequence, or for users that want to implement changes or new features in the pulse sequence overhead; sections that deal with Acode can be skipped by those who are not interested.

The material was originally written as documentation for the pulse sequence programming part of the Varian courses for users of Varian VXR-S, UNITY, and UNITY*plus* NMR spectrometers using VNMR software. Any input that would help improve this manual is very welcome.

Note that the specifics of the UNITY *INOVA* are *not* discussed—this will be covered in a future version of the manual. However, apart from some differences in the digital hardware, the basic concepts shown in this material are still valid also on the UNITY *INOVA* (see also the note below on software compatibility).

The present version of this manual specifically refers to VNMR version 5.1 and is partly inconsistent with earlier and later VNMR releases (in particular, the Acode structure for the UNITY *INOVA* is different, and the apdecode software mentioned in this manual will not work on the VNMR 5.2 software release). Most Acode printouts have been generated under VNMR 4.3 and are slightly different from Acode printouts obtained with apdecode under VNMR 5.1.

Rolf Kyburz Varian International AG, Zug/Switzerland February 1996

Acknowledgments

The author would like to thank his colleagues for their thorough proofreading, suggestions, and helpful discussions, as well as the Varian service personnel for withstanding the author's endless questions on the subject of this manual. In particular, Andy Myles (service supervisor in Darmstadt) has been extremely helpful in providing technical information.

Conventions in This Manual

The following notational conventions are used throughout all VNMR manuals:

• Typewriter-like characters are used to represent UNIX or VNMR commands, parameters, directories, and file names in the text of the manual; for example:

The shutdown command is in the /etc directory.

• Typewriter-like characters are also used for text displayed on the screen, including the text echoed on the screen as you enter commands; for example:

Self test completed successfully.

• Italicized typewriter-like characters are used for text displayed on the screen that is not the same every time; for example,

Abort at *some_address*

means the value of "*some_address*" depends upon when the abort command is made—what you might see on the screen is a message like this:

Abort at 47F82

• Special characters are used for keys on the keyboard and menu buttons on the screen; for example,

Press the Return key or select the Display button.

- Text shown between angled brackets in a syntax entry is optional. For example, if the syntax is seqgen s2pul<.c>, entering the ".c" suffix is optional, and typing seqgen s2pul.c or seqgen s2pul is functionally the same.
- Lines of text containing command syntax, examples of statements, source code, and similar material are often too long to fit the width of the page. To show that a line of text had to be broken to fit into the manual, the line is cut at a convenient point (such as at a comma near the right edge of the column), a backslash ("\") is inserted at the cut, and the line is continued as the next line of text. This notation will be familiar to C programmers. Note that the backslash is not part of the line and, except for C source code, should not be typed when entering the line.
- Because pressing the Return key is required at the end of almost every command or line of text you type on the keyboard, use of the Return key will be mentioned only in cases where it is *not* used. This convention avoids repeating the instruction "press the Return key" throughout most of this manual.

Chapter 1. **Overview**

This manual discusses the underlying functionality of the generation and execution of pulse sequences under VNMR. The idea is to provide the advanced pulse sequence programmer with a comprehensive insight into all aspects of pulse sequence generation. With better understanding, avoiding mistakes and writing more efficient pulse sequences should be possible, while at the same time moving towards more complex experiments.

1.1 Pulse Sequence Execution

Before going into details, let us review the various stages of pulse sequence generation and execution (see Figure 1):

- At the center of pulse sequence execution is the command go, which executes a compiled pulse sequence from seqlib, based on parameters from the current experiment and system configuration parameters (/vnmr/conpar). If external phase tables are used, these are retrieved from tablib, and pulse shapes or other waveform generator patterns are retrieved from shapelib. go stores its output (Acode, the data that result from executing the file in seqlib) in the directory /vnmr/acqqueue.
- Before go can be called, the parameters in the current experiment must be prepared. This is usually done by calling a macro (from maclib) with the name of the pulse sequence. This macro creates the necessary new parameters, usually from a library (parlib) with pulse sequence specific parameters, and it usually displays a text file with information about the pulse sequence.
- Before go can execute the pulse sequence file in seqlib, the pulse sequence file in psglib must be compiled using the command seqgen. This requires various files and libraries in a directory psg. seqgen itself internally is a complicated, multi-stage process, while for the user it looks like a simple command.
- The libraries of precompiled files required for seggen have to be created first. This is done using a command psggen, which again involves numerous source and header files in the directory psg. This step is normally not necessary for the user, because a complete set of psg libraries is part of the standard distribution tape.
- Once the pulse sequence file in seqlib has been successfully executed, the program Acqproc takes the Acode from /vnmr/acqqueue, transmits it via SCSI bus to the acquisition CPU (via the HAL, the Host-to-Acquisition Link), monitors the acquisition, stores the experimental data in the appropriate FID file within the experiment from which the acquisition has been started, and notifies VNMR at specific points in the acquisition (error condition, block size, and number of transients or entire experiment completed).
- In the acquisition CPU, the Acode is executed (interpreted), which results in information being fed into the pulse programmer or directly into the spectrometer that interacts with the magnet and the probe.

• The pulse programmer board controls almost all of the spectrometer (rf channels, waveform generators, frequency synthesis, receiver chain), including the analog-to-digital converter (ADC), which produces a new FID that is transferred to the sum-to-memory board (STM), which adds the current FID to the previous data in the HAL memory, from where Acqproc retrieves the final (or preliminary) FID onto the Sun memory (and disk).




Figure 1. Stages of pulse sequence execution

1.2 What to Expect in This Manual

In the following chapters, every single step in the above process is discussed in detail, starting with pulse sequence generation and the psggen step, continuing with go, then looking into the things that happen in the acquisition CPU and the spectrometer, while continuously adding knowledge about pulse sequence statements (sometimes also called elements or functions in VNMR documentation).

In the end, we will come back to the "origin" and discuss how the front end for the user (the user interface of a pulse sequence) is constructed.

Chapter 2. Sequence Generation: seqgen

The central step in the creation of a new pulse sequence is the compilation of the pulse sequence source code, which results in an executable file (with the name of the pulse sequence in seqlib, the library for compiled pulse sequences).

2.1 Modifying the File for dps


The VNMR command seqgen calls a UNIX shell script /vnmr/bin/seqgen with the same argument (the name of the pulse sequence with or without .c extension). seqgen can also be called from a UNIX environment:

seqgen sequencename<.c>

The UNIX file /vnmr/bin/seqgen is a lengthy shell script that first checks for all the files that are necessary for a successful compilation, then it checks for the presence of a program /vnmr/bin/dps_ps_gen. If this program is found, it is executed:

dps_ps_gen seqgencename.c \$HOME/vnmrsys/psglib 2

This results in a new file sequencenamedps.c (the name of the pulse sequence source file with dps inserted between the extension and the body of the name:



The new dps-modified file contains both the original text of the pulse sequence in the first part, plus a new, second part for the dps command, which graphically displays the pulse sequence (enclosed in #ifndef LINT... #endif, such that lint (the syntax checker) does not "look" at this part of the file):

#include <standard.h>

```
pulsesequence()
{
    /* equilibrium period */
    status(A);
    hsdelay(d1);
    /* tau delay */
    status(B);
    pulse(p1, zero);
    hsdelay(d2);
```

```
* observe period */
  status(C);
  pulse(pw,oph);
}
#ifndef LINT
extern FILE
             *dpsdata;
x_pulsesequence()
  fprintf(dpsdata, "\n status A %f ",pw);
  fprintf(dpsdata, " D_d1 %.9f ", d1);
  fprintf(dpsdata, "\n status B %f ",pw);
   fprintf(dpsdata, " P_p1 %.9f ", p1);
   fprintf(dpsdata, " D_d2 %.9f ", d2);
  fprintf(dpsdata, "\n status C %f ",pw);
  fprintf(dpsdata, " P_pw %.9f ", pw);
}
```

#endif

The second part of the dps-modified file contains all functions of the first part, prepended with "x_" (e.g., x_pulsesequence() instead of pulsesequence()). Comment is removed, and internal functions (pulse, delay, etc.) are replaced by functions that print some symbolic text into a temporary text file (curexp/dpsdata, referenced here as dpsdata). The VNMR dps command then decodes this text file. The internal structure and functionality of this text file will not be discussed any further here—this is beyond the scope of this manual.

If there was any C errors in the pulse sequence, such errors very likely would be propagated into the second part, leading to an increased number of error messages (despite the fact that the second part is "hidden from lint" through #ifndef LINT and #endif). For this reason, seqgen recompiles the pulse sequence *without* the dps modifications upon detection of errors with the compilation of the dps-modified file.

2.2 Running the make Program

All subsequent steps are performed using make, a UNIX program that facilitates the compilation of complex objects and also checks on the modification date of the files involved. make avoids an unnecessary compilation if an up-to-date executable already exists (see Chapter 3, "Object Library Generation: psggen," on page 29). make then calls cpp (the C preprocessor), lint (syntax checking), and cc (the C compiler and linker). If, at any of these stages, error messages are produced, the messages are stored in a separate file errmsg, which is displayed by seggen at the end of the operation.

The make program takes its instructions from /vnmr/acqbin/seqgenmake, a special "make file." If at the end of the entire process an executable is obtained, the dps part of the name is removed by the seqgen script and the executable pulse sequence is stored in seqlib. Figure 2 is a diagram of the make process.




Figure 2. Compiling from the make utility

2.3 Calling the C Preprocessor

Compiling a C program always involves first a call to the C preprocessor, which has the task of removing comments and resolving all the preprocessor statements in the C code. The C preprocessor cpp (or cc -P) is normally invoked implicitly as part of the cc command, and its output is directly fed into the next stage of the compilation. For the lint syntax check, however, the C preprocessor is called explicitly:

cc -P -DLINT -I\$HOME/vnmrsys/psg -I/vnmr/psg sequencename.c

This generates a file sequencename.i, which is then subjected to the lint check. The option -DLINT defines the preprocessor variable LINT (i.e., we are explicitly *including* source sections between #ifdef LINT and #endif, and explicitly *excluding* parts between #ifndef LINT and #endif). The -I options specify the path(s) to VNMR-specific include files.

The major task for the C preprocessor is the resolution of include statements, macros, and constants. In a first step, the include statements have to be identified and replaced with the contents of the file they indicate. Every C program has at least one include statement:

#include <stdio.h>

This statement incorporates the file /usr/include/stdio.h. VNMR pulse sequences usually also have a single include statement:

#include <standard.h>

This statement incorporates the file /vnmr/psg/standard.h (or the file \$HOME/vnmrsys/psg/standard.h, both paths are specified in the call to cc -P), and standard.h is mostly a collection of further include statements.

The program cpp recursively resolves all these lines, until no other include lines are found. Figure 3 is a diagram of the cpp program incorporating the include files.

The first include file from /vnmr/psg/standard.h is stdio.h (as for any C program). All other include files (oopc.h, acqparms.h, rfconst.h, aptable.h, power.h, macros.h, and apdelay.h) are from the psg directory,




Figure 3. Diagram of the cpp program and include files

either in /vnmr or in \$HOME/vnmrsys. The file oopc.h calls for another include file named /usr/kvm/sys/sys/types.h.

The purpose of all the include files is two-fold:

- To declare external definitions for functions and variables that are defined (or initialized) in another source module (which is precompiled and will simply be linked from an object library).
- To resolve constants and macros. These are pure preprocessor functions and must be resolved before the compiler can be called. These further modify the pulse sequence text for the C compiler (which doesn't understand constants or macros).

Most of the include files for pulse sequences deal with external declarations of functions and variables, or with the definition of constants. Macros are defined in the include file psg/macros.h (macros are just a more sophisticated type of constant). Because at least some of these include files contain information that can be very helpful to pulse sequence programmers, they are briefly discussed here:

- stdio.h defines the I/O functions for C. Nearly all C programs include stdio.h.
- psg/oopc.h defines some message structures and all the constants required in connection with general (object-oriented) routines like G_Pulse and G_Delay, which are described in the manual *VNMR User Programming*.
- psg/acqparms.h defines all (externally declared and initialized) variables related to acquisition parameters (pw, p1, d1, d2, d3, d4, etc.), externally defined code addresses (ct, oph, zero, one, two, three, v1, v2, v3, etc.) and objects (mostly structures), and the indices for the various rf channels.
- psg/rfconst.h defines constants such as TRUE and FALSE, fast bit offsets (hexadecimal constants) for the fast bits in VXR/UNIT- style output boards and acquisition control boards, device addresses (TODEV, DODEV, etc.), status field constants A (=0) to Z (=25), and a few others.

- psg/aptable.h defines constants and structures used in connection with tables, including the (externally defined) table addresses t1 to t60.
- psg/power.h defines a few constants used in connection with the G_Power general power statement, described in the manual *VNMR User Programming*.
- psg/macros.h defines constants used specifically with macros as well as all the macros used in pulse sequences. Most of the statements used in pulse sequences are actually macros that are converted into very few general statements such as G_Pulse and G_Delay.
- psg/apdelay.h defines constants (macro functions in reality) that specify the implicit delays involved with AP bus statements. Whenever possible, these delay constants (e.g., POWER_DELAY, SAPS_DELAY, OFFSET_DELAY) should be used instead of direct numeric values (e.g., 14.95e-6) in pulse sequences.

The reason why there are many different include files is that specific parts of the "overhead" defined in these files are used by various other source modules as well, and the "packaging" is defined by the needs of these other source modules. If there was only one source module (the pulse sequence itself), there would probably be a single include ("header") file.

Macros are a special kind of constant that allow for arguments. Let's take the definition of decpulse (from psg/macros.h) as an example:

```
#define decpulse(DECch,phaseptr) \
    G_Pulse(PULSE_DEVICE, DECch, \
    PULSE_WIDTH, length, \
    PULSE_PHASE, phaseptr, \
    PULSE_PRE_ROFF, 0.0, \
    PULSE_POST_ROFF, 0.0, \
    0)
```

The two arguments to the decpulse statement become arguments 4 and 6 of the G_Pulse statement; the other arguments in the general statement are determined by the macro name.

The C preprocessor dramaticallychanges the text of the original pulse sequence file, primarily by replacing include lines by the text of the include files (and recursively substituting text for include lines contained therein). The text of the pulse sequence function itself is different after the substitution of macros and constants (apart from the elimination of comment lines). Look at the simple pulsesequence function of a S2PUL (standard two-pulse) pulse sequence:

```
pulsesequence()
{
    /* equilibrium period */
    status(A);
    hsdelay(d1);
    /* tau delay */
    status(B);
    pulse(p1, zero);
    hsdelay(d2);
    /* observe period */
    status(C);
    pulse(pw,oph);
}
```

After the constants A, B, and C have been substituted, and the pulse macros have been replaced by the corresponding G_Pulse calls, the function consists of the following lines (the blank lines that result from the elimination of comments have been removed):

```
pulsesequence()
{
    status(0);
    hsdelay(d1);
    status(1);
    G_Pulse(1, p1, 5, zero, 0);
    hsdelay(d2);
    status(2);
    G_Pulse(1, pw, 5, oph, 0);
}
```

Note that the constants within the arguments to G_Pulse have been replaced by their numeric values: the result of the C preprocessor is not meant to be read by humans but rather adjusted for the C compiler. What is left are pure C control structures, function names, variables, and (numeric or string) constants. Remember that C knows symbolic constants only through preprocessor definitions.

The C preprocessor includes another useful feature that was briefly mentioned earlier in this chapter: conditional compilation. Based on the definition of a (preprocessor) variable, certain sections of a source file can be specifically excluded or included. This feature is used extensively in the "overhead" to pulse sequences (see Chapter 3, "Object Library Generation: psggen," on page 29). In pulse sequences, this feature is only used to hide things from lint, which would otherwise lead to an error message:

```
#ifndef LINT
static char SCCSid[] = "@(#)cosy.c 3.1...";
#endif
```

The second line would normally cause an error message from lint, because the array variable SCCSid is defined but not used (see below). Such preprocessor flags can either be defined in the file itself (#define LINT; note that this definition has a name, but no value) or through an argument (-Dflagname) with the preprocessor or the C compiler call, for example:

cc -P -DLINT -I\$HOME/vnmrsys/psg -I/vnmr/psg sequencename.c

2.4 Checking Syntax with lint

One of the drawbacks of the C programming language is the fact that its compiler has little error detection. The compiler only finds some very coarse syntactical errors, such as a mismatch in the parentheses or braces, and eventually also missing semicolons. Anything more subtle—such as a variable type mismatch, the use of uninitialized (or even undefined) variables, a mismatch in the number of arguments to a function, etc.—remains undetected by the compiler itself. With all the syntactical traps and pitfalls (e.g., using & instead of &&, or | instead of | | makes perfect sense in C), this weakness could make it a very difficult language for occasional pulse sequence programmers. Fortunately, there is lint, the C syntax checker.

The program lint can be called on any C source program. lint usually makes lots of output covering two aspects: definite errors and possible errors (the latter includes warnings about the possibility of exceeding the limits on arrays, if using an array). Error detection is not a problem with lint—it's rather that often the output is

overwhelming. For that reason, in checking pulse sequences, lint is called with some of its checking features turned off (e.g., options -n, -u, -v, -z; see man lint).

Like cc, lint first calls the C preprocessor cpp to resolve preprocessor statements (definitions, macros, and conditional sections). In the case of pulse sequences, this step would have the undesirable effect that the original macro syntax would be lost (and with it possibly errors in the arguments to macros) because macros would be resolved to object-oriented function calls, where argument type-checking is nearly impossible (at the lint level). For this reason, most of the file psg/macros.h is divided into two sections—one for the normal compilation (#ifndef LINT) and the other for the lint pass (#else...#endif). The lint version is generated in an explicit call to the C preprocessor:

cc -P -DLINT -I/vnmr/psg -I\$HOME/vnmrsys/psg sequencename.c

The lint part of psg/macros.h does not resolve the macros into their C function equivalents. It simply changes the macro name from lower to upper case (leaving the arguments intact), while constants are still replaced. The simplified S2PUL pulse sequence above would then look as follows in the corresponding s2pul.i file:

```
pulsesequence()
{
   status(0);
   hsdelay(d1);
   status(1);
   PULSE(p1, zero);
   hsdelay(d2);
   status(2);
   PULSE(pw, oph);
}
```

Subsequently, lint checks the syntax in this version of the file, and it "knows" the argument number and types for all the uppercase "pseudo functions."

How does lint know about syntax rules and, in particular, the syntax and types for externally defined functions and variables? The secret is a special library that is involved automatically with every lint call. For simple, standard C programs, this library is the file /usr/lib/lint/llib-lc.ln. This file is used automatically unless a specific rules file is specified with the lint call. For VNMR, a dedicated library file psg/llib-lpsg.ln has been created with a special lint call, using the -C option:

```
lint -a -n -z -DLINT -Cpsg lintfile.c
```

The file lintfile.c contains all the pulse sequence functions that are checked and all the externally defined variables and addresses. Different from the real functions, the functions in lintfile.c have a null body ({}). On the other hand, these functions have the same (complete) argument type declarations as their real equivalents.

As an example, let's take the lint definition of the rgpulse statement:

```
RGPULSE(pulsewidth, phaseptr, rx1, rx2)
double pulsewidth, rx1, rx2; codeint phasptr; { }
```

For the lint check, lint is called with all the preprocessor options (even though it actually operates on the preprocessed file sequencenamedps.i):

The output of the lint check is fed into the error message file errmsg (see Figure 4).




Figure 4. The lint check

Of course, VNMR is delivered with psg/llib-lpsg.ln, a complete lint library file. There seldom should be a need for a user to re-create this library file—perhaps only if somebody creates a new external function and wants to have it checked by lint. Note that unlike other files involved in pulse sequence generation, the lint library file *must* reside in /vnmr/psg unless the make file seqgenmake is modified for a local file.

Error messages and warnings from lint are added to the error messages file errmsg. The compilation does not stop upon error messages from the C preprocessor or the lint step; such messages are taken as warnings only and in the end lead to a message from seqgen: "Pulse sequence did compile, but may not function properly..."

2.5 Compiling and Linking

Once the syntax of the pulse sequence is checked, make proceeds to the compilation, which is done in two steps: first, the pulse sequence module is compiled and, second, the resulting object module is linked with other precompiled modules to form an executable program (see Figure 5).

The actual pulse sequence compilation is done by make with the command cc -c, which calls the C preprocessor implicitly. This time, of course, the -DLINT option is not specified, such that macros are resolved into the proper C functions. Also, sections that were skipped for the lint pass because they would lead to error messages, such as the definition of the SCCSid string (see above), are now included.




Figure 5. Compiling and linking to form an executable program

The compilation is performed with the following command string:

```
cc -O -I$HOME/vnmrsys/psg -I/vnmr/psg -c sequencenamedps.c
```

where -0 is the code optimizer option. Note that -0 is the same as -02, which means level 2 optimization (see man cc for more information). The -I options are for the preprocessor, which needs to be told where to find the VNMR-specific include files. With the -c option, the compiler does not produce an executable program, but rather an object module file sequencenamedps.o. In case of compilation errors, this file is missing, make aborts, and seqgen reports all errors that were detected up to that point¹. The resulting object module sequencenamedps.o is not executable, because it is the bare compiled pulse sequence function, with all the references to external variables and functions but without the actual code for them.

In order to arrive at an executable file, make calls cc again (in reality the call may differ slightly, depending on the software and hardware configuration; also, for various reasons, the syntax in the make file is somewhat different):

```
cc -O -Bstatic -s -L/vnmr/psg -o sequencenamedps \
    sequencenamedps.o -lpsglib -lparam -lm
```

The -0 option again specifies the code optimizer. All other options are not for cc itself, but for ld, the so-called link editor, which is called by cc:

The -s option causes ld to strip the so-called "symbol table and relocation bits" (used by debuggers) to save disk space. -Bstatic causes the externally referenced libraries to be linked (attached) statically (see below); the -L option permits specifying additional directories where ld searches for libraries; and the argument -o sequencenamedps specifies the name of the executable target file.

¹ The command /usr/ucb/cc (a shell script) works in two steps: the actual C compiler intermediately produces assembly language code (*.s), which is then converted to an object module using the as command (assembly language compiler). The first step can also be done explicitly using the -S option to the cc command.

The last three options (-lpsglib, -lparam, and -lm) specify the three libraries from where external modules are loaded: libpsglib.a and libparam.a (both normally in /vnmr/psg), and /usr/lib/libm.a. The latter contains the math libraries (extended floating point math functions), the former two contain all the VNMR-specific compiled external modules found in the directory /vnmr/psg. The link editor ld should normally be able to resolve all references to external functions and variables; otherwise, error messages result.

In the early releases of VNMR, the description above was the way the link loader worked. Meantime, pulse sequence overhead gained complexity, and each statically linked executable pulse sequence became 250 to 400 Kbytes (Sun-3 vs. Sun-4). With about 50 pulse sequences in the standard release, this takes about 10 to 20 Mbytes of disk space. Moreover, all of these sequences would have the *same* modules linked (i.e., a large portion of disk space is lost for 50 copies of almost the same software).

In this situation, it was very useful that Sun introduced the concept of dynamic binding and shared objects (libraries) with SunOS 4. Instead of binding (linking) all external modules at compile time, these objects are put into a special *shared* library: a library that is accessed and used by many object modules at run-time (when they are called)².

For pulse sequences, the run-time aspect of shared objects is not relevant, because the cases where two pulse sequences are called at the same time on a single system does not happen normally (the "call" in this case only means the seconds of execution time upon typing g_0). On the other hand, with run-time linking, every pulse sequence only occupies 16 or 24 Kbytes on the disk³, saving considerable disk space. By default, seqgen uses dynamic binding (run-time linking) when compiling pulse sequences.

How is dynamic run-time linking enabled and used in practice? To prepare for using run-time linking, a new file has to be created that combines the contents of the library files psg/libpsglib.a and psg/libparam.a in a suitable format. For VNMR, such files are included under the names psg/libpsglib.so.x.y and psg/libparam.so.x.y, where x and y indicate the major and minor revision number (in VNMR 4.3A the revision number is 5.0). The idea behind the revision numbers is that if several revisions of these files coexist in the same directory, only the last one (the one with the highest revision number) is used. The so stands for "shared objects."

During the compilation, static linking must be prevented. This outcome is achieved by omitting the -Bstatic option during the link loading (dynamic linking is the default). This way we get executables without the shared objects, and it turns out that for pulse sequences these files are only 16 or 24 Kbytes, again in multiples of 8 Kbytes (i.e., we save about 95% of the disk space that would be used with static binding). Invisibly for the user, at go time, the system now picks the executable from seqlib, links it with the psg/libpsglib.so.x.y library with the highest revision number, and executes the resulting code. The slow-down in go due to the run-time linking is negligible.

² For certain applications—not pulse sequences, though—this not only saves disk space, but also process/memory space because the shared objects (libraries) are only loaded into memory once and can the be accessed by several applications simultaneously.

³ The size of executables is always a multiple of 8 Kbytes, because programs are used (loaded into memory) in 8-Kbyte pages anyway.

One feature of run-time linking sometimes becomes apparent to the user: at link time (i.e., when executing the command), the revision dates of the current shared objects is compared with the date of the objects present at compile time (which is included in the executable file), and if there is a revision mismatch, this leads to error messages—it may even cause the program to crash. Normally, this should be no problem—but there are possible dangers when the shared objects are modified by the user, see Chapter 3, "Object Library Generation: psggen," on page 29.

One small point is left for seqgen: if the dps additions were included successfully, the pulse sequence name had the characters dps added all the way through the syntax check, compilation, and linking. In order for VNMR to find the executable pulse sequence, seqgen renames it in the end to the pulse sequence name without dps. Finally, all intermediate files are deleted, including the pulse sequence object file (sequencenamedps.o). This means that for every seqgen, the *entire* process is repeated, irrespective of whether the pulse sequence executable is already up-to-date or not.

Chapter 2. Sequence Generation: seqgen

Chapter 3. Object Library Generation: psggen

For creating the pulse sequence compilation libraries, over 80 C modules (mostly stored in /vnmr/psg) have to be compiled, using over 30 VNMR-specific include files. The compilation process is very similar to the compilation of a pulse sequence, except that no explicit C preprocessor call and lint syntax checking are involved.

3.1 How Are Object Libraries Generated?

The C preprocessor is called implicitly with the first compiler pass (cc -c). In the second C compiler pass (cc -o), no executables are created (the program, which is nothing but a single, big C program, is not complete because itt lacks the pulsesequence function), but instead, shared object libraries are generated.

Also different from the compilation of a pulse sequence, the intermediate object modules (*.o) are not discarded, but rather stored in object libraries (libpsglib.a for C modules supplied in /vnmr/psg, and libparam.a for C modules that are not distributed). This has a major advantage—upon recompilation, make only compiles C modules that have been updated because it compares the dates of the individual C modules and their corresponding object modules.

The process of creating the object libraries in the distributed software is so similar to the regeneration of these libraries upon user modifications that only the regeneration will be discussed in detail here.

3.2 Adding Changes to the Object Libraries

VNMR gives the user access to almost all of the pulse sequence-related source code and with that the ability to implement user modifications. In the past, this has also helped Varian to distribute bug fixes for the "pulse sequence overhead" as ASCII text —simple source code modifications for /vnmr/psg.

In accordance with the VNMR philosophy, such modifications are not done in /vnmr/psg directly, but rather in a local copy of that directory, which is generated with the setuserpsg shell script.

This script creates a local directory ~/vnmrsys/psg containing the two object libraries libpsglib.a and libparam.a, the library for run-time linking libpsglib.so.x.y(libparam.so.x.y is not modified by the user and is therefore still stored in /vnmr/psg), and the lint library file llib-lpsg.ln. The principle is that upon re-generation of the libraries, or upon compilation or execution of a pulse sequence, any file missing in the local library is taken from /vnmr/psg.

If a user ever needs to use modifications in this area, it most likely involves changes to a single source file. The user either can follow some given recipe (like to fix a bug in the pulse sequence overhead) or can go into /vnmr/psg and locate the file that needs modification (typically by using grep and an editor to locate the "critical spot"). The source file to be modified is then copied into the local psg directory and modified. The shell script command psggen is then called, for example:

```
cp /vnmr/psg/aptable.c ~/vnmrsys/psg
cd ~/vnmrsys/psg
vi aptable.c
psggen
```

As shown in Figure 6, the psggen shell script (in /vnmr/bin) invokes a rather complex series of events.




Figure 6. Events associated with the psggen shell script

First, psggen calls make with the file /vnmr/psg/makeuserpsg as the make file (make -fes makeuserpsg lib), and then almost all subsequent actions are done through make:

- For all header and source files (*.c, *.h, *.p) from /vnmr/psg that are *not* present in the local directory, a symbolic link is created. The make program expects to "see" a complete directory.
- All object files (*.o) are unpacked from the libpsglib.a library.
- The creation dates of all C source files and their corresponding object files are compared. Only those files in which the source file has a later modification date than the creation date of its object file are compiled. If you call psggen the first time in a newly generated local psg directory, all sources are compiled, because the original objects in libpsglib.a have been created at the Varian factory, whereas the source files carry a modification date from the time when they were stored in /vnmr/psg (i.e., for the first call of psggen all object files look "outdated"). Such a full compilation takes about five minutes on a SPARCstation 2 or equivalent. On subsequent calls, only the sources that have been modified since the last call to psggen are compiled. On a SPARCstation 2, psggen takes about one minute if only one module was updated.

- The make file (psg/makeuserpsg) contains all "object dependencies," which means that if a header file (*.h or *.p) is modified or updated, automatically all sources that include this file become part of the compilation.
- Depending on the above date comparisons, C modules are compiled (cc -c); the C preprocessor is called implicitly with that pass for each file.
- After the compilation, the object library (psg/libpsglib.a) is updated with the current object files.
- No executables are created in the second compiler (linker) pass, but instead, the shared object library is updated. For a local psg directory, this shared object library always has the name libpsglib.so.1.0.
- All object files are deleted from the directory. This seems similar to seqgen—but here, the objects are still available in archived form. The archive has the advantage that the directory listing can be kept small.
- In a last step, psggen removes all symbolic links (to source and header files in /vnmr/psg).

Now, pulse sequences can be compiled using the new libraries. Even if the external objects are not statically linked to the pulse sequence executable, pulse sequences should still be recompiled, to avoid internal inconsistencies.

The resulting pulse sequence generation libraries are active only for a single account, because they are in a local vnmrsys/psg directory. In principle, vnmr1 can make changes available to "the world" by copying the modified files into /vnmr/psg, although caution must be used with this operation, because errors can cause an inability to compile pulse sequences in general. You may have to reload VNMR to recover.

The main problem is that object archive files should not be copied using the cp command because this alters the date of the target file, which can lead to date inconsistencies, and subsequently make (seqgen!) refuses to work. Note also that make does not work if the system date is older than the modification date of source and object files (e.g., by mistake the system date is set backwards by several months or years). Another possible problem lies in the fact that the run-time linker takes the shared library with the highest revision date if several versions coexist; hence, you cannot copy libpsglib.so.l.0 from the local psg directory into /vnmr/psg. A procedure that should be safe is as follows:

1. Make a backup of /vnmr/psg:

```
cd /vnmr
mkdir psg.bk
cd psg
tar cf - * | (cd ../psg.bk; tar xvfBp -)
```

2. Transfer the local files:

```
cd; cd vnmrsys/psg
tar cf - * | (cd /vnmr/psg; tar xvfBp -)
cd ..
rm -r psg
cd /vnmr/psg
ls *.so.*
mv libpsglib.so.1.0 libpsglib.so.8.0
```

With the last command, the highest revision number is given to the updated copy (which was revision 1.0 in the local directory). The local directory can be deleted,

because all changes are now available from the system files. cp - p should also work instead of the tar pipe.

3.3 Adding a New Precompiled Object

Few users only need or want to make changes to the pulse sequence overhead; among those, even fewer come into a situation where they want to add a new file to the list of precompiled objects. This chapter was not written with the primary intention to teach how to do this, but rather to give some additional insight into the internal functionality of the psggen and seggen commands, and into pulse sequence-related software in general.

A situation where somebody *might* want to add a new file to the list of precompiled modules is the case of a frequently used, very long, custom-built include file. Such include files can lengthen the pulse sequence compilation time considerably; however, with a precompiled module, the same functions could be made available without slow-down in compilation speed.

First, we create a local psg directory and store the new module there. As an example, take the current include file /vnmr/psg/shape_pulse.c:

```
setuserpsg
cp /vnmr/psg/shape_pulse.c ~/vnmrsys/psg
cd ~/vnmrsys/psg
vi shape_pulse.c
```

We now need to change the C module in order for it to become a standalone module. Pulse sequence include files draw all their references to external definitions and the macro references from the standard include file standard.h. As a result, such files typically have no, or only specific, include lines. As a start we can insert all the include lines from /vnmr/psg/standard.h:

```
#include <stdio.h>
#include "oopc.h"
#include "acqparms.h"
#include "rfconst.h"
#include "aptable.h"
#include "power.h"
#include "macros.h"
#include "apdelay.h"
```

The include notation uses the double quotes (e.g., "oopc.h") to mean that local files (e.g., in psglib) are included. The angled bracket notation (e.g., <stdio.h>) causes cpp to only search the standard include libraries (e.g., ~/vnmrsys/psg, /vnmr/psg, /usr/include). If files are specified with absolute paths, the include notation doesn't matter.

Later, we may want to determine whether we need all these include lines. Most likely at least some of them can be deleted, but for the time being, too many include lines can't hurt. If include lines already exist in the new file, we insert them between the standard C include lines (such as #include <math.h>) and any VNMR include lines, avoiding duplication.

Our example file unnecessarily already contains a reference to stdio.h; therefore, we will not add that again. On the other hand, we may have to also include the additional definitions and references from standard.h:

For psggen to recognize the new source file, we need to modify the make file makeuserpsg. If psggen has not been called beforehand, we get the make file from the system psg directory:

```
cp /vnmr/psg/makeuserpsg .
vi makeuserpsg
```

We need to change the make file in two places: starting at around line 50 is an alphabetical list of source files. Here we insert a line for our new module at the appropriate place, with a backslash at the end of the line:

shape_pulse.c $\$

The other modification is in the section "Object dependencies," starting at about line 250, where individual compilation instructions and source code dependencies are listed alphabetically. Our additions are at approximately line 600:

Most of this is analogous to other entries; therefore, we don't really have to understand all the details of the make file language. A few basics, however, cannot be avoided:

- Constructs like \$(COMPILE.c) refer to macros defined elsewhere in the make file or alternatively in the arguments to the make command.
- The whole thing represents the rules and dependencies for the object file shape_pulse.o; all continuation lines *must* start with *tabs*, *not spaces*.

The lines with the header file names describe the object dependencies (i.e., the header file changes imply a recompilation of this particular module). As long as we are in the development stage for that particular module (and maybe don't know yet for sure what header files will finally be used), these lines aren't really required; only in the end we should make sure that all used include files are also listed in the object dependency list, to ensure proper updating in the case of header file changes. As a start, we could use two lines only:

```
shape_pulse.o : $(@:.o=.c)
    (umask 2; $(COMPILE.c) $(@:.o=.c) )
```

Our particular example includes math.h and makes calls to functions in the math library, like log10(). It is not necessary to specify \$(LIBS) in the last line (LIBS is defined as -lm, which would link /usr/lib/libm.a), because this flag is specified with the compilation of each pulse sequence.

Before we can call psggen, we need to make sure that a copy of the new source module exists in /vnmr/psg; otherwise, make results in error messages and refuses to work, then we can finally call psggen:

```
cp shape_pulse.c /vnmr/psg psggen
```

For the first time, we receive an error message, because make can't find a copy of shape_pulse.o in libpsglib.a. This message is harmless—the file is compiled anyway— and after psggen, libpsglib.a contains this object module as well.

You can check this with the entry:

```
ar t libpsglib.a | more
```

You can double-check this by recompiling the new module. To do this, first simulate a change in that module by altering its modification date:

```
touch shape_pulse.c psggen
```

You are now ready to debug the new source module (in the case of a new program) or to try and reduce the number of include files, in order to minimize the object dependencies. In the end, make sure the make file (makeuserpsg) contains object dependencies for all include files that are used. For complete debugging, we have to compile pulse sequences that use the new function. If previously the new source module was an include file, we need to delete that include line from the pulse sequence, because the file is not required any longer; it is now even incompatible as include file, because of its own include lines (which would lead to multiple error messages because of duplicate definitions).

The fact that the function now no longer is compiled together with the pulse sequence function leaves one open question: How does the pulse sequence (better: the compiler and linker) know that there is now a new function in some other object module? The answer is that external functions are found and incorporated properly as long *as they are of the default type* (int) *and do not return a value.* The vast majority of the functions used in pulse sequences are of the default type and do not return values— there is no need to declare them in the pulse sequence module. The same is true for the function that we just converted into a precompiled module.

In the standard include file, standard.h, we find a definition for the function getval that returns a double:

```
extern double getval();
```

This function is of the non-standard type (and returns a value), and hence it must be declared as an external function. Any module in psg that uses the function getval must also include the above line. The file standard.h also defines two more functions that are defined as type void:

```
extern void setprgmode(), prg_dec_off();
```

There is a potential disadvantage in making a file precompiled: lint no longer "sees" it, and therefore can no longer check the syntax in calls to functions that are now external. To compensate for this, you can add to the lint syntax file those functions in the new module that can be called in a pulse sequence. In our particular case, we would add the function definition (without contents, of course) to the file lintfile.c.

The lines to be added would be as follows in our case:

```
shape_pulse(shape,pws,phs,pwrtbl,phstbl,spwr,npulses,rx1,rx2)
char shape[MAXSTR];
codeint phs, pwrtbl, phstbl;
int npulses;
double pws, spwr, rx1, rx2;
{}
```

This "dummy" function should be added to the section with plain functions (starting at about line 350) of the file psg/lintfile.c; then the lint syntax file needs to be generated and copied into /vnmr/psg, see Section 2.4, "Checking Syntax with lint," on page 22:

```
cp /vnmr/psg/lintfile.c .
vi lintfile.c
lint -a -n -z -DLINT -Cpsg lintfile.c
cp llib-lpsg.ln /vnmr/psg
```

The last command allows seqgen to use the new lint library file; alternatively, we could modify the seqgen make file (seqgenmake, usually stored in

/vnmr/acqbin) and/or the seqgen shell script, to work with the local lint library. For that, we should not change the make file in /vnmr/acqbin, but rather create a local copy in ~/vnmrsys/psg (i.e., only modify the make file for the account that has a local lint library). This alternative is not described here. Chapter 3. Object Library Generation: psggen

Chapter 4. Time Events

In this chapter, we deduce more complex pulse sequence statements from simple delays, just to see what is their function. The syntax shown in this chapter is *not* the one used in the pulse sequence overhead; the idea here is not to explain how pulse sequence functions are defined exactly, but merely to understand their functionality in order to correlate the arguments of such a statement to what is happening in the spectrometer. At the same time, this is an opportunity to become familiar with the basic C pulse sequence syntax. Later we shall get a more detailed picture of what is actually happening in the acquisition CPU and the pulse programmer (see Chapter 8, "Acquisition CPU and Acode," on page 69, and subsequent chapters).

4.1 How Do Delays Work?

A delay is nothing but an exactly specified waiting time for the pulse programmer. In other words, a delay is a time during which the pulse programmer does nothing but hold the current status (we'll see later exactly what that means). Although even with a simple delay a fair amount of "internal software" is involved, we can think of a delay as a principal pulse sequence element. In fact, a delay or time event is one of the few very basic actions a pulse programmer can do.

Simple Delays

Even a function as simple as a delay has some logical decisions built in; namely, the property that if a length of zero is specified or if the specified duration is too short to be executed by the pulse programmer (see Chapter 9, "Pulse Programmers," on page 85), the entire pulse sequence statement is skipped. The idea is to avoid unnecessary dead times in a pulse sequence, as well as statements, that could lead to problems at execution time. A possible coding for the statement delay could be as follows:

```
delay(length)
double length;
{
    if (length >= MINDELAY) time_event(length);
    if ((length < 0.0) && (ix == 1))
        printf("A negative delay has been specified.\n");
}</pre>
```

where MINDELAY is a constant that defines the duration of the smallest executable delay (0.2 microseconds); and ix is the index of the current FID (starting at 1, up to arraydim). Most warnings are issued for the first FID only, in order to avoid flooding the VNMR text window with error messages in arrayed and multidimensional experiments.

The delay statement actually is a macro in /vnmr/psg that is resolved to the G_Delay statement, and the function time_event doesn't exist as such (instead, a statement delayer with a slightly different definition is used internally).

Delays With Homospoil Pulse

If we want to write or change the homospoil delay statement hsdelay, we first need a definition of what that statement is supposed to do. It is always better to first have an accurate definition before starting to code a program. Sometimes it is necessary to start with a simple definition that can be expanded during the implementation:

- The statement hsdelay should have a single argument, specifying the total duration of the homospoil delay: hsdelay(length).
- If the specified length is zero, the entire pulse sequence statement is skipped.
- If the homospoil flag hs is 'n' in the current status field, a simple delay should be performed. If hs does not cover the current status field, the last specified field of hs is taken instead (propagation of the last specified status field).
- Otherwise, execute a delay of length hst with the homospoil (Z1 gradient) turned on, followed by a normal delay, such that the total duration is the specified length.
- If the specified length is less than hst, issue a warning message and perform a simple delay instead.


A possible implementation of the hsdelay statement could be as follows:

```
hsdelay(length)
double length;
{
  int hspos = statusindex;
  int hslen = strlen(hs);
  if (hspos >= hslen) hspos = hslen - 1;
  if (length >= MINDELAY)
  {
    if (hs[hspos] == 'y')
    {
      if (length >= hst)
      ł
        HSgate(homospoil_bit,TRUE);
        delay(hst);
        HSgate(homospoil_bit,FALSE);
        delay(length - hst);
      }
      else
      {
        if (ix == 1)
          printf("delay shorter than hst - no homospoil \
          performed.");
        delay(length);
      }
    }
    else
      delay(length);
  }
  if ((length < 0.0) && (ix == 1))
    printf("A negative delay has been specified.\n");
}
```

where statusindex is a defined variable set by the status statement (to the value of its argument); hs and hst are standard parameters and do not have to be declared (the same as all parameters used in the s2pul pulse sequence). TRUE (1) and FALSE (0) are constants defined in /vnmr/psg (they are *not* part of the standard C definition).

The handling of the hs flag conforms to the general internal flag handling in VNMR, in that the last flag field is propagated throughout the rest of the string (e.g., hs='ny' is the same thing as hs='nyyyyyyy...'; flags can be up to 255 characters long). statusindex is a variable used and defined in /vnmr/psg, and HSgate is an undocumented internal statement¹ in /vnmr/psg. In principle, both can be used in pulse sequences.

The relatively complicated hsdelay statement above (for the "normal" case) leads to a fairly simple timing diagram:



In words, the homospoil (Z1 gradient) is turned on, a delay hst is executed, then the homospoil is turned off, and a second delay length-hst is executed. Although the function hsdelay has only one argument, it actually has two software-controlled variables: the specified length, and the hst duration of the homospoil pulse, which is a standard VNMR acquisition parameter. The strength (really a third parameter to hsdelay) of the homospoil pulse is adjustable in hardware only.

The hsdelay statement is a type of "automatic" or "high-level" statement in that users do not have to know how to turn on or off the homospoil gradient, or even what the standard parameter for the length of the homospoil pulse is, in order to use this statement in a pulse sequence. On the other hand, users often want *more control* in a pulse sequence, or perhaps for a specific reason they want to use a different parameter for the duration of the homospoil pulse. This is no problem, because the "low-level" statements can be taken from the above statement and used in a pulse sequence:

```
HSgate(homospoil_bit,TRUE);
delay(hst);
HSgate(homospoil_bit,FALSE);
delay(length - hst);
```

Other Delays

There are a few additional statements for delays: idelay for interactive parameter (length) adjustment during acquisition, incdelay and vdelay are delays where the actual length is calculated by the acquisition computer from a time base and a real-time counter. These other delays are no different than a simple delay in what they do to spins; they will not be discussed any further here.

¹ A price may be required to pay for the additional flexibility obtained in using undocumented lowerlevel statements like HSgate: The definition for such statements could change with a future VNMR release or the statement as such may cease to exist altogether, and such changes may not be documented in the standard VNMR manuals. It is up to the user to verify in /vnmr/psg that this statement still exists in its current form.

4.2 How Do Pulses Work?

Pulses on the Observe Channel

Several statements coexist for performing pulses on the observe channel; all of them are actually macros and are ultimately translated into calls to a generic function G_Pulse (see the manual *VNMR User Programming*). The most general of these statements is genpulse (not supported by dps at the time of this writing):

genpulse(length,phase,rx1,rx2,device);

where length is the duration of the pulse (see below) in seconds, phase is the pulse phase and is a reference to either a real-time variable (see below) or a phase table, rx1 and rx2 are "receiver gating times" (again see below), and finally device is the rf channel on which the pulse is to be performed (OBSch in the case of a pulse on the observe transmitter channel).² Note that within a pulse sequence all time events are in seconds, even if the parameter (of type "pulse") is in microseconds.

The statement genpulse is the most versatile, yet the most complex, of the statements for a pulse on a single channel at the level of the pulse sequence. It is not recommended except where the rf channel for a pulse should be kept under parameter control.

The statements normally used for a pulse on the observe transmitter channel are rgpulse, pulse, and obspulse. The function rgpulse has the following syntax:

rgpulse(length,phase,rx1,rx2);

This statement is equivalent to genpulse(length, phase, rx1, rx2, OBSch), i.e., the device name (number) is implicit. The next simpler statement is pulse:

pulse(length,phase);

The pulse statement is equivalent to rgpulse(length, phase, rof1, rof2) or genpulse(length, phase, rof1, rof2, OBSch), i.e., the device name (number) and the two receiver gating times are implicit; for the receiver gating times, parameters rof1 and rof2 are used. The simplest statement is obspulse:

obspulse();

This is equivalent to pulse(pw,oph), rgpulse(pw,oph,rof1,rof2), or genpulse(pw,oph,rof1,rof2,OBSch), i.e., in addition to the implicit parameters in the pulse statement, the pulse length pw and the pulse phase oph is implied.

All these statements are more or less "automatic" variants of calling a "pulse" statement in a pulse sequence. What actually is happening in the rf during the execution of such a pulse? Why do we need the "receiver gating times," what are their implicit

² Traditionally (and exclusively in VNMR releases up to 4.3), devices and channels were treated equivalently. They were named TODEV (observe transmitter channel), DODEV (decoupler), DO2DEV (second decoupler), and DO3DEV (third decoupler). These addresses can basically still be used in the current version of VNMR to specify any one of the rf channels. However, these device names hardcode the channels so the channels do not reflect the rf channel assignments made with the rfchannel parameter. If the rf channel assignment is to be made dynamically (as per rfchannel), the new device names OBSch (observe channel), DECch (decoupler channel), DEC2ch, and DEC3ch must be used. For compatibility, it is recommended to translate all occurrences of TODEV-type constants to OBSch and its equivalents, even if the rf channel argument.

and desirable values? To answer these questions, let us first have a look at how the rgpulse statement could be coded in C:

```
rgpulse(length,phase,rx1,rx2)
double length,rx1,rx2;
codeint phase;
{
   int rcvrflag = rcvroff_flag;
   if (length >= MINDELAY)
   {
      rcvroff();
      txphase(phase);
      delay(rx1);
      xmtron();
      delay(length);
      xmtroff();
      delay(rx2);
      if (!rcvrflag)
         rcvron();
   }
}
```

Like with the delay or hsdelay statements, unnecessary dead times in the pulse sequence should be avoided: if the specified length is zero, the *entire* pulse is skipped. This actually may have unintended implications, because in the past people have misadjusted the receiver gating times for functional delays (e.g., refocusing delays) in a pulse sequence: the fact that upon setting the pulse length to zero (for testing or calibration purposes), not just the pulse but also the two receiver gating delays disappear, has sometimes led to unexpected results.

The last if statement in the above function definition causes the receiver *not* to be turned back on if the receiver has been turned off globally (i.e., with a call to rcvroff() prior to calling rgpulse). The rcvroff_flag is used for UNITY-style and older rf only; for UNITY*plus* systems, a more complex construct is used within /vnmr/psg. (The value of rcvroff_flag needs to be stored in an intermediate variable, because the rcvroff() function alters that flag.)

What is the reason for using three time events in the above construction? Why can't we simply turn the amplifier on and off again? In order to discuss this, we need to take a detailed look at the timing diagram, which now (different from hsdelay) involves *several* devices: the transmitter, the receiver, and the transmitter phase setting.

As shown in Figure 7, it turns out that implicitly also the amplifier is involved in a particular way.

For the transmitter and the receiver (and the amplifier), simple gating (on/off) lines are involved. The 90-degree phase setting is slightly more complex, as it involves four possible values instead of a simple on/off information.

It turns out that we can deduce the four phase values from *two* on/off bits: an on/off bit (each bit represents the values 0 or 1, i.e., off or on) for the 180-degree phase shift and an on/off bit for the 90-degree phase shift.

In all, a pulse involves *five* on/off bits (fast bits, as we call them later): the transmitter, the receiver, the amplifier, the 180-degree phase shift, and the 90-degree phase shift.

Chapter 4. Time Events


| 180 | 90 | Phase shift |
|-----|-----|-------------|
| off | off | 0 degrees |
| off | on | 90 degrees |
| on | off | 180 degrees |
| on | on | 270 degrees |


In the table below, assume that we switch the transmitter phase from 90 to 180 degrees and that all other devices were in their default setting for a delay:

Obviously, a pulse is much more complex than a simple delay! Why all this? One of the points is that the receiver must be shut down while we are pulsing on the observe channel because both the pulse and the NMR signal are running partially through the same wire. Also, we want to prevent the NMR signal from going back into the amplifier, because this would cause a significant signal loss.

The diagram in Figure 8 shows the switching mode for receiving a signal ("receiver switched on"). In this mode, the T-switch (transmit switch) between the amplifier and the probe is open and the R-switch (read switch) is also open. This way, the NMR signal can only proceed into the preamplifier.

With the receiver switched off, both parts of the T/R-switch are closed, and the amplifier is now connected to the probe. The R-switch puts a specific point between the probe and the preamplifier to 0 volts (ground). At that point, the rf pulse is








Figure 8. Switching mode for receiving an NMR signal

reflected, building up a standing wave back in the direction of the amplifier (if there was no probe connected we would in fact reflect the entire power back into the amplifier!).

If now the cable between the zero voltage point and the probe connector has the right length (1/4 or 3/4 of the wavelength), the standing wave has maximum amplitude at the probe connector and all the power is directed into the probe, as shown in Figure 9.




Figure 9. Maximum power directed into the probe

This works very efficiently and avoids a switch *in* the line between the probe and the preamplifier (which would cause a significant signal loss). The only disadvantage is that the length of the cable needs to be adjusted with the frequency. If the length of the quarter-wavelength cable is not adjusted, a part of the amplitude is reflected back into the amplifier, and we would have less power available in the probe; this would affect the length of the 90-degree pulse.

In earlier rf schemes, both the T-switch and the R-switch were not active switches but rather a pair of passive, crossed diodes:



Such an arrangement would act as a closed switch for voltages above 1 volt (such as an rf pulse), but for voltages below around 0.5 volts, the diodes act as an open switch, so that low voltages (like an NMR signal) cannot pass. In reality, there is still some leakage across the diodes, but as far as the NMR signal is concerned, this is negligible.


With the advent of shaped pulses, this kind of rf switch became impossible because the low-level (voltage) part of a pulse shape would be heavily distorted. To avoid that, the passive diodes have been replaced by actively switched PIN diodes, a special type of

diode that are switched by applying a high dc voltage across them. The switching of the T/R-switch occurs anti-parallel to the receiver switching.

Even when these diode switches are opened (i.e., when the receiver is on), there is still some signal leakage across them. The amount is sufficient for the amplifier noise to enter the receiver chain and ruin the signal-to-noise ratio. Linear amplifiers (which are prerequisite for performing shaped pulses) produce relatively high levels of rf noise that can only be suppressed by turning off the last stage of the amplifier ("blanking the amplifier"). For that reason, the observe channel amplifier is automatically switched off ("blanked") whenever the receiver is on (anti-parallel switching).

All this explains the basic switching diagram shown in Figure 8, but does not explain why there should be *delays* before and after the actual pulse. In particular, we can't yet deduce any recommendations as to what length of delay is required before and after the pulse. Also, we would like to know whether there are any hidden additional delays not shown on this diagram.

The second question depends on the construction of the pulse programmer and cannot be answered at this point in time (in reality—as we shall see—there are no additional hidden delays involved in the rgpulse function).



To explain the first point we should check how the diagram in Figure 7 looks like in "real life"; This is shown in Figure 10.

Figure 10. Timing diagram in "real life"

The real reason for introducing additional delays around an rf pulse lies in the fact that we are dealing with real hardware. No matter how quickly a pulse programmer switches lines, it always takes a finite time for the various states to establish themselves. Such delays are already introduced by the pulse programmer (which has a finite switching time), but since this delays all states by the same amount, we don't have to consider this (as well as propagation delays through electrical cables) any further. In order to judge the time constants or the time to completely change status on the devices involved in an rgpulse statement, we have to consider the individual hardware involved. First, consider gate switching around the beginning of the pulse.

- The **receiver gate** (without the linked amplifier blanking) operates a number of TTL (transistor-transistor logic) gates and PIN-diodes (such as the T/R switch) with very rapid switching times (typically in the order of nanoseconds). For this reason, the receiver can be considered switching instantaneously relative to the time scale involved in pulse sequences—no particular delay is required to allow the receive and related hardware to switch on or off.
- The opposite is true for the **amplifier blanking**: when the blanking is removed, it takes a relatively long time for linear amplifiers to reach full output amplitude and phase stability. The problem is not so much the turn-on time, but rather the time it takes for the amplitude to stabilize. The time for blanking the output after a pulse is not so critical. UNITY and UNITY*plus* systems use linear (class A/B) amplifiers that are trimmed for short unblanking time: they reach full amplitude stability after 2 to 4 microseconds. Class A linear amplifiers in early high-field VXR spectrometers take up to 40 microseconds to stabilize after unblanking.

Class C amplifiers used in earlier equipment don't need to be blanked: they produce very little noise and are left on all the time. They are also much more efficient than linear amplifiers, but unfortunately they are unusable for pulse shaping. They are also not broadband by nature—an entire array of amplifiers is required to cover the full frequency range.

The amplifier appears to be the prime reason for having a delay prior to starting a pulse. If the amplifier was blanked (switched off) beforehand, pulses on the observe channel should be preceded by a delay of 10 microseconds, to make sure the pulse only starts when full output stability is reached. The unblanking time can be determined by performing a single pulse experiment (s2pul sequence with a sample of doped D₂O) using a very short pulse width (typically 0.2 to 0.5 microseconds) and arraying rof1 in the range of 0 to 10 microseconds.

On the other hand, if the amplifier was already unblanked beforehand (perhaps because another pulse is preceding the current rgpulse), the amplifier is still turned on (unblanked) and does not require a special delay. From an amplifier point-of-view, back-to-back pulses are no problem at all.

• To find out about the 90-degree ("quadrature") **phase shifting times**, we need to know how these phase shifts are generated. There are fundamental differences between a UNITY*plus* and prior generations of instruments in the way 90-degree and small-angle phase shifts are generated. In the UNITY*plus*, both kinds of phase shifts are generated on the same (transmitter) board, but by dedicated hardware, as shown in Figure 11 (the mechanism of programming small-angle phase shifts is addressed in Chapter 12, "AP Bus Traffic," on page 137). The 90-degree phase shifts are generated in a special circuitry (a "twisted ring shift register") that generates four 10.5 MHz frequencies (0, 90, 180, and 270 degrees) from 42 MHz, out of which one is selected by means of two digital input lines (90 and 180 degrees phase shift flags)³. The switching between the four phases is virtually immediate (TTL gates, switching times in the order of nanoseconds).

³ For wideline and multipulse spectroscopy, the 90-, 180-, and 270-degree phase angles can be fine adjusted in the order of millidegrees.




Figure 11. Phase shifting on UNITYplus systems

For the small-angle phase shifting, the 10.5-MHz frequency is split into 0- and 90degree components. Small-angle phase shifting takes 9 bits of input: values run from 0 to 512, with 360 steps for phase angles of 0 to 89.75 degrees, in 0.25-degree resolution. These are translated into digital amplitude values for the 0-degree and 90-degree components: a DAC converts the numeric values into dc voltages, which are then used in Gilbert multipliers to set the amplitude of the two orthogonal 10.5-MHz components. These two components are combined again, generating again a single 10.5-MHz frequency that has the (digital) small-angle phase shifting added in. This is the 10.5-MHz IF (intermediate frequency) that propagates the 90-degree and small-angle phase shifts into the transmitter frequency. The small-angle phase shifting circuitry has no effect on the timing of a pulse on the observe channel.

Figure 12 shows phase shifting on systems prior to the UNITY*plus*. On these systems, both 90-degree and small-angle phase shifts were generated using a 720-step sine look-up table. This generated a frequency of 2.5 MHz, which in turn was used to generate the 10.5 MHz IF.




Figure 12. Phase shifting on systems prior to UNITYplus

With this type of rf generation, a 2.5-MHz frequency is generated by reading the sine value from a 360-degree look-up table (720 steps, 0.5 degrees resolution) at a rate of 10 MHz. The resulting output (with a very rough digitization in the time axis) is then filtered in order to obtain a pure 2.5-MHz frequency that can then be mixed with 8 MHz to generate the 10.5 MHz IF. Small-angle phase shifts are nothing but an offset (up to 180 steps of 0.5 degrees) in the same look-up table.

The 90-degree phase shifts are selected by means of two digital input lines that cause look-up offsets in increments of 180 addresses.

What is now the timing of a 90-degree phase shift? The look-up offset is implemented immediately (TTL gates), but it only has an effect at the next 10 MHz trigger point, which can be up to 0.1 microseconds later. Also, the filter in the 2.5-MHz line further delays the response, such that with this type of rf a 90-degree phase shift is completed typically after 0.4 microseconds. Thus, delays of 0.5 to 1 microseconds are recommended to allow for the phase to settle.

• The actual **pulse switching** has several components: the rf gates and associated hardware on the transmitter board, the properties of the amplifier, and the inductance of the coil. Similar to the receiver gating, the rf gates on the transmitter board (TTL switches) can be thought of as switching instantaneously (nanoseconds range); therefore, the pulse turn-on time determined only by the amplitude rise time of the amplifier, which is usually determined using an oscilloscope. It turns out that even pulses of only 0.2 microseconds duration is fairly rectangular. Typically, within a few tens of nanoseconds, the output reaches full amplitude.

The pulse turn-on time is of no concern for the delay prior to the pulse, but it should at least be ensured that the rf is stable when the pulse is turned on. Also, we can do very little in terms of timing in a pulse sequence about turn-on effects from the inductance of the coil (see below for a discussion on probe ringdown).

The high-power amplifiers used with high-power solid-state NMR experiments $(1 \text{ kW }^{1}\text{H}/^{19}\text{F} \text{ or broadband})$ are at least as fast as the linear amplifiers used for liquids NMR (which feed the high-power amplifiers); thus, the amplifier considerations for solid-state systems are no different than for liquids experiments.

Taking all this into consideration, we can conclude that there are two cases for the delay preceding an rf pulse on the observe channel:

- If the amplifier was blanked beforehand, pulses on the observe channel should be preceded by a 10-microsecond delay, to allow for the amplifier to stabilize. At the same time, this delay is more than sufficient for the new phase to settle with all types of rf generation.
- If the amplifier was *not* blanked beforehand (either another pulse is preceding the current one or the receiver was switched off during the preceding [time] event), the delay prior to the pulse is mainly a phase-settling delay and can be reduced accordingly. On UNITY-type and earlier rf schemes, a delay of 1.0 microseconds should be more than sufficient for the phase to settle. With UNITY*plus*-type rf (with instantaneous phase switching), 0.2 microseconds is long enough—even back-to-back pulses do not cause problems.

What needs to be considered for determining the setting of the delay *following* a pulse on the observe channel? It turns out that this is simpler, because much less gating is involved and no phase change occurs.

- At the end of the pulse, the gates on the transmitter board close immediately; the **pulse turn-off time** should be in the order of nanoseconds. The amplifier turn-off time may be slightly longer, but this should not be a concern within most of the pulse sequence, except for the last pulse, where we should make sure that no "pulse break-through" occurs, which could affect preamplifier performance.
- The **receiver gates** operate at about the same speed as the transmitter gates.

Chapter 4. Time Events

• There are phenomena outside the actual rf that also need to be considered here, namely **probe ring-down** and related effects. The inductance (or better called the quality factor, *Q*) of the rf coil causes a finite delay in turning off the pulse (rf ring-down). With liquids probes, the rf ring-down time constant is typically 200 nanoseconds (high-field, ¹H) to several microseconds (low-gamma nuclei). Wideline and high-power CRAMPS probes have a lower *Q* and therefore much shorter ring-down time constants of 20 to 50 nanoseconds at proton frequency⁴.

Also (and often even more important), the coil experiences considerable mechanical forces during an rf pulse that shakes the coil mechanically (**acoustic ringing**), even if the coil is fixed firmly. These mechanical movements induce currents in the coil, and these currents in turn may affect the preamplifier and the receiver, distorting the beginning of the FID.

In conclusion, we have the following different situations for the delay after a pulse:

- If data are acquired immediately following the pulse (i.e., after an excitation pulse rather than after a refocusing pulse) or if pulses are performed in-between acquisition points, a delay is recommended after the last pulse in order to avoid probe ring-down saturating the preamplifier. For normal liquids experiments (¹H, ¹³C, ³¹P, ¹⁹F), 10 to 20 microseconds should be sufficient; for low-gamma nuclei, longer values are often appropriate (if the T_2 permits); and for typical solid-state NMR experiments, shorter values can be used because these probes have a much shorter ring-down time.
- For all other pulses in a pulse sequence, preamplifier saturation should not be a problem because it generally has sufficient time to recover. Hence, in a typical liquids experiments, the post-pulse delay can be set to zero for all pulses except the last pulse.
- If the receiver is switched off "globally" within a pulse sequence (e.g., switched off before the first pulse), all post-pulse delays *except for the last* can be set to zero.

Simple Pulses on Other RF Channels

Pulses on any rf channel can be performed using the genpulse statement the same way as used for the observe channel, except that the device OBSch is replaced by DECch, DEC2ch, or DEC3ch for the first, second, or third decoupler channel, respectively (see the footnote on page 40). As mentioned before, genpulse is not generally used (nor recommended) except for cases where the rf channel (the device) should be held under parameter control.

Typically, such pulses are performed using one of the following statements:

```
decrgpulse(length,phase,rx1,rx2);
decpulse(length,phase);
dec2rgpulse(length,phase,rx1,rx2);
dec3rgpulse(length,phase,rx1,rx2);
```

⁴ The rf ring-down time constant t_r is proportional to the quality factor Q, and inversely proportional to the rf frequency f: tr=Q/(3f). Liquids probes have a high Q of around 300 to 400, the Q of highpower wideline and CRAMPS probes is typically 10 times lower. After that time, the voltage across the coil has decayed to 1/e of its initial value. The total ring-down time t_{tr} for a pulse at power P (in watts), down to a residual level of 1 microvolt at an impedance of 50 ohms can be calculated as

$$t_{tr} = \left(-ln\left(\frac{1}{\sqrt{50P}}\right)\right) \times Q/(3f)$$

These statements are equivalents to rgpulse and pulse (for the second and third decoupler channel, no pulse equivalents have been defined). Unlike pulse, the preand post-pulse delays in decpulse are set to zero. All these statements are macros and are resolved to calls to the generic G_Pulse statement by the C preprocessor.

The gating scheme around decoupler pulses is basically the same as for pulses on the observe channel. A difference exists insofar as the decoupler amplifier(s) are usually set into continuous wave (CW) mode (i.e., they are constantly on⁵). Amplifier (noise) blanking is usually not required on these channels because their output is almost always be directed onto a different rf coil and can hardly reach the receiver⁶. For this reason a pre-pulse delay is only used for the phase to settle (particularly on UNITY and earlier systems; 0.5 to 1 microseconds are recommended), and a post-pulse delay is almost never necessary. On UNITY*plus* systems both delays can be left at zero length.

Simultaneous Pulses on Different RF Channels

The VNMR definition of simultaneous pulses is such that they should occur centered on top of each other. Based on this, a gating scheme can be deduced, for example, for simpulse (a simultaneous pulse on the observe channel and the first decoupler):

- 1. Gate receiver off, set phase on all channels involved.
- 2. Perform pre-pulse delay.
- 3. Switch on transmitter gate on the channel with the longer pulse.
- 4. Perform a delay; length: half the difference between the two pulse lengths.
- 5. Switch on transmitter gate on the channel with the shorter pulse.
- 6. Perform a delay; length: the shorter of the two pulse lengths.
- 7. Switch off transmitter gate on the channel with the shorter pulse.
- 8. Perform a delay; length: half the difference between the two pulse lengths.
- 9. Switch off transmitter gate on the channel with the longer pulse.
- 10. Perform post-pulse delay.
- 11. Gate receiver on, unless it has been switched off globally.

The C construct for such a statement uses the same basic calls as the rgpulse statement described above; if-else statements are used to decide which device to turn on first and which to turn off last. Similarly, statements for simultaneous pulses on three and four channels can be constructed. All these statements exist; they have the syntax:

```
simpulse(len1,len2,ph1,ph2,rx1,rx2);
sim3pulse(len1,len2,len3,ph1,ph2,ph3,rx1,rx2);
sim4pulse(len1,len2,len3,len4,ph1,ph2,ph3,ph4,rx1,rx2);
```

The indices in the arguments refer to the observe channel (1) and the first three decoupler channels (2, 3, 4). In all these statements, the observe channel is involved; therefore, (because in general we cannot, and do not want to, make assumptions about

 $^{^{5}}$ Except in situations where the decoupler nucleus is identical to the observe nucleus or falls into the same amplifier frequency band (3 H, 1 H, 19 F vs. the nuclei at 31 P frequency and below). In this case, decoupling occurs through the same amplifier as the observe pulses, which are in pulse mode and are therefore blanked.

⁶ Provided there is little or no crosstalk between the different coils in a probe. Exceptions would also be heteronuclear experiments on a single, double-tuned coil. In such cases the amplifier would have to be either put into pulsed mode or blanked explicitly, see Section 4.3, "Other State-Related Pulse Sequence Statements," on page 52.

Chapter 4. Time Events

the length of the individual pulses) for the pre-pulse delay the same considerations apply as for the pre-pulse delay to rgpulse, assuming that the amplifier for the longest pulse needs to be unblanked first, unless the receiver was already switched off during the previous time event.

The above statements are macros in which simpulse and sim3pulse are resolved into calls to the following C functions:

gensim2pulse(len1,len2,ph1,ph2,rx1,rx2,dev1,dev2); gensim3pulse(len1,len2,len3,ph1,ph2,ph3,rx1,rx2,dev1,dev2,dev3);

These functions take the corresponding number of device arguments (rf channel addresses, see the footnote on page 40); these devices can be specified in any order as long as the pulse lengths and the phases are specified in the same order. The macro sim4pulse is translated into a call to the function G_Simpulse that will not be discussed any further here. Similar to genrgpulse and genpulse, these functions are not supported by the dps command at this point in time, and their use is not recommended unless there is a need to keep the rf channel selection under parameter control (with the presence of the rfchannel parameter, those functions are not longer required).

Composite Pulses

From a pulse sequence or gating point of view, composite pulses are nothing but a sequence of ordinary pulses without explicit delays in-between. As outlined in the sections above, all pulses but the first can be preceded by a very short pre-pulse delay only, to allow for the phase shift to complete. On UNITY*plus* systems, even that delay can be set to zero (on other systems this may lead to a slight performance loss). For the pre-pulse delay of the first pulse, the same considerations apply as for a normal pulse (see above).

All post-pulse delays should be set to zero to avoid unnecessary gaps between the pulses. The maximum permissible gap between components of a composite pulse is determined by the question whether noticeable precession can occur during such delays. In liquids experiments, gaps of 1 microsecond should be harmless.

Considerations for the Delays Following the Last Pulse

Additional considerations apply to the post-pulse delay after the last rf pulse on the observe channel (better said as the pulse that generates the final, detected coherence), because that delay affects both the reliability of the autophasing (aph) as well as the first-order phase parameter (and with that the baseline flatness). Why is that?

In the last receiver stages before the ADC, the NMR signal (with IF $+/-\delta$) is mixed with the intermediate frequency in order to obtain audio frequencies (signals at -sw/2to +sw/2 Hz). At that level, any noise (and extra signals) outside the spectral range need to be filtered away in order to avoid folding noise and undesired signals into the spectrum (the extra noise would dramatically affect the signal-to-noise ratio). Such a filter delays the audio signal by a time that is inversely proportional to the filter bandwidth. This is taken care of by the software—the filter delay (shown in Figure 13) is calculated automatically and made part of the implicit acquisition, where prior to acquiring the first data point a delay of alfa+1/(beta*fb) is performed.




Figure 13. Filter delay

alfa is a standard parameter (in microseconds, defined as a pulse), fb is the filter bandwidth, and beta is a constant that depends on the type of audio filter in use. For 4-pole Butterworth filters (used in VXR and earlier spectrometers), beta is 2.0; for 8pole quasi-elliptical filters (the standard filter in UNITY and UNITY *plus* systems), beta is 1.29; and for 6-pole Bessel filters (used for spectral windows above 100 kHz), it is 2.33⁷. With this setup, the user can adjust the (last) post-pulse delay, while having an independent parameter for regulating the total delay between the last pulse and the first sampling point.

In 1983, D. Hoult et al.⁸ noted that by making the *total* delay between the last pulse and the first sampling point equal to the theoretical filter delay (for Butterworth filters, 1/(2*fb)), the best baseline flatness was obtained. This seems to disagree with the standard setup in Varian pulse sequences (with rof2 and alfa being 10 and 20 microseconds). What are the consequences of these additional delays? It turns out that with these delays, the first-order phase correction parameter lp has a negative value that depends on the spectral window sw. In fact, the lp parameter is decremented by 360 degrees with every dwell time by which we right-shift the first sampling point. If the condition described by Hoult is fulfilled, lp is close to zero—an approximation for the lp parameter can be calculated easily from rof2 and alfa⁹:

$$lp = -\frac{rof2 + alfa}{10^6 \times (1/sw)} \times 360 = -(rof2 + alfa) \times sw \times 360 \times 10^{-6}$$

The real problem is that due to the filter response, the first data point is never accurate. An offset in the first data point leads to a dc offset in the transformed spectrum $(l_{p=0})$. The first-order phase correction converts this dc offset into a sinoidal baseline distortion: at a correction of ± 360 degrees, a full sine wave is generated, and with every additional 360-degree deviation of l_{p} from zero, an additional sine wave is obtained. Also, as l_{p} increases, this baseline roll increases in amplitude (with very high values for l_{p} , it is not unusual that the baseline roll becomes larger than the NMR signals!).

Using the above formula, we can estimate that for proton spectra (spectral windows of 2000 to 10000 Hz) first-order phase corrections of -20 up to -100 degrees are obtained with the standard settings for rof2 and alfa. For typical carbon spectra (spectral windows of 10 to 50 kHz), the corrections are more severe with 100 to 500 degrees. In many cases—especially proton spectra on normal organic compounds—this baseline roll is noticeable but can easily be corrected by means of data processing software (bc command). In many other cases, such a baseline roll is undesirable—it can even ruin

 $^{^7}$ In the $^1\mathrm{H}/^{13}\mathrm{C}$ version of Gemini spectrometers, 3-pole Butterworth filters with a filter constant beta of 3.0 are used.

⁸ D.I. Hoult, C.-N. Chen, H. Eden and M. Eden, *J. Magn. Reson.* **51**, 110 (1983).

 $^{^9}$ The division by 10^6 is needed because within VNMR rof2 and alfa are defined as pulses and are entered in microseconds.

the quality of multidimensional spectra. To help cure this, the VNMR software contains the hoult command, which calculates rof2 to be one-third of the theoretical filter delay (1/(beta*fb)), as described in Hoult's paper, and then sets alfa to the *negative* value of rof2. This only works because the implicit delay that is executed is the *sum* of alfa and the filter delay.

The hoult command brings lp down to a value close to zero—but for many applications, especially biological NMR, this is still not good enough. Why? The reason for that is twofold: first, even the best filters are never 100.0% accurate (all electronic components have a limited accuracy) and hence the real filter delay is always slightly differ from the theoretical value. Second, it cannot be assumed that the excitation point (the point where all excited spins are in phase) is at the end of the rf pulse, because the precession *during* the pulse cannot be neglected.

Overall, it is impossible to calculate the exact filter delay. Consequently, an *empirical* timing correction can be used: a (1D) spectrum is first acquired, transformed, and phase-corrected, and then the above equation is modified to recalculate alfa from lp. This is done by calling the macro calfa:

$$alfa_{corr} = alfa + (lp \times 10^6)/(360 \times sw)$$

With this new value for alfa, the spectrum is then reacquired. It can be assumed, that this correction stays the same as long as the spectral window is not changed (although this probably is only true as long as the pulse width is also not altered).

Another implication from the delays between the last pulse and the first sampling point is the fact that the autophasing command aph first estimates a value for lp, assuming that both rof2 and alfa were used after the last pulse. It must therefore be ensured that the last pulse is in fact performed with a post-pulse delay rof2 (or an equivalent delay must be inserted after the last pulse); otherwise, the reliability of the autophasing can be severely affected.

4.3 Other State-Related Pulse Sequence Statements

Direct Gating

As shown in the sample coding from the previous sections, rf gates can also be operated directly. A number of simple, documented commands exist for the transmitter gates (xmtron/xmtroff, decon/decoff, dec2on/dec2off, dec3on/dec3off), receiver gate (rcvron/rcvroff), spare gates (splon/sploff, sp2on/sp2off), and decoupler amplifier blanking control (decblank/decumblank, dec2blank/dec2umblank, dec3blank/dec3umblank). Pulse statements operate these gates implicitly and automatically. If a user decides to operate gates explicitly, it is the user's responsibility to switch them back into "normal" mode at an appropriate time in the pulse sequence.

The statements for explicit 90-degree phase shifting (txphase, decphase, dec2phase, dec3phase) can be regarded as special (phase) gating statements (with the difference that they do not control rf gates or blank an amplifier). All these gating statements have no immediate effect on what is happening in the rf: they only influence the gate settings of subsequent time events.

Implicit Gating

The status statement can partly be regarded as implicit gating statement. Primarily, status changes the value of the variable statusindx to the specified flag field. At a secondary level, it also switches the decoupler gates depending on the corresponding status fields in the dm, dm2, and dm3 parameters. These only become active in subsequent time event. The definition is such that if a flag has less fields than addressed in the pulse sequence, the last flag field is taken instead (status takes an integer as argument; the constants A to Z [0 to 25] are defined in psg/rfconst.h).

In earlier rf schemes, the decoupler modulation mode (dmm) was also on direct status lines (two lines for up to four modulation modes) and was set by status the same way as the other gate lines; in UNITY and UNITY*plus* systems, the decoupler modulation modes (dmm, dmm2, dmm3) are set through a different mechanism, see Chapter 12, "AP Bus Traffic," on page 137.

4.4 Basic Purpose of a Pulse Sequence

The main purpose of a pulse sequence is to define a sequence of *times* (delays) and a series of associated *states* (on/off information). Pulses—even if composite or simultaneous on various rf channels—can be deduced from a sequence of simple delays with variable gating information. From that point-of-view, the heart of the pulse programmer can be diagrammed as shown in Figure 14.




Figure 14. Block diagram for a pulse programmer

This is not only a block diagram for what the pulse programmer is doing, but, as we see later, this diagram contains the "essence" of the pulse programmer hardware in all Varian NMR spectrometers from the XL, VXR, and UNITY, to the Gemini, UNITY*plus, GEMINI 2000,* and UNITY*INOVA.* Of course, there are additional features (we shall learn about rf devices that can not be addressed simply with a limited set of on/off states), but the central idea is already in the diagram—the pulse programmer is a state–times buffer with associated timing electronics that makes the code interpretation (the acquisition CPU) independent (asynchronous) of the pulse sequence timing.

Chapter 4. Time Events

Chapter 5. Submit to Acquisition: go

The go command (and related commands such as au and ga) has the central task of running the compiled (executable) sequence from the seqlib directory and submitting the experiment to the acquisition process Acqproc (go does not perform the acquisition itself). go then hands the control back to VNMR. It is up to Acqproc to manage the acquisition and notify VNMR about completed blocks, FIDs, and experiments.

Figure 15 outlines the processes and actions around the go command.

5.1 The Tasks for go

The command go directly executes and initiates the following steps:

- From the parameter seqfil, the go command finds out what pulse sequence is to be executed. If such a pulse sequence executable is not found in the user or the system seqlib directory, go aborts with an error message.
- If an executable pulse sequence is found, go strips all non-acquisition parameters off the current parameter tree and adds acquisition parameters from the global tree along with all the system configuration parameters. This parameter collection is then transferred to the pulse sequence.¹




Figure 15. Processes surrounding the go command

¹ Because only acquisition-, sample- and configuration-related parameters reach the pulse sequence, a getval statement on a display, processing, or spin simulation parameter *does not work* (a zero value is returned instead for numeric parameters)!

- The shim method from the directory shimmethods (as specified by the parameter method) is packed into an internal parameter comsstring and added to the set of transferred acquisition parameters.
- go then executes the sequence from seqlib. A compiled pulse sequence in the user's sequence library (~/vnmrsys/seqlib) takes precedence over a sequence (with the same name) from the system sequence library (/vnmr/seqlib).
- The execution of the pulse sequence first involves the run-time linker, which binds the other objects from psg/libpsglib.so.x.y (see Section 2.5, "Compiling and Linking," on page 24) to the pulse sequence. As part of the run-time linking, a revision check is made on the run-time link library psg/libpsglib.so.x.y. If the revision of that file (.x.y extension of the file name) is not correct, it results in a "run-time" error message or even in a crash. For such a case, check whether an old compiled pulse sequence (from a previous software release) was executed instead of an updated one in /vnmr/seqlib). If that happened, recompile the pulse sequence. It could also be that you used psggen (see Chapter 3, "Object Library Generation: psggen," on page 29) with an earlier software release, and the local psg directory (~/vnmrsys/psg) was not deleted when a new release was installed. In this situation, you must delete (or re-create and update) the local psg directory, and recompile the pulse sequence if it was a local file.¹
- If the pulse sequence does not return an error code upon execution, go notifies the acquisition process Acqproc of the new experiment in the acquisition queue. After that, the VNMR part of the go command is finished, and VNMR is ready for executing the next command.
- go usually appends a new experiment to the end of the acquisition queue, except if go (or au) was called with the argument 'next', in which case the experiment is queued immediately after the running experiment (insertion at the front of the queue) or immediately before the next acquisition is started, when called as part of conditional werr or wexp processing of an experiment that was started with the command au('wait').
- If VNMR cannot contact Acqproc after several attempts, an error message "maximum number of retries exceeded . . PSG aborted" or "The acquisition daemon 'Acqproc' is not active!" is issued, and no experiment is started. In this case, either Acqproc has died and needs to be started again, or for some reason VNMR is unable to contact Acqproc (it may be hanging), in which case Acqproc needs to be killed and restarted again. At this point, any running experiment needs to be restarted, and queued experiments need to be requeued (see Section 6.4, "Controlling Acqproc," on page 63). After this, call go again.
- The absence of Acqproc makes go non-functional on stand-alone data stations (except for go('acqi'), see also Section 5.3, "Using go('acqi')," on page 58).

¹ These problems are avoided if makeuser is called for every account and the local files are updated through this command. makeuser deletes old compiled pulse sequences (not the source file from psglib) and disables any local psg directory.

5.2 Tasks for the Pulse Sequence Executable

The primary task of the pulse sequence executable is to generate instructions and data that Acqproc later uploads into the acquisition computer plus data to be used by Acqproc to control that experiment. The resulting information is stored in several files in the acquisition queue directory /vnmr/acqqueue. The names for most of these files reflect the experiment name (expl to exp9), the name of the user (username), and the process-ID number of the pulse sequence executable (PID). Obviously there would be a problem if a user did not have write permission into /vnmr/acqqueue.

- Parameters are read from the set of acquisition parameters that is handed over by go; they are stored in local variables.
- Specified waveform generator (WFG) shapes and patterns are read from shapelib and transformed into a binary file expn.username.PID.RF in /vnmr/acqqueue. This file is later uploaded via HAL into the waveform generator(s). Together with expn.username.PID.RF, a second file named /vnmr/acqqueue/ldcontrol is generated or updated. This second file contains information about the size and structure of all waveforms and decoupling pattern to be used. If a system is not equipped with waveform generators, this file is not generated or updated; if waveform generators are installed, the file /vnmr/acqqueue/ldcontrol is updated even if no waveform generator is actively used in the experiment.
- The main task for the pulse sequence executable is the generation of instructions to be interpreted by the acquisition CPU called the *Acode*. The Acode is a sequence of 16- and 32-bit words, including instructions and a complete set of local variables used by the acquisition CPU, that are stored in a file in the directory /vnmr/acqqueue. The file name is expn.username.PID.Code. Typically, the Acode size is 1 to 3 Kbytes *per FID*. For multidimensional experiments, the Acode can easily fill many megabytes of disk space and can eventually create "disk full" problems.
- External phase tables are read from the specified file in tablib; they are built into the file with the instruction list for the acquisition CPU (see Chapter 9, "Pulse Programmers," on page 85).
- A file /vnmr/acqqueue/expn.username.PID is generated that contains acquisition control parameters the parameters that Acqproc needs to control the acquisition, store the FID in the proper place, and initiate the right actions when certain key conditions (error condition, block size completed, FID or experiment completed) are reached. This file is a standard ASCII parameter file, containing a set of parameters that are required for Acqproc to do its job (see also Chapter 6, "Acquisition Process," on page 59).

5.3 Using go('acqi')

Entering go('acqi'), or the equivalent macrogf, serves to generate the Acodes necessary to run acqi in interactive (FID or spectrum) mode. The same as entering go with no arguments, the pulse sequence executable is called, but in the case of go('acqi'), nothing is submitted to Acqproc, and the resulting files are read solely by the acqi program. For this reason, different naming conventions are used:

- acqi.Code replaces /vnmr/acqqueue/expn.username.PID.Code.
- acqi.RF replaces /vnmr/acqqueue/expn.username.PID.RF (only used when shaped pulses or rf patterns involving a waveform generator are involved).
- acqi.par replaces /vnmr/acqqueue/expn.username.PID; this file (also an ASCII parameter file) is meant for acqi directly; it tells acqi how to display the FID or the spectrum interactively (sp, wp, etc.).
- ldcontrol is the same as with the normal go command (only used if waveform generators are present).

Entering go('acqi') or gf is often also used to test a pulse sequence up to the point of Acode generation (without doing a real acquisition yet) or to permit inspection of the Acode and associated files. Note that the Acode generated by go('acqi') differs in several points from the code generated through a normal go:

- Only Acode for the first FID is generated,
- Parameter alock is set to 'n' (no automatic locking),
- Parameter load is set to 'n' (shim values are not loaded),
- Parameter wshim is set to 'n' (no automatic shimming occurs),
- Parameter ss is set to 0 (no steady-state pulses),
- Parameter dp is set to 'y' (32-bit acquisition),

The macro gf modifies even more parameters:

- Parameter cp is set to 'n' (the observe phase oph is not cycled automatically),
- Parameter gain is set to 'y' (use fixed gain),
- The lock to be set into the "fast loop" mode.

By definition, the parameters used for shimming on the FID or on the spectrum (as well as the parameters for interactive parameter adjustment) are those for the first increment of an arrayed experiment. Also, gf sets the parameters such that excessive startup delays are avoided (alock, wshim, ss, and gain parameters). For the FID mode in particular, cp='n' avoids jumping baselines and/or intensities by stopping the phase cycling (assuming the cp flag or oph are used in the pulse sequence).

Note also that the Acode segment for the first FID—as obtained via go('acqi') can *not* be taken as a model for subsequent code segments (both in size and contents), because it contains a number of rf initialization instructions that are not repeated for every other FID. The first Acode segment is always bigger than all following segments of an arrayed or multidimensional experiment.

Chapter 6. Acquisition Process

The acquisition process Acqproc operates in three phases: uploading and starting the acquisition operating system, uploading data for the current acquisition, and downloading acquired FIDs and associated information. Acqproc is the main software link from the host computer to the acquisition cabinet (the actual spectrometer)¹.

Acqproc is a background process that must be constantly running on a spectrometer. It is started at boot-up time through /vnmr/rc.vnmr (which is called by /etc/rc.local) if the file (flag) /etc/acqpresent is found. /vnmr/rc.vnmr, the UNIX kernel (/vmunix) with the Varian device driver, and the acqpresent flag are installed together with the VNMR software by the setacq installation script.

From a hardware point-of-view, the host computer is linked to the acquisition CPU through the HAL (Host-to-Acquisition Link) board, a board with a SCSI interface. Thus, the two computers are connected by a SCSI bus. Unless there is a second SCSI bus used exclusively for acquisition, the same SCSI bus is usually used for communication with hard disks, tape drives, and a CD-ROM drive (if present)². Acting as a special device on the SCSI bus, the HAL board requires a special device driver in the UNIX kernel (i.e., a special, proprietary communication protocol is used between the the host computer and the acquisition CPU).

6.1 Starting the Acquisition Operating System

Upon startup, Acqproc checks if the acquisition CPU has resident software. If not (e.g., after initialization), Acqproc uploads the acquisition operating system (i.e., the software that permanently resides in the acquisition CPU). The files that form the acquisition operating system are stored in the directory /vnmr/acq:

- rhmon.out is the part of acquisition CPU operating system that controls the communication with the host computer through the HAL board.
- autshm.out is the autoshimming part of the operating system.
- xrop.out is the central part of operating system (Acode decoder) for systems with output board (63-step FIFO).
- xrxrp.out is the central part of operating system (Acode decoder) for systems with acquisition control board (1024-step FIFO, mostly UNITY systems).
- xrxrh.out is the central part of operating system (Acode decoder) for systems with pulse sequence control board (2048-step FIFO, UNITY*plus* systems).

¹ The interactive acquisition program acqi talks directly to the acquisition CPU, be it for locking, for shimming, or for interactive acquisition.

² Because the electrical specification of the SCSI bus would impose a severe restriction with respect to the maximum cable length between the two computers (or any two SCSI devices), Varian uses a special "differential box" that adds an electrically driven branch to the SCSI bus. This box permits using a maximum cable length of over 10 meters between the host computer (actually the differential box) and the acquisition cabinet.

• xr.conf is the part of operating system that returns information about the rf configuration to Acqproc and the config program (/vnmr/bin/vconfig).

In earlier VNMR releases, a symbolic link named xr.out was created to one of the two files xrop.out and xrxrp.out when calling the setacq script, in which the user had to reply whether the system was equipped with 63-step or 1024-step loop FIFO. The task of selecting the right version of xrop.out, xrxrp.out, or xrxrh.out is taken over by Acqproc, which gets the necessary hardware information from the xr.conf part of the acquisition CPU operating system.

It is possible (although unlikely) that after the installation of a new version of VNMR, the acquisition CPU still has an old version of the acquisition operating system loaded. In such a situation, Acqproc may not "realize" the revision discrepancy, because the acquisition CPU is still running on uploaded software, which in turn leads to error messages or even a system hangup—at the very least when trying to communicate with the acquisition CPU. After loading VNMR (before starting the new Acqproc), the acquisition CPU *must* be reinitialized.

6.2 Queuing and Starting the Acquisition

Acqproc monitors the acquisition CPU. If no experiment is running or when the acquisition CPU is ready to accept the next Acode segment (instruction set and data for one FID), Acqproc uploads that Acode (and waveform generator data files if required) to the acquisition CPU (actually to the Host-to-Acquisition Link) and initiates the acquisition. Figure 16 shows the communications scheme controlled by Acqproc.

Acqproc controls not only the data acquisition, but also the queuing of acquisitions from different experiments and different users. The experiments in the queue are executed in the order of their submission³. "Acquisitions" here also include queue




Figure 16. Acquisition control by Acqproc

³ The only means of manipulating the queue is to abort the current acquisition (aa from the current experiment), to stop it temporarily and go to the next queued experiment (sa from the active experiment), to remove an acquisition from the queue (sa from the queued experiment before it becomes active), or to add a new experiment to the end of the queue.

items that do not generate NMR data as submitted by change, spin, lock, shim, sample, and su; these will not be discussed in detail here because their actions can also be made part of a "normal" experiment (as started with the go, ga, or au).

The process of starting up an experiment (after its submission by go) involves the following steps:

- If either the acquisition CPU is idle or has just terminated the previous experiment, the next item in the acquisition queue is launched.
- If waveform generator-shaped pulses, decoupling pattern, or shaped gradients are involved, the necessary waveform generator data are uploaded into the acquisition CPU memory from /vnmr/acqqueue/expn.username.PID.RF. The file /vnmr/acqqueue/ldcontrol contains the control information that permits Acqproc to instruct the acquisition CPU about the structures within the waveform generator data file, and where to load these structures (pattern, waveform).
- The Acode for the first FID is uploaded into the acquisition CPU as well, and its interpretation is initiated.
- In arrayed or multidimensional experiments, the Acode segment for the next FID is loaded ahead of time, such that the acquisition can switch from one FID to the next without having to upload Acode in between FIDs. This is called buffered acquisition; it dramatically reduces the dead-times between FIDs (no extra delays for data transfers or disk access).
- Acqproc also acts as a message handler and transmits interrupt messages (such as generated by the aa or sa commands) to the acquisition CPU.

6.3 Downloading the FID

Acqproc polls the acquisition CPU at regular (dynamic) intervals. If data are ready (after a block size or a completed FID or experiment), Acqproc downloads it onto the host computer memory and stores the data in the experiment directory from which the acquisition was started. This involves the following actions and mechanisms:

- The data transfer from the acquisition CPU (the HAL) to the host computer occurs through polling by Acqproc. Acqproc polls at regular intervals, but tries adjusting the polling rate to the rate at which data are ready in the acquisition CPU. This avoids unnecessary dead times.
- If data are ready (block size or FID completed, sa interrupt completed), Acqproc downloads the FID from the HAL memory to the host computer.
- Acqproc then checks whether a lock file ~/vnmrsys/lock_n.primary⁴ exists for the current experiment. If a lock file is found, its contents are inspected

⁴The lock file ~/vnmrsys/lock_n.primary has the experiment number coded in the file name (lock_n). The extension .primary indicates that it is a *primary lock file*. A *secondary lock file* is generated temporarily before generating the primary one, in order to avoid that two processes lock the same experiment at the same time (i.e., between detecting that "no lock file exists" and the creation of a lock file, another process could in theory create the same lock file; the secondary lock file avoids this). The lock file contains a code for the kind of lock (3 = foreground), the host name of the system on which the locking process runs, and the process-ID (PID) of the locking process. The PID permits detecting whether the locking process is still active.

in order to find out whether a foreground VNMR is active in the current experiment.

- The FID is stored in the appropriate experiment by over-writing or updating the applicable section of the FID file (~/vnmrsys/expn/acqfil/fid). The ownership of the FID file is transferred (chown) to the experiment owner. During this process, the active experiment is locked by Acqproc (acquisition lock).
- If a foreground VNMR is running in the current experiment, it is notified through send2Vnmr, and its (buffered) processed parameter tree is updated (ct and celem parameters). If conditional processing was specified (werr, wbs, wnt, or wexp parameters), this processing is initiated through the send2Vnmr call. The send2Vnmr call uses the following syntax:

```
ps -ax | /vnmr/acqbin/send2Vnmr "VNMR command string"
```

• If no foreground VNMR is found, or if no foreground VNMR process is running in the current experiment but conditional processing was specified, that processing is performed through a background VNMR call:

Vnmr -mback -nExp# "VNMR command string"

During this process, the active experiment is locked by VNMR (background lock).

- A log file /vnmr/acqqueue/Masterlog is maintained that contains information about submitted experiments, interrupts and completion points, as well as information about acquisition-related processing.
- When the experiment completes, the queue files expn.username.PID.Code, expn.username.PID.RF (if present), and expn.username.PID in /vnmr/acqqueue are deleted, and the next experiment is submitted to the acquisition CPU.

All the information that Acqproc needs for handling downloaded FIDs are contained in the parameter file /vnmr/acqqueue/expn.username.PID. This file includes the following parameters:

- Parameters dp, np, nf, and arraydim that contain data format information.
- Parameters bs, nt, ct, celem, and interleave (the latter corresponding to il) that allow Acqproc to determine whether a block or an FID are completed and which code segment to load next (if interleave='y', cycle block-wise through all FIDs; otherwise, switch to the next FID when the previous one is completed).
- Parameters werr, wbs, wnt, and wexp that specify what action needs to be initiated in VNMR in the case of an error condition, a completed block, a completed FID, or completed experiment.
- Parameter wait that specifies whether the following experiment can be started immediately following the current one (wait='n', the standard case) or whether the conditional (wexp or werr) processing should be completed first (wait='y'). Setting wait to 'y' allows for another experiment to be started on the same sample using au('next'), before switching to the next sample (i.e., at the front of the acquisition queue).
- Parameters gain and spin are used in connection with the Acquisition Status window. Unused are parameters priority, the queuing priority, and suflag.

Path and other file information:

- Directory fidpath stores the FID file.
- curexp is usually the parent of fidpath.
- userdir is usually the parent of curexp.
- The variable systemdir is usually /vnmr.
- id is the name of the parameter file itself; it is also the root part of the names for the Acode and waveform generator data files.
- date is a time stamp for the submission of the experiment

6.4 Controlling Acqproc

Because Acqproc needs to be able to store FIDs and modify parameter files in every user's experiment directories, it must be a process owned by root. Also, only root can use chown to transfer the ownership to the respective users. If Acqproc were owned by a standard VNMR user, only that user could acquire data⁵ (note that only one copy of Acqproc per system can run at a time).

The root ownership for the process Acqproc can create problems: in case of errors (in particular, experiencing communication failures on the SCSI bus) every user should be able to kill and restart Acqproc. And in the case of software troubles in the acquisition system, the acquisition CPU sometimes needs to be restarted, which should only be done while Acqproc is shut down (otherwise very likely Acqproc will not be able to communicate with the HAL board). UNIX processes can only be killed by their respective owners, or by root; this implies that only root can kill Acqproc, and root must restart it, unless the SUID protection bit (4000) is set. A consequence could be that many people need to know the root password, which of course is totally unacceptable from a security point of view.

The solution for this situation is a special account acqproc that is equipped with root privileges (UID and GID are the same as for root) and which may or may not be equipped with a password. In networked environments, a password (different from the root password, of course) is certainly recommended; however, this precludes resetting Acqproc from the .rootmenu to generate such an account. A shell script /vnmr/bin/makesuacqproc is provided, which must be executed as root.

When logging into the account acqproc, no shell is obtained (as with usual UNIX accounts), but a shell script /vnmr/bin/execkillacqproc (only executable by root) is executed instead, after which the user is automatically logged out again immediately. /vnmr/bin/execkillacqproc kills Acqproc if it is found in the process table; otherwise, it starts Acqproc (as root). This means that for restarting (killing and starting) Acqproc, the user must call su acqproc twice. Broadcast messages indicate whether Acqproc has been aborted or started (this message is displayed in every active UNIX window and the VNMR master window).

⁵ The other users could open up the permissions in their directories and files, but this would not solve all problems; also, it is undesirable to loosen the file system security by making so many files writable by everybody (or at least the group nmr).

For SunView, by adding the line

"Fix Acquisition" su acqproc

to a user's ~/.rootmenu file, this function becomes available in the SunView.rootmenu. This can only work if the account acqproc has no password, because there is no keyboard input to commands that are executed from the .rootmenu.

For OpenWindows, the menu file is ~/.openwin-menu, and the line to be added is as follows:

"Fix Acquisition" exec su acqproc

Chapter 7. Digital Components

For the remaining sections of this manual, having at least a rough idea about the structure of the digital part of the spectrometer is essential. Figure 17 gives on overview of the relevant components that are involved—more details are provided in the following chapters.




Figure 17. Digital components of a spectrometer

Chapter 7. Digital Components

It turns out that there are two CPU boards in the acquisition computer: the Host-to-Acquisition Link (HAL) and the acquisition CPU. Both use a Motorola MC68000 CPU chip and have a serial port for a diagnostics terminal. The two computers fulfil different tasks (the HAL is dedicated to communications with the host computer and it holds the actual NMR data), but they operate on the same bus, share the same address space, and work on the same data. Therefore, in subsequent chapters we treat them both as one single computer (the fact that there are two CPUs is irrelevant to the user).

7.1 Main Boards

The tasks for the main components in the digital part of the spectrometer can be summarized as follows:

• The main purpose of the **HAL board** is to act as a link between the host computer and the acquisition CPU. Through its SCSI interface, the host computer (Acqproc) can upload not only Acode (see Chapter 8, "Acquisition CPU and Acode," on page 69), but also the acquisition operating system through DMA (direct-to-memory access) into the acquisition CPU. After an experiment, the NMR data are downloaded through the same interface (see Chapter 6, "Acquisition Process," on page 59).

The other purpose of the HAL board is to hold the current NMR data in its RAM. For arrayed and multidimensional experiments, there is sufficient RAM to hold both the current and the next FID. After completion of the scans for an FID, this permits continuing with the next data set without delay ("buffered acquisition").

• The main task of the **acquisition CPU** is the interpretation of the Acode, generating the information that is fed into the pulse programmer. With the aid of the programmer, it also perform the tasks of autolocking and autoshimming (the spinner speed regulation is done by dedicated peripheral circuitry).

The acquisition CPU is equipped with boot PROMs, containing a primitive operating system that is running after switching on or resetting the CPU, up until the real acquisition operating system is loaded (see Section 6.1, "Starting the Acquisition Operating System," on page 59). The acquisition operating system fills most of its RAM, which during acquisition also holds the Acode, the data, and instructions needed for performing an experiment (see Chapter 8, "Acquisition CPU and Acode," on page 69).

The acquisition CPU sends information to the pulse programmer and exchanges information with the automation board through an on-board parallel I/O channel (see Chapter 13, "Acquisition CPU Communication," on page 145).

The acquisition CPU has two serial ports, one of which can be used for a diagnostics monitor (either to display diagnostic information during NMR experiments or to run PROM-based on-board diagnostics software). In UNITY and earlier spectrometers, the other serial port was used to exchange information with the VT controller (the Oxford VTC-4); in UNITY*plus* spectrometers, this is done via automation board.

• In a UNITY *plus* spectrometer, the **automation board** drives five serial ports, through which information is exchanged with the sample changer (ASM-100 or SMS), the magnet leg pneumatics control circuitry (eject/insert, slow drop, and bearing air flows), the spinner control circuitry, and the VT controller (one serial port is currently unused). In UNITY and earlier spectrometers, the second serial

port of the acquisition CPU was used for the VT controller. Instead, on these systems the automation board also drove the lock power, phase and gain controls, and the receiver gain; all these functions are now addressed via the AP bus (see Chapter 12, "AP Bus Traffic," on page 137).

In a UNITY*plus,* the automation board also contains the bootup mode selector and some battery-buffered RAM ("zero-power RAM") to store the rf configuration.

- The main task of the **pulse programmer** is to accurately control the timing during an NMR experiment. It also acts as a buffer for timing and control information; it directly controls most of the rf gates and has a dedicated output channel (the AP bus, see Chapter 12, "AP Bus Traffic," on page 137) through which it can send information to most parts of the spectrometer. It controls the shim gradients and all devices that set frequency offsets, attenuations and pulse amplitudes, phase shifts and pulsed field gradient amplitudes (see also Chapter 9, "Pulse Programmers," on page 85).
- The ADC (Analog-to-Digital Converter) board receives two audio signals (0 and 90 degrees) from the receiver board (via an audio filter). Each of the channels is fed into sample-and-hold and ADC circuitry. The data sampling is triggered by the pulse programmer. If necessary, the data are also scaled down on the ADC board. At its output, there is a FIFO (first-in-first-out) buffer, from which the data are fed directly into the sum-to-memory board.
- The STM (Sum-to-Memory) board reads in the current FID from the HAL board through DMA and adds in the new FID from the FIFO buffer on the ADC board, checks for mathematical overflow (at which point the FID will be scaled by another bit), and stores the result back in the HAL memory. The STM board does complex additions according to the current receiver phase (see Chapter 18, "Acquiring Data," on page 205).

7.2 Bus Structures

The main components in the digital part of the spectrometer are linked with different bus structures that carry the data traffic:

- SCSI (Small Computer System Interface) bus is an 8-bit parallel bus that links the HAL board, the host computer, and most host computer peripherals (disks, tapes, CD-ROM, etc.).
- Acquisition CPU bus is a 16-bit parallel bus structure that is integrated in the cardcage backplane (VERSAbus). It carries the data traffic between the acquisition CPU, the HAL board, and the STM board. It also provides dc power to the other boards in the same card cage (e.g., the automation board).
- **I/O bus** is an 8-bit parallel bus that is driven by the acquisition CPU. This bus connects the acquisition CPU, automation board, and pulse programmer.
- AP (analog port) bus is the most important link between the pulse programmer and most of the rf devices. Apart from the gating lines, the AP bus carries information on frequency offsets, attenuation and power modulation levels, pulsed field gradient amplitudes, small-angle phase shifts—all the numeric information for the rf part in general. The AP bus is also used to transfer shapes and patterns to the waveform generators. On the UNITY*plus*, lock power, lock phase, lock gain, and the receiver gain are set over the AP bus (on UNITY and earlier spectrometers,

this was done with the automation board). The AP bus is 16-bit parallel (see also Chapter 12, "AP Bus Traffic," on page 137).

 On UNITY plus systems, the gating information ("fast lines") is routed in its own bus, the HS (high-speed) bus (see Chapter 9, "Pulse Programmers," on page 85).

The components shown in Figure 17 do not cover the entire digital part of the spectrometer, but mainly the computing part. Especially on more recent instruments, such as the UNITY*plus*, digital components reach much further into the spectrometer:

- The waveform generators (see Chapter 16, "Waveform Generators," on page 163) are completely digital.
- Each transmitter board consists of two parts: the rf board and a digital control board (see Section 16.1, "How Does a Waveform Generator Fit Into the System?," on page 163 and Section 4.2, "How Do Pulses Work?," on page 40).

Because virtually everything in the spectrometer is digitally controlled, digital components are found on every subunit of the spectrometer (see also Section 12.2, "What Devices are Driven by the AP Bus?," on page 139).

Chapter 8. Acquisition CPU and Acode

After the acquisition CPU is "running" (i.e., once the acquisition operating system has been uploaded and started) and "active" (performing an experiment), it in essence holds the following data blocks:

- Acquisition operating system (see Section 6.1, "Starting the Acquisition Operating System," on page 59), including the Acode interpreter.
- Current FID.
- FID-specific data and instructions (Acode).

8.1 CPU Address Space

The operating system and the Acode are located in the RAM of the acquisition CPU; the FID is stored in the RAM of the HAL board. For the rest of this chapter, we will discuss the Acode alone. With respect to the acquisition CPU and the HAL board, you don't have to be concerned with anything but the data in the address space of the acquisition CPU and the way these data are interpreted and used.

The CPU address space also includes status registers on the STM (sum-to-memory) board. This enables communicating with a board driven by firmware ("hard-coded software") and permits transferring acquisition information (number of points, observe phase) to the STM board (see also Chapter 18, "Acquiring Data," on page 205).

8.2 Looking at Acode

Acode (acqi.Code from go('acqi') or expn.username.PID.Code for normal acquisitions, both in /vnmr/acqqueue, see also the previous chapters) is primarily a binary file, mostly consisting of 16-bit integers, interspersed with bit patterns (organized in 16-bit binary words) and some 32-bit long integers. As 16-bit units dominate (and everything is organized in 2- or 4-byte units), the od -i command (or od -s under Solaris 2.3) is an almost adequate way to look at this file:

| UNITY400:vnmr1 - 1> od -i /vnmr/acqqueue/acqi.Code | | | | | | | | |
|--|--------|--------|--------|--------|--------|------|--------|--------|
| 0000000 | 0 | 1 | 0 | б | 0 | 421 | 0 | -28096 |
| 0000020 | 0 | 16 | 0 | 0 | 0 | 0 | 0 | 0 |
| 0000040 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| 0000060 | 0 | 0 | 0 | 1 | 0 | 0 | 268 | 94 |
| 0000100 | 1 | 293 | 3 | 106 | 106 | 85 | 4 | 0 |
| 0000120 | 0 | 0 | 0 | 0 | 64 | 0 | 16 | 16 |
| 0000140 | 0 | 0 | 0 | 0 | 512 | 1024 | 1536 | 0 |
| 0000160 | 4096 | 8192 | 12288 | 0 | 0 | 1 | 0 | 0 |
| 0000200 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| 0000220 | 0 | 0 | 0 | 0 | 63 | 0 | 0 | 0 |
| 0000240 | 0 | 1 | 2 | 3 | 0 | 0 | 0 | 0 |
| 0000260 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| 0000300 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 1542 |
| 0000320 | 0 | 0 | 0 | 0 | 20 | 31 | 59 | 1 |
| 0000340 | 15 | 1 | 0 | 0 | 0 | 8 | 0 | 16 |
| 0000360 | 0 | 24 | 0 | 24 | 11 | 148 | 15 | 2 |
| 0000400 | 0 | 0 | 0 | 256 | 0 | 512 | 0 | 768 |
| 0000420 | 0 | 768 | 11 | 164 | 53 | 0 | 1 | 19 |
| 0000440 | 5 | 157 | 1 | 4 | 6 | б | -21696 | -17651 |
| 0000460 | -25856 | -21676 | -17428 | -25809 | -25713 | 6 | 8 | -22752 |

| 0000500 | -18688 | -26690 | -18449 | -18514 | -18437 | -18498 | -26879 | -18688 |
|---------|--------|--------|--------|--------|--------|--------|--------|--------|
| 0000520 | 159 | 1 | 8 | -22752 | -18688 | -26690 | -18449 | -18514 |
| 0000540 | -18437 | -18498 | -26879 | -18688 | 0 | 16 | 257 | 0 |
| 0000560 | 68 | 1 | 360 | 65 | 769 | 0 | 59 | 2867 |
| 0000600 | 1025 | 79 | 16 | 55 | 59 | 2966 | 514 | 4095 |
| 0000620 | 0 | 4095 | 6 | 3 | -21608 | -17454 | -25856 | -25856 |
| 0000640 | 6 | 1 | -21614 | -17504 | б | 1 | -21605 | -17664 |
| 0000660 | 6 | 1 | -21615 | -17536 | б | 1 | -21616 | -17647 |
| 0000700 | 6 | 7 | -32255 | -32236 | -32224 | -32207 | -32191 | -32171 |
| 0000720 | -32156 | -32143 | 16 | 258 | 0 | 59 | 2866 | 1025 |
| 0000740 | 79 | 16 | 30 | 59 | 2982 | 514 | 4095 | 0 |
| 0000760 | 4095 | 6 | 3 | -21592 | -17622 | -25856 | -25856 | 6 |
| 0001000 | 1 | -21598 | -17472 | 6 | 1 | -21589 | -17664 | 6 |
| 0001020 | 1 | -21599 | -17536 | 6 | 1 | -21600 | -17646 | 6 |
| 0001040 | 1 | -21688 | -17663 | 6 | 1 | -21708 | -17579 | 6 |
| 0001060 | 1 | -21707 | -17664 | 6 | 1 | -21683 | -17631 | 6 |
| 0001100 | 1 | -21693 | -17663 | 6 | 1 | -21687 | -17664 | 6 |
| 0001120 | 1 | -21706 | -17658 | 6 | 1 | -21605 | -17528 | 6 |
| 0001140 | 1 | -21589 | -17528 | 151 | 4129 | 71 | 63 | 2 |
| 0001160 | 440 | 30000 | 250 | 0 | 151 | 12297 | 7 | 9 |
| 0001200 | 95 | 4 | 151 | 12787 | 7 | 9 | 75 | 1 |
| 0001220 | 7 | 1792 | 15363 | 0 | 0 | 5229 | 27904 | 0 |
| 0001240 | 43 | 1 | 0 | 256 | 8390 | 4123 | 19 | 5 |
| 0001260 | 98 | 150 | 0 | 0 | 6 | 1 | -21605 | -17528 |
| 0001300 | 150 | 0 | 0 | 6 | 1 | -21589 | -17528 | 150 |
| 0001320 | 0 | 0 | 150 | 0 | 0 | 150 | 0 | 0 |
| 0001340 | 150 | 0 | 0 | 16 | 1 | 74 | 150 | 0 |
| 0001360 | 1 | 150 | 0 | 1 | 151 | 4489 | 150 | 0 |
| 0001400 | 3 | 151 | 4369 | 150 | 0 | 1 | 151 | 4489 |
| 0001420 | 150 | 0 | 0 | 150 | 0 | 0 | 39 | 74 |
| 0001440 | 39 | 16 | 1 | 74 | 16 | 2 | 74 | 152 |
| 0001460 | 8487 | 4117 | 90 | 99 | 1 | 0 | -28096 | 8390 |
| 0001500 | 4123 | 7 | 97 | 20 | 232 | | | |
| 0001512 | | | | | | | | |

This Acode is for a simple experiment (s2pul, p1=0, d1=0, d2=0) and obtained with the gf macro ("real" Acode for arrayed or multidimensional experiments can be *very* long, see Section 5.2, "Tasks for the Pulse Sequence Executable," on page 57). Looking at the above code, we can already (roughly) distinguish two parts: the top part that consists mainly of "simple numbers" and the bottom part that has a lot of "funny numbers" (indicating that there are things other than 16- or 32-bit integers).

The -i option causes od to interpret the file given in the argument as 16-bit signed integers; bit patterns that have the most significant bit set, large unsigned 16-bit integers, negative or very large 32-bit integers result in large negative numbers being shown. od -1 would interpret most 32-bit integers correctly (those which start at an odd address), but would misinterpret most of the file by taking two 16-bit numbers for a single 32-bit number; bit pattern would again mostly be shown as very large numbers. Other options for od (for octal or hexadecimal output) make the result even less interpretable. Clearly, a tool is needed that properly takes all the bits and pieces of the Acode apart.

Methods of Interpreting the Contents of Acode Files

The information on the contents of the Acode is found (somehow) in /vnmr/psg, although not in a single place or file, but rather convoluted and spread over many source and header files. In general it is *not* necessary to know about the contents of the Acode, although for people who do low-level pulse sequence programming (or modify files in /vnmr/psg), an Acode decoder can sometimes provide extremely valuable debugging information—especially because the dps command does not interpret certain low-level commands and does not show phase tables, real-time phase math, and even some of the "normal" statements, such as offset.

• In the remaining parts of this manual, Acode contents are shown as output by a program that has been submitted to the user library (userlib/bin/apdecode). This program shows Acode structures in detail (all contents); the instruction part is decoded in that only the instruction codes themselves are shown, whereas arguments to code functions are (partly) decoded.

apdecode is certainly not a perfect solution because it doesn't try decoding everything; it will also not work with all possible rf configurations. Also, because the Acode undergoes extensive changes with every VNMR release, apdecode is undergoing constant adaptations, and it, in general, only works with the VNMR release for which it has been written. The main purpose for showing apdecode output in this manual is to illustrate specific contents of the Acode and to show the effect that specific functions (phase calculations, phase tables) have on Acode.

An alternative option for looking at Acode (and related things) is the 'debug' option to the go command, which will display extensive (debugging) information in the window from which VNMR was called (the parent window). For anyone who is not an expert, go('debug') cannot be recommended because its output is much too detailed to be of real value. The output of apdecode is much more compact (but may be incomplete at times).

• The output of the Acode decoder apdecode on the same Acode as shown above looks as follows:

Number of traces: 1 Code Start: 6 Code Ends: 421

| aaddr f | addr 1 | addr v | value | comment |
|---------|--------|--------|-------|--|
| | | | | ############# LC (Low Core) Structure #################################### |
| 6 | 0 | 03 | 37440 | LC->np (long) |
| 8 | 2 | 2 | | LC->nt (long) |
| 10 | 4 | 4 | 0 | LC->ct (long) |
| 12 | б | б | 0 | |
| 14 | 8 | 8 | 0 | LC->rsum (long) |
| 16 | 10 | 10 | 0 | LC->dpts (long) |
| 18 | 12 | 12 | 0 | LC->autop (long) |
| 20 | 14 | 14 | 0 | LC->stmar (long) |
| 22 | 16 | 16 | 0 | LC->stmcr (long) |
| 24 | 18 | 18 | 0 | LC->rtvptr (long) |
| 26 | 20 | 20 | 1 | LC->elemid (long) |
| 28 | 22 | 22 | 0 | LC->squi (long) |
| 30 | 24 | 24 | 268 | LC->idver |
| 31 | 25 | 25 | 94 | LC->o2auto (Offset to AUTOD) |
| 32 | 26 | 26 | 1 | LC->ctctr |
| 33 | 27 | 27 | 293 | LC->dsize (blocks) |
| 34 | 28 | 28 | 3 | LC->asize (blocks) |
| 35 | 29 | 29 | 106 | LC->codeb (Offset to Code) |
| 36 | 30 | 30 | 106 | LC->codep |
| 37 | 31 | 31 | 85 | LC->status |
| 38 | 32 | 32 | 4 | LC->dpf (0=int, 4=long) |
| 39 | 33 | 33 | 0 | LC->maxscale |
| 40 | 34 | 34 | 0 | LC->icmode |
| 41 | 35 | 35 | 0 | LC->stmchk |
| 42 | 36 | 36 | 0 | LC->nflag |
| 43 | 37 | 37 | 0 | LC->scale |
| 44 | 38 | 38 | 64 | LC->check |
| 45 | 39 | 39 | 0 | LC->oph |
| 46 | 40 | 40 | 16 | LC->bsval |
| 47 | 41 | 41 | 16 | LC->bsctr |
| 48 | 42 | 42 | 0 | LC->ssval |
| 49 | 43 | 43 | 0 | LC->ssctr |
| 50 | 44 | 44 | 0 | LC->ctcom |
| 51 | 45 | 45 | | LC->dptab 0 0x200 0x400 0x600 |
| 55 | 49 | 49 | | LC->obsptb 0 0x1000 0x2000 0x3000 |

| 59 | 53 | 53 | 0 | LC->rfphpt |
|----------|-------|----|------|---|
| 60 | 54 | 54 | 0 | LC->curdec (unused) |
| 61 | 55 | 55 | 1 | LC->cpf (cycle phase flag) |
| 62 | 56 | 56 | 0 | LC->maxconst |
| 63 | 57 | 57 | 0 | LC->tablert |
| 64 | 58 | 58 | 0 | LC->output card status reg |
| 65 | 59 | 59 | 0 | LC->analog port status reg |
| 66 | 60 | 60 | 0 | LC->analog port data reg |
| 67 | 61 | 61 | 0 | LC->analog port address/control reg |
| 68 | 62 | 62 | 0 | LC->input card status reg |
| 69 | 63 | 63 | 0 | LC->input card data req |
| 70 | 64 | 64 | 0 | LC->input card ocsr reg |
| 71 | 65 | 65 | 0 | LC->stm card status req |
| 72 | 66 | 66 | 0 | LC->tpwrr (xmtr power) |
| 73 | 67 | 67 | 0 | LC->dpwrr (dec power) |
| 73 | 68 | 68 | 0 | |
| | | | | LC->tphsr (xmtr phase shift) |
| 75 | 69 | 69 | 0 | LC->dphsr (dec phase shift) |
| 76 | 70 | 70 | 63 | LC->dlvlr (dec level) |
| 77 | 71 | 71 | 0 | LC->srate (High Speed Rotor Freq) |
| 78 | 72 | 72 | 0 | LC->rttmp (temp real time->interlock) |
| 79 | 73 | 73 | 0 | LC->spare1 (unused) |
| 80 | 74 | 74 | 0 | LC->zero |
| 81 | 75 | 75 | 1 | LC->one |
| 82 | 76 | 76 | 2 | LC->two |
| 83 | 77 | 77 | 3 | LC->three |
| 84 | 78 | 78 | 0 | LC->v1 |
| 85 | 79 | 79 | 0 | LC->v2 |
| 86 | 80 | 80 | 0 | LC->v3 |
| 87 | 81 | 81 | 0 | LC->v4 |
| 88 | 82 | 82 | 0 | LC->v5 |
| 89 | 83 | 83 | 0 | LC->v6 |
| 90 | 84 | 84 | 0 | LC->v7 |
| 91 | 85 | 85 | 0 | LC->v8 |
| 91 92 | 86 | 86 | 0 | LC->v9 |
| | | | | |
| 93 | 87 | 87 | 0 | LC->v10 |
| 94 | 88 | 88 | 0 | LC->v11 |
| 95 | 89 | 89 | 0 | LC->v12 |
| 96 | 90 | 90 | 0 | LC->v13 |
| 97 | 91 | 91 | 0 | LC->v14 |
| 98 | 92 | 92 | 0 | LC->(v15) |
| 99 | 93 | 93 | 0 | LC->(v16) |
| | | | | |
| | | | | ############# AUTOD: Automation Data ################################## |
| 100 | 94 | 0 | 0 | AUTOD->checkmask (long) |
| 102 | 96 | 2 | 0 | AUTOD->when_mask: load='n' wshim='n' |
| 103 | 97 | 3 | 1542 | AUTOD->control_mask |
| 104 | 98 | 4 | 0 | AUTOD->best |
| 105 | 99 | 5 | 0 | AUTOD->loops |
| 106 | 100 | 6 | 0 | AUTOD->sample_mask (tray location) |
| 107 | 101 | 7 | 0 | AUTOD->sample_error |
| 108 | 102 | 8 | 20 | AUTOD->recgain |
| 109 | 103 | 9 | 31 | AUTOD->lockpower |
| 110 | 104 | 10 | 59 | AUTOD->lockgain |
| 111 | 105 | 11 | 1 | AUTOD->lockphase |
| | 105 | | - | Norob / rockphase |
| | | | | ############## Instruction Section #################################### |
| 112 | 106 | 0 | 1 5 | |
| 112 | 106 | 0 | 15 | |
| 100 | 1.0.0 | | | allBits = 0x18, addr = 0xb0094 |
| 126 | 120 | 14 | 15 | SETPHATTRibutes CH2 0 0x100 0x200 0x300 |
| | | | | allBits = 0x300, addr = 0xb00a4 |
| 140 | 134 | 28 | 53 | ACQBITMASK 0 |
| 142 | 136 | 30 | 1 | CBEGIN |
| 143 | 137 | 31 | 19 | INITialize acq |
| 144 | 138 | 32 | 5 | CLEAR |
| 145 | 139 | 33 | 157 | LOCKFILTER fast = 1, slow = 4 |
| 148 | 142 | 36 | 6 | APBOUT 7 items 0xab40 0xbb0d 0x9b00 0xab54 0xbbec |
| | | | | 0x9b2f 0x9b8f |
| 157 | 151 | 45 | 6 | APBOUT 9 items 0xa720 0xb700 0x97be 0xb7ef 0xb7ae |
| | | | č | 0xb7fb 0xb7be 0x9701 0xb700 |
| 168 | 162 | 56 | 159 | TUNE_FREQ CH1 9 words |
| 200 | | 50 | | 0xa720 0xb700 0x97be 0xb7ef 0xb7ae 0xb7fb 0xb7be |
| | | | | 0x9701 0xb700 |
| | | | | 0127,01 0AD100 |

| 180 | 174 | 68 | 0 | NO_OP |
|-------|------------|-----|-----|--|
| 181 | 175 | 69 | 16 | SETPHAS90 CH1c 0 |
| 184 | 178 | 72 | 68 | PHASESTEP CH1 360 units (90.00 degrees) |
| 187 | 181 | 75 | 65 | SETPHASE CH1f 0 |
| 190 | 184 | 78 | 59 | APChipOUT APaddr 11, reg 51, +logic, 1 byte |
| | | | | max 79, offset 16, value 55 |
| 196 | 190 | 84 | 59 | APChipOUT APaddr 11, reg 150, -logic, 2 bytes |
| 100 | 100 | 04 | 55 | max 4095, offset 0, value 4095 |
| 202 | 106 | 0.0 | c | APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 |
| 202 | 196 | 90 | 6 | |
| 208 | 202 | 96 | 6 | APBOUT 2 items 0xab92 0xbba0 |
| 212 | 206 | 100 | 6 | APBOUT 2 items 0xab9b 0xbb00 |
| 216 | 210 | 104 | 6 | |
| 220 | 214 | 108 | 6 | APBOUT 2 items 0xab90 0xbb11 |
| 224 | 218 | 112 | б | APBOUT 8 items 0x8201 0x8214 0x8220 0x8231 0x8241 |
| | | | | 0x8255 0x8264 0x8271 |
| 234 | 228 | 122 | 16 | SETPHAS90 CH2c 0 |
| 237 | 231 | 125 | 59 | APChipOUT APaddr 11, reg 50, +logic, 1 byte |
| | | | | max 79, offset 16, value 30 |
| 243 | 237 | 131 | 59 | |
| 245 | 257 | 191 | 55 | max 4095, offset 0, value 4095 |
| 0.4.0 | 0.4.2 | 100 | ~ | |
| 249 | 243 | 137 | 6 | APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 |
| 255 | 249 | 143 | 6 | APBOUT 2 items 0xaba2 0xbbc0 |
| 259 | 253 | 147 | 6 | APBOUT 2 items 0xabab 0xbb00 |
| 263 | 257 | 151 | 6 | APBOUT 2 items 0xabal 0xbb80 |
| 267 | 261 | 155 | 6 | APBOUT 2 items 0xaba0 0xbb12 |
| 271 | 265 | 159 | б | APBOUT 2 items 0xab48 0xbb01 |
| 275 | 269 | 163 | 6 | APBOUT 2 items 0xab34 0xbb55 |
| 279 | 273 | 167 | 6 | |
| 283 | 277 | 171 | | APBOUT 2 items 0xab4d 0xbb21 |
| | 281 | 175 | | APBOUT 2 items 0xab43 0xbb01 |
| 287 | | | | APBOUT 2 items 0xab43 0xbb01 APBOUT 2 items 0xab49 0xbb00 |
| | 285 | 179 | 6 | |
| | | 183 | 6 | |
| 299 | 293 | 187 | | APBOUT 2 items 0xab9b 0xbb88 |
| 303 | 297 | 191 | 6 | |
| 307 | 301 | 195 | 151 | EVENT1_TWRD 1.000 usec |
| 309 | 303 | 197 | 71 | GAINAutomation |
| 310 | 304 | 198 | 63 | SETVT Oxford PID 440, temp 3000.0, vtc 25.0 |
| 315 | 309 | 203 | 0 | NO_OP |
| 316 | 310 | 204 | 151 | |
| 318 | 312 | 206 | 7 | STartFIFO |
| 319 | 313 | 207 | 9 | StopFIFO |
| | 314 | | | - |
| | | 208 | 95 | PADelay 4 words |
| 322 | 316 | | 151 | EVENT1_TWRD 500 msec |
| | 318 | 212 | | STartFIFO |
| 325 | 319 | 213 | 9 | StopFIFO |
| 326 | 320 | 214 | 75 | SHIMAutomation mode = 1, 7 words |
| | | | | 0x700 0x3c03 0 0 0x146d 0x6d00 0 |
| 341 | 335 | 229 | 43 | NOISE loop 256 pts, dwell 199 usec + 1.000 usec |
| 342 | 336 | 230 | 19 | INITialize acq |
| 343 | 337 | 231 | 5 | CLEAR |
| 344 | 338 | 232 | 98 | NextSCan |
| 345 | 339 | 233 | 150 | HighSpeedLINES (void) |
| 348 | 342 | 236 | 6 | APBOUT 2 items 0xab9b 0xbb88 |
| 352 | 346 | | 150 | HighSpeedLINES (void) |
| | | 240 | | |
| 355 | 349 | 243 | 6 | |
| 359 | 353 | 247 | 150 | HighSpeedLINES (void) |
| 362 | 356 | 250 | 150 | HighSpeedLINES (void) |
| 365 | 359 | 253 | 150 | HighSpeedLINES (void) |
| 368 | 362 | 256 | 150 | HighSpeedLINES (void) |
| 371 | 365 | 259 | 16 | SETPHAS90 CH1 zero |
| 374 | 368 | 262 | 150 | HighSpeedLINES RXOFF |
| 377 | 371 | 265 | 150 | HighSpeedLINES RXOFF |
| 380 | 374 | 268 | 151 | EVENT1_TWRD 10.000 usec |
| 382 | 376 | 270 | 150 | HighSpeedLINES RXOFF TXON |
| 385 | 379 | 273 | 151 | EVENT1_TWRD 7.000 usec |
| | | | | — |
| 387 | 381 | 275 | 150 | HighSpeedLINES RXOFF |
| 390 | 384 | 278 | 151 | EVENT1_TWRD 10.000 usec |
| 392 | 386 | 280 | 150 | HighSpeedLINES (void) |
| 395 | 389 | 283 | 150 | HighSpeedLINES (void) |
| 398 | 392 | 286 | 39 | ASSIGNFUNC zero oph |
| 401 | 205 | 200 | 16 | SETPHAS90 CH1 zero |
| | 395 | 289 | | |
| 404 | 395 398 | 289 | 16 | SETPHAS90 CH2 zero |

407 401 295 152 EVENT2_TWRD 296 usec + 850 nsec 410 404 298 90 SETInputCardMode 414 408 302 99 ACQXX loop np=37440, dwell 199 usec + 1.000 usec 417 411 305 7 STartFIFO 418 412 306 97 HouseKEEPing 419 413 307 20 BRANCH Offset 232 Total code size = 421 words / 842 Bytes / 0.8 KB

8.3 Structure of Acode Files

The Acode consists of four parts:

- File header
- Data structure LC (low-core)
- Second data structure AUTOD (automation data)
- Instruction section

In the case of arrayed and multidimensional experiments, the Acode contains one file header plus three other parts: the LC data structure, AUTOD data structure, and an instruction section *per FID*. The LC and AUTOD data structures are always the same (within one software release at least): rigid assemblies of 16- and 32-bit numbers according to a predefined scheme. The instruction section is variable in length and contents, and varies considerably with the hardware configuration. The last part of the instruction sections depends on the current pulse sequence. As mentioned before, the first instruction segment is always much longer than all others, because of all the initialization instructions.

Let us now have a more detailed look at the different parts of the Acode. No attempt will be made to explain all parts of the Acode—in particular, only those parts of the data structures are explained that are relevant to the following chapters or which are needed to explain the Acode by itself.

Acode File Header

The Acode has a variable-size file header, consisting of 32-bit (long) integers. The first number describes the *number of code segments* (corresponds to the value of the arraydim, except for go('acqi') or gf, which always produces one code segment only). In the above example (obtained with gf), there is one code segment only.

The next (32-bit integer) number is the *offset to the first code segment* as 16-bit (2-byte) address. In this case (one code segment), the header is three 32-bit words long, resulting in an offset of 6 for the first segment (the first number of the header is located at address 0).

The following 32-bit numbers are the number of 16-bit words up to (and including) the last code word for every code segment, or the *addresses of the next code segments*, if any more segments follow. For a standard acquisition (not in Acode generated with go('acqi') or gf), there are arraydim code segments; therefore, the header is 2*(arraydim + 2) words (16 bits or 2 bytes each) long, or 4*(arraydim + 2) bytes.

The header information is used by Acqproc to extract individual code segments, because in arrayed or multidimensional experiments the code is transferred to the acquisition CPU segment by segment during the acquisition.

LC Data Structure

For the rest of the Acode, the program apdecode lists three different offsets (in 16bit words): an absolute word count throughout the entire Acode file (aaddr), a word count per code segment (faddr), and a local word count or address (laddr)—this kind of code offset is used by the Acode interpreter itself: addresses are either offsets within the LC or AUTOD structures, or within the instruction segment (see below).

As explained previously, the term LC stands for "low core." This terminology is historic and has little meaning with respect the memory organization in the acquisition CPU. Both structures (AUTOD and LC) are stored in the acquisition memory the way they are stored in /vnmr/acqqueue, and the same is done with the instruction section, but the three sections are *not* stored in consecutive memory locations (in particular, the instruction section is *not* stored after the two data structures). After they are loaded into their proper memory locations, the two data structures form the "workspace" for the instruction section. Such a workspace could also be "built" by the acquisition software, but there is of course a good reason for creating an image of that workspace together with the instruction section as part of the Acode: the two data structures also serve as container for numeric information that is transferred between the two computers. This becomes obvious when we look at some of the contents of LC (only selected parts are explained here).

LC contains many items that sound very familiar to the pulse sequence programmer: it turns out that all real-time variables (ct, v1, one, etc., see below) are actually (16-bit) *addresses* to some specific word (16-bit integer) within the LC structure. Because they are addresses and not normal integers, a special type codeint was created, describing C variables that contain addresses (code offsets) within the LC structure. As we see later, by using those codeint variables, the pulse sequence program is able to generate Acode that contains these addresses. Based on that, during the acquisition the Acode interpreter can do real-time calculations and operations using the addressed memory locations.

LC starts with a series of long (32-bit and two 16-bit word) integers, containing:

- LC->np, the number of points to acquire, corresponding to the np parameter in VNMR; LC->np is the C syntax for the nt element of a structure named LC—we use this syntax to differentiate the structure element from the np variable in C and the np parameter in VNMR.
- LC->nt, the number of transients to acquire, corresponding to the nt parameter.
- The number of acquired transients, corresponding to the ct parameter in VNMR; as we will see below, all so-called real-time variables in a pulse sequence are 16-bit integers, and all real-time math operations are 16-bit operations; the codeint variable ct in a pulse sequence therefore is defined as the "low-order" half of the LC->ct structure element: ct is first set to the offset to LC->ct and then incremented by one (this is found in /vnmr/psg/psg.c); all code offsets are in (16-bit) words, *not in bytes*. LC->ct is incremented *after* every scan (i.e., during the first scan ct and LC->ct is 0, and the acquisition for one FID is finished when the two long integers LC->ct and LC->nt are the same).

- The index of the current FID and Acode segment (in C pulse sequence code corresponding to the variable ix in the C code of a pulse sequence; the index for the first FID is 1, not 0). LC->ix is an unsigned, 32-bit integer; therefore, the theoretical maximum number of elements is 4,294,967,296. This may sound excessive, but earlier releases of VNMR had LC->ix defined as a 16-bit integer, which limited the number of elements to 32767, a value that is easily exceeded in 3D or 4D experiments.
- "squi", the current set of "quiescent" states: in this location the acquisition software stores (at real time) the set of gating information (the "fast bits") to be used with the next pulse programmer event, and to (or from) which new gate setting is added and subtracted.

The remainder of LC consists of 16-bit words and integers, among others including:

- LC->o2auto, the offset to the AUTOD structure (within the code segment, in 16bit words), in other words: the size of the LC structure in 16-bit words.
- LC->codeb, the offset to the instruction section within the code segment (faddr in the apdecode output) or the size of the sum of the LC and AUTOD structures.
- LC->codep, a pointer to the current code word.
- LC->dpf, the "double-precision flag"—an integer that contains either 0 (for 16bit acquisitions) or 4 (for 32-bit acquisitions, dp='y').
- LC->dsize, the data size in 512-byte blocks (for dp='y': 4*np/512, rounded up to the next full block).
- LC->asize, the Acode size in 512-byte blocks (minimum: 3 blocks).
- LC->scale, the number of binary scaling operations during an acquisition: whenever the maximum or minimum number in an FID would exceed the numeric range at the selected precision, the ADC output and the current FID are scaled down by a factor of 2 (this corresponds to a right-shift in the binary numbers). For 16-bit acquisitions (dp='n'), the range is from -32768 to +32767; for 32-bit acquisitions (dp='y'), it is from -2,147,483,648 to +2,147,483,647.
- LC->maxscale, the maximum number of permissible scaling operations (rightshifts): before starting any acquisition, the system measures the amount of noise by acquiring 256 data points at the conditions of the experiment. From that noise the maximum number of scaling operations is determined, such that after the scaling the noise is still properly digitized. If at that point more scaling would be required, the acquisition is stopped with the message "maximum number of transients accumulated".
- LC->cpf, the "cycle phase flag": 0 stands for cp='n', 1 for cp='y' in VNMR; if cp='n', the observe phase oph remains constant; otherwise, it is incremented with the ct counter.
- LC->ssval, the number of steady-state pulses, corresponding to the ss parameter in VNMR. A codeint (real-time) variable ssval contains the *address* (the offset) to LC->ssval, such that in a pulse sequence we can instruct the acquisition CPU to perform (mathematical) operations based on the value stored in LC->ssval.
- LC->ssctr, the counter location, in which the acquisition CPU counts the steady-state pulses; this location is initially set to the value of LC->ssval and the

decremented after each transient. During "real" transients, LC->ssctr remains set to zero. Also here (equivalent to ssval), a codeint variable ssctr exists that allows performing real-time calculations based on the value of LC->ssctr.

- LC->bsval, the block size, corresponding to the bs parameter in VNMR. A codeint (real-time) constant bsval contains the offset to LC->bsval, such that in a pulse sequence we can instruct the acquisition CPU to perform (mathematical) operations based on the value stored in LC->bsval.
- LC->bsctr, the counter location, in which the acquisition CPU counts through the block size; at the beginning of the acquisition this location is set to the value of LC->ssval and is decremented after each transient; whenever the location LC->bsctr contains zero, the FID is stored, and LC->bsctr is reset to the value of LC->bsval. Again, a codeint variable (bsctr) exists that allows performing real-time calculations based on the value of LC->bsctr, or even recalculating the value of LC->bsctr (e.g., for implementing dynamic block sizes).
- LC->oph, the observe phase real time variable: oph in a pulse sequence is nothing but the address (the local offset, laddr in the apdecode output) to that structure element. If cp='n' (LC->cpf is 0), LC->oph contains 0, otherwise (cp='y', and LC->cpf is 1) LC->oph is the same as ct (the low-order half of LC->ct). As the observe phase can only assume four values (0, 1, 2, 3), it doesn't matter if oph contains the values from ct that are incremented up to nt-1: the software simply looks at the last two (the two least-significant) bits of that 16-bit number, which is identical to a *modulo 4* function, resulting in a sequence 0, 1, 2, 3, 0, 1, 2, 3, etc. Of course, oph can also be calculated (or set from a phase table, see the following chapters); in any case, it is *not* necessary to perform a modulo function to ensure values in the proper range of 0 to 3—for the modulo (mod 4) function is implicit!
- LC->v1, LC->v2,...LC->v14, 14 storage locations for results of phase calculations, flags, etc.: in a pulse sequence these locations can be addressed via the codeint ("real-time") variables v1, v2, ... v14. Since real-time math happens detached from the C compiler and the execution of the C pulse sequence, the only way to instruct the acquisition CPU to perform math with specified variables is to the addresses on which mathematical and logical operations are to be performed. There are two kinds of such addresses: addresses to the variables v1, v2, ... v14, oph, ssctr, and bsctr (the latter two are rarely used as variables), and addresses to constants (see below) and "system variables" like ct (as well as ssval, ssctr, bsval, bsctr, tablert). The integers LC->v1 up to LC->v14 are initialized with zero value, unless the user specifies a different value using the initval statement. This also is the only way to set the contents of these locations directly from within C; once the experiment is running, the values can and will only be changed by real-time math.
- LC->zero, LC->one, LC->two, LC->three, four locations with the numeric values 0, 1, 2, and 3: these numbers are accessible via the real-time (codeint) variables zero, one, two, and three. As numeric C constants (or variables) cannot exist in the Acode translation, these four frequently used numeric values are stored in real-time locations and made accessible by the address variables zero, one, two, three. They all could theoretically constructed from a single value (zero or one), but it is more convenient to have at least these three numbers available for real-time math.

• LC->tablert, this is the location in which phase values from a table are stored in table calls (e.g., the pulse(pw,t1) statement), unless the table value is extracted into a real-time variable (v1, v2, ... v14, using the getelem statement, or into oph with the setreceiver statement, see also Chapter 11, "Phase Tables," on page 115).

The other elements (addresses) in the LC structure are used internally within the acquisition CPU and are rarely of interest to the user, not even for debugging.

The LC structure (together with AUTOD) is defined in /vnmr/psg/lc.h. It is initialized in /vnmr/psg/psg.c, which also defines and sets the corresponding codeint (real-time) variables. The LC and AUTOD structures are to be regarded as non-user modifiable, because there are matching counterparts within the VNMR module and in the acquisition operating system.

The AUTOD Data Structure

In the current software, the AUTOD structure serves two purposes: it provides parameters and intermediate storage locations for the autoshimming routine, and it contains the current values for the operations that are performed via the automation control board (see Chapter 7, "Digital Components," on page 65). Compared to the LC structure, AUTOD is small, containing one long (32-bit) integer and ten 16-bit integers:

- AUTOD->when_mask is a 16-bit integer that contains the information from the VNMR load and wshim parameters.
- AUTOD->checkmask, AUTOD->control_mask, AUTOD->best, AUTOD->loops are used for the autoshimming.
- AUTOD->samplemask is the sample changer location (0 = not used).
- AUTOD->sample_error contains sample changer error codes.
- AUTOD->recgain is the receiver gain (value as set by the parameter, or as set by the autogain (if gain='n').
- AUTOD->lockpower, AUTOD->lockgain, AUTOD->lockphase are the lock parameters—again, either as set by the parameters or (in case the autolock is selected) as set by the autolock function.

In earlier VNMR releases, AUTOD was larger and contained many more functions:

- 32 shim gradient coil values. With the increasing number of gradients and the many different possible gradient sets, this became obsolete and was moved into the instruction section. No fixed space is allocated for the gradient settings.
- The shim method as ASCII text (up to 128 characters). This was moved into the instruction section, where the method text appears in parsed (semi-interpreted) form. No fixed amount of space is allocated for the shim method any longer.
- Knobs information—memory space that was allocated for processes internal to the acquisition CPU. This was removed from AUTOD as well.

The Instruction Section

The instruction section of each Acode segment is structured by itself and consists of two subsections:

- The first part does all the initializations and is executed only once per FID.
- The second part contains the code that is generated by the pulse sequence function itself and is looped over nt.

We will now take a more detailed look at the two parts of the instruction section, without discussing individual functions (these will be discussed in later chapters).

The Initialization Part

The initialization part of the instruction section contains definitions and initializations. To a very large extent, this part of the code depends on the rf configuration. The more channel and the more devices that exist, the longer this section becomes. Of course, it also depends on the parameter settings (gradient settings, frequencies, power levels, etc.), but these mainly change the numeric (and binary) values used in this section and not so much its length.

In arrayed and multidimensional experiments, this section also depends on the FID index. Certain segments and instructions (as indicated below) only show up in the instruction segment for the first FID and, in general, the instruction segments for subsequent FIDs are slightly shorter. The difference is about 50 to 150 Acode words (about 100 to 300 bytes), depending on the spectrometer configuration.

At the beginning of the initialization part, the binary values for the fast bits that determine the 90-degree phase shifts are defined for each channel, together with an internal address used in connection with phase shifting. This part only shows up in the Acode for the first FID.

| 112 | 106 | 0 | 15 | SETPHATTRibutes CH1 | 0 | 0x8 | 0x10 | 0x18 |
|-----|-----|----|----|---------------------|-------|--------|-----------------|-------|
| | | | | allBits = | 0x18 | , addr | $= 0 \times b0$ | 094 |
| 126 | 120 | 14 | 15 | SETPHATTRibutes CH2 | 0 | 0x100 | 0x200 | 0x300 |
| | | | | allBits = | 0x300 | , addr | $= 0 \times b0$ | 0a4 |

The next, large section includes initialization output to virtually every digital and rf device in the spectrometer (including magnet leg control and shims), ensuring that they are in their proper state, as defined by configuration and non-pulse sequence dependent acquisition parameters. Most of this section consists of AP output (see Chapter 12, "AP Bus Traffic," on page 137). The lock filter response times are set, all shim gradient settings, frequencies, and power levels are set, amplifiers, rf switches, and relays are set into a proper state. In this section, the TUNE_FREQ instructions (setting the PTS frequency synthesizers), the instructions that initialize the waveform generators, and a few other functions are only present as part of the first Acode segment.

| 140 | 134 | 28 | 53 | ACQBITMASK 0 |
|-----|-----|----|-----|---|
| 142 | 136 | 30 | 1 | CBEGIN |
| 143 | 137 | 31 | 19 | INITialize acq |
| 144 | 138 | 32 | 5 | CLEAR |
| 145 | 139 | 33 | 157 | LOCKFILTER fast = 1, slow = 4 |
| 148 | 142 | 36 | 6 | APBOUT 7 items 0xab40 0xbb0d 0x9b00 0xab54 0xbbec |
| | | | | 0x9b2f 0x9b8f |
| 157 | 151 | 45 | 6 | APBOUT 9 items 0xa720 0xb700 0x97be 0xb7ef 0xb7ae |
| | | | | 0xb7fb 0xb7be 0x9701 0xb700 |
| 168 | 162 | 56 | 159 | TUNE_FREQ CH1 9 words |
| | | | | 0xa720 0xb700 0x97be 0xb7ef 0xb7ae 0xb7fb 0xb7be |
| | | | | 0x9701 0xb700 |
| 180 | 174 | 68 | 0 | NO_OP |
| | | | | |

| 181 175 69 16 SETPHAS90 CH1c 0 184 178 72 68 PHASESTEP CH1 360 units (90.00 degrees) 187 181 75 65 SETPHASE CH1f 0 190 184 78 59 APChipOUT APaddr 11, reg 51, +logic, 1 byte max 79, offset 16, value 55 196 190 84 59 APChipOUT APaddr 11, reg 150, -logic, 2 bytes max 4095, offset 0, value 4095 202 196 90 6 APBOUT 2 items 0xab92 0xbba0 028 212 206 100 6 APBOUT 2 items 0xab91 0xbb80 0202 216 100 6 APBOUT 2 items 0xab90 0xbb10 028 022 214 108 6 APBOUT 2 items 0xab91 0xb20 0x8221 0x8231 0x824 220 214 108 6 APBOUT 2 items 0x8261 0x8214 0x8220 0x8231 0x824 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 0, value 4095 243 237 | | | | | | |
|--|-----|-----|-----|-----|------------|--|
| 187 181 75 65 SETPHASE CH1f 0 190 184 78 59 APChipOUT APaddr 11, reg 51, +logic, 1 byte max 79, offset 16, value 55 196 190 84 59 APChipOUT APaddr 11, reg 150, -logic, 2 bytes max 4095, offset 0, value 4095 202 196 90 6 APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 208 202 96 6 APBOUT 2 items 0xab91 0xbba0 212 206 100 6 APBOUT 2 items 0xab91 0xbb00 216 210 104 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 2 items 0xab20 0x8250 233 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xab40 0xbb20 | 181 | 175 | 69 | 16 | SETPHAS90 | CH1c 0 |
| 190 184 78 59 APChipOUT APaddr 11, reg 51, +logic, 1 byte max 79, offset 16, value 55 196 190 84 59 APChipOUT APaddr 11, reg 150, -logic, 2 bytes max 4095, offset 0, value 4095 202 196 90 6 APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 208 202 96 6 APBOUT 2 items 0xab92 0xbba0 212 206 100 6 APBOUT 2 items 0xab91 0xbb80 216 210 104 6 APBOUT 2 items 0xab90 0xbb10 216 210 104 6 APBOUT 2 items 0xab91 0xb80 220 214 108 6 APBOUT 2 items 0xab91 0xb80 224 218 112 6 APBOUT 2 items 0xab90 0xb21 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 409 10x420 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba0 0xbb00 | 184 | 178 | 72 | 68 | PHASESTEP | |
| max 79, offset 16, value 55 196 190 84 59 APChipOUT APaddr 11, reg 150, -logic, 2 bytes max 4095, offset 0, value 4095 202 196 90 6 APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 208 202 96 6 APBOUT 2 items 0xab92 0xbba0 212 206 100 6 APBOUT 2 items 0xab91 0xbb00 216 210 104 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 4 items 0xab90 0xbb11 224 218 112 6 APBOUT 4 items 0xab90 0xbb11 224 218 112 6 APBOUT 4 items 0xab90 0xbb11 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba2 0xbbc0 255 249 143 6 APBOUT 2 items 0xaba1 0xbb80 | 187 | 181 | 75 | 65 | SETPHASE | CH1f 0 |
| 196 190 84 59 APChipOUT APaddr 11, reg 150, -logic, 2 bytes max 4095, offset 0, value 4095 202 196 90 6 APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 208 202 96 6 APBOUT 2 items 0xab91 0xbbd2 0x9b00 0x9b00 212 206 100 6 APBOUT 2 items 0xab91 0xbb80 212 206 104 6 APBOUT 2 items 0xab91 0xb80 220 214 108 6 APBOUT 2 items 0xab91 0xb80 220 214 108 6 APBOUT 2 items 0xab91 0x8214 0x8220 0x8231 0x824: 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba 0xbb2a 0x9b00 0x9b00 255 253 147 6 APBOUT 2 items 0xaba 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba 0xbb12 271 265 159 6 APBOUT 2 items 0xab3 0xbb12 271 265 | 190 | 184 | 78 | 59 | APChipOUT | APaddr 11, reg 51, +logic, 1 byte |
| max 4095, offset 0, value 4095 202 196 90 6 APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 208 202 96 6 APBOUT 2 items 0xab98 0xbbd2 0x9b00 0x9b00 212 206 100 6 APBOUT 2 items 0xab90 0xbb00 216 210 104 6 APBOUT 2 items 0xab90 0xbb11 220 214 108 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 8 items 0xab201 0x8214 0x8220 0x8231 0x824 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 166, -logic, 1 byte <max 0,="" 4095<="" 79,="" offset="" td="" value=""> 243 237 131 59 APChipOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba8 0xbb12 249 243 137 6<!--</td--><td></td><td></td><td></td><td></td><td></td><td>max 79, offset 16, value 55</td></max> | | | | | | max 79, offset 16, value 55 |
| 202 196 90 6 APBOUT 4 items 0xab98 0xbbd2 0x9b00 0x9b00 208 202 96 6 APBOUT 2 items 0xab92 0xbba0 212 206 100 6 APBOUT 2 items 0xab91 0xbb00 216 210 104 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 8 items 0xab90 0xb214 0x8220 0x8231 0x824 0x8255 0x8264 0x8211 0x8220 0x8231 0x824 0x8255 0x8264 0x8271 0x8231 0x824 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095 249 243 137 6 APBOUT 2 items 0xaba8 0xbb00 0x9b00 0x9b00 255 249 143 | 196 | 190 | 84 | 59 | APChipOUT | APaddr 11, reg 150, -logic, 2 bytes |
| 208 202 96 6 APBOUT 2 items 0xab92 0xbba0 212 206 100 6 APBOUT 2 items 0xab9b 0xbb00 216 210 104 6 APBOUT 2 items 0xab91 0xbb80 220 214 108 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 8 items 0x8201 0x8214 0x8220 0x8231 0x824 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 243 137 6 APBOUT 4Paddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba8 0xbb2a 0x9b00 0x9b00 263 257 151 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xab34 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb01 275 269 <t< td=""><td></td><td></td><td></td><td></td><td></td><td>max 4095, offset 0, value 4095</td></t<> | | | | | | max 4095, offset 0, value 4095 |
| 212 206 100 6 APBOUT 2 items 0xab9b 0xbb00 216 210 104 6 APBOUT 2 items 0xab91 0xbb80 220 214 108 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 8 items 0x8201 0x8214 0x8220 0x8231 0x824: 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba8 0xbbc0 263 257 151 6 APBOUT 2 items 0xaba1 0xbb80 264 257 155 6 APBOUT 2 items 0xab34 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb01 275 269 163 6 | 202 | 196 | 90 | 6 | APBOUT 4 | items 0xab98 0xbbd2 0x9b00 0x9b00 |
| 216 210 104 6 APBOUT 2 items 0xab91 0xbb80 220 214 108 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 8 items 0x8201 0x8214 0x8220 0x8231 0x824 0x8255 0x8264 0x8271 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba1 0xbb00 263 257 151 6 APBOUT 2 items 0xaba1 0xbb00 263 257 151 6 APBOUT 2 items 0xaba4 0xbb01 271 265 159 6 APBOUT 2 items 0xab4 0xbb1 271 265 159 6 APBOUT 2 items 0xab4 0xbb1 275 269 163 6 APBOUT 2 items 0xab | 208 | 202 | 96 | 6 | APBOUT 2 | items 0xab92 0xbba0 |
| 220 214 108 6 APBOUT 2 items 0xab90 0xbb11 224 218 112 6 APBOUT 8 items 0x8201 0x8214 0x8220 0x8231 0x824 0x8255 0x8264 0x8271 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 2 items 0xaba2 0xbbc0 255 249 143 6 APBOUT 2 items 0xaba1 0xbb80 263 257 151 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba4 0xbb01 275 269 163 6 APBOUT 2 items 0xab34 0xbb5 279 273 167 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab44 0xbb55 279 273 167 6 APBOUT 2 items 0xab43 0xbb01 | 212 | 206 | 100 | 6 | APBOUT 2 | items 0xab9b 0xbb00 |
| 224 218 112 6 APBOUT 8 items 0x8201 0x8214 0x8220 0x8231 0x824 0x8255 0x8264 0x8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba0 0xbb12 271 265 159 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab40 0xbb21 | 216 | 210 | 104 | 6 | APBOUT 2 | items 0xab91 0xbb80 |
| Dx8255 Dx8264 Dx8271 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba1 0xbb12 271 265 159 6 APBOUT 2 items 0xab34 0xbb12 275 269 163 6 APBOUT 2 items 0xab40 0xbb12 275 269 163 6 APBOUT 2 items 0xab40 0xbb12 276 273 </td <td>220</td> <td>214</td> <td>108</td> <td>6</td> <td>APBOUT 2</td> <td>items 0xab90 0xbb11</td> | 220 | 214 | 108 | 6 | APBOUT 2 | items 0xab90 0xbb11 |
| 234 228 122 16 SETPHAS90 CH2c 0 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba4 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb1 275 269 163 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb01 276 273 167 6 APBOUT 2 items 0xab43 0xbb01 276 289 183 6 APBOU | 224 | 218 | 112 | 6 | APBOUT 8 | items 0x8201 0x8214 0x8220 0x8231 0x8241 |
| 237 231 125 59 APChipOUT APaddr 11, reg 50, +logic, 1 byte max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba4 0xbb12 271 265 159 6 APBOUT 2 items 0xab30 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab35 0xbb00 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 284 175 6 APBOUT 2 items 0xab49 0xbb01 285 179 6 APBOUT 2 items 0xab49 0xbb01 <t< td=""><td></td><td></td><td></td><td></td><td></td><td>0x8255 0x8264 0x8271</td></t<> | | | | | | 0x8255 0x8264 0x8271 |
| max 79, offset 16, value 30 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba1 0xbb00 263 257 151 6 APBOUT 2 items 0xaba0 0xbb12 271 265 159 6 APBOUT 2 items 0xaba4 0xbb01 275 269 163 6 APBOUT 2 items 0xab30 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab35 0xbb00 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 2923 187 6 APBOUT 2 items 0xab40 0xbb01 291 285 179 6 APBOUT 2 items 0xab40 0xbb06 | 234 | 228 | 122 | 16 | SETPHAS90 | CH2c 0 |
| 243 237 131 59 APChipOUT APaddr 11, reg 166, -logic, 2 bytes max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xaba1 0xbb00 263 257 151 6 APBOUT 2 items 0xaba1 0xbb10 267 261 155 6 APBOUT 2 items 0xab48 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb01 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 284 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab40 | 237 | 231 | 125 | 59 | - | |
| max 4095, offset 0, value 4095 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xabab 0xbbc0 263 257 151 6 APBOUT 2 items 0xabab 0xbb10 263 257 151 6 APBOUT 2 items 0xabab 0xbb12 267 261 155 6 APBOUT 2 items 0xab40 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb55 279 273 167 6 APBOUT 2 items 0xab40 0xbb21 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT | | | | | | max 79, offset 16, value 30 |
| 249 243 137 6 APBOUT 4 items 0xaba8 0xbb2a 0x9b00 0x9b00 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xabab 0xbb00 263 257 151 6 APBOUT 2 items 0xabab 0xbb10 263 257 151 6 APBOUT 2 items 0xabab 0xbb12 267 261 155 6 APBOUT 2 items 0xab48 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab43 0xbb12 277 273 167 6 APBOUT 2 items 0xab40 0xbb21 283 277 171 6 APBOUT 2 items 0xab40 0xbb01 285 179 6 APBOUT 2 items 0xab40 0xbb00 <tr< td=""><td>243</td><td>237</td><td>131</td><td>59</td><td>-</td><td></td></tr<> | 243 | 237 | 131 | 59 | - | |
| 255 249 143 6 APBOUT 2 items 0xaba2 0xbbc0 259 253 147 6 APBOUT 2 items 0xabab 0xbb00 263 257 151 6 APBOUT 2 items 0xaba1 0xbb80 267 261 155 6 APBOUT 2 items 0xaba0 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab40 0xbb01 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 284 175 6 APBOUT 2 items 0xab40 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab40 0xbb06 299 293 187 6 APBOUT 2 items 0xab49 0xbb06 299 293 | | | | | | |
| 259 253 147 6 APBOUT 2 items 0xabab 0xbb00 263 257 151 6 APBOUT 2 items 0xabal 0xbb10 267 261 155 6 APBOUT 2 items 0xabal 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab40 0xbb01 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 284 175 6 APBOUT 2 items 0xab43 0xbb01 297 281 175 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab46 0xbb06 299 293 187 6 APBOUT 2 items 0xab49 0xbb88 303 297 | | | | | APBOUT 4 | items 0xaba8 0xbb2a 0x9b00 0x9b00 |
| 263 257 151 6 APBOUT 2 items 0xabal 0xbb80 267 261 155 6 APBOUT 2 items 0xaba0 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab40 0xbb00 283 277 171 6 APBOUT 2 items 0xab43 0xbb01 297 281 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb01 295 289 183 6 APBOUT 2 items 0xab46 0xbb06 299 293 187 6 APBOUT 2 items 0xab49 0xbb88 303 297 191 6 APBOUT 2 items 0xab40 0xbb88 | 255 | 249 | 143 | 6 | | |
| 267 261 155 6 APBOUT 2 items 0xaba0 0xbb12 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab35 0xbb00 283 277 171 6 APBOUT 2 items 0xab40 0xbb01 291 285 179 6 APBOUT 2 items 0xab40 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | 259 | 253 | 147 | 6 | APBOUT 2 | items 0xabab 0xbb00 |
| 271 265 159 6 APBOUT 2 items 0xab48 0xbb01 275 269 163 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab35 0xbb00 283 277 171 6 APBOUT 2 items 0xab4d 0xbb21 287 281 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab49 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | 263 | 257 | 151 | 6 | APBOUT 2 | items 0xabal 0xbb80 |
| 275 269 163 6 APBOUT 2 items 0xab34 0xbb55 279 273 167 6 APBOUT 2 items 0xab35 0xbb00 283 277 171 6 APBOUT 2 items 0xab4d 0xbb21 287 281 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | 267 | | | | | |
| 279 273 167 6 APBOUT 2 items 0xab35 0xbb00 283 277 171 6 APBOUT 2 items 0xab4d 0xbb21 287 281 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | | 265 | | | APBOUT 2 | items 0xab48 0xbb01 |
| 283 277 171 6 APBOUT 2 items 0xab4d 0xbb21 287 281 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | | | | | | |
| 287 281 175 6 APBOUT 2 items 0xab43 0xbb01 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | | | | | | |
| 291 285 179 6 APBOUT 2 items 0xab49 0xbb00 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | | | | | | |
| 295 289 183 6 APBOUT 2 items 0xab36 0xbb06 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | 287 | 281 | | | | |
| 299 293 187 6 APBOUT 2 items 0xab9b 0xbb88 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | | | | | | |
| 303 297 191 6 APBOUT 2 items 0xabab 0xbb88 | | | | | | |
| | | | | | | |
| 307 301 195 151 EVENT1_TWRD 1.000 usec | | | | | | |
| | 307 | 301 | 195 | 151 | EVENT1_TWR | D 1.000 usec |

The receiver gain is set to the value specified in the corresponding VNMR parameter, and the VT controller is initialized (a temperature setting of 3000 degrees switches off the active VT regulation, the PID regulation parameters are set with every acquisition, and the VT cut-over value is set at the same time). Unless a temperature array is specified, the SETVT instruction is only used in the first FID.

| 309 | 303 | 197 | 71 | GAINAutomation | | | | | |
|-----|-----|-----|-----|------------------|------|------|---------|-----|------|
| 310 | 304 | 198 | 63 | SETVT Oxford PID | 440, | temp | 3000.0, | vtc | 25.0 |
| 315 | 309 | 203 | 0 | NO_OP | | | | | |
| 316 | 310 | 204 | 151 | EVENT1_TWRD | 10 | msec | | | |
| 318 | 312 | 206 | 7 | STartFIFO | | | | | |
| 319 | 313 | 207 | 9 | StopFIFO | | | | | |
| | | | | | | | | | |

Next, the preacquisition delay (pad) is executed; in Acode segments other than for the first FID, this is omitted.

| 320 | 314 | 208 | 95 | PADelay 4 words | | |
|-----|-----|-----|-----|-----------------|-----|------|
| 322 | 316 | 210 | 151 | EVENT1_TWRD | 500 | msec |
| 324 | 318 | 212 | 7 | STartFIFO | | |
| 325 | 319 | 213 | 9 | StopFIFO | | |

Depending on the wshim parameter, autoshimming is now performed on the first FID only, or for each increment. The shim method has become part of the Acode instructions (it does *not* show up in ASCII text, but rather in parsed form).

326 320 214 75 SHIMAutomation mode = 1, 7 words 0x700 0x3c03 0 0 0x146d 0x6d00 0

The last element of the initialization part is the noise measurement: for the first FID, 256 noise points are acquired using the spectral window (as well as the gain and filter settings, of course) of the "real" FID. This noise measurement is used to calculate the maximum number of down-scalings (right-shifts) that can be performed under the

current conditions (in case the maximum or minimum integer is reached in the FID), such that the noise is still sufficiently digitized.

 341
 335
 229
 43
 NOISE loop 256 pts, dwell 199 usec + 1.000 usec

 342
 336
 230
 19
 INITialize acq

 343
 337
 231
 5
 CLEAR

With these instructions, the pulse sequence independent part of the instruction section terminates. The keyword for the pulse-sequence-related section is NextSCan, or the Acode instruction 92.

Pulse-Sequence-Related Part

The pulse-sequence-related part of the instruction section can be as short as 80 Acode words, as in the example (s2pul) below. In other cases it can be up to several thousand Acode words. The maximum Acode size for a single FID (including the initialization section) currently is 10,000 words.

The starting point for the pulse sequence section is always the instruction 92 (NextSCan). Typically, the pulse sequence function itself is a sequence of math functions (not in this example), AP output (APBOUT) statements, instructions that set fast bits (HighSpeedLINES, SETPHAS90), and time events. Any pulse sequence statement call that can possibly change the fast bits (state information) generates a HighSpeedLINES call. In this particular example, we recognize only three time events: pw was set to 7 microseconds, and rof1 and rof2 both set to 10 microseconds. The other pulse (p1) was set to zero (thus the time events associated with p1 do not show up), and the two delays are set to zero and are therefore skipped.

| 344 | 338 | 232 | 98 | NextSCan | |
|-----|-----|-----|-----|----------------|---------------|
| 345 | 339 | 233 | 150 | HighSpeedLINES | (void) |
| 348 | 342 | 236 | б | APBOUT 2 items | 0xab9b 0xbb88 |
| 352 | 346 | 240 | 150 | HighSpeedLINES | (void) |
| 355 | 349 | 243 | 6 | APBOUT 2 items | 0xabab 0xbb88 |
| 359 | 353 | 247 | 150 | HighSpeedLINES | (void) |
| 362 | 356 | 250 | 150 | HighSpeedLINES | (void) |
| 365 | 359 | 253 | 150 | HighSpeedLINES | (void) |
| 368 | 362 | 256 | 150 | HighSpeedLINES | (void) |
| 371 | 365 | 259 | 16 | SETPHAS90 | CH1 zero |
| 374 | 368 | 262 | 150 | HighSpeedLINES | RXOFF |
| 377 | 371 | 265 | 150 | HighSpeedLINES | RXOFF |
| 380 | 374 | 268 | 151 | EVENT1_TWRD | 10.000 usec |
| 382 | 376 | 270 | 150 | HighSpeedLINES | RXOFF TXON |
| 385 | 379 | 273 | 151 | EVENT1_TWRD | 7.000 usec |
| 387 | 381 | 275 | 150 | HighSpeedLINES | RXOFF |
| 390 | 384 | 278 | 151 | EVENT1_TWRD | 10.000 usec |
| 392 | 386 | 280 | 150 | HighSpeedLINES | (void) |
| 395 | 389 | 283 | 150 | HighSpeedLINES | (void) |
| | | | | | |

The remainder of the code is related to the implicit acquisition. The observe phase (oph) is set to zero, because this code was generated using the gf macro that disables the observe phase cycling (to avoid jumping dc levels in the real-time FID display). All rf channels are then reset to zero phase to avoid a center glitch due to rf leakage. Of course, this is only relevant in real acquisitions with cp='y' (i.e., *with* observe phase cycling). After that, the filter delay, alfa+1/(beta*fb), is executed.

The instruction 90 (SETInputCardMode) transmits the observe phase to the sum-tomemory board (see Chapter 18, "Acquiring Data," on page 205). This is done automatically before acquiring the first data point (after which the receiver phase can't be changed any longer); here, the entire FID is acquired in a single instruction.

 398
 392
 286
 39
 ASSIGNFUNC
 zero oph

 401
 395
 289
 16
 SETPHAS90
 CH1 zero

 404
 398
 292
 16
 SETPHAS90
 CH2 zero

 407
 401
 295
 152
 EVENT2_TWRD
 296 usec + 850 nsec

 410
 404
 298
 90
 SETInputCardMode

 414
 408
 302
 99
 ACQXX loop np=37440, dwell 199 usec + 1.000 usec

 417
 411
 305
 7
 STartFIFO

The last instruction before "branching" (jumping) back to the instruction NextSCan is for housekeeping. Housekeeping is necessary for the sum-to-memory card to return to its proper state, finishing storing the FID in the RAM. It also increments the ct counter (LC->ct), decrements block size or steady-state counters (LC->bsctr or LC->ssctr), checks for the end of the an acquisition block size or for the termination of the steady-state scans, and starts the necessary actions (see also Chapter 18, "Acquiring Data," on page 205).

418 412 306 97 HouseKEEPing 419 413 307 20 BRANCH Offset 232

8.4 Acode Interpretation

The acquisition operating system interprets the Acode instruction section. This process involves a number of different tasks:

- Setting internal status registers.
- Setting or reading status registers on other boards, such as the sum-to-memory board or the HAL.
- Math operations on internal registers (see also Chapter 10, "Phase Calculations," on page 95).
- Table manipulations (in essence, extraction of single values from tables, see also Chapter 11, "Phase Tables," on page 115) and the main task.
- Constructing FIFO words that are then suitably "packed" and sent to the pulse programmer through the host I/O bus (see Chapter 7, "Digital Components," on page 65).

FIFO Flow

The last acquisition operating system task, constructing FIFO words, needs to be coordinated with the pulse programmer; otherwise, the pulse programmer might obtain a few short time events or AP bus words, and execute them, and then run out of FIFO words before the acquisition CPU is able to deliver the continuation of the pulse sequence. With a 0.2 microsecond minimum time event, or 1.15 microseconds (2.15 microseconds for older systems) per AP bus word, the execution of these events can easily be faster than the delivery of information from the acquisition CPU. If the FIFO (see Chapter 9, "Pulse Programmers," on page 85) runs empty, this generates an error message "FIFO underflow". On the other hand, if the FIFO is full, the pulse programmer sets a status register that indicates to the acquisition CPU that the FIFO is full, and no more FIFO words are generated.

Initially, the FIFO is empty, but it is stopped (no information is released at the "other" end); this way the acquisition CPU can feed the FIFO words from the initialization part without danger that the FIFO runs empty. Actually, it is exactly the initialization part that consists of only fast (AP bus) events. Only after most of the initialization data has been fed into the pulse programmer and the VT controller has been initialized, the FIFO output is started with Acode instruction 7 (STartFIFO, at local offset 206 in

the above example). The following block (preacquisition delay) is surrounded by Acode instruction 9 (StopFIFO) and another STartFIFO, and also the pulse sequence code itself is surrounded by these function.

In theory, it is possible to have only very few, fast events in the pulse sequence statement (e.g., no delay, a very short pulse, and a very short acquisition at a very large spectral window). In such a case, it must be ensured that the sequence of events during the execution of the pulse sequence is not disrupted by a FIFO underflow.

It can be assumed that during the execution of a pulse sequence, the pulse programmer normally does not run empty or spend extra waiting times in a "stopped" state, because this would also disrupt steady states.

Acode Size Limitations, Acode Buffering

The acquisition CPU has somewhat less than 22 Kbyte of memory space reserved for storing Acodes (including the LC and AUTOD structures and the instruction segment). That defines the absolute maximum for the Acode size on these systems: about 10900 words in the instruction segment. There hasn't been an experiment yet that exceeded this limit (but we certainly could construct one, for the sake of the argument), so this should not be a point to be concerned about for pulse sequence programmers.

The topic of Acode size is more complex than that, however. If we were to fill the Acode space in the acquisition CPU to more than 50%, then only one Acode segment can be held in the acquisition CPU at any given time during an experiment. This means, that after finishing one increment of a multi-FID experiment the system would first have to upload the next Acode segment from the host (and this again can only happen when the acquisition CPU is *polled* by the host) before it can continue with the acquisition of the next trace. This leads to inter-increment delays of over one second (the actual length of that delay is unpredictable) and totally disrupts the steady state in arrayed and *n*D experiments (apart from lengthening the overall experiment time, making any calculation/prediction of the experiment duration very unreliable).

If, however, the Acode segment size is less than 50% of the available memory, more than one Acode segment will be stored in the acquisition CPU, and when one increment is finished, the interpretation of the next Acode segment can start immediately (the system will store as many Acode segments in the acquisition CPU as it can, given the current memory limitations). So, to maintain any steady-state across multiple increments, we should not exceed the limit of approximately 5400 words (10800 bytes) per Acode segment¹. To estimate the size of an Acode segment, enter go('acqi') in VNMR and check the size of file /vnmr/acqqueue/acqi. Code².

A limit of 5400 words still is fairly large for an Acode segment, but it is definitely not impossible to exceed that limit, in particular in some complex sequences with a long, explicitly coded spin locking sequence, or when using either extremely long phase tables and/or trying to use apshaped_pulse with complex pulse shapes (over 1000 to 2000 slices).

¹ We haven't confirmed whether 5400 words is the real and accurate limit beyond which the multiple buffering is lost. The pulse sequence generation software has a built-in limit of 11000 words per Acode segment (or about 10900 words in the instruction section).

² This size is approximate because go('acqi') calculates the Acode only for the first increment (which always has some extra instruction overhead), and also because go('acqi') temporarily alters some of the parameters (see also Section 5.3, "Using go('acqi')," on page 58).

For most cases there is a solution that avoids the problem:

- In tables use the division return factor (see also Chapter 11, "Phase Tables," on page 115).
- The number of shape slices is definitely limited with apshaped_pulse (see also Section 16.5, "What If a Waveform Generator Is Not Available," on page 189).
- As for DANTE-type pulses, explicitly coded programmed spin locking sequences, etc., it is very inefficient to code such elements using the rgpulse statement. rgpulse generates several extra HighSpeedLINES Acode instructions that can be avoided by using constructs like

```
txphase(phase);
delay(delta);
xmtron(); delay(length); xmtroff();
```

which is much more economic in terms of Acode space and execution (see also Section 20.2, "Sideband Suppression in MAS Experiments," on page 230).

• If you have a waveform generator, use that for modulated spin locking instead of programming it explicitly.

Chapter 9. Pulse Programmers

In earlier chapters, we have shown out how a pulse programmer could be constructed in the context of a Varian spectrometer. Of course, this was a rough sketch, showing only the general working aspect, rather than a detailed picture of its internal functionality. In this chapter, we look at the *real* functionality of this board, which is the most central piece of hardware for the execution of a pulse sequence.

9.1 Layout of the Pulse Programmer

Figure 18 shows a functional diagram of the pulse programmer. The pulse programmer gets its input via the CPU I/O bus from the acquisition CPU (see also Chapter 7, "Digital Components," on page 65). This information is fed into a FIFO buffer, a memory buffer into which information is fed sequentially and which releases the




Figure 18. Structure of the pulse programmer

Chapter 9. Pulse Programmers

information in the same sequence through the second port¹. It is a special kind of dualport memory containing words of a certain width. The number of words in a FIFO buffer is often called the *depth* or the *length* of the FIFO.

Since the width of the CPU I/O bus (16-bits parallel) is not the same as the width of the FIFO (28, 36, or 54 bits), the input logic (bus decoder) composes FIFO words from several I/O bus words.

There are actually two FIFO buffers in sequence. Normally, the information (FIFO words) "falls through" to the end of the second FIFO, or (if there is already information in the buffer) down to the position behind the last word that was fed into the buffer previously—until both buffers fill up. At that point, the pulse programmer signals "FIFO full" to the acquisition CPU, which then temporarily stops producing FIFO words until free space is again available in the FIFO.

The information that comes out at the end of the buffer can be fed back into the second FIFO: a multiplexer (a parallel "Y switch" with two inputs and one output) determines whether the information going into the second FIFO is taken from the first FIFO (the normal case) or from the output. The latter is used for the *hardware looping*, in which the information in the second FIFO (the loop FIFO) is circulating for a predefined number of cycles. At the beginning of the last loop pass, the loop control circuitry (a 16-bit loop counter) switches the input (the multiplexer) back to the preloop FIFO.

A part of every FIFO word is state-related information, often called *fast bits* (because they can be switched instantaneously with any FIFO word. At the output, these bits directly drive rf gates and switches, set phases in 90-degree increments, and blank (or unblank) amplifiers.

From the remaining bits in the FIFO word, 16 are used to define the time; out of these, 12 define the *time count* (1 to 2^{12} ; i.e., 1 to 4096 units), the remaining four are used to define the *time base* (seconds, milliseconds, microseconds, "nanoseconds"). The timing information is fed into the timing circuitry, which activates the corresponding timer. When the timer has counted down to zero, a trigger signal is given to the FIFO, causing it to release the next FIFO word. Alternatively, the newer versions of the pulse programmer board can disable the timers for specific FIFO words and cause the following FIFO word to be released upon sensing an external trigger signal. This allows synchronizing a pulse sequence with external events (such as a rotor period for solids experiments, or the heart beat or a respiration cycle for imaging experiments).

Because there are only four different time bases, only two bits are needed to select the time base. The four control bits at the same time can *redefine* the timing part of a FIFO word to be information that is fed into the AP bus by which 16 devices can be addressed directly (i.e., a part of the 16 timing bits serves as address). The rest is numeric information that is transferred to the addressed device (e.g., the small-angle phase shifters).

Alternatively, the AP bus can be used in indirect mode, in which the address is first transmitted in a separate AP word, followed by AP words containing the numeric information. This way, an unlimited number of devices (several thousand) can be addressed, and numeric information of arbitrary precision can be transmitted (e.g.,

¹ This differs greatly from the way information is stored and recalled randomly in RAM. In a FIFO, no addressing is involved: only one word at a time can be read out, and the read-out sequence is given by the order of the input.

frequency information for the PTS synthesizer). During AP bus traffic, the timing circuitry is still activated, ensuring that every AP bus word remains on the bus for a specific, well-defined time (1.15 microseconds for UNITY*plus* systems, 2.15 microseconds on earlier systems), such that there is sufficient time for the addressed devices to decode the address and read the information off the bus.

The last four bits in the FIFO word include the "command to convert" (CTC) bit that triggers the ADC. Unlike all other bits, this bit is not held up during an entire time event, but is electronically reformed to a short trigger pulse. The remaining three bits are used for loop and FIFO control.

All Varian pulse programmers have essentially the same layout—be it the output board used in XLs and VXRs, the acquisition control board used in VXRs and UNITY spectrometers, the pulse sequence control board used in UNITY*plus* spectrometers, or similar boards in other Varian spectrometers. With few exceptions, the differences are only quantitative: both the depth and the width of the FIFO has changed over the years (adding more fast bits), the timing resolution has been improved, and the external trigger has been added.

In summary, the various pulse programmers can be characterized as shown in Table 1. The differences can be summarized as follows:

- The FIFO width varies from 28 bits (Gemini) and 36 bits (XL, VXR, and UNITY) to 54 bits (UNITY*plus*); the number of fast bits in essence is equal to the total number of bits per FIFO word minus 20.
- The length of the loop FIFO has been increased from 63 words (Gemini, XL, and VXR) to 1024 words (VXR-S and UNITY) and 2048 words (UNITY*plus*).
- The preloop FIFO does not exist on the Gemini output board. It was half the size of the loop FIFO in the output board and is of the same depth as the loop FIFO in newer systems.
- On the acquisition control board, the timing resolution has been improved to 25 nanoseconds (100 nanoseconds in earlier boards).
- In the UNITY*plus*, the AP bus has been speeded up by almost a factor of two (see Chapter 12, "AP Bus Traffic," on page 137).

| | Output Board (GEMINI) | Output Board | Acquisition Control Board | Pulse Sequence Control Board |
|------------------------|--------------------------|-------------------------|------------------------------|---------------------------------|
| Use | Gemini | XL, VXR, early VXR-S | late VXR-S, UNITY | UNITYplus |
| FIFO width (fast bits) | 28 (8) bits | 36 (16) bits | 36 (16) bits | 54 (34) bits |
| Pre-loop FIFO | none | 32 words | 1024 words | 2048 words |
| Loop FIFO | 63 words | 63 words | 1024 words | 2048 words |
| Timing Resolution | 100 nsec | 100 nsec | 25 nsec | 25 nsec |
| AP Bus Speed | 2.15 µsec/word | 2.15 µsec/word | 2.15 µsec/word | 1.15 µsec/word |
| External Trigger | no | no | yes | yes |

Table 1. Pulse programmer characteristics

• The external trigger only became available with the acquisition control board.

Fast bits are covered in the next section, and timing characteristics are considered in more detail at the end of this chapter. The consequences of various preloop and loop FIFO lengths are discussed more in Section 14.3, "Hardware Loops," on page 150.

9.2 Fast Bits

Fast bits have been provided for those states that possibly need to be changed with every time event, even as short as 0.2 microseconds, such as transmitter gates, 90-degree phase shifts, receiver gates, and amplifier blanking, plus eventually external devices, such as a laser for CIDNP experiments; also waveform generators need to be triggered with a fast bit.

The demand for fast bits has been growing over the years. The Gemini had only very modest needs in terms of fast bits (using only two rf channels), a UNITY*plus* spectrometer with up to six rf channels is much more demanding in terms of the number of fast-switching lines. For a long time, systems (all except Gemini and UNITY*plus*) were equipped with a fixed number of 16 fast bits—with the consequence that in systems with third rf channel and waveform generators (particularly UNITY systems) some of the fast bits had to be reassigned by shifting some of the fast bit functionality (the decoupler modulation mode dmm in particular) to the AP bus. This cut-over happened with the transition from AP interface board type 2 to type 3 of the same board (only type 3 permits controlling dmm via the AP bus).

Table 2 lists the assignment for the fast bits in pulse programmers with 16 fast bits. The hexadecimal number in the first row describes the numeric value of the corresponding

| HEX Value | AP Interface Type 2 | AP Interface Type 3 |
|-----------|---------------------------------------|---------------------------------------|
| 0x1 | VAR1 (Varian-reserved line 1) | WFG1 (fast line for WFG #1) |
| 0x2 | VAR2 (Varian-reserved line 2) | WFG2 (fast line for WFG #2) |
| 0x4 | SP1 (spare bit 1) | SP1 (spare bit 1) |
| 0x8 | SP2 (spare bit 2) | SP2 (spare bit 2) |
| 0x10 | DECLVL (decoupler level switching) | DECUPLR2 (decoupler 2, gate) |
| 0x20 | MODMA (dmm, bit 1) | DC2_90 (decoupler 2, 90 deg. phase) |
| 0x40 | MODMB (dmm, bit 2) | DC2_180 (decoupler 2, 180 deg. phase) |
| 0x80 | HomoSpoilON (homospoil gate) | HomoSpoilON (homospoil gate) |
| 0x100 | DECPP (homodecoupler gating) | WFG3 (fast line for WFG #3) |
| 0x200 | DC90 (decoupler, 90 deg. phase) | DC90 (decoupler, 90 deg. phase) |
| 0x400 | DC180 (decoupler, 180 deg. phase) | DC180 (decoupler, 180 deg. phase) |
| 0x800 | DECUPLR (decoupler gate) | DECUPLR (decoupler gate) |
| 0x1000 | RFP90 (observe xmtr. 90 deg. phase) | RFP90 (observe xmtr. 90 deg. phase) |
| 0x2000 | RFP180 (observe xmtr. 180 deg. phase) | RFP180 (observe xmtr. 180 deg. phase) |
| 0x4000 | TXON (observer xmtr. gate) | TXON (observer xmtr. gate) |
| 0x8000 | RXOFF (receiver off gate) | RXOFF (receiver off gate) |

Table 2. Fast bit assignments, output boards and acquisition control boards

bit within the 16-bit "fast-bit word" or within the "quiescent states" word in the LC structure. This information is taken from the header file rfconst.h found in /vnmr/psg. The CTC (command-to-convert) bit does not count as a regular fast bit, because it behaves differently (see above).

Evidently, it is a bad idea to try addressing the ambiguous bits (DECLVL/DECUPLR2, MODMA/DC2_90, MODMB/DC2_180, and DECPP/WFG3) directly (using low-level statements taken from /vnmr/psg), because this would generate code that is only correct for one particular type of hardware and could create havoc when executed on a different system. Actually, *no* fast bit should be addressed directly (in the UNITY*plus*, the assignments changed altogether, and in fact, the gate statement for manipulating fast bits directly no longer exists in VNMR 4.3 or later.)

The pulse sequence control board in the UNITY*plus* has a totally different fast bit assignment, as can be seen from Table 3. The fast bits are organized in groups of five bits per rf channel; each channel has its own "receiver gate" (the amplifier blanking), a transmitter gate, a gate for the waveform generator, and two gates for the 90-degree phase shifts. The observe is switched together with the blanking line for the observe channel amplifier. Four additional bits control the homospoil pulse, the rotor synchronization hardware, and the two spare gating lines (for triggering external devices). The latter two are no longer part of LC->squi (the 32-bit "quiescent states" word in the LC structure), but are handled separately.

9.3 Timers and Timer Words

All pulse programmers use four different timers, and the time count (the number of time units to be counted down by the timer) is a 12-bit binary number, resulting in values between 1 and 4096. The four timers include:

- Seconds timer (1 to 4096 seconds)
- Milliseconds timer (1 to 4096 milliseconds)
- Microseconds timer (1 to 4096 microseconds)
- "Nanosecond" timer, which is in units of 100 nanoseconds for output boards and on the Gemini (resulting in a theoretical range of 0.1 to 409.6 microseconds), or in units of 25 nanoseconds for acquisition control boards and pulse sequence control boards (resulting in a theoretical range of 0.025 to 102.4 microseconds).

The minimum timer word is 0.2 microseconds on all pulse programmers. Timer words shorter than that cannot be executed and are suppressed already at C level in the pulse sequence software (they would result in a pulse programmer error message).

A problem seems to exist with the timers listed above, in that it seems that time events above 4 seconds can only be executed in steps of 1 second, time events between 4 milliseconds and 4 seconds can be executed with a precision of 1 milliseconds, and time events between 102 microseconds (409 microseconds on older boards) and 4 milliseconds can be performed with a precision of 1 microsecond only. This seems rather limited and would in fact have a severe impact on multidimensional experiments, because most likely the 2D increments could probably not be spaced properly due to the "granularity" of the duration of time events. The same limitation would apply to the 1D dwell time, which would impose severe restrictions to the settability of the sw parameter. All this is clearly unacceptable.

| HEX Value | RF channel | Function |
|------------|-------------|--------------------------------------|
| | KI' channei | |
| 0x1 | channel 1 | amplifier blanking / receiver gating |
| 0x2 | channel 1 | transmitter gate |
| 0x4 | channel 1 | waveform generator gate |
| 0x8 | channel 1 | 90 degrees phase shift |
| 0x10 | channel 1 | 180 degrees phase shift |
| 0x20 | channel 2 | amplifier blanking ("receiver") |
| 0x40 | channel 2 | transmitter gate |
| 0x80 | channel 2 | waveform generator gate |
| 0x100 | channel 2 | 90 degrees phase shift |
| 0x200 | channel 2 | 180 degrees phase shift |
| 0x400 | channel 3 | amplifier blanking ("receiver") |
| 0x800 | channel 3 | transmitter gate |
| 0x1000 | channel 3 | waveform generator gate |
| 0x2000 | channel 3 | 90 degrees phase shift |
| 0x4000 | channel 3 | 180 degrees phase shift |
| 0x8000 | channel 4 | amplifier blanking ("receiver") |
| 0x10000 | channel 4 | transmitter gate |
| 0x20000 | channel 4 | waveform generator gate |
| 0x40000 | channel 4 | 90 degrees phase shift |
| 0x80000 | channel 4 | 180 degrees phase shift |
| 0x100000 | channel 5 | amplifier blanking ("receiver") |
| 0x200000 | channel 5 | transmitter gate |
| 0x400000 | channel 5 | waveform generator gate |
| 0x800000 | channel 5 | 90 degrees phase shift |
| 0x1000000 | channel 5 | 180 degrees phase shift |
| 0x2000000 | channel 6 | amplifier blanking ("receiver") |
| 0x4000000 | channel 6 | transmitter gate |
| 0x8000000 | channel 6 | waveform generator gate |
| 0x10000000 | channel 6 | 90 degrees phase shift |
| 0x20000000 | channel 6 | 180 degrees phase shift |
| 0x40000000 | | homospoil gate |
| 0x80000000 | | rotor synchronization |
| [0x10] | | spare gate 1 |
| [0x20] | | spare gate 2 |

 Table 3. Fast-bit assignments on pulse sequence control boards

The solution is to perform *double-precision timer words*, where necessary. If there is a remainder after performing the major part of a time event in one unit (e.g., seconds for time events above 4.096 seconds), a *second* time event is included that performs the remainder of the time event in the next smaller unit. Single-precision timer words are used only if there is no remainder from a single time unit. With these two options for time events, we get the following possibilities:

- Time events above 4100 seconds are performed as a double timer word with the seconds timer; the maximum time event is 8192 seconds. The round-off error (up to 0.5 seconds) is less than 0.012% for this range of durations.
- Time events above 4.1 seconds can be performed with millisecond precision, with a round-off error of less than 0.012% (up to 0.5 milliseconds absolute).
- Time events between 4.2 milliseconds (4.5 milliseconds on boards with 100 nanoseconds resolution) and 4.1 seconds are performed with microsecond precision, the maximum round-off error being again below 0.012% (up to 0.5 microseconds absolute).
- Time events smaller than 4.2 milliseconds (4.5 milliseconds on boards with 100 nanoseconds resolution) are performed with the timing resolution of the pulse programmer (100 or 25 nanoseconds). By nature, the round-off error can be considerable in this range: for the output board it is up to 50 nanoseconds absolute, or 12 ppm at the upper end, and 25% at the lower end (0.2 microseconds). With newer boards, it is up to 12.5 nanoseconds in absolute terms, or up to 3 ppm at the upper limit and below 6% for the shortest time intervals.

The reality is somewhat more complicated. On all pulse programmers it takes 150 nanoseconds for the hardware to decode the time base, store the time count in the time counter, and start the corresponding timer. Consequently, there is a 150 nanoseconds dead-time involved with every timer word (300 nanoseconds with double-precision time events). For small time events, the software corrects for this error by adjusting the nanosecond count, for larger durations this error remains uncorrected, but is negligible (0.007% maximum).

One more thing was neglected up to now: the numeric range of a 12-bit binary number is actually from 0 to 4095, but the timers, on the other hand, always perform at least one time count to begin with (0 input results in a count of 1). Thus, the 12-bit input range translates to a count range of 1 to 4096 units.

Table 4 summarizes the various possibilities with single- and double-precision time events. Single precision (except for those with the smallest units) is, of course, only performed if there is no remainder. The somewhat "odd" lower limits are due to the fact that smaller delays will be performed with the next smaller time base (and because there is some overlap in the ranges of the individual time bases). Whenever the "nanoseconds" time base is involved on boards with 25-nanosecond resolution, 6 or 12 counts are subtracted from the "nanoseconds" part for single- and double-precision time events, respectively; these time events will be accurate to 25 nanoseconds.

On the output board (63-word loop FIFO), short time events have slightly different characteristics, because the timing resolution is 100 nanoseconds only and timing errors are not corrected for this hardware. Table 5 covers output boards.

It is not necessary to fully understand all of the tables shown here, but this information should help in understanding how certain events are translated into single- or doubleprecision timer words. The question whether double-precision timer words are used or not isn't relevant for most cases either, but may play an important role in the execution of hardware loops, especially on systems with 63-word loop FIFO (see also Section 14.3, "Hardware Loops," on page 150).

9.4 Problems with Timer Word Errors

Compared with the total duration of a timer word, the possible timing errors seem small and, overall, the timer word precision seems more than sufficient for the common NMR experiments. There is one important exception, however: in nD experiments the evolution delay is performed as a single time event that is calculated for every increment (see also Section 19.1, "Indirect Time Domain Incrementation," on page

| Туре | Time Base(s) | [Counts (min/max)] Duration range | Resolution | Timing Error |
|--------|----------------|--|-----------------|-----------------|
| double | sec / sec | [4096 + 1 up to 4096 + 4096] 4097 sec - 8192 sec | 1 sec | 300 nsec longer |
| single | sec | 5 sec - 4096 sec | 1 sec | 150 nsec longer |
| double | sec / msec | [4 + 97 up to 4096 + 4096] 4.097 sec - 4100.096 sec | 1 msec | 300 nsec longer |
| single | msec | 5 msec - 4096 msec | 1 msec | 150 nsec longer |
| double | msec / µsec | [4 + 97 up to 4096 + 4096] 4.097 msec - 4.100096 sec | 1 µsec | 300 nsec longer |
| single | μsec | a "nsec" time event is always added to co | prrect for dead | times |
| double | µsec / 25 nsec | [102 + 11 (+ 12) / 4096 + 4096 (+ 12)] 102.575 µsec - 4198.7 µsec | 0.025 µsec | accurate |
| single | 25 nsec | [2 (+ 6) - 4096 (+ 6)] 0.2 µsec - 102.55 µsec | 0.025 µsec | accurate |

Table 4. Single- and double-precision timer word characteristics

Table 5. Single- and double-precision timer word characteristics, output boards

| Туре | Time Base(s) | [Counts (min/max)] Duration range | Resolution | Timing Error |
|--------|-----------------|---|------------|-----------------|
| double | sec / sec | [4096 + 1 up to 4096 + 4096] 4097 sec - 8192 sec | 1 sec | 300 nsec longer |
| single | sec | 5 sec - 4096 sec | 1 sec | 150 nsec longer |
| double | sec / msec | [4 + 97 up to 4096 + 4096] 4.097 sec - 4100.096 sec | 1 msec | 300 nsec longer |
| single | msec | 5 msec - 4096 msec | 1 msec | 150 nsec longer |
| double | msec / µsec | [4 + 97 up to 4096 + 4096] 4.097 msec - 4.100096 sec | 1 µsec | 300 nsec longer |
| single | μsec | 411 µsec - 4096 µsec | 1 µsec | 150 nsec longer |
| double | µsec / 100 nsec | [409 + 7 up to 4096 + 4096] 409.7 μsec - 4505.6 μsec | 0.1 µsec | 300 nsec longer |
| single | 100 nsec | [2 up to 4096] 0.2 μsec - 409.6 μsec | 0.1 µsec | 150 nsec longer |

215). What counts in the end is the *difference* between the increments, which should be as equal as the dwell times during the standard acquisition. Any random variation in the "evolution dwell time" will translate to noise in the indirect dimension of the final spectrum and, much worse, any periodic variation will lead to extra signals (sidebands) in the indirect domain.

What is the nature of these timing errors in nD experiments? Typically (with spectral windows of a few kHz) the evolution time for first increments on an nD experiment is below 4 msec (i.e., it is performed accurately with the full timing precision of the pulse programmer, 25 nsec on most systems). But for those increments that have an evolution time *greater than 4 msec* (see the above tables) two errors occur:

- The delays are rounded off to full µsec time events. This will lead to a periodic error in the course of the evolution time incrementation.
- The timer word overhead is no longer corrected for (i.e., evolution times larger than 4.097 msec are 300 nsec too long, or 150 nsec for single-precision millisecond time events).

We could set the spectral window in the indirect dimensions such that the evolution time is always a multiple of a full μ sec (which certainly would be inconvenient), but then still there would be a single incrementation discontinuity when crossing the 4.097 msec limit, beyond which the timing overhead is not corrected. This alone probably wouldn't be too bad, but worse is that when the evolution time happens to "fall onto" a full msec duration, which will occur periodically, the timing overhead would only be half of what it is for the other time events above 4.097 msec. This can translate to visible artifacts (sidebands) in nD spectra with a large dynamic range.

To correct for these problems, as of VNMR 5.1 the module /vnmr/psg/delay.c has been modified such that any time event above 4 msec and below 4 sec is performed in *three timer words*:

- a (single) msec timer word, with the msec portion of the duration *minus 2 msec*.
- a double (usec/nsec) timer word of 2 msec, *plus* the sub-msec fraction, *minus the timing overhead of the first timer word*.

This change has the effect of producing triple-precision timer words for durations between 4 msec and 4 sec, and with this *all time events up to 4 sec are performed at the full precision of the pulse programmer* (25 nsec on most systems). These "triple-precision time events" will, of course, lead to extra FIFO words in the pulse programmer; however, this only affects delays above 4 msec, and hence should not lead to problems such as having too many events in a hardloop, etc.

9.5 Timer Words and Fast Bits in the Acode

For acquisition control boards and pulse sequence control boards, time events are encoded with the instructions 151 and 152, EVENT1_TWRD and EVENT2_TWRD. EVENT1_TWRD instructions are followed by a 16-bit timer word (the time base is in the 4 most-significant bits), EVENT2_TWRD are followed by two 16-bit timer words.

Chapter 9. Pulse Programmers

Fast bits are set by separate instructions: gate-related information (excluding the 90degree phase shifts) is set by instruction 150 (HighSpeedLINES, followed by a 32bit pattern), 90-degree phase shifts are set with the instruction 16 (SETPHAS90, followed by a channel identifier and the address of a location in the LC structure):

| 362 | 356 | 250 | 150 | HighSpeedLINES | (void) |
|-----|-----|-----|-----|----------------|-------------|
| 365 | 359 | 253 | 150 | HighSpeedLINES | (void) |
| 368 | 362 | 256 | 150 | HighSpeedLINES | (void) |
| 371 | 365 | 259 | 16 | SETPHAS90 | CH1 zero |
| 374 | 368 | 262 | 150 | HighSpeedLINES | RXOFF |
| 377 | 371 | 265 | 150 | HighSpeedLINES | RXOFF |
| 380 | 374 | 268 | 151 | EVENT1_TWRD | 10.000 usec |
| 382 | 376 | 270 | 150 | HighSpeedLINES | RXOFF TXON |
| 385 | 379 | 273 | 151 | EVENT1_TWRD | 7.000 usec |
| 387 | 381 | 275 | 150 | HighSpeedLINES | RXOFF |
| 390 | 384 | 278 | 151 | EVENT1_TWRD | 10.000 usec |
| 392 | 386 | 280 | 150 | HighSpeedLINES | (void) |

Any pulse sequence statement that can possibly alter the fast bits (like status or the various gating instructions within the G_pulse statement, even if the pulse duration is zero) produces a HighSpeedLINES instruction in the Acode. Because many time events in a pulse sequence may be set to zero (like in the above example using the s2pul pulse sequence), it is not unusual to find a series of HighSpeedLINES with identical arguments and no time events in-between. Of course, only instructions 151 and 152 (EVENT1_TWRD and EVENT2_TWRD) actually produce FIFO words (1 and 2, respectively); the HighSpeedLINES and SETPHAS90 only change specific bits in the "quiescent states" register LC->squi.

For pulse programmers with 16 fast bits (output boards and acquisition control boards), the two spare lines (set and unset in the pulse sequence code by the statements splon(), sploff(), sp2on(), and sp2off()) are part of the standard fast bit pattern and are set the same as the other gating bits, through the HighSpeedLINES instruction. For the pulse sequence control board, the two spare lines are set with a special instruction 156 (SPARE12), followed by a 16-bit word containing the numeric equivalent to the two fast bits (0x10 for spare line 1, 0x20 for the spare line 2), because these two bits are not part of the standard 32-bit fast bit word LC->squi.

Earlier software releases (prior to UNITY*plus* and the pulse sequence controller board) did not have the Acode instructions 150 to 156; at that time, the high-speed lines (including the two spare lines) were set together with the time event instructions EVENT1 and EVENT2 (46 and 47). These instructions were followed by one or two timer words, plus an extra 16-bit word with the fast bits. A HighSpeedLINES instruction was not required; only the 90-degree phase shifts were set by the same instruction SETPHAS90 (16) as today.

Chapter 10. Phase Calculations

Phase cycling in pulsed NMR experiments is the alteration of rf phases as a function of the scan number, thereby co-adding experiments with different rf phases, by altering the phases in a pulse sequence as a function of the transient counter variable ct. There are many reasons to do phase cycling, including:

- Coherence pathway selection, such as multiple-quantum filtering, f₁ or f₂ quadrature selection, etc.
- Cancellation of artifacts arising from spectrometer imperfections, such as channel imbalance (quadrature images), dc offsets (center glitch), phase errors, pulse imperfections (off-resonance effects).
- Cancellation of phase errors due to J-coupling.
- Cancellation of phase errors due to (sometimes deliberately) improper refocusing.

10.1 How Do Phase Calculations Work?

The most natural thing seems to be to calculate the phases directly from the ct counter via some mathematical algorithm. Of course, this presumes that such an algorithm exists: we need a mathematical prescription that generates a given vector (a one-dimensional array of scalar numbers; for 90-degree phase shifting typically integer numbers between 0 and 3) from the vector of natural numbers (i.e., the ct counter). It can be shown that for any *repetitive* sequence of numbers, such an algorithm exists.

Phase cycles by nature are repetitive; therefore, the calculation is feasible, supposing we have the necessary mathematical tools. The complexity of the mathematical procedure (the number of mathematical operations and the storage requirements for intermediate results) heavily depends on the (repetitive) sequence of numbers (the target vector) that is to be created; the means that VNMR provides should actually be sufficient for creating *any* phase cycle used in NMR experiments.

True random sequences of numbers cannot be generated this way, but pseudo-random sequences of sufficient length and quality can be generated using integer math (using modulo functions, see Section 10.8, "Real-Time Random Numbers," on page 112).

The Tools

The basis for phase mathematics is formed by two main ingredients: *operators* (the math operations themselves) and *operands* (the objects—numbers and storage locations—that are used and manipulated by the operators). A third ingredient, *real-time logical decisions*, is sometimes also used for the construction of more complex phase cycles. Also, *C constructs* are often used to make phase calculations dependent on certain VNMR parameters like phase.

C constructs can only happen on a *per-FID* basis, because after the Acode generation C constructs can have no effect on the execution of real-time math. Also, C math operators cannot be used for real-time math: phase math operators are interpreted by the acquisition CPU in "real time"; operands are addresses (offsets in the LC structure), on which the real-time math operations are performed.

Real-Time Operands

The operands for real-time math can be divided into constants and variables:

- Simple numeric constants: zero, one. two, three (LC elements with numeric values 0, 1, 2 and 3, respectively).
- Constants ssval and bsval, holding the number of steady-state and blocksize transients.
- Predefined variables ct, ssctr, and bsctr, holding the transient counter (least-significant half only, see "LC Data Structure" on page 75), the steady-state transient counter, and the block size transient counter, respectively.
- True real-time variables v1, v2, ... v14.

It obviously doesn't make sense to change the values of the numeric constants and, with very few exceptions, to modify other constants and predefined variables. Therefore, results are normally placed only in the true real-time variables v1 to v14.

Real-Time Operators

A full set of basic math operators exists for the construction of phase cycles. These include unary operators (one argument only), and operators with two and three arguments. Always the last argument is the *target operand* and that contains the result of the calculation (it is therefore *modified* by the operator), preceding operands are not modified by the calculation. Table 6 shows the set of real-time math operators (ph1, ph2, etc. are real-time operands, see above.

With all these operators, the last argument is the target operand and is modified by the real-time calculation. Math results can, of course, be used in subsequent calculations. Complex algorithms are built by chaining a series of simple operations.

All calculations are integer operations. Divisions always result in an integer, and fractional numbers are truncated to the next lower integer number. All operations are

| Туре | Syntax | Meaning |
|---------------------------|---|---|
| assignment | assign(ph1,ph2); | ph2 = ph1 |
| increment / decrement | <pre>incr(ph1); decr(ph1);</pre> | ph1 = ph1 + 1 ph1 = ph1 - 1 |
| addition / subtraction | add(ph1,ph2,ph3); sub(ph1,ph2,ph3); | ph3 = ph1 + ph2 ph3 = ph1 - ph2 |
| multiplication | dbl(ph1,ph2); mult(ph1,ph2,ph3); | ph2 = 2 * ph1 ph3 = ph1 * ph2 |
| division | hlv(ph1,ph2); divn(ph1,ph2,ph3); | ph2 = ph1 / 2 ph3 = ph1 / ph2 |
| modulo | <pre>mod2(ph1,ph2); mod4(ph1,ph2); modn(ph1,ph2,ph3);</pre> | ph2 = ph1 % 2 ph2 = ph1 % 4 ph3 = ph1 % ph2 |

Table 6. Real-time math operators

performed in 16-bit integer format with a range of -32768 up to 32767. Integer overflow is suppressed—results are subjected to an implicit (mod 32768) operation.

Variable operands can contain positive as well as negative numbers. In most of the typical applications for phase calculations (e.g., the calculation of 90-degree or small-angle phase shift steps, the use of a real-time variable is as a table index) the function for which the result is used has a very limited numeric range, such as 0 to 3 for 90-degree phase shifts, or 0 to 7 for 45-degree phase steps. In these cases, the value in the variable undergoes an implicit modulo function (without altering the variable itself): for 90-degree phase shifts, the two least-significant bits are extracted (-1 mod 4 is 3), for 45-degree phase steps, the three least-significant bits are taken (-1 mod 8 is 7). In other words, as long as the number of possible (phase) values is a power of two (2, 4, 8, ...), a negative number *n* in phase cycling behaves like 32768 - |n|.

As mentioned above, in most cases the result of a real-time calculation is subject to some modulo function in the end (either implicitly or explicitly). For instance, with 90-degree phase shifting it normally does not matter whether a variable contains the numbers 3, 7, 11, or -1. In fact, many people tend to think that in such situations "4 is the same as 0" and mentally apply a modulo function to every intermediate result in a phase calculation. This thinking is certainly not correct, but it still leads to the proper result, *as long as no division* (divn, hlv) *is involved*.

New Real-Time Numeric Constants

The real-time numeric constants zero, one, two, and three exist for calculations. In principle, three of those could be calculated from a single value, for example:

| assign(zero,v1); | /* | 0 | */ |
|----------------------|----|---|----|
| <pre>incr(v1);</pre> | /* | 1 | */ |
| dbl(v1,v2); | /* | 2 | */ |
| add(v1,v2,v3); | /* | 3 | */ |

In other words, three of the four constants are not really a necessity, but are a pure convenience for the programmer (actually, we save some Acode space because saving three more constants costs three 16-bit words, while the above calculations take up 12 AP words). Sometimes we would like to make calculations with much higher numeric values, the value 127 perhaps (for an example see Section 10.8, "Real-Time Random Numbers," on page 112). Again, such numbers can be calculated:

| dbl(two,v1); | /* | 4 */ |
|----------------------------|----|--------|
| dbl(v1,v1); | /* | 8 */ |
| <pre>mult(v1,v1,v1);</pre> | /* | 64 */ |
| dbl(v1,v1); | /* | 128 */ |
| decr(v1); | /* | 127 */ |

On the other hand, it would seem convenient if we could just create a new numeric constant. This is in fact possible, using the initval statement:

```
initval(127.0,v1);
```

This initializes v1 (LC->v1) with the numeric value 127 (instead of 0). This way we haven't wasted a single word of Acode! However, there is a big *caveat* in this method: v1 in fact contains the value 127 when the above statement is used, but this happens only once: upon typing go, when the Acode is generated (*better*: when the LC structure is created and initialized). If *any* real-time calculation overwrites the value of v1 (i.e., if v1 is used as target operand), the initial value is lost for the current transient and any

subsequent transients for that FID! The consequence is that the variable that is initialized using initval is lost as a real-time variable for that pulse sequence.

There is also a potential danger that the programmer uses initval at the beginning of a pulse sequence and then by mistake overwrites the value. This leads to very subtle pulse sequence errors that can be hard to find and debug (certainly, the syntax checker does not find such errors!).

The conclusion is that initval allows simplifying the pulse sequence code, but in general it is probably better to avoid using it, thus saving real-time variables and avoiding pulse sequence errors. We should be careful when using initval to store power levels in real-time variables (this practice was common in Varian pulse sequences for a long time). If by mistake power levels are altered by real-time calculations, excessive rf load and heating may damage the probe and precious samples. It is better to avoid setting power levels through real-time variables (see also Chapter 12, "AP Bus Traffic," on page 137).

The fact that initval happens at go time and not at real time means that it doesn't matter where the initval statement is placed. It can be the last statement in the pulse sequence, but the value can be used in the first real-time calculation.

Phase Calculations in the Acode

Phase calculations are the part of the pulse sequence program that is transferred into Acode almost directly. Every statement creates one Acode instruction, and every realtime math argument translates to an Acode word (LC address pointer). The following example is generated from a slightly modified version of the standard inadqt.cpulse sequence in VNMR 4.3:

| | v3 | ct | MOD4FUNC | 37 | 272 | 378 | 384 |
|-----|-----|------|------------|----|-----|-----|-----|
| | v9 | ct | HLVFUNC | 34 | 275 | 381 | 387 |
| | v9 | v9 | HLVFUNC | 34 | 278 | 384 | 390 |
| | v10 | zero | ASSIGNFUNC | 39 | 281 | 387 | 393 |
| | v1 | v9 | ASSIGNFUNC | 39 | 284 | 390 | 396 |
| | v9 | v9 | HLVFUNC | 34 | 287 | 393 | 399 |
| | v2 | v9 | ASSIGNFUNC | 39 | 290 | 396 | 402 |
| | v9 | v9 | HLVFUNC | 34 | 293 | 399 | 405 |
| | v9 | v9 | HLVFUNC | 34 | 296 | 402 | 408 |
| | v9 | v9 | MOD2FUNC | 36 | 299 | 405 | 411 |
| | v1 | vl | DBLFUNC | 33 | 302 | 408 | 414 |
| v1 | v1 | v9 | ADDFUNC | 29 | 305 | 411 | 417 |
| | oph | vl | ASSIGNFUNC | 39 | 309 | 415 | 421 |
| | v8 | v2 | DBLFUNC | 33 | 312 | 418 | 424 |
| v2 | v2 | v9 | ADDFUNC | 29 | 315 | 421 | 427 |
| oph | oph | v8 | ADDFUNC | 29 | 319 | 425 | 431 |
| | | v3 | DBLFUNC | 33 | 323 | 429 | 435 |
| v4 | v4 | v3 | ADDFUNC | 29 | 326 | 432 | 438 |
| v3 | v9 | v3 | ADDFUNC | 29 | 330 | 436 | 442 |
| oph | oph | v4 | ADDFUNC | 29 | 334 | 440 | 446 |
| oph | oph | v10 | ADDFUNC | 29 | 338 | 444 | 450 |
| | | | | | | | |

Considering the fact that this code generates several phase cycles that are 64 steps long, the above code can be regarded as being quite efficient: only 70 Acode words (21 instructions) are used to generate these phase cycles (and this isn't even the shortest possible coding!).

10.2 Case 1: Decoding Phase Calculations

Anyone programming sequences with real-time math has two hurdles to overcome in connection with the real-time math: to understand how a sequence with real-time math works, and then to find an appropriate real-time math algorithm that generates the desired phase cycle. The understanding part is the easier of the two tasks and is therefore discussed first.

If we take the above example, we might find a phase cycling (calculation) section in the pulse sequence that looks as follows (assuming that there are no comments):

```
mod4(ct,v3);
hlv(ct,v9); hlv(v9,v9);
assign(v9,v1);
hlv(v9,v9); assign(v9,v2);
hlv(v9,v9); hlv(v9,v9); mod2(v9,v9);
dbl(v1,v1); add(v9,v1,v1);
assign(v1,oph);
dbl(v2,v8);
add(v9,v2,v2);
add(v8,oph,oph);
dbl(v3,v4); add(v3,v4,v4);
add(v3,v9,v3);
add(v4,oph,oph);
```

The only safe way (apart from using a computer program) to evaluate the resulting phase cycles is to write down the result of every single calculation step. For the beginner, an easy method is to write the vectors in columns (at least for shorter and less complex phase cycles):

| | mod4 (ct,v3) | hlv (ct,v9) | hlv (v9,v9) | assign (v9,v1) | hlv (v9,v9) | assign (v9,v2) | hlv (v9,v9) |
|----|-----------------|----------------|----------------|-------------------|----------------|-------------------|----------------|
| ct | v3 | v9 | v9 | v1 | v9 | v2 | v9 |
| 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| 1 | 1 | 0 | 0 | 0 | 0 | 0 | 0 |
| 2 | 2 | 1 | 0 | 0 | 0 | 0 | 0 |
| 3 | 3 | 1 | 0 | 0 | 0 | 0 | 0 |
| 4 | 0 | 2 | 1 | 1 | 0 | 0 | 0 |
| 5 | 1 | 2 | 1 | 1 | 0 | 0 | 0 |
| 6 | 2 | 3 | 1 | 1 | 0 | 0 | 0 |
| 7 | 3 | 3 | 1 | 1 | 0 | 0 | 0 |
| 8 | 0 | 4 | 2 | 2 | 1 | 1 | 0 |
| 9 | 1 | 4 | 2 | 2 | 1 | 1 | 0 |
| 10 | 2 | 5 | 2 | 2 | 1 | 1 | 0 |
| 11 | 3 | 5 | 2 | 2 | 1 | 1 | 0 |
| 12 | 0 | 6 | 3 | 3 | 1 | 1 | 0 |
| 13 | 1 | 6 | 3 | 3 | 1 | 1 | 0 |
| 14 | 2 | 7 | 3 | 3 | 1 | 1 | 0 |
| 15 | 3 | 7 | 3 | 3 | 1 | 1 | 0 |
| 16 | 0 | 8 | 4 | 4 | 2 | 2 | 1 |
| 17 | 1 | 8 | 4 | 4 | 2 | 2 | 1 |
| 18 | 2 | 9 | 4 | 4 | 2 | 2 | 1 |

Chapter 10. Phase Calculations

| | mod4 (ct,v3) | hlv (ct,v9) | hlv (v9,v9) | assign (v9,v1) | hlv (v9,v9) | assign (v9,v2) | hlv (v9,v9) |
|----|-----------------|----------------|----------------|-------------------|----------------|-------------------|----------------|
| ct | v3 | v9 | v9 | v1 | v9 | v2 | v9 |
| 19 | 3 | 9 | 4 | 4 | 2 | 2 | 1 |
| 20 | 0 | 10 | 5 | 5 | 2 | 2 | 1 |
| 21 | 1 | 10 | 5 | 5 | 2 | 2 | 1 |
| 22 | 2 | 11 | 5 | 5 | 2 | 2 | 1 |
| 23 | 3 | 11 | 5 | 5 | 2 | 2 | 1 |
| 24 | 0 | 12 | 6 | 6 | 3 3 | | 1 |
| 25 | 1 | 12 | 6 | 6 | 3 | 3 | 1 |
| 26 | 2 | 13 | 6 | 6 | 3 | 3 3 | |
| 27 | 3 | 13 | 6 | 6 | 3 | 3 | 1 |
| 28 | 0 | 14 | 7 | 7 | 3 | 3 | 1 |
| 29 | 1 | 14 | 7 | 7 | 3 | 3 | 1 |
| 30 | 2 | 15 | 7 | 7 | 3 | 3 | 1 |
| 31 | 3 | 15 | 7 | 7 | 3 | 3 | 1 |
| 32 | 0 | 16 | 8 | 8 | 4 | 4 | 2 |
| 33 | 1 | 16 | 8 | 8 | 4 | 4 | 2 |
| 34 | 1 | 17 | 8 | 8 | 4 | 4 | 2 |
| 35 | 2 | 17 | 8 | 8 | 4 | 4 | 2 |

For more complex phase cycles, the vectors are better written in rows instead of columns, either on a piece of paper or, even better, directly in the pulse sequence as a comment behind every single function call:

```
mod4(ct,v3);
                 /* v3 = 01230123 */
hlv(ct,v9); /* v9 = 001122334455 ... */
hlv(v9,v9); /* v9 = 00001112222333344 ... = ct/4 */
assign(v9,v1); /* v1 = ct/4 */
hlv(v9,v9); /* v9 = ct/8 */
assign(v9,v2); /* v2 = ct/8 */
hlv(v9,v9); /* v9 = ct/16 */

hlv(v9,v9); /* v9 = ct/16 */

hlv(v9,v9); /* v9 = ct/32 */

mod2(v9,v9); /* v9 = (32*0) (3

dbl(v1,v1); /* v1 = 00002222
                 /* v9 = (32*0) (32*1) */
                  /* v1 = 00002222444466666 ... */
                  /* = 0000222200002222 (quadrature phases) */
add(v9,v1,v1); /* v1 = 000022220000222200002222
                            11113333111133331111333311113333 */
assign(v1,oph); /* oph = v1 */
dbl(v2,v8); /* v8 = 000000022222222 */
add(v9,v2,v2); /* v9 = 000000022222222000000022222222
                            111111133333331111111133333333 */
add(v8,oph,oph); /* oph = 0000222222220000000222222220000
                            1111333333311111111333333331111 */
                 /* v4 = 02460246 */
dbl(v3,v4);
add(v3,v4,v4)
                  /* v4 = 03690369 = 03210321 */
add(v3,v9,v3); /* v3 = 01230123230123010123012323012301
                            12301230301230121230123030123012 */
add(v4,oph,oph); /* oph = 03212103210303210321210321030321
                            10323210321010321032321032101032 */
```

Phase cycles can be rather long. To avoid excessive typing, a shorthand syntax is very helpful (although it makes addition of phase vectors more difficult). It is recommended to use the shorthand syntax that is also used in the definition of phase tables (covered in the section "Shorthand Notation" on page 116):

(a b c d)n repeat the entire sequence within the parentheses n times (e.g., for n=2: a b c d a b c d)
[a b c d]n repeat each individual element within the brackets n times (e.g., for n=2: a a b b c c d d)

Such a simplified comment would look as follows:

```
mod4(ct,v3);
                 /* v3 = 0 1 2 3 */
hlv(ct,v9); /* v9 = [ 0 1 2 3 ]2 */
hlv(v9,v9); /* v9 = [ 0 1 2 3 ]4 */
assign(v9,v1); /* v1 = [ 0 1 2 3 ]4 */
hlv(v9,v9); /* v9 = [ 0 1 2 3 ]8 */
assign(v9,v2); /* v2 = [ 0 1 2 3 ]8 */
hlv(v9,v9); /* v9 = [ 0 1 2 3 ]16 */
hlv(v9,v9); /* v9 = [ 0 1 2 3 ]32 */
mod2(v9,v9); /* v9 = [ 0 1 ]32 */
dbl(v1,v1); /* v1 = [ 0 2 ]4 */
add(v9,v1,v1); /* v1 = [ 0 2 0 2 0 2 0 2 1 3 1 3 1 3 1 3 ]4 */
assign(v1,oph); /* oph = [ 0 2 0 2 0 2 0 2 1 3 1 3 1 3 1 3 ]4 */
dbl(v2,v8); /* v8 = [ 0 2 ]8 */
add(v9,v2,v2); /* v9 = [ 0 2 0 2 1 3 1 3 ]8 */
add(v8,oph,oph); /* oph = [ 0 2 2 0 0 2 2 0 1 3 3 1 1 3 3 1 ]4 */
dbl(v3,v4); /* v4 = 0 2 0 2 */
add(v3,v4,v4); /* v4 = 0 3 2 1 */
add(v3,v9,v3); /* v3 = ( 0 1 2 3 0 1 2 3 2 3 0 1 2 3 0 1 )2
                      (12301230 30123012)2*/
add(v4,oph,oph); /* oph = ( 0 3 2 1 2 1 0 3 2 1 0 3 0 3 2 1 )2
                           (10323210 32101032)2*/
```

10.3 Case 2: Creating Phase Math for Given Phase Tables

While the reconstruction of phase cycles from real-time math statements is not very difficult and can be solved by applying pure logic (e.g., in a computer program), the opposite task is quite different: nobody has yet come up with a simple logical approach for the construction of real-time math statements for a given phase cycle. The problem of creating a real-time math algorithm is largely a question of intuition, experience, and (mental) pattern recognition. Generating phase cycling math is also not an "exact science" in the sense that for any phase cycle there are many algorithms that lead to the correct result. The general goal is to create short and understandable algorithms, but what finally counts is the resulting phase cycle. For the execution of the pulse sequence in the pulse programmer, it normally doesn't matter *how* a phase cycle was generated.

Simple Phase Cycles

As we shall see later, most phase cycles (except perhaps the most simple ones) can be constructed by co-adding simpler (sub-)cycles: simple phase cycles can form the tools from which we can construct complex cycles. We should, therefore, first learn how to construct some simple phase cycles from ct.

Phase Incrementation

Incrementing phases are generated either from ct directly or by using division operation to slow ct down (the final mod4 operations are optional for quadrature phases):

```
assign(ct,v1); /* 0 1 2 3 */
mod4(ct,v1); /* 0 1 2 3 */
hlv(ct,v1); /* 0 0 1 1 2 2 3 3 4 4 5 5 ... */
mod4(v1,v1); /* 0 0 1 1 2 2 3 3 4 4 5 5 ... */
hlv(ct,v1); /* 0 0 1 1 2 2 3 3 */
hlv(ct,v1); /* [ 0 1 2 3 ]2 */
hlv(v1,v1); /* [ 0 1 2 3 ]4 */
hlv(v1,v1); /* [ 0 1 2 3 ]8 */
divn(ct,three,v1); /* [ 0 1 2 3 ]3 */
```

Alternating Phases

Alternating phases are obtained by doubling an incrementing sequence (of 90-degree phase shifts). Again, for quadrature phases, the final mod4 operations are optional:

| dbl(ct,v1); | /* | 0 | 2 | 0 | 2 */ |
|--|----------------|-------------|-------------|-------------|---|
| <pre>add(ct,ct,v1);</pre> | /* | 0 | 2 | 0 | 2 */ |
| <pre>hlv(ct,v1); dbl(v1,v1); mod4(v1,v1);</pre> | /* | 0 | 0 | 2 | 1 2 2 3 3 4 4 */ 2 4 4 6 6 8 8 */ 2 */ |
| <pre>hlv(ct,v1); mod2(v1,v1); dbl(v1,v1);</pre> | /* | 0 | 0 | 1 | 1 2 2 3 3 4 4 */ 1 0 0 1 1 */ 2 0 0 2 2 */ |
| hlv(v1,v1); dbl(v1,v2); | /* /* /* |] [] | 0 0 0 | 1 2 1 | 2 3]2 */ 2 3]4 */]4 */ 2 3]8 */]8 */ |
| <pre>divn(ct,three,v1); dbl(v1,v1); mod4(v1,v1);</pre> | /* | [| 0 | 2 | 4 6]3 */ |

Decrementing Phases

Decrementing phase cycles can be obtained by adding three copies of an incrementing phase cycle (multiplying an incrementing cycle by the constant three) or by subtracting an incrementing phase from a fixed phase:

```
add(ct,ct,vl); /* 0 2 4 6 ... */
add(ct,vl,vl); /* 0 3 6 9 ... */
mod4(vl,vl); /* 0 3 2 1 0 3 2 1 */
mult(ct,three,vl); /* 0 3 6 9 ... = 0 3 2 1 */
hlv(ct,vl); /* 0 0 1 1 2 2 3 3 ... */
```

```
add(three,one,v2); /* 4 (= 0) */

sub(v2,v1,v1); /* 4 4 3 3 2 2 1 1 0 0 -1 -1 */

mod4(v1,v1); /* [ 0 3 2 1 ]2 */

hlv(ct,v1); /* 0 0 1 1 2 2 3 3 4 ... */

sub(zero,v1,v1); /* 0 0 -1 -1 -2 -2 -3 -3 -4 ... */

mod4(v1,v1); /* 0 0 3 3 2 2 1 1 0 ... */
```

Shifted Pattern

```
add(ct,two,v1); /* 2 3 4 5 6 7 8 9 ... */

hlv(v1,v1); /* 1 1 2 2 3 3 4 4 ... */

hlv(v1,v1); /* 0 0 1 1 1 1 2 2 ... */

mod2(v1,v1); /* 0 0 1 1 1 1 2 2 ... */

sub(one,v1,v1); /* 1 1 0 0 0 0 1 1 */

mult(v1,three,v3); /* 0 0 3 3 3 3 0 0 */

add(ct,one,v1); /* 1 2 3 4 5 6 7 8 ... */

hlv(v1,v1); /* 0 1 1 2 2 3 3 4 ... */

dbl(v1,v1); /* 0 2 2 0 0 2 2 0 */
```

Complex Phase Cycles

The trick for generating complex phase cycles with real-time math is to recognize "internal periodicities," or repeating patterns *within* a phase cycle, and (more difficult) a shifted pattern in a phase cycle. The target is to "decompose" a complex phase cycle into statements that can either be generated directly or through one of the short algorithms presented above ("reverse synthesis"). In the pulse sequence, we can then compose the complex phase cycle simply by adding up the elements. We now look at a few examples (all used for quadrature phase shifting with a range of 0 to 3).

0 0 2 2 2 2 2 0 0 Pattern

This phase cycle can either be regarded as a shifted pattern (see above) or it can be thought of as being identical to $(0 \ 0 \ 2 \ 2 \ 2 \ 4 \ 4)$. We can split that into two groups of four phases, within which the phase is alternated: $(0 \ 0 \ 2 \ 2)$, $(2 \ 2 \ 4 \ 4)$. The second of these groups is shifted by two units compared to the first group. This leads to the following decomposition:

 $0 \ 0 \ 2 \ 2 \ 2 \ 2 \ 4 \ 4 \ = \ 0 \ 0 \ 2 \ 2 \ 0 \ 0 \ 2 \ 2 \ + \ 0 \ 0 \ 0 \ 0 \ 2 \ 2 \ 2 \ 2$

We can, therefore, generate this phase cycle from two simple cycles. Note that an intermediate result from the second part can be used for the first part:

hlv(ct,v1); /* 0 0 1 1 2 2 3 3 */ hlv(v1,v2); /* 0 0 0 1 1 1 1 1 2 2 2 3 3 */ dbl(v1,v1); /* 0 0 2 2 */ dbl(v2,v2); /* 0 0 0 0 2 2 2 2 2 2 */ add(v1,v2,v1); /* 0 0 2 2 2 2 4 4 */ mod4(v1,v1); /* 0 0 2 2 2 2 0 0 */

Chapter 10. Phase Calculations

The treatment as shifted pattern (as shown previously) is slightly more efficient:

| <pre>add(ct,two,v1);</pre> | /* | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | */ |
|----------------------------|----|---|---|---|---|---|---|---|---|--------|
| hlv(v1,v1); | /* | 1 | 1 | 2 | 2 | 3 | 3 | 4 | 4 | */ |
| hlv(v1,v1); | /* | 0 | 0 | 1 | 1 | 1 | 1 | 2 | 2 | */ |
| dbl(v1,v1); | /* | 0 | 0 | 2 | 2 | 2 | 2 | 4 | 4 | */ |

0213Pattern

This phase cycle can be interpreted as being composed from two subcycles $(0 \ 2)$ and $(0 \ 0 \ 1 \ 1)$, with the second two pulses shifted by one unit (90 degrees) against the first two. This leads to the following algorithm:

| dbl(ct,v2); | /* | 0 | 2 | 0 | 2 | */ | | | |
|----------------|----|---|---|---|---|-----|---|---|----|
| hlv(ct,v1); | /* | 0 | 0 | 1 | 1 | 2 2 | 3 | 3 | */ |
| mod2(v1,v1); | /* | 0 | 0 | 1 | 1 | */ | | | |
| add(v1,v2,v1); | /* | 0 | 2 | 1 | 3 | */ | | | |

0231 Pattern

Exchanging the last two phases in the previous pattern changes things completely. We get a cycle composed of the simple cycle $(0 \ 0 \ 1 \ 1)$ and a shifted pattern $(0 \ 2 \ 2 \ 0)$:

| add(ct,one,v2); | /* | 1 | 2 | 3 | 4 | */ | | | |
|-----------------|----|---|---|---|---|-----|---|---|----|
| hlv(v2,v2); | /* | 0 | 1 | 1 | 0 | */ | | | |
| dbl(v2,v2); | /* | 0 | 2 | 2 | 0 | */ | | | |
| hlv(ct,v1); | /* | 0 | 0 | 1 | 1 | 2 2 | 3 | 3 | */ |
| mod2(v1,v1); | /* | 0 | 0 | 1 | 1 | */ | | | |
| add(v1,v2,v1); | /* | 0 | 2 | 3 | 1 | */ | | | |

02132031 Pattern

Splitting this cycle, which is equivalent to $(0 \ 2 \ 1 \ 3 \ 2 \ 4 \ 3 \ 5)$, into two groups of four leads to the two subcycles, $(0 \ 0 \ 0 \ 2 \ 2 \ 2 \ 2)$ and $(0 \ 2 \ 1 \ 3 \ 0 \ 2 \ 1 \ 3)$, with the $(0 \ 2 \ 1 \ 3)$ part of the second subcycle being composed from two more subcycles, $(0 \ 0 \ 1 \ 1)$.

| dbl(ct,v2); | /* | 0 | 2 | 0 | 2 */ |
|----------------|----|---|---|---|--------------|
| hlv(ct,v1); | /* | 0 | 0 | 1 | 1 2 2 3 3 */ |
| hlv(v1,v3); | /* | [| 0 | 1 | 2 3]4 */ |
| dbl(v3,v3); | /* | [| 0 | 2 |]4 */ |
| mod2(v1,v1); | /* | 0 | 0 | 1 | 1 */ |
| add(v1,v2,v1); | /* | 0 | 2 | 1 | 3 */ |
| add(v1,v3,v1); | /* | 0 | 2 | 1 | 3 2 0 3 1 */ |

It turns out that there is a much simpler solution if we take four groups of two phases, in which case the elements are $(0 \ 2)$ and $(0 \ 0 \ 1 \ 1 \ 2 \ 2 \ 3 \ 3)$:

| hlv(ct,v1); | /* | 0 | 0 | 1 | 1 | 2 | 2 | 3 | 3 | */ |
|----------------|----|---|---|----|---|---|---|---|---|----|
| dbl(ct,v2); | /* | 0 | 2 | */ | / | | | | | |
| add(v1,v2,v1); | /* | 0 | 2 | 1 | 3 | 2 | 0 | 3 | 1 | */ |

0231201320130231 Pattern

Splitting this phase cycle into four groups of four phases we get an A-B-B-A pattern with $A=(0\ 2\ 3\ 1)$ and $B=(2\ 0\ 1\ 3)$. B is phase-shifted by 180 degrees (2 units); this becomes obvious if we think about the equivalent of 4 and 0 or 5 and 1 for quadrature phase shifting: B is (2 4 5 3). The decomposition of phase cycle (0 2 3 1 2 0 1 3 2 0 1 3 0 2 3 1) is equal to (0 2 3 1) plus (0 0 0 0 2 2 2 2 2 2 2 2 0 0 0 0), a simple shifted pattern, [0 2]8 left-shifted by four steps.

| add(ct,one,v2); | /* | v2 | = | 1 | 2 | 3 | 4 | */ | |
|-----------------|----|----|---|---|----|---|---|----|--------------------------|
| hlv(v2,v2); | /* | v2 | = | 0 | 1 | 1 | 0 | */ | |
| dbl(v2,v2); | /* | v2 | = | 0 | 2 | 2 | 0 | */ | |
| hlv(ct,v1); | /* | v1 | = | 0 | 0 | 1 | 1 | 2 | 2 3 3 */ |
| mod2(v1,v1); | /* | v1 | = | 0 | 0 | 1 | 1 | */ | |
| add(v1,v2,v1); | /* | v1 | = | 0 | 2 | 3 | 1 | */ | |
| dbl(two,v3); | /* | v3 | = | 4 | */ | | | | |
| add(ct,v3,v4); | /* | v4 | = | 4 | 5 | б | 7 | 8 | 9 10 11 12 13 */ |
| dbl(v3,v3); | /* | v3 | = | 8 | */ | | | | |
| divn(v4,v3,v4); | /* | v4 | = | 0 | 0 | 0 | 0 | 1 | 1 1 1 1 1 1 1 2 */ |
| dbl(v4,v4); | /* | v4 | = | [| 0 | 2 | 2 | 4 | 4 6 6 8 8 10]4 |
| | | v4 | = | [| 0 | 2 | 2 | 0 |]4 */ |
| add(v1,v4,v1); | /* | v1 | = | 0 | 2 | 3 | 1 | 2 | 4 5 3 2 4 5 3 0 2 3 1 |
| | | v1 | = | 0 | 2 | 3 | 1 | 2 | 0 1 3 2 0 1 3 0 2 3 1 */ |

Note that even for quadrature phases the "implicit mod4 function" can only be applied if no division is involved (after a division it is all right again)!

03212103210303211032321032101032Pattern

The best approach for this phase cycle is to start by splitting into two groups, and then splitting the resulting 16-step cycle into four groups of four steps each:

These are all subcycles we have discussed previously: $(0 \ 3 \ 2 \ 1)$ is a simple decrementing cycle, $[0 \ 2 \ 2 \ 0]4$ was used in the previous phase calculation, and the subcycle $[0 \ 1]16$ can be calculated as $(ct/16) \mod 2$:

```
dbl(two,v1);
                                 /* v1 = 4 */
add(ct,v1,v2);
                                /* v2 = 4 5 6 7 8 9 10 11 12 13 .. */
                                /* v1 = 8 */
dbl(v1,v1);
divn(v2,v1,v2); /* v2 = 0 0 0 0 1 1 1 1 1 1 1 1 2 .. */
dbl(v2,v2); /* v2 = [ 0 2 2 4 4 6 6 8 8 10 ... ]4
                                    v2 = [ 0 2 2 0 ]4 */
dbl(v1,v1);
                                /* v1 = 16 */
divn(ct,v1,v3);
mod2(v3,v3);
                                /* v3 = ct/16 = [ 0 1 2 3 ]16
                                /* v3 = [ 0 1 ]16

      sub(zero,ct,v4);
      /* v4 = 0 -1 -2 -3 = 0 3 2 1 */

      add(v4,v2,v4);
      /* v4 = 0 3 2 1 2 1 0 3 2 1 0 3 0 3 2 1

      add(v4,v3,v4);
      /* v4 = 0 3 2 1 2 1 0 3 2 1 0 3 0 3 2 1

                                /* v4 = 0 3 2 1 2 1 0 3 2 1 0 3 0 3 2 1 */
                                             1 0 3 2 3 2 1 0 3 2 1 0 1 0 3 2 */
```

This phase cycle is longer than the previous one, but it requires fewer math statements. If we skip the mod2 step, we would even obtain a phase cycle of 64 steps, which uses less Acode (39 words) than if we would store the resulting phase table in the Acode (as 16-bit numbers)! Divisions normally multiply the phase cycle length by the divisor, multiplications and modulo steps usually shorten the length of a phase cycle.

Phase Cycles for Many Pulses

By tradition, Varian pulse sequences were mostly written with a code section containing the real-time phase math for all pulses (and the receiver), followed by the actual pulse sequence. This has the advantage of keeping things apart for more clarity, and also it often permits using intermediate results or the final phase cycle for one pulse also for the calculation of other phase cycles, reducing the total number of math statements for the entire pulse sequence. On the other hand, if there are many pulses with many different phase cycles, there may be a problem with the limited number of real-time variables (14 total). There may be more phase cycles in a single pulse sequence than real-time variables. Also, such a convoluted, complex phase cycling section is often difficult to decode.

The other extreme would be to calculate each phase *just before using it*. This keeps the individual phase calculations simpler and avoids conflicts with the number of real-time variables (after having calculated and used one phase cycle, all variables can be reused for the next calculation). On the other hand, this method may dramatically increase the total number of math statements in a pulse sequence, can make it difficult to read (pulses and delays are interspersed with phase calculations), and blows up the Acode. In particular, this method also neglects the phase *relations* within a pulse sequence—the fact that intermediate results and entire phase cycles can often be reused for other phase calculations is not incidental, but due to the close relationship between all phases in a pulse sequence! More about that in Chapter 11, "Phase Tables," on page 115.

10.4 Real-Time Logical Decisions

The Acode interpreter not only has the ability of performing a linear interpretation, but it can also make decisions and perform branching based on the contents of real-time variables (a looping capability is also built in, see Chapter 14, "Repeating Events," on page 147). The use of real-time decisions for the conditional creation of FIFO words (conditional pulses or pulse sequence fragments) is discussed in more detail in Section 15.2, "Real-Time Decisions," on page 160. Here, we just discuss the application of real-time decisions to the construction of phase cycles using real-time math.

Let's look at a phase cycle used in a phase-cycling-based implementation of the E.COSY experiment (basically a linear combination of different double-quantum filtered COSY experiments). In this experiment, the first two pulses undergo a complicated scheme of 45-degree phase shifts:

0 1 0 7 0 1 2 7 0 1 4 7 0 1 6 7 0 1 0 7 0 1 2 7 0 1 4 7 0 3 6 5

There must be a direct mathematical way to construct this kind of phase cycle, but that algorithm would be excessively complicated. Instead, let us take a different approach by separating the phases for the even and odd scans (the phases are in 45-degree increments, the range of values therefore is 0 to 7) :

| 0 | 1 | 0 | 7 | 0 | 1 | 2 | 7 | 0 | 1 | 4 | 7 | 0 | 1 | 6 | 7 | 0 | 1 | 0 | 7 | 0 | 1 | 2 | 7 | 0 | 1 | 4 | 7 | 0 | 3 | 6 | 5 | = |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 0 | | 0 | | 0 | | 2 | | 0 | | 4 | | 0 | | 6 | | 0 | | 0 | | 0 | | 2 | | 0 | | 4 | | 0 | | б | | |
| | 1 | | 7 | | 1 | | 7 | | 1 | | 7 | | 1 | | 7 | | 1 | | 7 | | 1 | | 7 | | 1 | | 7 | | 3 | | 5 | |

The new (sub-)phase cycles look much more manageable than the original one, at least the subcycle for the odd pulses (for the even scans we need some exception handling, because the last two phases behave differently from the others).

What we now need is some construct that can be described as follows:

```
if (odd scan)
   (use phase cycle a); /* 0 0 0 2 0 4 0 6 */
else
   (use phase cycle b); /* 1 7 1 7 1 7 1 7 */
```

Fortunately, this is very easy to do. We just use a real-time variable that is zero for all the odd transients by applying mod2 to ct, and then we can use real-time branching to calculate specific phases for both the odd and even transients (for more information in real-time decisions, see Section 15.2, "Real-Time Decisions," on page 160):

```
mod2(ct,v1); /* 0101 */
ifzero(v1);
    ... /* phase selection for odd transients */
elsenz(v1);
    ... /* phase selection for even transients */
endif(v1);
```

Using this mechanism, we can switch between two phase cycles with every scan. Unlike shown above, the partial phase cycles do not have gaps, but the "empty positions" are filled with suitable numbers (such that the partial phase cycles are easy to calculate).

 Full cycle:
 0
 1
 0
 7
 0
 1
 4
 7
 0
 1
 6
 7
 0
 1
 0
 7
 0
 1
 4
 7
 0
 1
 6
 7
 0
 1
 0
 7
 0
 1
 4
 7
 0
 1
 6
 7
 0
 1
 0
 1
 6
 5

 Cycle
 A:
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1
 1
 7
 7
 1

The two resulting subcycles are definitely easier, but still not trivial to construct using standard math—but we can use exactly the same algorithm to construct the two subcycles:

Chapter 10. Phase Calculations

Subcycle A can be derived from two trivial phase cycles by switching between the two secondary subcycles with every pair of scans. Subcycle B can be constructed in a similar way (x stands for a non-zero value in the flag variable):

Cycle B: 1 1 7 7 1 1 7 7 1 1 7 7 1 1 7 7 1 1 7 7 1 1 7 7 1 1 7 7 1 1 7 7 1 1 7 7 3 3 5 5 Cycle E: 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 3 3 5 5 7 Cycle F: 1 1 7 7 1 1 1 7 7 1 1 7 7 1

In the actual coding, we first construct the four secondary subcycles (one of them subcycle C—is trivial) and the three flag variables, and then construct the final phase cycle with a series of nested real-time if statements (this allows bypassing the intermediate construction of the subcycles A and B). The flag variable "E/F" is best created as a shifted incrementing cycle:

Flag E/F: 1 1 1 1 2 2 2 2 3 3 3 3 4 4 4 4 5 5 5 5 6 6 6 6 7 7 7 7 0 0 0 0

This is constructed using modulo 8 of ct/4+1. The entire phase "calculation" can be coded as follows:

```
mod2(ct,v10);
                             /* v10 = 0 1 = flag A/B */
hlv(ct,v1); /* v1 = ct/2 ^/
mod2(v1,v11); /* v11 = 0 0 1 1 = flag C/D */
hlv(v1,v12); /* v12 = ct/4 */
' v12 = ct/4 + 1 */
                            /* v1 = ct/2 */
incr(v12); /* v12 = ct/4 + 1 */
dbl(two,v13); /* 4 */
dbl(v13,v13); /* 8 */
modn(v12,v13,v12) /* [ 1 2 3 4 5 6 7 0 ]2 = flag E/F */
v2 = 0 0 1 1 */
/* v2 = 0 0 2 2 */
mult(v2,three,v3);
add(v2,three,v2);
incr(v3);
/* v3 = 0 0 6 6 */
/* v2 = 3 3 5 5 = 0
/*
                            /* v2 = 0 0 1 1 */
                            /* v2 = 3 3 5 5 = subcvcle E */
incr(v3);
hlv(v1,v1);
dbl(v1,v1);
modn(v1,v13,v1);
                            /* v3 = 1 1 7 7 = subcycle F */
                            /* v1 = ct/4 */
                             /* v1 = [ 0 2 4 6 8 10 12 ... ]4 */
                             /* v1 = [ 0 2 4 6 ]4 = subcycle D */
    zero(v10); /* odd scans */
ifzero(v11); /* scans 1(,2),5(,6),etc. */
ifzero(v10);
       assign(zero,v1); /* subcycle C */
    endif(v11); /* (subcycle D is default in v1) */
    senz(v10); /* even scans */
ifzero(v12); /* scans 30,32 */
elsenz(v10);
       assign(v2,v1); /* select subcycle E */
    endif(v12);
endif(v10);
stepsize(45.0, TODEV);
xmtrphase(v1);
 . . .
```

The full E.COSY phase cycle is six times longer than the one shown above, with yet another "exception" in the sixth loop. Although it is obviously simpler to generate such phase cycles using tables, the main intention here was to demonstrate the possibility of generating "arbitrary" phase cycles using real-time math in combination with real-time decisions.

10.5 Steady-State Phase Cycling

All phase calculations up to now were based on the ct real-time variable. This works fine for all standard pulses up to nt=32768. After that many transients, the phase cycling continues as if ct would restart at 0, because the real-time variable ct is only the low-order half of LC->ct. For scan 32768, LC->ct is set to 32768, which is 1 in the high-order half (word) and 0 in the low-order word. This should be fine, because the number of steps in most phase cycling schemes is in powers of two anyway; even if this is not the case. After that many scans, a disruption in the phase cycling has no noticeable effect on the final spectrum.

Nevertheless it is a *bad idea* to use ct as exclusive basis for phase cycling. During steady-state pulses, the ct counter remains at 0, and therefore no phase cycling occurred during steady-state transients! This means that no matter how big a number of steady-state scans is selected, the steady state during the ss pulses is *different* from the steady state reached after the first few real scans. In all but the most primitive experiments (s2pul with cp='n'), the true steady state is only reached after the first "real" scans! This is particularly bad with cancellation experiments, unless some presaturation scheme (homospoil-90-homospoil, or two long orthogonal spin locking pulses at the beginning of each scan) is used to erase all residual coherence from previous scans.

The proper solution is *not* to use ct directly for phase calculations, but to construct a counter that varies though the steady-state pulses. This is possible by using the ssctr real-time variable. ssctr is equal to the value of ssval for the first steady-state scan and is decremented after each ss scan. When the ssctr variable is zero, the real scans start. Various schemes have been proposed to combine the ct and ssctr counters, like the following:

After this, v9 would be used instead of ct for all phase calculations. This scheme takes into account the fact that ssctr "counts backwards" (it is decremented, instead of incremented like ct) so phase cycling starts the same way as with the "real" scans.

This scheme is certainly better than using ct only, but it still has some deficiencies. After the steady-state scans, the phase cycling counter jumps back to zero and may, therefore, disrupt the steady state. This can again be particularly bad with cancellation experiments (e.g., using an odd number of steady-state scans with double-quantum filtered experiments). In principle, this scheme requires setting ss to the total length of the phase cycle! The idea behind steady-state scans is to "make the spins believe that we have been performing an infinite number of scans in the past," before the "real" scans start. This can be achieved by performing *the last elements* of the phase cycle rather than the beginning elements during the steady-state scans. This can be achieved in a very simple way:

sub(ct,ssctr,v9); /* -ss,-ss+1,.. -1,0,1,.. nt-1 */

This works fine, because the signed binary value -1 is equivalent to the unsigned binary value 1,111,111,111,111,111 (decimal 65535) or 2^{n} -1. As long as the length of the phase cycle is a power of two, this will execute the last elements of a phase cycle. As there are only very few examples of pulse sequences where the phase cycle length is

not a power of two (COSY-3 is one of them), this should be made the default way to generate a base counter for phase cycling.

10.6 C Constructs and Phase Calculations

C constructs are used to implement parameter-dependent phase cycling. The most frequently used parameter-dependent phase cycling is the type for phase-sensitive multidimensional NMR, where a coding similar to the type below is applied. A pseudo-variable phase is used to differentiate between the TPPI (time-proportional phase incrementation) and hypercomplex (States-Haberkorn-Ruben) types of phase-sensitive 2D (*n*D) spectra. The Varian convention is to use phase=3 for TPPI spectra (adding 90 degrees to the phase shift of the first pulse with every time increment) and phase=1, 2 for hypercomplex experiments, where phase=1 is unshifted, phase=2 has 90 degrees added to the pulses prior to the evolution phase (and relevant to the phase of the observed signal), as well as to the observe phase. In the coding below, v1 is supposed to be the phase of the pulse prior to the evolution (e.g., in a NOESY experiment):

```
int t1_counter = (int) (d2 * swl + 0.1);
int phase1 = (int) getval("phase");
if (phase1 == 2)
    incr(v1);
else if (phase1 == 3);
{
    initval((double) t1_counter), v10);
    add(v1,v10,v1);
}
```

This is only the most basic implementation for phase-sensitive NMR. Modern sequences use refinements for the hypercomplex method that are discussed in more detail in Chapter 19, "Multidimensional Experiments," on page 215.

There is another nice example of the use of C constructs for implementing variable phase cycles: in the relayh sequence's relayed COSY variant (relay greater than 0), in which the number of refocusing elements is determined by the parameter relay. The refocusing pulses cannot (or at least should not) be cycled only synchronously, because this would lead to an accumulation of the errors due to imperfections in the refocusing pulses. In the Varian implementation of this sequence, C constructs are used in two places to adjust the phase cycling to the actual setting of the relay parameter:

- The length of the 00112233 phase cycling is calculated in a C for loop.
- The phase cycling for the refocusing periods is calculated "on the fly", while coding the relay (refocusing) intervals such that the 90- degree pulse after the first relay has the fastest phase alternation and the phase inversion of 90-degree pulses following subsequent relay periods is slowed down progressively (by a factor of two per relay period).

This is a slightly simplified version of this pulse sequence:

```
/* relayh - relayed cosy, including regular cosy */
pulsesequence()
{
    int i, relay = (int) (getval("relay") + 0.5);
    double tau = getval("tau");
```

```
/* v1 = 00112233 */
hlv(ct, v1);
for (i = 0; i < relay + 1; i++)
hlv(v1,v1);
dbl(ct,oph);
add(ct,v1,v2);
add(oph,v1,oph)
                                /* [ 0 1 2 3 ]2**(relay+2) */
                                /* \text{ oph} = 0202 \ 0202 \ */
                                /* v2 = 0123 0123 + 00112233 */
                                /* oph = 0202 0202 + 00112233 */
                                 /* v3 = 0011 2233 4455 ... */
hlv(ct,v3);
status(A);
   hsdelay(d1);
                                 /* preparation period */
status(B);
   pulse(pw,v1);
                               /* evolution period */
   delay(d2);
   pulse(pw,v2);
                               /* start of mixing period */
   for (i = 0; i<relay; i++) /* relay coherence */</pre>
   {
      delay(tau/2);
      pulse(2.0*pw,v2);
      delay(tau/2);
      hlv(v3,v3);
dbl(v3,v4);
                              /* v3 = [ 0 1 2 3 ]2**(relay+1) */
                              /* v4 = [ 0 2 ]2**(relay+1) */
                              /* v5 = v4+v2 (with 00112233) */
      add(v2,v4,v5);
      pulse(pw,v5);
   }
status(C);
```

10.7 Why Phase Calculations?

}

There is no doubt that real-time phase calculations are not the simplest way to generate phase cycling in pulse sequences, because it may be difficult (some users may even call it painful) for inexperienced users to decode phase calculations, or even more to generate a real time calculation algorithm for a given phase table if a pulse sequence is published with the phase cycles in the form of numeric tables. On the other hand, there are definitely also arguments for phase calculations (as opposed to using tables):

- For very long phase cycles, the calculations are more efficient in terms of Acode space (see also Chapter 11, "Phase Tables," on page 115).
- Scientists who generate new pulse sequences do *not* think in phase tables but rather in algorithms (like phase cycling individual pulses and/or the observe phase for a specific coherence level selection, multiple-quantum filtering, etc.). Real-time phase calculations offer a way to directly implement such algorithms.
- Understanding the internal phase cycling mechanisms from long numeric tables may be as difficult (if not more difficult) than trying to understand real-time math. Changing the phase cycling order with tables may be more difficult than with real-time math.

10.8 Real-Time Random Numbers

There are many examples of a systematic variation of a variable during the execution of a pulse sequence, such as:

- Incrementation of the evolution time as function of the increment number.
- Systematic incrementation or decrementation of the mixing time or associated time intervals in some NOESY-type experiments.
- Time-proportional phase incrementation for phase-sensitive *n*D experiments.

There are also examples of random variation of variables, like the variation of the mixing time in some (other) types of NOESY experiments. In this case, the randomization happens from increment to increment, and a C construct can be used to generate individual random numbers for every code segment:

```
#define CONSTANT 1073741824
double rv;
if (ix == 1)
    srandom(getpid());
rv = ((double) (random() - CONSTANT))/((double) CONSTANT);
mix = mix * (1.0 + rv * getval("mixvar"))
```

The randomizer must be initialized once (with the first increment) by calling the srandom function (calling it with the current process-ID ensures true randomization, also between subsequent calls of the same pulse sequence). The function random then generates a random number between 0 and 2,147,483,647 (the maximum positive 32-bit integer). mixvar is a scaling factor (between 0 and 1.0 in this case, some implementation may use a scaling between 0 and 100.0) that allows controlling the degree of randomization (the proportion of the variation range and the total duration of the delay. The implementation above varies the delay in both directions.

In some rare cases—like some Z-filters, or randomization of d1 for avoiding phasing problems in the case of long T_2 -relaxation (sharp lines) and short repetition rates—such a randomization should occur within a single increment, and systematic variation is not desired or not possible. This is a tricky issue, because the acquisition operating system (the Acode interpreter) does not provide a random number generator.

This problem has two solutions: randomization using a table of random numbers (see also Section 11.7, "Using Tables as Source for Random Numbers," on page 135) and randomization using a real-time pseudo-random number generator.

An algorithm for the second solution has been presented in *Magnetic Moments*¹; the coding below is adapted from that article:

```
double random = getval("random");
if ((d1 - random/2.0) > 0.2e-6)
  delay(d1 - random/2.0);
if (random > 0.0)
{
  initval((double) ((ix*13) % 511),v14);
                                                     /* seed */
                     /* v13 = 4 = phasecycle */
  initval(4.0,v13);
                                                 /* v7 = 8 */
  dbl(two,v7); dbl(v7,v7);
  mult(v7,v7,v7);
                                                 /* v7 = 64 */
                                                /* v8 = 128 */
  dbl(v7,v8);
  decr(v7);
                                                 /* v7 = 63 */
```

¹ R. Boyko, B. Sykes & G. Gray, *Magnetic Moments*, III.3, p.4.

This method generates random numbers in the range of 0 to 126. The authors of the algorithm verified that this mechanism generates nearly ideally distributed numbers within the specified range. The initialization of v14 (the random seed) with a function of the increment number ix also ensures true randomization between different traces. Unlike the variation using a C randomization based on the process-ID, this type of randomizer generates the *same random numbers with every* go *command*. Ultimately, the real-time random number is used to determine the (real-time) loop cycles over a small delay (1/126 of the total variation range). Loops are discussed in detail in Chapter 14, "Repeating Events," on page 147.

Note that in the above code, the real-time variable v14 (the random seed) is initialized using the initval statement. This number is then modified *by purpose* (in the vast majority of the cases this is "forbidden"). The random variation occurs at the beginning of every phase cycle. It is a good idea to complete every phase cycle under identical conditions in order to ensure proper subtraction and cancellation.

}

Chapter 10. Phase Calculations

Chapter 11. Phase Tables

The basic syntax and the mechanisms of using phase tables are simple and easy to understand, construct, and use.

11.1 Basic Syntax

A simple, external ASCII file, such as the file /vnmr/tablib/sample of ~/vnmrsys/tablib/sample, takes up the actual phase tables:

```
/* sample phase table file */
t1 = 0 0 0 0 1 1 1 1 2 2 2 2 3 3 3 3
    /* 1st pulse */
t2 = 0 2 1 3 1 3 2 0 2 0 3 1 3 1 0 2
    /* 2nd pulse, oph */
```

The table file syntax is nearly trivial. Comments are defined the same way as in a C program:

- Tables have names between t1 and t60.
- Table values are positive integer numbers (16-bit, range between 0 and 32767), separated by spaces.
- Table values are separated from the table name by an *equal sign* (=) *that must be surrounded by spaces*.
- Table length is arbitrary. If necessary, tables can be extended over several lines.

In the pulse sequence, the table file must be loaded explicitly (specifying its name); thereafter, the specified tables (t1, t2 in the above example) can be used directly:

```
#include "standard.h"
pulsesequence()
{
    loadtable("sample");
    status(A);
    hsdelay(d1);
    status(B);
    pulse(p1,t1);
    hsdelay(d2);
    status(C);
    pulse(pw,t2);
    setreceiver(t2);
}
```

The table file name is relative. In this example ~/vnmrsys/tablib/sample has preference over /vnmr/tablib/sample. The table names t1 and t2 are used like real-time variables. They are of type codeint, but unlike real-time variables do *not* point to some element in the LC structure. As we will see later, the table names are used in the C code only to differentiate between the various tables (the table names are lost in the Acode, see below). The receiver phase is defined in the real-time variable oph. If the receiver phase should be taken from a phase table, oph must be set from a table. This is usually done with the setreceiver statement, but the same can also be achieved by the getelem statement:

```
getelem(t2,ct,oph);
```

This explicitly extracts element number ct from the table t2. getelem can also be used to extract phases from tables into other real-time variables, for example:

getelem(t1,ct,v10);

This way, real-time math can be performed on phases from phase tables (real-time math does not work with table names). Applications for this feature will be discussed in detail in Section 11.6, "Combining the Best of the Two Worlds," on page 129.

Shorthand Notation

Two types of shorthand notation exist for external phase table files: repeated table sequences and repeated table elements.

Repeated Table Sequences

By definition, any table is a repeated sequence. If the table index (the pointer to the current table element) is larger than the number of elements in the table, the table lookup automatically restarts at the beginning of the table. Thus, the table index is automatically taken *modulo of the table length*. This occurs individually for each table. For example, even if the total phase cycle length (in a pulse sequence) is 16, it is not necessary to define:

This would be a waste of Acode space. The much shorter table

t1 = 1 3

defines exactly the same table, but requires much less space. We don't need a shorthand notation for this case¹. On the other hand, a shorthand notation is useful in the case of a table like:

 t1
 =
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0
 2
 0

In shorthand notation this would be written as

t1 = (02)8 (13)8 (20)8 (31)8

The sequence between the parentheses is repeated as many times as indicated by the repetition number that must follow the closing parenthesis *without space* (the spaces between parentheses and the numbers inside are optional).

Repeated Table Elements

Very often tables consist of repeated elements, such as

t1 = 0 0 0 0 0 0 0 0 0 2 2 2 2 2 2 2 2 1 1 1 1 1 1 1 1 1 3 3 3 3 3 3 3 3 3

In shorthand notation, this can be dramatically simplified and written as

t1 = [0 2 1 3]8

¹ There is an exception: table-to-table math functions (ttadd, ttsub, ttmult, ttdiv) require that the target table (which is also one of the two operand tables) is at least as long as the second operand table. For this reason, we may have to write down a table sequence more than once (usually done using shorthand notation). See the manual *VNMR User Programming* for more information.

The syntax around the brackets is the same as with the parentheses, except that the index after the brackets indicates the number of repetitions for each table element between the brackets.

Brackets and parentheses can be used sequentially in within the same table but cannot be nested. Note that the shorthand notation is just an easy and simple way of *defining and entering* long phase tables—*the tables are internally generated at full length* (i.e., no Acode space is saved by using shorthand notation, see also below).

Advanced Features

Several advanced features, including the division factor and an automatically incrementing index, make working with tables easier.

Division Factor

The last table described above can also be written differently:

 $t1 = \{ 0 2 1 3 \} 8$

Again, the same syntax as with the parentheses and the brackets applies to the braces ("curly brackets"). This looks like exactly the same thing as the shorthand syntax using the brackets and, from a superficial point-of-view, this is true. The index after the braces defines how many times each element is repeated before passing to the next element in the table. However, internally the braces work differently. In the case of the braces, the table generated internally is only the table elements between the braces.

The index after the braces defines what is called the *division factor* (sometimes also called the *division return factor*). Before reading a table element, the table index (the pointer to the current table element) is *divided* by the division factor, before the modulo function (modulo the table length) is performed and the table element is extracted. Every table has a division factor. The default setting for the division factor is 1, and the braces simply alter that value.

In other words, the braces not only are a shorthand notation for tables with repeated elements, they also define a "shorthand table" internally (i.e., different from the brackets, the braces also save Acode space, see below). *Therefore, it is recommended to use braces instead of brackets wherever possible*².

Brackets and parentheses cannot be nested within themselves, but they can be nested within braces, as seen in the following examples:

 $t1 = \{ (0 2) 4 (1 3) 4 \} 4$ $t2 = \{ [0 2] 4 (1 3) 4 \} 4$

In conventional syntax, these tables would be written as follows:

| t1 | = | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 |
|----|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| | | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 |
| | | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 |
| | | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 |

² There is an exception to this as well. Table-to-table math operations (ttadd, ttsub, ttmult, and ttdiv) directly operate on the vectors as they are defined internally, without taking into consideration any division factors. This can lead to unexpected results if division factors different from 1 are used, apart from complications due to the restrictions with respect to the table lengths in such operations (see the previous footnote).

| t2 | = | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
|----|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| | | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| | | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 |
| | | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 |

Note that the division factor always applies to the entire table; whereas, brackets (and parentheses) can also be used on parts of a table only:

t1 = [0 2]4 2 0 2 0 2 0 2 0

This table is perfectly acceptable, but the same line with braces instead of brackets *is not*.

The maximum division factor is 64. For tables with a higher number of repetitions on the individual table elements, the following definitions can be used:

 $t1 = \{ 0 0 0 0 1 1 1 1 2 2 2 2 3 3 3 3 \} 64$ $t1 = \{ [0 2 1 3] 4 \} 64$

The division factor can not only be defined in the table file, but it can also be defined or altered within the pulse sequence, using the setdivnfactor statement, for example:

```
setdivnfactor(t1,16);
```

This feature is used frequently in the type of coding described in Chapter 19, "Multidimensional Experiments," on page 215.

Tables with Automatically Incrementing Index

The table index (the pointer to the current table element) is normally either specified explicitly (using the getelem statement), or the transient counter ct is automatically (implicitly) taken as index (when specifying a table address t1...t60 as phase argument to statements like pulse, txphase, etc.). This way, the table index is incremented automatically with every transient (implicit use) or set as specified by the index argument to getelem.

However, there are applications where we would like to "scan" through a entire table (or parts of it) during a single transient, for instance, when the table defines the amplitudes of a pulse shape or the phases of an explicitly programmed decoupling pattern. Using table features discussed up to now, this could be programmed in the following manner:

```
assign(one,v9);  /* table index */
initval(32.0,v10);  /* loop cycles (no. of table elements) */
loop(v10,v11);
getelem(t1,v9,v1);  /* extract table value */
rgpulse(pw,v1,0.0,0.0);  /* perform one pulse element */
incr(v9);  /* increment table index */
endloop(v10);
```

Much of the code above is for the pointer generation and incrementation after reading an using each table element in turn (note that the number of loop cycles is *not necessarily* equal to the number of table elements). To make it easier to use tables for reading out successive amplitude and phase elements from tables within a single pulse sequence pass, phase tables have been equipped with an *autoincrement attribute*. This attribute is unset by default. It can be set in the table file by using "+=" instead of the simple equals sign "=" between the table name and the table elements:

 $t1 += \{ 0 2 2 0 2 0 0 2 \} 4$

Alternatively, the autoincrement attribute can be set with the setautoincrement statement (there is no function to unset the autoincrement attribute):

```
setautoincrement(t1);
```

With this feature, a "scan through a phase table" can be realized in a much simpler way:

Note that with the autoincrement attribute the table index is reset *once per FID* (when creating the code). If a table is not read out completely during one scan, the table index is incremented from *somewhere within the table* with the next scan.

Autoincrementing tables can be used not only directly (e.g., as phase argument to a pulse or txphase statement) but also through the setreceiver and getelem statements; however, in the case of setreceiver and getelem, the specified table index is *disregarded*.

How Does a Table Work?

To the user, a table consists of a table *name*, table *values*, and the *division factor* and *autoincrement* attributes (*autoincrement* is a flag). At the C level, each table has two more attributes: the current table *index* and the table *length*. The table index is only used in connection with the autoincrement attribute. Both the table length and the division factor are used to extract a table element either from an implicit or a specified table index:

```
extract_index=((index / division_factor) % table_length)
```

If the extraction index is *not* specified ("direct use" of phase tables as argument to statements like pulse, txphase, etc.), ct is used as table index (unless the table has the autoincrement attribute set). The variable ct is the default table index.

With the getelem statement, the table value is extracted into the specified real-time variable. With the setreceiver function, it is extracted into the oph real-time variable. Where tables are used "directly," the table value is extracted into a dedicated real-time variable tablert (the address of LC->tablert):

```
txphase(t3);
```

is equivalent to

```
getelem(t3,ct,tablert);
txphase(tablert);
```

In other words, statements that take a real-time variable as an argument check whether the real-time variable is a table name or a normal variable, and then take the corresponding actions.

For example, the txphase statement can be written as follows³:

```
(void) txphase(ph)
codeint ph;
{
    if ((ph >= t1) && (ph <= t60))</pre>
```

³ The actual coding is quite different (a setquadphase statement does not actually exist), but the equivalent to this construct is still used in the real software (the statement setphase90 in /vnmr/psg/rfchan_device.c).

```
{
   getelem(ph,ct,tablert);
   setquadphase(tablert, TODEV);
}
else
   setquadphase(ph, TODEV);
}
```

This also works with autoincrementing tables, because this simply causes the ct index to be disregarded at a lower level (in the Acode interpreter).

Note that with the standard VNMR software the default table index is ct (e.g., whenever a table that is not autoincrementing is used directly). Therefore, there is *no phase cycling during steady-state pulses*. It is strongly recommended to not use the table variables directly, but to use getelem with a pseudo ct counter:

```
sub(ct,ssctr,v13);
getelem(t1,v13,v1);
getelem(t2,v13,t2);
```

Extracting table values into the *corresponding* real-time variables ($t1 \rightarrow v1$, $t2 \rightarrow v2$, etc.) makes it easier to read and understand the sequence.

11.2 Inline Phase Tables

It is possible to avoid a separate file for tables by defining the tables as (static) arrays in the pulse sequence source file itself between the header (include lines) and the function pulsesequence:

```
static int table1[8] = {0,2,1,3,3,1,2,0};
static int table2[16] = {0,0,0,0,1,1,1,1,2,2,2,2,2,3,3,3,3};
```

The name of the static variable is arbitrary. This static array of integers is then converted to a table using the settable statement:

```
settable(t1,8,table1);
settable(t4,16,table2);
```

In contrast to the table definition using external files, tables defined through settable can also contain negative numbers (i.e., numbers in the full range of signed 16-bit integers, from -32768 to 32767). Division factors and the autoincrement attribute have to be set separately (after the call to settable), using the setautoincrement and setdivnfactor statements.

With the exception of cases where negative numbers are required in a table (e.g., if the table is used to set the amplitude of pulsed field gradients), using inline tables in general is not recommended because of the following reasons:

- The definition syntax is more complicated (there is lint checking on it, though),
- No shorthand syntax is available.
- Division factors and autoincrement attribute cannot be set together with the table definition (it is less obvious what the tables actually are if division factors are involved).
- Changing phase tables requires recompiling the pulse sequence.

11.3 Table Math

Tables cannot be created or defined only as static objects. Scalar and vector math are possible with tables. Four statements are provided for scalar operations on tables:

```
tsadd(tablename,scalarvalue,modulovGalue);
tssub(tablename,scalarvalue,modulovalue);
tsmult(tablename,scalarvalue,modulovalue);
tsdiv(tablename,scalarvalue,modulovalue);
```

The argument tablename is one of the table addresses, t1 to t60. The operation (addition, subtraction, multiplication, division) is applied to every value in the table. At the end, a modulo function is applied to the result, unless the last argument (moduloval) is set to zero, in which case the numeric range of the resulting table often exceeds the numeric range of the table prior to the execution of the scalar operation (in the vast majority of the cases, this has no adverse effects due to the implicit modulo operation that occurs when using the table). Obviously, the divisor (second argument) to tsdiv cannot be zero, because division by zero would lead to mathematical overflow.

There is also a set of table-to-table (vector-to-vector) statements available for more complex math with tables:

```
ttadd(target_table,operand_table,modulovalue);
ttsub(target_table,operand_table,modulovalue);
ttmult(target_table,operand_table,modulovalue);
ttdiv(target_table,operand_table,modulovalue);
```

The mathematical operation in these statements is applied to each "equivalent" pair of numbers in two tables, and the result is stored in one of the two specified tables. It is not possible to generate new tables using table math. Also, the operand table cannot be longer than the target table, because table math cannot alter the length of tables.

The following example should illustrate the functionality of table-to-table math. Let's assume we have defined the following two tables:

t1 = 0 0 0 0 2 2 2 2t2 = 0 1 2 3

In the pulse sequence, we now apply the following function:

```
ttadd(t1, t2, 0);
```

This causes each element in table t2 to be added to the corresponding element in table t1. If the operand table is shorter than the target table, the shorter table is expanded to the same length before performing the operation. The result of the above function is the following table t1:

t1 = 0 1 2 3 2 3 4 5

(no modulo value was specified; therefore, values above 3 are obtained). The other table-to-table operations use the same working principle. For ttdiv, the operand table must not contain zeroes, of course, because this would lead to math overflow.

There is a fundamental difference between table math and real-time math using variables. Different from the real-time variables, table math occurs in the host computer, when executing the pulse sequence code. That means that a table may not be changed *within* a pulse sequence. Once a table has been built into the Acode (see below), it cannot be modified. In fact, a table cannot be modified using table math after

it has been used the first time (be it as argument to a pulse or phasing function, through getelem or setreceiver calls)—this will lead to a run-time error message.

This precludes the alteration of tables within composite pulses and similar pulse sequence elements (like BIRD pulses) that use multiple closely related ("parallel") phase cycles. For such cases it is either necessary to define separate tables or to use real-time math to derive the related phases (see also Section 11.6, "Combining the Best of the Two Worlds," on page 129).

11.4 Phase Tables in the Acode

Tables are incorporated into the instruction segment of the Acode *at the point where they are used for the first time*. Unused tables are discarded. The table values are stored as 16-bit integers. Evidently, long tables can occupy a lot of Acode space.

Besides the table values, four more 16-bit integers per table define the table length, the autoincrement flag, the division factor, and the table index. The table name itself is not at Acode level. Tables are identified by their address in the instruction segment (actually the offset of the first attribute). The table size attribute also permits the Acode interpreter to calculate the offset to the next instruction after the table.

Let's first have a look at the code generated by a getelem statement (extracted from the Acode for a noesy experiment):

```
36325739 ASSIGNFUNC
369
                                         zero v5
372
        366260105TABLE 261size 2, autoinc 0, divn_ret 1, ptr 0
                        0 2
379 373267106TASSIGN table 261 ct v1
383 377271105TABLE 272size 2, autoinc 0, divn_ret 16, ptr 0
                        0 2
390 384278106TASSIGN table 272 ct v2
394 388282105TABLE 283size 4, autoinc 0, divn_ret 2, ptr 0
                       0 1 2 3
403 397291106TASSIGN table 283 ct v3

        407
        40129529
        ADDFUNC
        v1 v5 v1

        411
        40529929
        ADDFUNC
        v1 v2 oph

        415
        40930329
        ADDFUNC
        v3 oph oph

419 413307105TABLE 308size 2, autoinc 0, divn_ret 8, ptr 0
                      0 1
426 420314106TASSIGN table 308 ct v4
430 42431829 ADDFUNCo ph v4 oph
434 42832229 ADDFUNC
                                       v1 v4 v1

        438
        43232629
        ADDFUNC
        v2
        v4
        v2

        442
        43633029
        ADDFUNC
        v3
        v4
        v3
```

Each of the four getelem calls in this part of the pulse sequence generates a TABLE instruction (code 105), which inserts the table into the Acode as described above, followed by a TASSIGN instruction (code 106), which extracts a table value into a real-time variable. TASSIGN is followed by the table address, the table index (ct in this case), and the target address (a real-time variable in LC). For each of the tables, the TABLE instruction is inserted at the point where a table is first referred to.

The setreceiver statement generates the same kind of Acode (except that the target address is LC->oph):

 846
 840
 734
 105
 TABLE
 735
 size
 4, autoinc
 0, divn_ret
 1, ptr
 0

 0
 1
 2
 3
 3
 3
 3
 55
 849
 743
 106
 TASSIGN
 table
 735
 ct oph

Even using tables directly (here as phase arguments to a simpulse statement) does not require a new Acode instruction. The target address is LC->tablert in this case.

603 597 491 150 HighSpeedLINES (void) 606 600 494 105 TABLE 495 size 4, autoinc 0, divn_ret 8, ptr 0 1 2 3 0 615 609 503 106 TASSIGN table 495 ct tblrt 619 613 507 16 SETPHAS90 CH1 tblrt 622 616 510 105 TABLE 511 size 4, autoinc 0, divn_ret 1, ptr 0 0 1 2 3
 631
 625
 519
 106
 TASSIGN
 table
 511 ct tblrt

 635
 629
 523
 16
 SETPHAS90
 CH2 tblrt
 638 632 526 150 HighSpeedLINES (void) 641 635 529 150 HighSpeedLINES XOFF 644 638 532 150 HighSpeedLINES RXOFF DECRG 647 641 535 150 HighSpeedLINES RXOFF DECRG 650 644 538 151 EVENT1_TWRD 1.000 usec 652 646 540 150 HighSpeedLINES RXOFF DECRG DEC 655 649 543 151 EVENT1_TWRD 7.925 usec 657 651 545 150 HighSpeedLINES RXOFF TXON DECRG DEC 660 654 548 151 EVENT1_TWRD 27.725 usec 662 656 550 150 HighSpeedLINES RXOFF DECRG DEC 665 659 553 151 EVENT1_TWRD 7.925 usec 667 661 555 150 HighSpeedLINES RXOFF DECRG 670 664 558 150 HighSpeedLINES DECRG 673 667 561 150 HighSpeedLINES (void)

11.5 Tables vs. Real-Time Calculations

In this section, we make a point-by-point comparison of tables and real-time calculations. The comparison is illustrated by examples.

Point-to-Point Comparison

Tables and real-time calculations can be compared on length, complexity, Acode space required, understandability, constructionally, and changeability.

Possible Length, Complexity

Real-time calculations permit calculating phase cycles that are virtually unlimited in length (a phase cycle with 1024*1024 steps is no problem!), but limited in complexity ("arbitrary" phase cycles may require an enormous amount of real-time calculations that are difficult to construct and understand). *Tables* are limited in length (up to 8192 steps per table, or about 9500 steps maximum per pulse sequence in total, excluding division factors), but of unlimited complexity.

Acode Space Consumption

Real-time calculations are inefficient for relatively short and simple phase cycles like 0 2 1 3 3 1 2 0, but long phase cycles can often be constructed with relatively few math statements. *Tables* are efficient as long as they are short. Long tables can consume a lot of Acode space. Acode space efficiency can be crucial in 3D and 4D experiments (although these normally have very short phase cycles).

Chapter 11. Phase Tables

Easy to Understand? Easy to Construct?

Finding out the phase cycle generated by *real-time calculations* is often not easy. Constructing a real-time math algorithm for a given, non-trivial phase table is definitely a difficult task for non-specialists. Given a complete (phase) table in literature, a *table file* is trivial to construct, and it also is trivial to see what the actual numeric sequence is. On the other hand, with longer tables it can be *more difficult* to understand what the *underlying phase cycling algorithms* are than with real-time math, because with the latter such algorithms can often be directly coded. For the designer of a new pulse sequence, real-time math may be the more adequate choice. NMR scientists often think more in terms of phase-cycling algorithms, rather than final phase tables (of course, in publications it is easier to simply print out phase tables rather than having to explain the various phase cycling elements and algorithms).

Easy to Change?

Changing a phase cycle generated by *real-time math* normally requires analyzing and rebuilding the entire algorithm and can be time-consuming. Changing a phase cycle from full *tables* first requires understanding and analyzing the underlying algorithms, and then changing (all) the tables. With long tables, this can be as complicated as changing real-time calculations, maybe even more difficult sometimes.

Comparison by Examples

To emphasize the differences, and in particular the advantages and disadvantages of the two methods for generating phase cycles, a pulse sequence with relatively long phase cycles (not excessive, though), inadqt.c, has been selected for the following comparison (*the sequence was simplified for the purpose of this chapter*).

Using Calculations Only

Let's first have a look at a sequence that is written with phase calculations only:

```
pulsesequence()
{
  double tau = 1.0 / (4.0 * getval("jcc"));
  int phase = (int) (getval("phase") + 0.5);
  mod4(ct, v3);
  hlv(ct, v9);
  hlv(v9, v9);
  assign(v9, v1);
  hlv(v9, v9);
  assign(v9, v2);
  hlv(v9, v9);
  hlv(v9, v9);
  mod2(v9, v9); /* v9 = F2 guad. image suppression */
  dbl(v1, v1);/* v1 = suppresses artifacts from
                          imperfections in 1st 90 deg. pulse */
  add(v9, v1, v1);
  assign(v1, oph);
  dbl(v2, v8);
  add(v9, v2, v2)i/* v2 = suppresses artifacts from
                         imperfect 180 refocusing pulse */
  add(v8, oph, oph);
  dbl(v3, v4);
```

```
add(v3, v4, v4);
add(v3, v9, v3);
add(v4, oph, oph);/* v3 = selects DQC during t1
                       evolution period */
assign(zero, v10);/* v10 = F1 guadrature */
if (phase == 2)
   incr(v10);
add(v10, oph, oph);
status(A);
   hsdelay(d1);
status(B);
   stepsize(45.0, TODEV);
   xmtrphase(v10);
   rgpulse(pw, v1, rof1, 0.0);
   delay(tau);
   rgpulse(2.0*pw, v2, rof1, 0.0);
   delay(tau);
status(C);
   rgpulse(pw, v9, rof1, 0.0);
   xmtrphase(zero);
   delay(d2);
   pulse(pw, v3);
status(D);
```

The phase calculation section is rather long. There are five different phases, and the phase cycle length is 64 steps. Even an experienced spectroscopist could spend many minutes figuring out what the phase cycle is for the various pulses. In the Acode, the phase calculations occupy 70 words, which isn't too bad for phase cycles of this length.

```
275
     26916398 NextSCan
276
     27016437 MOD4FUNCct v3
    27316734 HLVFUNCct v9
279
282 27617034 HLVFUNCv9 v9
285 27917339 ASSIGNFUNCv9 v1
288 28217634 HLVFUNCv9 v9
291 28517939 ASSIGNFUNCv9 v2
294 28818234 HLVFUNCv9 v9
297
     29118534 HLVFUNCv9 v9
300
     29418836 MOD2FUNCv9 v9
     29719133 DBLFUNCv1 v1
303
306 30019429 ADDFUNCv9 v1 v1
310 30419839 ASSIGNFUNCv1 oph
313 30720133 DBLFUNCv2 v8
316 31020429 ADDFUNCv9 v2 v2
320 31420829 ADDFUNCv8 oph oph
324 31821233 DBLFUNCv3 v4
327 32121529 ADDFUNCv3 v4 v4
331 32521929 ADDFUNCv3 v9 v3
     32922329 ADDFUNCv4 oph oph
335
339 33322739 ASSIGNFUNCzero v10
342 33623029 ADDFUNCv10 oph oph
346 340234 6 APBOUT2 items 0xa511 0xb57c
350 344238150 HighSpeedLINESDECUP
353 347241151EVENT1_TWRD1500 msec
355 349243150 HighSpeedLINESDECUP
358 35224668 PHASESTEP CH190 units (45.00 degrees)
361 35524965 SETPHASE CH1v10
364
     35825216 SETPHAS90 CH1v1
364 35825210 SEIFINASSO CHIVI
367 361255150HighSpeedLINESRXOFF DECUP
370 364258151EVENT1_TWRD40.000 usec
```

}

Chapter 11. Phase Tables

372 366260150HighSpeedLINESRXOFF TXON DECUP 375 369263151EVENT1 TWRD10.400 usec 377 371265150HighSpeedLINESRXOFF DECUP 380 374268150HighSpeedLINESDECUP 377271152EVENT2_TWRD6 msec + 250 usec 383 386 38027416 SETPHAS90CH1 v2 389 383277150 HighSpeedLINERXOFF DECUP 392 386280151EVENT1_TWRD40.000 usec 394 388282150 HighSpeedLINESRXOFF TXON DECUP 397 391285151EVENT1_TWRD20.800 usec 399 393287150HighSpeedLINESRXOFF DECUP 396290150HighSpeedLINESDECUP 402 405 399293152EVENT2_TWRD6 msec + 250 usec 408 402296150HighSpeedLINESDECUP 411 40529916 SETPHAS90CH1 v9 414 408302150 HighSpeedLINESRXOFF DECUP 417 411305151EVENT1_TWRD40.000 usec 419 413307150 HighSpeedLINESRXOFF TXON DECUP 422 416310151EVENT1_TWRD10.400 usec 424 418312150 HighSpeedLINESRXOFF DECUP 427 421315150 HighSpeedLINESDECUP 42431865 SETPHASECH1 zero 430 433 42732116 SETPHAS90CH1 v3 436 430324150 HighSpeedLINESRXOFF DECUP 439 433327151EVENT1_TWRD10.000 usec 441 435329150 HighSpeedLINESRXOFF TXON DECUP 444 438332151EVENT1_TWRD10.400 usec 446 440334150HighSpeedLINESRXOFF DECUP 449 443337150HighSpeedLINESDECUP 452 446340150HighSpeedLINESDECUP 44934316 SETPHAS90CH1 zero 455 458 45234616 SETPHAS90CH2 zero 461 455349150HighSpeedLINESRXOFF DECUP 464 458352151EVENT1_TWRD10.000 usec 466 460354150HighSpeedLINESDECUP 469 463357152EVENT2_TWRD122 usec + 350 nsec 472 46636090 SETInputCardMode 476 47036499 ACQXX loop np=4096, dwell 80.000 usec 479 473367 7 STartFIFO 480 47436897 HouseKEEPing 481 47536920 BRANCH Offset 159 ______ 483 words / 966 Bytes / Total code size 0.9 KB _____

The Acode segment for the first FID is 477 words or 954 bytes. We now want to compare this with another version of the same pulse sequence, where all the math statements have been replaced by phase tables.

Using "Pure Tables"

If we remove all the phase calculations and write a table file instead, the pulse sequence text is dramatically simplified. The only phase-related statements left are those that deal with the f_1 quadrature detection (see also Chapter 19, "Multidimensional Experiments," on page 215). Apart from the fact that tables are used instead of phase calculations, the functionality of the following version is exactly the same as with the version above (at least as far as the pulse programmer is concerned).

```
pulsesequence()
{
   double tau = 1.0 / (4.0 * getval("jcc"));
   int phase = (int) (getval("phase") + 0.5);
   loadtable("inadqt");
   setreceiver(t5);
   assign(zero,v6);
   if (phase == 2)
   {
      incr(v6);
      incr(oph);
   }
   status(A);
      hsdelay(d1);
   status(B);
      stepsize(45.0, TODEV);
      xmtrphase(v6);
      rgpulse(pw, t1, rof1, 0.0);
      delay(tau);
      rgpulse(2.0*pw, t2, rof1, 0.0);
      delay(tau);
   status(C);
      rgpulse(pw, t3, rof1, 0.0);
      xmtrphase(zero);
      delay(d2);
      pulse(pw, t4);
   status(D);
}
```

Shorthand syntax allows us to simplify the coding of all the phase tables involved. Unfortunately, only three of the five tables can be shortened using division factors two tables must be left at full length. These tables also reveal a potential disadvantage of the shorthand syntax: in order to correlate the different phase tables, we have to (mentally) translate the tables back to full length!

Very rarely phase cycling occurs isolated in one table only. Changing the phase of a pulse most likely causes the observe phase, if not also the phase of other pulses, to change as well. It would, therefore, help if the tables were written at full length, on a single line per table, such that the phase cycles could be vertically correlated. This would of course prohibit any shorthand syntax and division factors and would, therefore, make the tables even less efficient (in terms of Acode space consumption). Also, long (full) tables aren't necessarily more readable due to their length (each of them would fill two or four lines if the line length is limited to less than 80 characters).

In the Acode, the phase tables (at least those with a division factor of 1) occupy a considerable number of words:

| 275 | 269 | 163 | 98 | NextSCan |
|------------|------------|------------|------------|---|
| 276 | 270 | 164 | 105 | |
| | | | | |
| | | | | 2 1 0 3 0 3 2 1 |
| | | | | 0 3 2 1 2 1 0 3 |
| | | | | 2 1 0 3 0 3 2 1 |
| | | | | 1 0 3 2 3 2 1 0 |
| | | | | 3 2 1 0 1 0 3 2 |
| | | | | 1 0 3 2 3 2 1 0 |
| | | | | 3 2 1 0 1 0 3 2 |
| 345 | 339 | 233 | 106 | TASSIGN table 165 ct oph |
| 349 | 343 | 237 | 39 | ASSIGNFUNC zero v6 |
| 352 | 346 | 240 | б | APBOUT 2 items 0xa511 0xb57c |
| 356 | 350 | 244 | 150 | HighSpeedLINES DECUP |
| 359 | 353 | 247 | 151 | EVENT1_TWRD 1500 msec |
| 361 | 355 | 249 | 150 | HighSpeedLINES DECUP |
| 364 | 358 | 252 | 68 | PHASESTEP CH1 90 units (45.00 degrees) |
| 367 | 361 | 255 | 65 | SETPHASE CH1 v6 |
| 370 | 364 | 258 | 105 | TABLE 259 size 16,autoinc 0,divn_ret 4,ptr 0 |
| | | | | 0 2 0 2 0 2 0 2 |
| | | | | 1 3 1 3 1 3 1 3 |
| 391 | 385 | 279 | 106 | TASSIGN table 259 ct tblrt |
| 395 | 389 | 283 | 16 | SETPHAS90 CH1 tblrt |
| 398 | 392 | 286 | 150 | HighSpeedLINES RXOFF DECUP |
| 401 | 395 | 289 | 151 | EVENT1_TWRD 40.000 usec |
| 403 | 397 | 291 | 150 | HighSpeedLINES RXOFF TXON DECUP |
| 406 | 400 | 294 | 151 | EVENT1_TWRD 10.400 usec |
| 408 | 402 | 296 | 150 | HighSpeedLINES RXOFF DECUP |
| 411 | 405 | 299 | 150 | HighSpeedLINES DECUP |
| 414 | 408 | 302 | 152 | EVENT2_TWRD 6 msec + 250 usec |
| 417 | 411 | 305 | 105 | TABLE 306 size 8, autoinc 0,divn_ret 8, ptr 0 |
| | | | | 0 1 2 3 1 2 3 0 |
| 430 | 424 | 318 | 106 | TASSIGN table 306 ct tblrt |
| 434 | 428 | 322 | 16 | SETPHAS90 CH1 tblrt |
| 437 | 431 | 325 | 150 | HighSpeedLINES RXOFF DECUP |
| 440 | 434 | 328 | 151 | EVENT1_TWRD 40.000 usec |
| 442 | 436 | 330 | 150 | HighSpeedLINES RXOFF TXON DECUP |
| 445 | 439 | 333 | 151 | EVENT1_TWRD 20.800 usec HighSpeedLINES RXOFF DECUP |
| 447 450 | 441 444 | 335 | 150 150 | |
| 450 453 | 444 | 338 341 | 150 | HighSpeedLINES DECUP EVENT2_TWRD 6 msec + 250 usec |
| 455 | 450 | 344 | 152 | HighSpeedLINES DECUP |
| 459 | 453 | 347 | 105 | TABLE 348 size 2, autoinc 0, divn_ret 32, ptr 0 |
| 155 | 155 | 517 | 105 | 0 1 |
| 466 | 460 | 354 | 106 | TASSIGN table 348 ct tblrt |
| 470 | 464 | 358 | 16 | SETPHAS90 CH1 tblrt |
| 473 | 467 | 361 | 150 | HighSpeedLINES RXOFF DECUP |
| 476 | 470 | 364 | 151 | EVENT1_TWRD 40.000 usec |
| 478 | 472 | 366 | 150 | HighSpeedLINES RXOFF TXON DECUP |
| 481 | 475 | 369 | 151 | EVENT1_TWRD 10.400 usec |
| 483 | 477 | 371 | 150 | HighSpeedLINES RXOFF DECUP |
| 486 | 480 | 374 | 150 | HighSpeedLINES DECUP |
| 489 | 483 | 377 | 65 | SETPHASE CH1 zero |
| 492 | 486 | 380 | 105 | TABLE 381 size 64, autoinc 0,divn_ret 1,ptr 0 |
| | | | | 0 3 2 1 0 3 2 1 |
| | | | | 0 3 2 1 0 3 2 1 |
| | | | | 0 3 2 1 0 3 2 1 |
| | | | | 0 3 2 1 0 3 2 1 |
| | | | | 1 0 3 2 1 0 3 2 |
| | | | | 1 0 3 2 1 0 3 2 |
| | | | | 1 0 3 2 1 0 3 2 |
| | | | | 1 0 3 2 1 0 3 2 |

| 561 | 555 | 449 | 106 | TASSIGN table 38 | 1 ct tblrt |
|-------|--------|------|---------|------------------|-----------------------|
| 565 | 559 | 453 | 16 | SETPHAS90 | CH1 tblrt |
| 568 | 562 | 456 | 150 | HighSpeedLINES | RXOFF DECUP |
| 571 | 565 | 459 | 151 | EVENT1_TWRD | 10.000 usec |
| 573 | 567 | 461 | 150 | HighSpeedLINES | RXOFF TXON DECUP |
| 576 | 570 | 464 | 151 | EVENT1_TWRD | 10.400 usec |
| 578 | 572 | 466 | 150 | HighSpeedLINES | RXOFF DECUP |
| 581 | 575 | 469 | 150 | HighSpeedLINES | DECUP |
| 584 | 578 | 472 | 150 | HighSpeedLINES | DECUP |
| 587 | 581 | 475 | 16 | SETPHAS90 | CH1 zero |
| 590 | 584 | 478 | 16 | SETPHAS90 | CH2 zero |
| 593 | 587 | 481 | 150 | HighSpeedLINES | RXOFF DECUP |
| 596 | 590 | 484 | 151 | EVENT1_TWRD | 10.000 usec |
| 598 | 592 | 486 | 150 | HighSpeedLINES | DECUP |
| 601 | 595 | 489 | 152 | EVENT2_TWRD 22 u | usec + 350 nsec |
| 604 | 598 | 492 | 90 | SETInputCardMode | 2 |
| 608 | 602 | 496 | 99 | ACQXX loop np=40 | 96, dwell 80.000 usec |
| 611 | 605 | 499 | 7 | STartFIFO | |
| 612 | 606 | 500 | 97 | HouseKEEPing | |
| 613 | 607 | 501 | 20 | BRANCH | Offset 159 |
| | | | | | |
| ===== | | | ======= | | |
| Total | L code | size | 615 | words / 123 | 0 Bytes / 1.2 KB |
| ===== | | | | | |

In this case, the Acode segment for the first FID is 609 words or 1218 bytes (compared to 477 words or 954 bytes with phase calculations). The Acode size is increased by 254 bytes per FID (28%). This may seem rather marginal, however, we should also consider 2D and 3D experiments.

Furthermore, there are certainly cases of complex pulse sequences that generate a lot of Acode by themselves. In such cases, it is not impossible that phase tables can cause the Acode size to exceed the upper limit (10000 words), in which case a sequence would not execute at all. There are also 1D pulse sequences such as cccc.c (composite pulse inadequate) with a large number of different phase cycles of 256 or even up to 1024 steps each. We can certainly imagine cases where the tables themselves would not fit into the available Acode space unless they can be abbreviated.

11.6 Combining the Best of the Two Worlds

Obviously, neither pure calculations nor pure phase tables provide an optimum solution for programming phase cycles. What we are looking for is a way of defining the phase cycling that is at the same time all of the following:

- Easy to understand and analyze.
- Easy to create and compose.
- Easy to change the phase cycling scheme (this is mostly a feature for pulse sequence developers).
- Visualizes the internal phase cycling algorithms.
- Does not use excessive Acode space.

A possible solution to this problem is to use a combination of phase tables and realtime calculations. Why not use phase tables to define the (simple) *phase cycling elements* and real-time math to realize *phase cycling algorithms* (i.e., phase dependencies, synchronous phase variation, etc.)?

Chapter 11. Phase Tables

It turns out that most phase cycling elements are simple—the most common elements are the following:⁴

{ 0123}n
{ 02}n
{ 0321}n
{ 0131}n
{ 0131}n
{ 0131}n

Furthermore, most phase cycling algorithms (in the sense of phase cycling *relations*) are almost trivial. In most sequences, it simply involves adding up certain phases in order to obtain the observe phase and to add some "macroscopic" phase cycling to all phase cycles (at least for pulses on the observe channel).

The most striking case for using this solution is in the area of composite pulses. These pulses use phases that are normally simply orthogonal or anti-phase relative to each other, and it would be a real waste to define phase table for each of the pulses. It is much simpler to extract one phase and to calculate the others from that one.

How can we realize this concept? Let's take the same pulse sequence as used above, inadqt.c, analyze it, and see whether the new concept is feasible. We first start with the full phase table file from the last version:

In order to find out about phase relations, it is best (at least for beginners) to first write out the full phase cycles:

| t1 = | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 | 0 | 0 | 0 | 0 | 2 | 2 | 2 | 2 |
|------|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 | 1 | 1 | 1 | 1 | 3 | 3 | 3 | 3 |
| t2 = | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 3 | 3 | 3 | 3 | 3 | 3 | 3 | 3 |
| | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 3 | 3 | 3 | 3 | 3 | 3 | 3 | 3 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| t3 = | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| t4 = | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 |
| | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 |
| t5 = | 0 | 3 | 2 | 1 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 0 | 3 | 2 | 1 |
| | 1 | 0 | 3 | 2 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 1 | 0 | 3 | 2 | 1 | 0 | 3 | 2 | 3 | 2 | 1 | 0 | 3 | 2 | 1 | 0 | 1 | 0 | 3 | 2 |

The one fact common to *all* these phase tables is that after 32 steps all phase cycles, including the observe phase ± 5 , are repeated with a phase shift of 90 degrees (all tables are incremented by 1). This is the quad image suppression in f_2 (the observed dimension)⁵ and is the slowest of the phase cycles in this sequence. It turns out that ± 3

⁴ One could even further simplify this list, saying that the most frequently used phase cycling elements are just { 0 1 }n and { 0 2 }n; { 0 1 2 3 }n and { 0 2 1 3 }n are just combinations of the two simpler ones: { 0 1 2 3 }n = { 0 1 }n + { 0 2 }(2n), and { 0 2 1 3 }n = { 0 2 }n + { 0 1 }(2n). The cycle { 0 3 2 1 }n is just the reversal of { 0 1 2 3 }n.

⁵ This forms part of a CYCLOPS phase cycling (A.D. Bain, *J. Magn. Reson.* **56**, 418 (1984); D.I. Hoult, R.E. Richards, *Proc. Roy. Soc. London Ser.* A **344**, 311 (1975)), which is present only partially in this pulse sequence. The "other half" of the full cycle ($\{ 0 \ 2 \ \} 64$), which would reduce any axial signal ("center glitch") in f₂, is not present (nor required) in this case.

consists of only this phase cycle. It can, therefore, be regarded as constant phase (0), with quad image suppression added. If we "subtract" (isolate) that phase cycle from all tables, we end up with a much simpler list of tables. The subtracted phase cycle is stored in a new table that we can add to all the other tables using real-time math:

You might think that we have already reached a sufficient degree of simplification (t1 and t2 can be shortened dramatically using division factors), but this is not the end. The observe phase often is simply a function of some (or all) of the pulse phases. In this case, it is obvious that both t1 and t4 (and t3?) are contained in the observe phase cycle. If we subtract these phases from the observe phase, we obtain

If we then write this cycle as

we recognize that it is in fact twice the values of t2! This mechanism has been called EXORCYCLE⁶. Each time the phase of a refocusing pulse is altered by 90 degrees, the phase of the observed signal (and hence the receiver phase) changes by 180 degrees. Therefore, the observe phase in the inadgt experiment can be expressed as

t5 = t1 + 2*t2 + t4

In fact, we don't need to define the long observe phase cycle as a table at all. This phase will be calculated in the pulse sequence. The phase table can, therefore, be reduced to the following list:

/* inadqt phase tables */

```
/* 1st 90 */
t1 = \{ 0 2 \} 4
t2 = \{ 0 1 2 3 \} 8
                             /* 180 */
                             /* 2nd 90 */
t3 =
      0
t4 = 0321
                             /* 3rd 90 */
t6 = \{ 0 1 \} 32
                            /* macroscopic phase cycling */
/* phase cycling is performed in the following order:
1) t4 (nt=4): double-quantum coherence selection (4-step cycle)
2) t1 (nt=8): suppression of artifacts due to imperfections in
               the first 90 (2-step cycle)
3) t2 (nt=32): EXORCYCLE phase cycling to remove artifacts from
               imperfect spin refocusing (4-step cycle)
4) t6 (nt=64): f2 quad image suppression through macroscopic phase
                cycling (2-step cycle) */
```

In the pulse sequence, we first create a pseudo ct counter to also obtain steady-state phase cycling. With this counter as index, we then extract all phases from the tables, calculate the observe phase using the formula we have obtained from the analysis, add the "macroscopic" { 0 1 }32 phase cycle (t6), and finally adjust the phases for f_1 quadrature detection.⁷

⁶G. Bodenhausen, R. Freeman, and D.L. Turner, J. Magn. Reson. 27, 511 (1977).

Chapter 11. Phase Tables

{

}

```
pulsesequence()
   double tau = 1.0 / (4.0 * getval("jcc"));
   int phase = (int) (getval("phase") + 0.5);
   loadtable("inadqt");
   sub(ct, ssctr, v10);
                               /* pseudo ct counter */
                               /* extract phases from tables */
   getelem(t1, v10, v1);
   getelem(t2, v10, v2);
   getelem(t3, v10, v3);
   getelem(t4, v10, v4);
   getelem(t6, v10, v6);
   dbl(v2, oph);
                               /* calculate observe phase */
   add(oph, v1, oph);
   add(oph, v4, oph);
   add(v1, v6, v1);
                               /* add macroscopic phase cycling */
   add(v2, v6, v2);
   add(v3, v6, v3);
   add(v4, v6, v4);
   add(oph, v6, oph);
   assign(zero, v5);
                               /* f1 quadrature (phase=1,2) */
   if (phase == 2) incr(v5);
   add(oph, v5, oph);
   status(A);
     hsdelay(d1);
   status(B);
      stepsize(45.0, TODEV);
      xmtrphase(v5);
      rgpulse(pw, v1, rof1, 0.0);
      delay(tau);
      rgpulse(2.0*pw, v2, rof1, 0.0);
      delay(tau);
   status(C);
      rgpulse(pw, v3, rof1, 0.0);
      xmtrphase(zero);
      delay(d2);
      pulse(pw, v4);
   status(D);
```

You might argue that the new calculation section complicates the sequence again (compared to the version with "pure tables"). Also, we again are using what you might call "real-time calculations that are difficult to understand." However, a closer inspection of the above code shows that there is no black magic in the real-time math statements used here. We simply add up the various phases-the math operations involved are almost trivial. Also, you have to realize that (with the exception of some pulsed field gradient and some trivial 1D experiments) the phase cycling to a pulse sequence is usually as important as the pulses and the delays themselves. A phase cycling section as long as the example above (especially if it also nicely visualizes the internal phase cycling algorithms and relations involved) is certainly justified.

⁷ Remember that this pulse sequence has been simplified for this chapter. More explanation on how to manipulate the phase cycles in order to achieve quadrature detection in various modes in nD spectra is given in Chapter 19, "Multidimensional Experiments," on page 215.

Finally, we can double-check how much Acode we are using by combining phase tables and real-time calculations:

| 275 | 269 | 163 | 98 | NextSCan | |
|-----|-----|-----|-----|-----------------|-----------------------------------|
| 276 | 270 | 164 | 30 | SUBFUNC | ct ssctr v10 |
| 280 | 274 | 168 | 105 | | 2, autoinc 0,divn_ret 4, ptr 0 |
| | | | | 0 2 | |
| 287 | 281 | 175 | 106 | TASSIGN table 1 | 69 v10 v1 |
| 291 | 285 | 179 | 105 | TABLE 180 size | 4, autoinc 0,divn_ret 8,ptr 0 |
| | | | | 0 1 2 | 3 |
| 300 | 294 | 188 | 106 | TASSIGN table 1 | 80 v10 v2 |
| 304 | 298 | 192 | 105 | TABLE 193 size | 1, autoinc 0, divn_ret 1, ptr 0 0 |
| 310 | 304 | 198 | 106 | TASSIGN table 1 | 93 v10 v3 |
| 314 | 308 | 202 | 105 | TABLE 203 size | 4, autoinc 0, divn_ret 1, ptr 0 |
| | | | | 032 | |
| 323 | 317 | 211 | 106 | TASSIGN table 2 | 03 v10 v4 |
| 327 | 321 | 215 | 105 | | 2, autoinc 0,divn_ret 32,ptr 0 |
| | | | | 0 1 | , |
| 334 | 328 | 222 | 106 | TASSIGN table 2 | 16 v10 v6 |
| 338 | 332 | 226 | 33 | DBLFUNC | v2 oph |
| 341 | 335 | 229 | 29 | ADDFUNC | oph v1 oph |
| | | | | | |
| 345 | 339 | 233 | 29 | ADDFUNC | oph v4 oph |
| 349 | 343 | 237 | 29 | ADDFUNC | v1 v6 v1 |
| 353 | 347 | 241 | 29 | ADDFUNC | v2 v6 v2 |
| 357 | 351 | 245 | 29 | ADDFUNC | v3 v6 v3 |
| 361 | 355 | 249 | 29 | ADDFUNC | v4 v6 v4 |
| 365 | 359 | 253 | 29 | ADDFUNC | oph v6 oph |
| 369 | 363 | 257 | 39 | ASSIGNFUNC | zero v5 |
| 372 | 366 | 260 | 29 | ADDFUNC | oph v5 oph |
| 376 | 370 | 264 | 6 | APBOUT 2 items | 0xa511 0xb57c |
| 380 | 374 | 268 | 150 | HighSpeedLINES | DECUP |
| 383 | 377 | 271 | 151 | EVENT1_TWRD | 1500 msec |
| 385 | 379 | 273 | 150 | HighSpeedLINES | DECUP |
| 388 | 382 | 276 | 68 | PHASESTEP CH1 | 90 units (45.00 degrees) |
| 391 | 385 | 279 | 65 | SETPHASE | CH1 v5 |
| 394 | 388 | 282 | 16 | SETPHAS90 | CH1 v1 |
| 397 | 391 | 285 | 150 | HighSpeedLINES | RXOFF DECUP |
| 400 | 394 | 288 | 151 | | 40.000 usec |
| 402 | 396 | 290 | 150 | _ | RXOFF TXON DECUP |
| 405 | 399 | 293 | 151 | EVENT1_TWRD | 10.400 usec |
| 407 | 401 | 295 | 150 | — | RXOFF DECUP |
| 410 | 404 | 298 | 150 | HighSpeedLINES | DECUP |
| 413 | 407 | 301 | 152 | EVENT2 TWRD | 6 msec + 250 usec |
| 416 | 410 | 304 | 16 | SETPHAS90 | CH1 v2 |
| | | | | | |
| 419 | 413 | 307 | 150 | | RXOFF DECUP |
| 422 | 416 | 310 | 151 | EVENT1_TWRD | 40.000 usec |
| 424 | 418 | 312 | 150 | HighSpeedLINES | RXOFF TXON DECUP |
| 427 | 421 | 315 | 151 | EVENT1_TWRD | 20.800 usec |
| 429 | 423 | 317 | 150 | HighSpeedLINES | RXOFF DECUP |
| 432 | 426 | 320 | 150 | HighSpeedLINES | DECUP |
| 435 | 429 | 323 | 152 | EVENT2_TWRD | 6 msec + 250 usec |
| 438 | 432 | 326 | 150 | HighSpeedLINES | DECUP |
| 441 | 435 | 329 | 16 | SETPHAS90 | CH1 v3 |
| 444 | 438 | 332 | 150 | HighSpeedLINES | RXOFF DECUP |
| 447 | 441 | 335 | 151 | EVENT1_TWRD | 40.000 usec |
| 449 | 443 | 337 | 150 | HighSpeedLINES | RXOFF TXON DECUP |
| 452 | 446 | 340 | 151 | EVENT1_TWRD | 10.400 usec |
| 454 | 448 | 342 | 150 | HighSpeedLINES | RXOFF DECUP |
| 457 | 451 | 345 | 150 | HighSpeedLINES | DECUP |
| 460 | 454 | 348 | 65 | SETPHASE | CH1 zero |
| 463 | 457 | 351 | 16 | SETPHAS90 | CH1 v4 |
| 466 | 460 | 354 | 150 | HighSpeedLINES | RXOFF DECUP |
| 469 | 463 | 357 | 151 | EVENT1_TWRD | 1.000 usec |
| 471 | 465 | 359 | 150 | HighSpeedLINES | RXOFF TXON DECUP |
| 474 | 468 | 362 | 151 | EVENT1_TWRD | 10.400 usec |
| 1/1 | 100 | 502 | 191 | 2,0011_100 | 10.100 abec |

| 476 | 470 | 364 | 150 | U. abcood INEC | RXOFF DECUP |
|-------|--------|--------|-----|-----------------|------------------------|
| - • • | - • • | | | HighSpeedLINES | |
| 479 | 473 | 367 | 150 | HighSpeedLINES | DECUP |
| 482 | 476 | 370 | 150 | HighSpeedLINES | DECUP |
| 485 | 479 | 373 | 16 | SETPHAS90 | CH1 zero |
| 488 | 482 | 376 | 16 | SETPHAS90 | CH2 zero |
| 491 | 485 | 379 | 150 | HighSpeedLINES | RXOFF DECUP |
| 494 | 488 | 382 | 151 | EVENT1_TWRD | 10.000 usec |
| 496 | 490 | 384 | 150 | HighSpeedLINES | DECUP |
| 499 | 493 | 387 | 152 | EVENT2_TWRD | 122 usec + 350 nsec |
| 502 | 496 | 390 | 90 | SETInputCardMod | e |
| 506 | 500 | 394 | 99 | ACQXX loop np=4 | 096, dwell 80.000 usec |
| 509 | 503 | 397 | 7 | STartFIFO | |
| 510 | 504 | 398 | 97 | HouseKEEPing | |
| 511 | 505 | 399 | 20 | BRANCH | Offset 159 |
| | | | | | |
| ===== | | ====== | | | |
| Total | L code | size | 513 | words / 10 | 26 Bytes / 1.0 KB |
| ===== | | | | | |

Here, the Acode segment for the first FID is 507 words (1014 bytes), compared to 477 words, 954 bytes, with pure real-time calculations, and 609 words, 1218 bytes, with "pure tables": we have almost reached the Acode space efficiency of real-time math, while providing a C code that is superior to both other approaches in terms of clarity and simplicity, while making it easier to change things and experiment with different versions of the phase cycling.

Supposing we decided that the phase cycling on the first pulse is less important than the EXORCYCLE (on the refocusing pulse); therefore, we wanted to cycle the phase of the refocusing pulse before alternating the phase of the first pulse. Instead of changing complex phase tables or rebuilding complicated real-time math, we could simply supply an alternate phase table file:

| t1 | = | { | 0 | 2 | }16 | | /* | 1st | 90 | */ | | | |
|----|---|---|---|---|-----|----|----|------|------|-------|-------|---------|----|
| t2 | = | { | 0 | 1 | 23 | }4 | /* | 180 | */ | | | | |
| t3 | = | | 0 | | | | /* | 2nd | 90 | */ | | | |
| t4 | = | | 0 | 3 | 2 1 | | /* | 3rd | 90 | */ | | | |
| tб | = | { | 0 | 1 | }32 | | /* | macr | cosc | copic | phase | cycling | */ |

We wouldn't even have to change the text of the pulse sequence.

This is, of course, mostly a point for pulse sequence designers; "standard" users seldom alter the phase cycling of a pulse sequence. Also, the sequence of importance and relative order of phase cycling certainly is a non-trivial issue. In cancellation experiments (like double-quantum filtering or double-quantum spectroscopy like the inadqt experiment shown above), we always do the main cancellation first, then we would like to start cancelling the biggest artifacts (like the artifacts from off-resonance effects, pulse missettings and imperfections, rf inhomogeneity) and gradually proceed to the minor artifacts (like f₂ quadrature images and axial signals).

The reason for this order is that in the case of changing conditions (such as environmental or instrumental instabilities, or a decaying sample), the cancellation quality suffers if excessive time passes between the two scans that are supposed to subtract the artifact from itself (by alternating its phase). The best chance for a good cancellation is if the two scans follow each other immediately; hence, double-quantum filters and the like are performed with first priority. Minor artifacts leave only very small residual signals, even if they are imperfectly cancelled⁸.

Obviously, changing existing pulse sequences (whether they are written with real-time math or with full phase tables) to the new style would require first analyzing the phase cycling, which (as mentioned before) may be a non-trivial task for many users; however, we do not suggest that people rewrite all existing (working) pulse sequences. It's more the idea to provide a new, powerful concept for people that generate *new* pulse sequences.

11.7 Using Tables as Source for Random Numbers

As an alternative to generating random numbers within a single FID using real-time calculations (e.g., for Z-filtering, or for d1-randomization, see also Section 10.8, "Real-Time Random Numbers," on page 112), such random numbers can, of course, also be taken from a table. The table solution provides the option of using a fairly large range of random numbers (up to the full or positive range of 16-bit integers, depending on the method of definition). The length of the sequence of random numbers from tables is almost unlimited (up to 8192 numbers), but this length has to be "paid" in Acode space (2 bytes per number). For small numeric values, the method using real-time calculations is superior, because long sequences can be generated with a limited number of math statements, but this solution is probably limited to a linear distribution within a given range.

For non-linear (e.g., Gaussian) distributions or for random numbers with large numeric ranges (e.g., values between 0 and 4095), the table approach is the preferable solution. Linearly distributed random numbers are best generated using a suitable C construct, such as the following for 1024 values between 0 and 4095:

```
static int randomtable[1024];
pulsesequence()
{
    int i;
    if (ix == 1) srandom(getpid());
    for (i = 0; i < 1024; i++)
        randomtable[i] = \
            (int) (4095.999*((double) random())/2147483647.0);
    settable(t1,1024,randomtable);
    ...</pre>
```

For non-linear distributions (like a Gaussian distribution of values), a suitable C construct has to be found. For an approximate Gaussian distribution, the data can also be taken from an FID with noise only (no signal). The corresponding procedure would be as follows (1024 numbers, Gaussian distribution between 0 and 4095):

 In VNMR, set the following parameters and then acquire an FID: pw=0 (avoid signal) sw=100000 (the noise increases with sw) np=n (where n is the desired number of random numbers) nt≥1 (for 16-bit ADC) or nt≥4 (for 15-bit ADC) or nt≥256 (for 12-bit ADC)

⁸ The question of the optimum order in the phase cycling depends on the relative intensity of the artifacts (which may be a question of parameter selection, or a question of rf and probe quality), and on the environmental and instrumental stability. There is no global recipe for the phase cycling order. It can be assumed that the VNMR pulse sequences have been optimized in this respect; however, depending on the conditions and the instrument configuration, an alternate phase cycling order may sometimes improve the results.

ss=0 increase the gain such that the ADC is almost filled

• In UNIX (Bourne shell script):

```
#!/bin/sh
cd
cd
cd vnmrsys/expn/acqfil
echo -n "t50 = " > $HOME/vnmrsys/tablib/randomtable
od -iw2 fid +74 | awk `
    {    if ((NR > 0) && (((NR - 1) % 8) == 0)) printf("\n\t")
        printf(" %7d",($2+32768.0)/65535.0*4095.999)
    }' >> $HOME/vnmrsys/tablib/randomtable
```

This generates a table file (table ±50) with eight numbers per line.

Chapter 12. AP Bus Traffic

States that are driven by the fast bits can change instantaneously (at least as far as the programming goes), they can be set and unset within time intervals as short as 200 nanoseconds, with a timing precision of 25 (or 100) nanoseconds. It would be nice if all spectrometer states could be changed that quickly! Unfortunately, with instruments and experiments becoming more and more complex, for a fully equipped modern spectrometer this would require over a thousand fast lines. Also, pulse programming would be rather inefficient considering the fact that with most time events only one or two of the fast bits would be altered, and yet we would have to feed the pulse programmer with over 100 bytes for every FIFO word (and mostly the same bytes over and over again). The system designers, therefore, divided the devices into three classes:

- Devices that require very fast switching (down to 200 nanoseconds for switching in both directions): typically rf gates and 90-degree phase shifts (see Section 9.2, "Fast Bits," on page 88).
- Devices that need to be switched with accurate timing, but not necessarily at the speed of the fast bits, devices that would take longer anyway to switch to another state, and devices that are not switched on and off in rapid succession, but rather keep their state for many time events (or the entire sequence). These are the units for which the AP (analog port) bus was built.
- Devices that do not need to be synchronized with the pulse sequence, such as the sample changer, spinner control, sample insert and eject, and the like (these aren't even handled by the pulse programmer, see Chapter 13, "Acquisition CPU Communication," on page 145).

It was decided to limit the fast bits to those devices that really need the speed of the state bits. Most other devices are addressed by the AP bus.

12.1 What Is the AP Bus

As the name indicates, the AP bus is a true bus structure. In computer terminology, a *bus* involves a linear structure (a cable, sometimes a computer backplane) with several parallel lines (16 in the case of the AP bus). The lines are used to transmit address and control information as well as data to a (potentially large) number of devices that are all hooked up to the AP bus with special interface chips (the AP bus chip). Some devices (like attenuators and the linear amplifiers on UNITY systems) share a common interface to the AP bus, the AP interface card.

The information sent over the AP bus is mostly numeric; a maximum of 8 bits (1 byte) of information can be transmitted per AP word. Depending on the device, two different transmission modes can be used:

- Direct binary information, as used by attenuators, phase and amplitude modulators, shim or pulsed field gradient DAC values, etc. Depending on the way the hardware works (positive or negative logic), this can be normal or bit-inverted (negated) binary information.
- BCD (binary coded decimal) information. Here, every digit of a decimal number is transmitted in a separate binary AP word. This may sound inefficient (only 4 out

of 8 data bits are used per AP word), but it is adequate for some devices like the PTS frequency synthesizers and other similar devices (decoupler modulator [dmf], offset synthesizers, etc.). These devices operate in (frequency) *decades* and would otherwise have to generate the decimal information internally.

Apart from numeric values, the information transmitted through the AP bus can also be of simple binary nature, like the "command" or flag that sets a linear amplifier into CW mode or back into pulsed mode, or a "strobe", a command that causes devices to carry out previously transmitted numeric information. The strobe is used for smallangle phase shifting on the UNITY*plus* (where the phase value is transmitted in two words but only carried out when the second word is received), or in the case of PTS frequency synthesizers with latching: With these, all decades are first transmitted in BCD mode, but the frequency only changes when the strobe is received in a separate AP word¹.

For the vast majority of the cases, the AP bus operates in "indirect" mode (i.e., the first word sent over the AP bus contains the address and the following words contain the information to be transferred to the device). In this mode, the address fills 12 bits of the first AP word. Thanks to this long address, the number of devices that can be driven by the AP bus is virtually unlimited (4096 addresses are available). In very few cases (e.g., small-angle phase shifting on UNITY and earlier systems), the AP bus is used in direct addressing mode, where address and data are transmitted with the same AP word. This mode is limited to 16 devices (4 address bits).²

The AP words pass the pulse programmer within normal FIFO words. The difference to a time event is that the timing part of the FIFO word in this case contains the AP bus bits instead of time and time base. This part of the FIFO word is fed either into the timer(s) or into the AP bus, depending on a special bit within the FIFO word that differentiates timer words and AP words. The definition of the AP bus is such that every word must be held on the bus lines for 1 microsecond (UNITY*plus*) or 2 microseconds (UNITY and earlier systems) to allow for the devices to decode the address and read the information off the bus. AP words are, therefore, timed (accurately) by the pulse programmer. It takes 150 nanoseconds to decode the timer or AP word and initialize the timer circuitry; the total time per AP word is therefore 1.15 microseconds on a UNITY*plus* and 2.15 microseconds on all earlier systems. This is the AP bus cycle time.

The time to transmit information to any device is always a multiple of the bus cycle time. For instance, it takes 9 AP words to change the frequency on a PTS synthesizer on a UNITY*plus*. Therefore, it takes 10.35 microseconds to change a frequency with the offset statement. Complete tables for the number of FIFO and AP words involved, as well as for time requirements are given in the manual *VNMR User Programming*.

With respect to the fast bits, AP words behave like normal (fixed length) time events. The fast bit part of the FIFO words is filled with the current settings of the fast bits.

¹ On PTS synthesizers without latching, each decade changes (or keeps) its frequency value.

 $^{^2}$ The addressing is actually more complex than that. There are 16 "board addresses" (4 bits) that are shared between devices operating in direct mode (not used on UNITY*plus*) and devices that are addressed via AP bus chip. For each of the board addresses (AP bus chips), there are 256 register subaddresses, which allows for a variety of functions to be performed. Each AP bus chip can address many devices.

Fast bits are slightly different on older systems with an *output board* (63-word loop FIFO). On these systems, the fast bits maintain the settings of the *last executed time event*. If an AP word would immediately follow a pulse without post-pulse delay, that pulse would be prolonged by the duration of the AP bus event(s).

To avoid this, AP functions that generate AP words also generate a 200-nanosecond delay preceding the AP words. The total duration of AP events on these systems (early VXR systems) is, therefore, 0.2 microseconds longer than on systems with the acquisition control board (UNITY). The programmer can use the apovrride statement to turn off this additional 0.2-microsecond delay for the *next* AP event (e.g., for the case of multiple AP events in sequence or if the fast bits were already set correctly from the previous delay).

The AP bus overall is a fast transport medium: even in BCD mode, the equivalent of over 3.5 million bits (UNITY*plus*) or 1.8 million bits (on earlier systems) can be transmitted per second (in binary mode the transfer rate is twice as great).

12.2 What Devices are Driven by the AP Bus?

As instruments become more complex, the number of devices addressed by AP bus has been growing constantly over the years. In VXR spectrometers, the following devices were driven through the AP bus:

- PTS frequency synthesizers (not on fixed-frequency rf channels).
- Offset synthesizers with fixed frequency and broadband (as opposed to direct synthesis) rf channels.
- Decoupler high-power level (class C decoupler amplifier).
- Decoupler low-power attenuator (class C decoupler amplifier).
- Decoupler modulator frequency.
- Audio filter bandwidth.
- Small-angle phase shifts (with direct synthesis rf only).
- Shim DAC values.

The UNITY spectrometer added a fair number of additional devices; including the following:

- RF power attenuators (one per channel with linear amplifiers), 63 or 79 dB.
- Fine attenuators (optional, mostly used with systems having a solid-state module); available for the first two channels only; the same port or AP words can be used to drive an additional ("third") 63-dB attenuator to facilitate pulse shaping without waveform generators.
- Linear amplifier status (cw and pulse modes).
- Waveform generator shapes and pattern (see also Chapter 16, "Waveform Generators," on page 163).
- Waveform generator shape selection and time scaling.
- Amplitude of a pulsed field gradient (PFG amplifier).
- Imaging gradient amplitudes.

The AP interface board (of which several versions exist) interfaces rf coarse and fine attenuators and AMT linear amplifiers with the AP bus. The last version of this board (type 3) allow setting the decoupler modulation mode (dmm parameter) by the AP bus.

Finally, in the UNITY plus the following devices were added to the list:

- Lock power, lock gain, and lock phase.
- Receiver gain.
- Linear amplitude modulator.

The AP interface board is not present in UNITYplus spectrometers.

12.3 AP Bus Words in the Acode

AP bus traffic is predominant in the initialization section of the instruction section in the Acode, as can be seen from the following part of the (interpreted) Acode (note that with the exception of rf attenuator commands, the contents of the AP bus words are *not* decoded with the apdecode command used here):

| 148 | 142 | 36 | б | APBOUT | 7 | items 0xab40 0xbb0d 0x9b00 0xab54 0xbbec 0x9b2f 0x9b8f |
|-----|-----|-----|-----|-----------|---|---|
| 157 | 151 | 45 | б | APBOUT | 9 | items 0xa720 0xb700 0x97be 0xb7ef |
| | | | | | | 0xb7ae 0xb7fb 0xb7be 0x9701 |
| | | | | | | 0xb700 |
| 168 | 162 | 56 | 159 | TUNE_FREQ |) | CH1 9 words 0xa720 0xb700 0x97be 0xb7ef |
| | | | | | | 0xb7ae 0xb7fb 0xb7be 0x9701 0xb700 |
| 180 | 174 | 68 | 0 | NO_OP | | |
| 181 | 175 | 69 | 16 | SETPHAS90 |) | CH1c 0 |
| 184 | 178 | 72 | 68 | PHASESTEP | þ | CH1 360 units (90.00 degrees) |
| 187 | 181 | 75 | 65 | SETPHASE | | CH1f 0 |
| 190 | 184 | 78 | 59 | APChipOUI | | APaddr 11, reg 51, +logic, 1 byte |
| | | | | | | max 79, offset 16, value 55 |
| 196 | 190 | 84 | 59 | APChipOUT | | APaddr 11, reg 150, -logic, 2 bytes |
| | | | | | | max 4095, offset 0, value 4095 |
| 202 | 196 | 90 | б | APBOUT | 4 | items 0xab98 0xbbd2 0x9b00 0x9b00 |
| 208 | 202 | 96 | б | APBOUT | 2 | items 0xab92 0xbba0 |
| 212 | 206 | 100 | б | APBOUT | 2 | items 0xab9b 0xbb00 |
| 216 | 210 | 104 | б | APBOUT | 2 | items 0xab91 0xbb80 |
| 220 | 214 | 108 | 6 | APBOUT | 2 | items 0xab90 0xbb11 |
| 224 | 218 | 112 | 6 | APBOUT | 8 | items 0x8201 0x8214 0x8220 0x8231 |
| | | | | | | 0x8241 0x8255 0x8264 0x8271 |
| 234 | | 122 | | SETPHAS90 | | CH2c 0 |
| 237 | 231 | 125 | 59 | APChipOUI | | APaddr 11, reg 50, +logic, 1 byte |
| | | | | | | max 79, offset 16, value 30 |
| 243 | 237 | 131 | 59 | APChipOUI | | |
| | | | | | | max 4095, offset 0, value 4095 |
| 249 | | 137 | | APBOUT | | items 0xaba8 0xbb2a 0x9b00 0x9b00 |
| 255 | | 143 | | APBOUT | | items 0xaba2 0xbbc0 |
| 259 | 253 | 147 | | APBOUT | | items 0xabab 0xbb00 |
| 263 | 257 | 151 | | APBOUT | | items 0xabal 0xbb80 |
| 267 | 261 | 155 | | APBOUT | | items 0xaba0 0xbb12 |
| 271 | 265 | 159 | | APBOUT | | items 0xab48 0xbb01 |
| 275 | 269 | | 6 | APBOUT | | items 0xab34 0xbb55 |
| 279 | 273 | 167 | | APBOUT | | items 0xab35 0xbb00 |
| 283 | 277 | 171 | 6 | APBOUT | | items 0xab4d 0xbb21 |
| 287 | 281 | 175 | | APBOUT | | items 0xab43 0xbb01 |
| 291 | 285 | 179 | | APBOUT | | items 0xab49 0xbb00 |
| 295 | 289 | 183 | | APBOUT | | items 0xab36 0xbb06 |
| 299 | 293 | 187 | | APBOUT | | items 0xab9b 0xbb88 |
| 303 | 297 | 191 | 6 | APBOUT | 2 | items 0xabab 0xbb88 |

Mostly for AP bus traffic, the code 6 (APBOUT) is used. This Acode instruction is grouped with an Acode word describing the number of AP words that follow (actually, the number of words minus one) and the AP bus part (shown here as hexadecimal value) of each FIFO word that will be generated. The code 59 (APCOUT) is used specifically for rf attenuators and power modulators. The waveform generator instructions (using the AP bus) have a special format.

12.4 Timing Considerations

The timing on the AP bus (or of AP bus FIFO words) is inherent and implicit. No delay needs to be specified for the AP bus traffic to occur. Code like

```
delay(tau);
offset(getval("offset1"),TODEV);
pulse(pw,oph);
```

is acceptable and complete. In general, inherent (hidden) AP delays are negligible compared to most delays (*n*D evolution, J-evolution, refocusing), which are mostly are an order of milliseconds or longer. However, there is a danger of losing coherence due to chemical shift evolution (precession) in (hypothetical) constructs like the following:

```
pulse(pw,v1);
delay(tau);
pulse(2.0*pw,v2);
delay(tau);
offset(getval("offset1"),TODEV);
rlpower(tpwr,TODEV);
pulse(pw,v3);
```

The AP bus delays after the second refocusing delay are over 26 microseconds on a UNITY with latching in the frequency synthesis (8.05 microseconds on a UNITY *plus*). This makes the second refocusing interval longer by an amount that can cause considerable precession due to chemical shift evolution with large spectral windows. With any refocusing, make sure the two refocusing delays are *absolutely equal*. Therefore, the hidden AP delays in the above construct should be added to the first delay or (more correctly) be subtracted from the delay adjacent to the AP events. At the same time, symmetry with respect to the pre- and post-pulse delays should be ensured:

```
rgpulse(pw,v1,rof1,0.0);
delay(tau-rof1);
rgpulse(2.0*pw,v2,rof1,0.0);
delay(tau-rof1-15.05e-6-2.15e-6);
offset(getval("offset1"),TODEV);
rlpower(tpwr,TODEV);
rgpulse(pw,v3,rof1,0.0);
```

In this example, the programmer has looked up the length of the hidden AP delays in the tables in the manual *VNMR User Programming* (note that at the level of the pulse sequence, all time events are defined in seconds, irrespective of the VNMR parameter definition involved). Because these delays are specific to certain types of hardware (a UNITY without latching on the observe PTS synthesizer, in this case), this construct (and with it the pulse sequence) becomes specific to an instrument and cannot be ported to other systems without rechecking (and maybe adjusting) the delay corrections.

Chapter 12. AP Bus Traffic

Therefore, even if the example of code above is correct in every aspect, it should be regarded as bad programming practice! VNMR software provides the tools for making this *correct code for every architecture:*

```
rgpulse(pw,v1,rof1,0.0);
delay(tau-rof1);
rgpulse(2.0*pw,v2,rof1,0.0);
delay(tau-rof1-OFFSET_DELAY-POWER_DELAY);
offset(getval("offset1"), TODEV);
rlpower(tpwr,TODEV);
rgpulse(pw,v3,rof1,0.0);
```

The file /vnmr/psg/apdelay.h defines a series of such (pseudo) constants (macros that select the right constant for the current architecture) to be used instead of fixed delays. Table 7 gives a list of the predefined AP delay constants. The AP delays

| Constant Name | Comments |
|---------------------|---|
| POWER_DELAY | coarse and fine attenuators with linear amplifiers (power, rlpower) |
| SAPS_DELAY | small-angle phase shifting (xmtrphase, dcplrphase, dcplr2phase, dcplr3phase) |
| OFFSET_DELAY | offset, obsoffset, decoffset, dec2offset, dec3offset (PTS synthesizer without latching) |
| OFFSET_LTCH_DELAY | offset, obsoffset, decoffset, dec2offset, dec3offset (PTS synthesizer with latching) |
| WFG_START_DELAY | starting one waveform generator (start of shaped_pulse, decshaped_pulse, dec2shaped_pulse, dec3shaped_pulse) |
| WFG_STOP_DELAY | <pre>stopping one waveform generator (at end of shaped_pulse, decshaped_pulse, dec2shaped_pulse, dec3shaped_pulse (0.0 on UNITY<i>plus</i>)</pre> |
| WFG2_START_DELAY | starting two waveform generators (e.g., start of simshaped_pulse) |
| WFG2_STOP_DELAY | stopping two waveform generators (e.g., at end of simshaped_pulse) (0.0 on UNITY <i>plus</i>) |
| WFG3_START_DELAY | starting three waveform generators (e.g., start of sim3shaped_pulse) |
| WFG3_STOP_DELAY | stopping three waveform generators (e.g., at end of sim3shaped_pulse) (0.0 on UNITYplus) |
| PRG_START_DELAY | obsprgon, decprgon, dec2prgon, prg_dec_on |
| PRG_STOP_DELAY | obsprgoff, decprgoff, dec2prgoff, prg_dec_off (0.0 on UNITY <i>plus</i>) |
| SPNLCK_START_DELAY | start of spinlock, decspinlock, dec2spinlock, genspinlock |
| SPNLCK_STOP_DELAY | end of spinlock, decspinlock, dec2spinlock, genspinlock (0.0 on UNITY <i>plus</i>) |
| SPN2LCK_START_DELAY | start of gen2spinlock |
| SPN2LCK_STOP_DELAY | end of gen2spinlock (0.0 on UNITYplus) |
| GRADIENT_DELAY | rgradient, vgradient; start and end of zgradpulse |

 Table 7. Predefined AP bus delay constants.

associated with the setting of the fine attenuators using the pwrf and rlpwrf statements (6.45 microseconds on UNITY, 4.6 microseconds on UNITY plus systems) are *not* covered by constants defined in /vnmr/psg/apdelay.h.

Another place it is strongly recommended to compensate for "hidden" AP bus delays is in the evolution time in nD experiments, where constructs like the following one can be used (at the same time we compensate for the precession during two adjacent 90-degree pulses):

Note that UNITY and VXR systems were normally equipped with PTS synthesizers *without* latching, and on UNITY*plus* and UNITY*INOVA* systems, OFFSET_DELAY and OFFSET_LTCH_DELAY have identical values; therefore, using OFFSET_DELAY covers the vast majority of the hardware configurations.

Note also that status changes can cause hidden AP delays to occur (see the manual *VNMR User Programming*).

Chapter 12. AP Bus Traffic

Chapter 13. Acquisition CPU Communication

Apart from the pulse programmer, the acquisition CPU has yet an other channel to communicate with devices in the spectrometer. It can use one of its RS-232 ports, and it can use the host I/O bus that is also used to transfer information to and from the pulse programmer. For the latter, the CPU communicates with the automation control board (see Chapter 7, "Digital Components," on page 65).

Because it is not passing the pulse programmer, such communication with acquisition devices can only be loosely coordinated with the timing of the acquisition. It has to happen when the FIFO (the pulse programmer) is stopped (or hasn't been started yet). This is appropriate for devices that typically have to be set or regulated once per sample (sample changer, magnet leg pneumatics, etc.), where the reaction time is seconds or minutes rather than microseconds (VT controller) or that would not possibly change during the pulse sequence (lock power, gain, phase).

13.1 Regular Pulse Sequence Communication

What exactly are the tasks that are achieved this way? Let's first list the devices for VXR and UNITY spectrometers:

- Communication to and from the VT controller is done from one of the two serial (RS-232) ports of the acquisition CPU (the other port being used for acquisition diagnostics purposes).
- The automation control board is the communication link to the sample changer (ASM-100 or SMS), the magnet leg pneumatics (eject/insert, slow drop period, bearing control), the spinner control circuitry, the lock parameters (power, gain, phase), and finally the receiver gain setting. This board is connected with the acquisition CPU via the same bus (the host I/O bus) as the pulse programmer.

On the UNITY *plus* spectrometer, the tasks that were achieved via these channels have changed slightly:

- Communication with the VT controller has been moved to the new automation control board.
- Lock parameters (power, gain, phase) and receiver gain are now set via the AP bus.
- The tasks left for the automation control board include the magnet leg pneumatics control, and interaction with slow devices that involve two-way communication, like the sample changer (ASM-100 or SMS) and the VT controller¹.

¹ The interfaces to the sample changer and the VT controller are standard serial (RS-232) ports. For setup and diagnostics purposes, these devices can also be operated off-line, with a "dumb" terminal. Both devices should not be disconnected as long as the software is configured for them; otherwise, because of the two-way communication scheme, acquisitions can abort with error messages. Before disconnecting or switching off either of these two devices, always reconfigure VNMR by either using config, or by setting vttype=0 (to disable VT control) or traymax='n' (to disable the sample changer).

13.2 Diagnostics and Error Output

One of the two serial (RS-232) ports on the acquisition CPU can be used for diagnostics purposes by hooking up a terminal. Different types of information can be obtained this way, depending on the acquisition bootup selector switch (a thumbwheel beside the acquisition CPU cardcage on VXR and UNITY spectrometers, or a toggle switch on the UNITY *plus* automation control board):

- If the switch is in position "0", error output is obtained only for serious incidents (CPU hangup and the like). This information may be valuable to software and hardware engineers for debugging purposes, but is seldom used in the field.
- If the switch is in position "1" or "2", continuous diagnostics output is made during an experiment. Position "1" creates output for terminals with cursor addressing (typically Televideo 915 or 925 terminals), position "2" is for "dumb" terminals. With dumb terminals, the output appears line by line; with addressable terminals, the output is organized as full-page display and is, therefore, easier to read and follow.

It is important to know that diagnostics (as opposed to error) output is generated during the housekeeping delay at the end of each scan. If the bootup selector switch is in non-zero position, the diagnostics output is generated, regardless whether a terminal is connected or not. This makes the length of the housekeeping delays *strictly non-deterministic*!

Hooking up a terminal may cause a further slow-down if the communication speed or the terminal are slow (low baud rate or handshaking, due to a slow terminal). Typically, with non-zero setting on the bootup selector, only about 5 scans can be performed per second, because of long housekeeping delays. Position "0" should, therefore, be the normal setting².

Even with the acquisition bootup selector switch in position "0", the residual housekeeping delay is still in the order of milliseconds. If an acquisition should be faster than that, this is only possible by doing multi-FID scans (using multiple acquire calls in the pulse sequence in connection with the nf and cf parameters). This way, we can acquire several FIDs with a single housekeeping delay. This technique is used extensively in multiecho imaging experiments.

 $^{^{2}}$ The position of the bootup switch is checked with every go; therefore, it is *not required* to reboot the acquisition CPU when changing the bootup selector setting.

Chapter 14. Repeating Events

Repetitious events in a pulse sequence may be defined better by C loops, real-time loops, or hardware loops. Each type of loop is discussed in this chapter.

14.1 C Loops

Any C construct can basically be used in a pulse sequence; therefore it is *formally* correct to generate successive, repetitive real-time events in a pulse sequence using a C loop:

```
int i;
for (i = 0; i < 1000; i++)
{
    rgpulse(pw,zero,rof1,0.0);
    rgpulse(pw,one,rof1,0.0);
}</pre>
```

This is not only correct C code, but it also creates correct Acodes that should execute properly—maybe! What is the problem?

- This construct is extremely inefficient in that it creates huge amounts of Acode (1000 copies of the same sequence of four time events with HSLINES instructions in-between), such that there can easily be an Acode overflow (10000 Acode words per FID maximum).
- Depending on the length of the time events involved (pw and rof1 in the above example), the Acode interpretation could become the rate-determining step, in which case the sequence would abort with a "FIFO underflow" message.

This does not mean that C loops are "forbidden" in pulse sequences, but in general they should not be used in sections that generate Acodes. A good example for an exception to this rule is the relayh pulse sequence that adds a variable number of relay periods (depending on the parameter relay) to a standard COSY pulse sequence using a C loop (see Section 10.6, "C Constructs and Phase Calculations," on page 110).

The following code from /vnmr/psg/shape_pulse.c is maybe an interesting exception. Here, a while loop construct is used to set a real-time variable (v12) to the value of a C (integer) variable (npulses), without using initval. This is particularly interesting for external procedures: initval statements reserve a variable completely and also restrict its use. With this construct, the variable can still be used elsewhere in the pulse sequence:

```
assign(zero,v12);
                                 /* v12 = 0 */
assign(2010,.11,
mult(three,three,v13);
                                /* v13 = 9 (RT increment) */
i = 9;
                                /* i = v13 (C increment) */
                                /* npulses > 0? */
while (npulses)
ł
   if (npulses >= i)
                               /* difference > increment? */
   {
      add(v12,v13,v12);
npulses -= i;
                               /* add increment */
                                /* calculate remainder */
   }
                                 /* difference<increment */</pre>
   else
```

```
{
    divn(v13,three,v13); /* divide increment (RT) */
    i /= 3; /* divide increment (C) */
}
}
```

This construct is not optimized for large numbers, but it shows a way to circumvent initval calls even in cases where a real-time variable should be initialized to a value that is parameter-based (i.e., not known prior to the function call).

14.2 Real-Time Loops

A much more adequate way to define repeating events is to use *real-time loops* (sometimes also called "software loops"). Different from C, where there are three different looping mechanisms (the for loop, the while loop, and the do . . . while loop), there is only one mechanisms for generating loops in the Acode interpretation, as shown in the following example of an explicitly encoded decoupling sequence on the transmitter channel (WALTZ-16, the text shown includes the table file and a part of a hypothetical pulse sequence):

```
t1 += 2 2 0 0 2 0 2 0 2 0 2 0 2
     0 0 2 2 0 2 0 2 0 2 0
     0 0 2 2 0 2 0 2 0 2 0
     2 2 0 0 2 0 2 0 2 2 0 2
/* 4 */
dbl(two,v12);
mult(v12,v12,v12); /* 16 = composite pulses per cycle */
xmtron();
loop(v10,v11);
 loop(v12,v13);
   txphase(t1);
   delay(pw);
   txphase(t1);
   delay(2.0*pw);
   txphase(t1);
   delay(3.0*pw);
 endloop(v13);
endloop(v11);
xmtroff();
```

In this example, the number of WALTZ loop cycles (v10) is calculated first. The WALTZ sequence itself is constructed using a nested loop of composite pulses, whereby the pulse phases are taken from an autoincrementing table. Real-time loops can easily be nested over to several levels. The corresponding Acode instruction segment looks as follows:

| 314 | 308 | 202 | 33 | DBLFUN | С | | t | wo | v12 | | | | |
|-----|-----|-----|-----|--------|------|------|-----|------|-------|--------|-----|-------|---|
| 317 | 311 | 205 | 31 | MULtFU | NC | | v | 12 | v12 | v12 | | | |
| 321 | 315 | 209 | 150 | HighSp | eedI | INES | ΤX | ON | | | | | |
| 324 | 318 | 212 | 39 | ASSIGN | FUNC | 2 | ze | ro | v11 | | | | |
| 327 | 321 | 215 | 42 | IFMInu | sFUN | 1C | v | 10 | one | Offset | = | 315 | |
| 331 | 325 | 219 | 39 | ASSIGN | FUNC | 2 | ze | ro | v13 | | | | |
| 334 | 328 | 222 | 42 | IFMInu | sFUN | 1C | v | 12 | one | Offset | = | 309 | |
| 338 | 332 | 226 | 105 | TABLE | 227 | size | 48, | auto | inc 1 | , divn | ret | 1,ptr | 0 |
| | | | | 2 | | 2 | 0 | 0 | 2 | 0 | 2 | 0 | |
| | | | | 2 | | 2 | 0 | 2 | 0 | 0 | 2 | 2 | |
| | | | | 0 | | 2 | 0 | 2 | 0 | 0 | 2 | 0 | |

| | | | | 0 | 0 | 2 | 2 | 0 | 2 | 0 | 2 |
|-----|-----|-----|-----|-----------|-------|-------|-------|------|--------|---|-----|
| | | | | 0 | 0 | 2 | 0 | 2 | 2 | 0 | 0 |
| | | | | 2 | 0 | 2 | 0 | 2 | 2 | 0 | 2 |
| 391 | 385 | 279 | 106 | TASSIGN | table | 227 | tb | lrt | | | |
| 394 | 388 | 282 | 16 | SETPHAS9 | 0 | CH1 | l tb | lrt | | | |
| 397 | 391 | 285 | 151 | EVENT1_T | WRD | 8.0 |)00 u | sec | | | |
| 399 | 393 | 287 | 106 | TASSIGN | table | 227 | tb | lrt | | | |
| 402 | 396 | 290 | 16 | SETPHAS9 | 0 | CH1 | l tb | lrt | | | |
| 405 | 399 | 293 | 151 | EVENT1_T | WRD | 16. | 000 | used | 2 | | |
| 407 | 401 | 295 | 106 | TASSIGN | table | 227 | tb | lrt | | | |
| 410 | 404 | 298 | 16 | SETPHAS9 | 0 | CH1 | l tb | lrt | | | |
| 413 | 407 | 301 | 151 | EVENT1_T | WRD | 24. | 000 | used | 2 | | |
| 415 | 409 | 303 | 27 | INCRFUNC | | v1 | 3 | | | | |
| 417 | 411 | 305 | 42 | IFMInusF | UNC | v1 | 3 | v12 | Offset | = | 226 |
| 421 | 415 | 309 | 27 | INCRFUNC | | v1 | 1 | | | | |
| 423 | 417 | 311 | 42 | IFMInusF | UNC | v1 | 1 | v10 | Offset | = | 219 |
| 427 | 421 | 315 | 150 | HighSpeed | dLINE | s (vo | oid) | | | | |
| | | | | | | | | | | | |

The two arguments specified to the loop statement specify the loop *count* (v10 and v12 in this example) and the loop *counter* (the variable that is used to count the loop cycles). The loop count variable itself is not altered by the loop construct. The loop counter variable is also supplied to the endloop statement—mainly to define which loop statement the endloop call refers to (it also makes it easier to read the code and helps avoiding mistakes). Strictly speaking, the argument to endloop isn't really necessary, because loops always must be *either sequential or hierarchically nested*: it doesn't make sense to start a loop inside an other loop, but to terminate it outside.

The Acode interpreter works on a more primitive level than the pulse sequence language itself. First, the loop counter (v11 for the outer loop, v13 for the inner loop) is initialized to zero. Then, the system checks whether "loop counts - one" is negative (i.e., the loop count is zero), in which case the entire loop segment is skipped (jump to offset 315 for the outer loop, to offset 309 for the inner loop). At the end of each loop cycle (supposing a non-zero loop count was specified), the loop counter is incremented, and if "loop counter - loop count" is negative (i.e., "loop counter < loop count"), a *branch* to the first instruction inside the loop is performed (jump to offset 226 for the inner loop, to offset 219 for the outer loop).

There are some principal limitations with the above implementation of a delay with WALTZ decoupling. The length of the delay has a distinct "granularity", in that it is a multiple of the duration of a WALTZ cycle¹. Also, the WALTZ modulation can only run synchronously, and the only way to program events (e.g.: pulses) simultaneous to this delay is to break up the delay into several segments, with the additional events inbetween. Of course, this would even increase problems with the "duration granularity".

In general, the above implementation of an explicit modulation scheme is quite elegant and efficient (in terms of Acode space usage); however, it may suffer from some further limitation—through the software loop we feed the pulse programmer with a large number of FIFO words. If the events inside the loop are now all very short (a few microseconds only), the Acode interpretation may become the rate-determining step, and the system may run into a FIFO underflow problem (the sequence would abort at that point). The real-time loop is a very efficient, flexible and elegant tool, but it may not be suited for very fast loops as they are needed in many cases like CRAMPS and

¹ One could in principle always round down the WALTZ loop cycles and perform an additional delay to fill the remainder such that the total delay is accurate, but then that remainder would not be properly modulated (unless more complex coding is used).

similar experiments. Note also that the loop count is a 16-bit real-time variable, and the number of loop cycles is therefore limited to 32767.

As shown in the above example, a real-time loop cycle may include phase changes (and phase calculations). It could even include real-time *if* statements (see Section 15.2, "Real-Time Decisions," on page 160). The loop cycles can therefore be variable not only in the phase, but even in the number and kind of events inside the loop.

14.3 Hardware Loops

The hardware loop (i.e., the possibility to cycle words in the loop FIFO of the pulse programmer) was created to compensate for the basic deficiency of real-time loops: the lack of ultimate speed (i.e., to cover the difference between the speed of Acode interpretation and the possible speed of propagation of FIFO words in the pulse programmer). Typical examples for experiments that require the hardware looping capability of the pulse programmer are multipulse and CRAMPS type of experiments such as MREV-8 or BR-24 (a part of the MREV-8 sequence is shown here):

```
pulsesequence()
{
  double tau = getval("tau"),
         dtau = tau - pw - rof1 - rof2;
   initval(np/2.0, v9);
   . . .
  delay(d1);
  rgpulse(pw,v4,rof1,rof2); /* prep pulse */
  starthardloop(v9);
     delay(dtau);
     rgpulse(pw,v4,rof1,rof2);
                                      x */
     delay(dtau);
     rgpulse(pw,v3,rof1,rof2);
                                   /* -y */
     delay(tau+dtau);
     rgpulse(pw,v1,rof1,rof2);
                                   /* y */
     delay(dtau);
     rgpulse(pw,v2,rof1,rof2);
                                   /* -x */
     delay(tau+dtau-2.0e-7);
     acquire(2.0,2.0e-7);
                                   /* acquire */
     rgpulse(pw,v2,rof1,rof2);
                                   /* -x */
     delay(dtau);
     rgpulse(pw,v3,rof1,rof2);
                                   /* -y */
     delay(tau+dtau);
     rqpulse(pw,v1,rof1,rof2);
                                   /* y */
     delay(dtau);
     rgpulse(pw,v4,rof1,rof2);
                                   /* x */
     delay(tau);
   endhardloop();
```

}

In the case of the MREV-8 sequence, a series of 8 pulses and 9 delays (typically 25 single-precision time events in total, with rof2 set to zero) is inserted *into each sampling interval* (the dwell time being 10 to 100 microseconds). The BR-24 multipulse experiment even asks for 24 pulses and 25 delays (73 time events in total) to be squeezed into the same time interval. All these time events typically are on the order of 1 to a few microseconds. There would be no chance for this experiment to work with a soft loop, but for the pulse programmer with its hardware looping capability this is no

problem, as long as the number of events in the loop doesn't exceed the number of words in the loop FIFO. Older systems with output boards (63-word FIFO) cannot perform BR-24 experiments, but all newer systems (using acquisition control boards and pulse sequence control boards) should have no limitation in the area of multipulse experiments.

Hardware loops have some inherent limitations:

- The length of a hardware loop is limited to the size of the loop FIFO: 63 FIFO words for output boards, 1024 FIFO words for acquisition control boards, and 2048 FIFO words for pulse sequence control boards (real-time or software loops can be as big as an entire code segment). When calculating the number of FIFO words in a hardware loop, be aware of double-precision time events. The manual *VNMR User Programming* contains detailed lists on the number of FIFO words involved in the statements that use the AP bus.
- All loop cycles are identical (in real-time or software loops, real-time decisions can be made inside the loop, and phase can be altered either by real-time calculations, by using auto-incrementing tables, or by recalculating table indices).
- Hardware loops can not be nested, because only one loop is implemented in real hardware (the loop FIFO of the pulse programmer).
- The limitation in the number of loop cycles is the same as for real-time loops, because the (hardware) loop counter is also a 16-bit number (i.e., *the maximum number of loop cycles is 32767*).

Multiple hardware loops can be used sequentially. If hardware loops are placed back-to-back (no time event in-between), there is a restriction in the duration of all but the last (back-to-back) hardware loops, in that the total duration of all FIFO words in a loop must be *at least* 0.4 microseconds for each FIFO word in the loop that follows.

With output boards, there are additional restrictions and limitations in that there must be at least a single time event (e.g., a delay of 0.2 microseconds) between any two hardware loops. A hardware loop cannot be shorter than 6.3 microseconds per cycle (the "fall-through time" of the loop FIFO), and with multiple hardware loops following each other, the length of all but the first loop cycles must be at least 80 to 100 microseconds. Also, try to limit the number of events between sequential hardware loops (but keep at least one event in-between); otherwise, it is possible that the loop FIFO runs empty when starting the second loop. Due to the limitation in pre-loop FIFO size, there may not be enough space to pre-load the second loop. Explicit hardware looping is not available on Gemini spectrometers².

Note that also the implicit acquisition and an explicit acquisition over the full FID³ is performed using hardware looping (see Chapter 18, "Acquiring Data," on page 205).

There are also some *programming restrictions* with respect to hardware loops: the use of *real-time math* and of the use of autoincrementing *tables* is not permitted inside hardware loops. This is mainly because it would make programmers believe that real-time math and the incrementation of table pointers also continues during the execution of hardware loops. In addition to that, *real-time decisions* (see Section 15.2, "Real-

² The Gemini pulse programmer is not equipped with a pre-loop FIFO; therefore, no events can be pre-loaded after a hardware loop. The consequence is that only one hardware loop per pulse sequence can be performed at the very end of the sequence. This almost by definition is the acquisition loop. ³ More exactly, the acquire statement for more than two data points (per call) outside a hardware loop.

Time Decisions," on page 160) and *real-time loops* are not permitted inside hardware loops. This mostly has to do with the way hardware loops are programmed in Acode⁴.

Therefore, C constructs (loops and calls to statements) are the only way to simplify the coding of hardware loops at C level. This also implies that hardware loops may not be very efficient in terms of Acode space. Let's see how the WALTZ-16 decoupling sequence would be coded in a pulse sequence—the same or similar types of constructs have been used in explicitly programmed MLEV-16, MLEV-17, and similar modulation schemes.

```
#include <standard.h>
waltza()
{
  txphase(two); delay(3.0*pw); txphase(zero); delay(4.0*pw);
  txphase(two); delay(2.0*pw); txphase(zero); delay(3.0*pw);
  txphase(two); delay(1.0*pw); txphase(zero); delay(2.0*pw);
  txphase(two); delay(4.0*pw); txphase(zero); delay(2.0*pw);
  txphase(two); delay(3.0*pw);
}
waltzb()
{
  txphase(zero); delay(3.0*pw); txphase(two); delay(4.0*pw);
  txphase(zero); delay(2.0*pw); txphase(two); delay(3.0*pw);
  txphase(zero); delay(1.0*pw); txphase(two); delay(2.0*pw);
  txphase(zero); delay(4.0*pw); txphase(two); delay(2.0*pw);
  txphase(zero); delay(3.0*pw);
}
pulsesequence()
ł
  status(A);
   delay(d1);
  status(B);
   rqpulse(pw,zero,rof1,0.0);
   initval(getval("tau")/(16.0*6.0*pw),v10);
   xmtron();
   starthardloop(v10);
     waltza(); waltzb(); waltzb(); waltza();
   endhardloop();
   xmtroff();
  status(C);
   pulse(pw,oph);
}
```

We could of course have followed the coding scheme of the real-time loop solution and split up the sequence in simple composite (90-180-270) pulses. As mentioned before, the hardware loop solution is not very efficient in Acode space usage, as you can see

⁴ The number of FIFO words inside a hardware loop (which is an argument to the HDLOOP Acode instruction) is calculated by the time when the Acode is built (i.e., when typing go), and not during execution time.

from the printout below⁵ (what really counts, of course, is the experimental flexibility, not the ultimate level of Acode economy):

| | | | | • | | | |
|------|-----|-----|-----------|-------------|---------------|----------|-----|
| 317 | 311 | 205 | 0 | NO_OP | | | |
| 321 | 315 | 206 | 67 | HWLOOP with | acq interrupt | 36 words | v10 |
| 322 | 316 | 210 | 16 151 | SETPHAS90 | CH1 two | | |
| 325 | 319 | 213 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 327 | 321 | 215 | 16 | SETPHAS90 | CH1 zero | | |
| 330 | 324 | 218 | 151 | EVENT1_TWRD | 32.000 usec | | |
| 332 | 326 | 220 | 16 | SETPHAS90 | CH1 two | | |
| 335 | 329 | 223 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 337 | 331 | 225 | 16 | SETPHAS90 | CH1 zero | | |
| 340 | 334 | 228 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 342 | 336 | 230 | 16 | SETPHAS90 | CH1 two | | |
| 345 | 339 | 233 | 151 | EVENT1_TWRD | 8.000 usec | | |
| 347 | 341 | 235 | 16 | SETPHAS90 | CH1 zero | | |
| 350 | 344 | 238 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 352 | 346 | 240 | 16 | SETPHAS90 | CH1 two | | |
| 355 | 349 | 243 | 151 | EVENT1_TWRD | 32.000 usec | | |
| 357 | 351 | 245 | 16 | SETPHAS90 | CH1 zero | | |
| 360 | 354 | 248 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 362 | 356 | 250 | 16 | SETPHAS90 | CH1 two | | |
| 365 | 359 | 253 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 367 | 361 | 255 | 16 | SETPHAS90 | CH1 zero | | |
| 370 | 364 | 258 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 372 | 366 | 260 | 16 | SETPHAS90 | CH1 two | | |
| 375 | 369 | 263 | 151 | EVENT1_TWRD | 32.000 usec | | |
| 377 | 371 | 265 | 16 | SETPHAS90 | CH1 zero | | |
| 380 | 374 | 268 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 382 | 376 | 270 | 16 | SETPHAS90 | CH1 two | | |
| 385 | 379 | 273 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 387 | 381 | 275 | 16 | SETPHAS90 | CH1 zero | | |
| 390 | 384 | 278 | 151 | EVENT1_TWRD | 8.000 usec | | |
| 392 | 386 | 280 | 16 | SETPHAS90 | CH1 two | | |
| 395 | 389 | 283 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 397 | 391 | 285 | 16 | SETPHAS90 | CH1 zero | | |
| 400 | 394 | 288 | 151 | EVENT1_TWRD | 32.000 usec | | |
| 402 | 396 | 290 | 16 | SETPHAS90 | CH1 two | | |
| 405 | 399 | 293 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 407 | 401 | 295 | 16 | SETPHAS90 | CH1 zero | | |
| 410 | 404 | 298 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 412 | 406 | 300 | 16 | SETPHAS90 | CH1 zero | | |
| 415 | 409 | 303 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 417 | 411 | 305 | 16 | SETPHAS90 | CH1 two | | |
| 420 | 414 | 308 | 151 | EVENT1_TWRD | 32.000 usec | | |
| 422 | 416 | 310 | 16 | SETPHAS90 | CH1 zero | | |
| 425 | 419 | 313 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 427 | 421 | 315 | 16 | SETPHAS90 | CH1 two | | |
| 430 | 424 | 318 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 432 | 426 | 320 | 16 | SETPHAS90 | CH1 zero | | |
| 435 | 429 | 323 | 151 | EVENT1_TWRD | 8.000 usec | | |
| 437 | 431 | 325 | 16 | SETPHAS90 | CH1 two | | |
| 440 | 434 | 328 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 442 | 436 | 330 | 16 | SETPHAS90 | CH1 zero | | |
| 445 | 439 | 333 | 151 | EVENT1_TWRD | 32.000 usec | | |
| 447 | 441 | 335 | 16 | SETPHAS90 | CH1 two | | |
| 450 | 444 | 338 | 151 | EVENT1_TWRD | 16.000 usec | | |
| 452 | 446 | 340 | 16 | SETPHAS90 | CH1 zero | | |
| 455 | 449 | 343 | 151 | EVENT1_TWRD | 24.000 usec | | |
| 457 | 451 | 345 | 16 | SETPHAS90 | CH1 two | | |
| | | | | | | | |

⁵ The WALTZ-16 solution using real-time loops shown previously (see Section 14.2, "Real-Time Loops," on page 148) required 113 Acode words (226 bytes). The coding with a hardware loop shown here takes up 184 Acode words or 368 bytes. This is *not* necessarily a typical example.

Chapter 14. Repeating Events

| 460 | 454 | 348 | 151 | EVENT1_TWRD | 24.000 usec |
|-----|-----|-----|-----|----------------|-------------|
| 462 | 456 | 350 | 16 | SETPHAS90 | CH1 zero |
| 465 | 459 | 353 | 151 | EVENT1_TWRD | 32.000 usec |
| 467 | 461 | 355 | 16 | SETPHAS90 | CH1 two |
| 470 | 464 | 358 | 151 | EVENT1_TWRD | 16.000 usec |
| 472 | 466 | 360 | 16 | SETPHAS90 | CH1 zero |
| 475 | 469 | 363 | 151 | EVENT1_TWRD | 24.000 usec |
| 477 | 471 | 365 | 16 | SETPHAS90 | CH1 two |
| 480 | 474 | 368 | 151 | EVENT1_TWRD | 8.000 usec |
| 482 | 476 | 370 | 16 | SETPHAS90 | CH1 zero |
| 485 | 479 | 373 | 151 | EVENT1_TWRD | 16.000 usec |
| 487 | 481 | 375 | 16 | SETPHAS90 | CH1 two |
| 490 | 484 | 378 | 151 | EVENT1_TWRD | 32.000 usec |
| 492 | 486 | 380 | 16 | SETPHAS90 | CH1 zero |
| 495 | 489 | 383 | 151 | EVENT1_TWRD | 16.000 usec |
| 497 | 491 | 385 | 16 | SETPHAS90 | CH1 two |
| 500 | 494 | 388 | 151 | EVENT1_TWRD | 24.000 usec |
| 502 | 496 | 390 | 150 | HighSpeedLINES | (void) |
| | | | | | |

The HDLOOP instruction that starts the hardware loop in the Acode is followed by the number of FIFO words inside the loop (*not* the number of Acode words!), and the loop count (a real-time variable). The endhardloop statement does not generate an Acode instruction by itself, but it causes the number of loop FIFO words to be written into the HDLOOP instruction.

Chapter 15. Decisions

The use of C-based (or previously Pascal-based) decisions in pulse sequences has a long tradition in Varian pulse sequences.

15.1 Decisions and Branchings in C

C-based decisions are based on both qualitative as well as quantitative checks. They are used to increase the reliability of many pulse sequences, to make them easier to use, and to increase their flexibility.

Aborting a Sequence in Case of a Improperly Set Parameter

Aborting a sequence can involve duty cycle calculations and checking, checks for parameter settings that could possibly damage the spectrometer hardware (like rf coils or amplifiers in the case of excessive rf power), and avoiding "impossible" parameter settings (like trying to decouple a 2D spectrum with antiphase magnetization). The following example is from the HMQC pulse sequence:

```
if ((dm[A] == 'y') || (dm[B] == 'y'))
{
    printf("DM must be set to either 'nny' or 'nnn'.\n");
    abort(1);
}
```

Standard output created by the printf function shows up in the VNMR text window. The abort function terminates the go command and the Acode generation through the compiled pulse sequence, and Acode generated up to that point is discarded.

Ensuring Compatibility with Spectrometer Hardware

Many pulse sequences contain statements that can only be used on certain types of spectrometer hardware. Therefore, many sequences include constructs such as:

```
if (newtrans)
{
   stepsize(base, TODEV);
   xmtrphase(v1);
}
else
   phaseshift(base,v1,TODEV);
```

Although programmers prefer to write "universal" pulse sequences, it would often make a pulse sequence very complicated if we tried to cover the entire range of spectrometers that can be driven by VNMR software. Thus, some sequences use constructs such as the following:

As shown in Table 8, there are a number of predefined flag variables that can be used in such constructs. Their definitions can be found in the header file /vnmr/psg/acqparms.h.

| Flag Variable | Meaning / Function |
|---------------|---|
| newtrans | True if system equipped with direct synthesis (as opposed to broadband or fixed frequency) rf on observe channel. |
| newdec | True if system equipped with direct synthesis (as opposed to broadband or fixed frequency) rf on (first) decoupler channel. |
| newtransamp | True if system equipped with linear (as opposed to class C) amplifiers on observe channel. |
| newdecamp | True if system equipped with linear (as opposed to class C) amplifiers on (first) decoupler channel. |
| vttype | 0 = none, 1 = Varian, 2 = Oxford |
| Hlfreq | Proton frequency of instrument: 200, 300, 400, 500, 600, 750 |
| automated | True if system has computer-controlled lock, gain, and decoupler power. This flag is left over from old software; all spectrometers running VNMR are automatic. |
| fifolpsize | Size of loop FIFO: 63 (output board), 1024 (acquisition control board), or 2048 (pulse sequence control board). |
| NUMch | Number of rf channels that are configured: 2, 3, or 4. |

| Table 8. | PSG hardware | flag and | configuration | variables |
|----------|--------------|----------|---------------|-----------|
|----------|--------------|----------|---------------|-----------|

Issuing Warning Messages in the Case of Questionable Parameter Settings

Sometimes it may not be necessary to abort the pulse sequence, but it is perhaps adequate to have the sequence display a warning message:

```
if ((newtransamp) && (rof1 < 9.9e-6) && (ix == 1))
printf("Warning: ROF1 is less than 10 usec.\n");</pre>
```

In the case of arrayed and multidimensional experiments, it is strongly recommended to only print a warning for the first increment; otherwise, the VNMR text window may be flooded with error messages. Of course, there may be cases where a check should be performed for each increment. The above example should just remind the operator that rof1 is below its normal value for linear amplifiers (which could result in pulse amplitude instabilities).

Turning On or Off Pulse Sequence Features

A very popular and powerful feature in VNMR pulse sequences is the use of *flags* for enabling and disabling individual pulse sequence statements. This allows combining many pulse sequences (in the sense of a defined sequence of pulses and delays) into a single pulse sequence program. There are numerous examples of pulse sequences with additional experimental flexibility through C decisions, like the HMQC sequence with its optional nulling of protonated signals through a BIRD inversion pulse (simplified coding), such as the following:

if ((null > 0.0) && (mbond[A] == 'n'))
{

```
rgpulse(pw,v1,rof1, 0.0);
delay(bird-rof1);
decrgpulse(pwx,v9,rof1,0.0);
simpulse(2.0*pw,2.0*pwx,v1,v1,1.0e-6,0.0);
decrgpulse(pwx,v9,1.0e-6,0.0);
delay(bird-rof1);
rgpulse(pw,v2,rof1,0.0);
hsdelay(null);
```

The second conditional in this example ensures than BIRD nulling is not used when the long-range option (mbond='y') is used.

Other examples include sequences like hetcor (with a built-in long-range option and the possibility to suppress proton multiplets using the hmult='n', which causes an inversion pulse to be replaced by a BIRD pulse), and many others. Sometimes, numeric parameters are used to make decisions (like checking whether the delay null is nonzero in the above example), but mostly such switches are based on flag parameters (such as mbond and hmult in the above examples).

Dynamically Arranging the Sequence of Events

Sometimes it is necessary to rearrange the sequence of events depending on the value of a parameter. One example is found in the HMQC pulse sequence, where there is a 180-degree proton pulse in the middle of the evolution time. The principal sequence of events is as follows:

```
decpulse(pwx,v3);
delay(d2/2.0-rof1);
pulse(2.0*pw,v4);
delay(d2/2.0-rof2);
decpulse(pwx,v5);
```

This construct *looks* correct, but it has several deficiencies. The duration of the central pulse is not compensated for in the evolution time, and the precession during the 90-degree pulses surrounding the evolution time is not taken into account (see Chapter 19, "Multidimensional Experiments," on page 215). Additionally, the first increment the delay d2 is zero, giving an error message about negative delays.

Apart from that (which is just a programming problem), the first increment has a considerable gap between the two 90-degree X-pulses, which causes problems in the first trace and baseline distortions in the final 2D spectrum. In order to do it properly, the pulses on both rf channels should be performed on top of each other for the first increment. Unfortunately, the two X-pulses don't have the same phase; therefore, a single simpulse statement cannot be used. One can minimize interpulse dead times by keeping the amplifiers unblanked for parts of a pulse sequence:

```
corr = 2.0*pwx/3.1416 + pw + 1.0e-6;
if (d2/2.0 > corr)
{
    rcvroff();
    decrgpulse(pwx,v3,rof1,0.0);
    delay(d2/2.0-corr);
    rgpulse(2.0*pw,v4,1.0e-6,0.0);
    delay(d2/2.0- corr);
    decrgpulse(pwx,v5,1.0e-6,0.0);
    rcvron();
}
else
```

```
{
    simpulse(pw,pwx,v4,v3,rof1,0.0);
    simpulse(pw,pwx,v4,v5,1.0e-6,0.0);
}
```

This is a simplified solution: Strictly speaking, this is still not perfect. If pwx is less than pw, there is still a gap between the X-pulses, but with most configurations used for this type of experiment, the X-pulses are longer than the proton pulses, and in this case the above solution is accurate for the first increment. Also, increments other than the first one can fall into a domain where the X pulses overlap with the proton 180, which may require some modifications to the above algorithm.

An interesting example is the combination of two pulse sequences into a single one for the same experiment. The example shown here is a constant-time heteronuclear correlation experiment, where a simultaneous inversion pulse moves within a fixed time interval¹. Such experiments often suffer from limited resolution because the moving pulse at some point (increment) reaches the end of the fixed delay. In the example shown here, the fixed interval contains a BIRD pulse, which would normally set a limit to the number of increments. Thanks to use of C decisions (and to the fact that an Acode segment is calculated individually for each 2D increment), we can apply the trick of letting the simultaneous inversion "jump over the BIRD pulse" and continue behind it up to the very end of the fixed delay, which doubles the achievable resolution in f₁ (only relevant parts of the pulse sequence are shown here):

```
decpulse(pp,t1);
                                     /* excitation pulse */
if (d2/2.0 < d3)
{
   if (d2/2.0 > rof1)
                                     /* start of evolution time */
     delay(d2/2.0-rof1);
  rgpulse(pw,t1,rof1,0.0);
                                     /* 180 H + 180 X */
  simpulse(2.0*pw,2.0*pp,t2,t3,1.0e-6,0.0);
  rgpulse(pw,t1,1.0e-6,0.0);
   if (d3 - d2/2.0 > rof1)
     delay(d3-d2/2.0-rof1);
  decrgpulse(pp,t1,rof1,0.0); /* BIRD-pulse */
  delay(tau-rof1);
  rgpulse(pw,t1,rof1,0.0);
  simpulse(2.0*pw,2.0*pp,t2,t4,1.0e-6,0.0);
  rgpulse(pw,t1,1.0e-6,0.0);
  delay(tau-rof1);
  decrgpulse(pp,t1,rof1,0.0);
  delay(d3 -tau -rof1);
}
else
{
  delay(d3-rof1);
                                   /* BIRD-Pulse */
  decrgpulse(pp,t1,rof1,0.0);
  delay(tau-rof1);
  rgpulse(pw,t1,rof1,0.0);
  simpulse(2.0*pw,2.0*pp,t2,t4,1.0e-6,0.0);
  rgpulse(pw,t1,1.0e-6,0.0);
  delay(tau-rof1);
  decrgpulse(pp,t1,rof1,0.0);
```

¹ M. Perpick-Dumont, W.F. Reynolds & R. Enriquez, *Magn. Reson. Chem.* **26**, 358 (1988); W.F. Reynolds, S. McLean, M. Perpick-Dumont & R. Enriquez, *Magn. Reson. Chem.* **26**, 1068 (1988).

Implicit Decisions

. . .

}

Of course, decisions are not limited to the pulse sequence itself. Decisions are involved internally in most pulse sequence statements, whether to do the right thing for the given hardware or to avoid unnecessary Acode and FIFO words (like time events with zero duration).

Decisions Set by the status Statement

Another class of decisions is hidden in the status statement (see also "Implicit Gating" on page 53). Most pulse sequences are divided into several basic sections by inserting status calls at suitable points in the pulse sequence:

```
status(A);
hsdelay(d1);
status(B);
pulse(p1,zero);
hsdelay(d2);
pulse(pw,oph);
status(C);
```

This is a simplified version of the s2pul pulse sequence. It is split into three basic sections: status A, which includes the relaxation period (d1); status B, which includes the two pulses and the evolution time (d2); and status C, which is valid during the acquisition time. The argument to the status function is a simple integer: A is defined as 0, B is defined as 1, etc., up to Z, which is defined as 25. status (A) simply means that during the following time events, the first ("0th" in C syntax) character in predefined multifield flag variables (dm, dm2, dm3, dmm, dmm2, dmm3, homo, homo2, homo3, and hs) is active up to the next status call with a different value in the argument. These flag variables can be up to 255 characters long (beyond the predefined constants A to Z), but in practice only 3 to 5 fields are used.

These multifield flags are a very easy way of controlling the gating of transmitters and modulation modes, and enabling or disabling homospoil pulses in hsdelay calls, which significantly enhances the flexibility of most pulse sequences.

The definition is that if more fields are addressed in the pulse sequence than there are characters in a flag string, the last character is propagated as much as necessary (i.e., dm = 'ny' is the equivalent to $dm = 'nyyyyyyyyyy \dots$ '). At the end of a scan (and at the end of the experiment), the system *implicitly returns to* status(A). This is most relevant to experiments that use high-power decoupling, or decoupling during the acquisition in general. For example, when using s2pul with decoupling, it is better not to use dm = 'y', but rather dm = 'nyy', even if delay d1 is not used. This switches off the decoupler at the end of the experiment and avoids sample heating due to the rf power fed into the probe continuously. If a relaxation delay is to be used with s2pul

(and p1=0), we can use d2 instead of d1, and dm='nyy' for full decoupling or dm='nyn' for gated decoupling.

Note that with some configurations (UNITY *plus* and UNITY spectrometers), the status statement may not just change some fast bits (status lines), but may involve AP bus FIFO words (depending on which flags change at that point) that take a finite time (see Chapter 12, "AP Bus Traffic," on page 137).

Checking Flag Parameters

Simple flag tests are easy. Constructs like if (dm[A] == 'y') are adequate in most cases and work also with flags that are newly created for a specific pulse sequence (like if (mbond[A] == 'n') in a previous example in this chapter). Problems occur when flag fields *other than the first* should be tested: if (dm[C] == 'y') sounds like a correct test—but what if the parameter dm has been set to 'y'? In this case, dm[C] returns a random value (dm[B] would return a null character, but still not 'y'!), and the test would fail, because the test returns FALSE instead of TRUE. A proper construct for testing flag fields, other than the first one, is not that simple:

```
int flagsize,
    index = statusindx;
char xflag[MAXSTR];
getstr("xflag",xflag);
...
flagsize = strlen(xflag);
if (index >= flagsize)
    index = flagsize - 1;
if (xflag[index] == 'y')
    ...
```

The variable statusindx is the currently active status, as set by the status statement. Another solution would be to ensure (by means of a macro or parameter entry restrictions) that this particular flag has the required length. Fortunately, the vast majority of user-created flags use only one field.

15.2 Real-Time Decisions

Decisions that need to be made from transient to transient cannot be programmed in C, because only one Acode segment is generated per FID, which performs all scans. In such a case, we need real-time decisions—decisions made by the Acode interpreter.

Programming Real-Time Decisions

Suppose that in a particular sequence, with every odd scan we wanted to add a refocusing period to a pulse sequence. This could be achieved as follows:

There is only one logical test for real-time decisions: the test whether a real-time variable (v10 in this case) is zero. To fulfil this condition for every odd scan, we

| calculate ct mod | a 2. The above construct g | generates the following | ng Acode (pw=7, |
|------------------|----------------------------|-------------------------|-----------------|
| tau=1/140, ro: | f1=10): | | |

| 298 | 292 | 186 | 36 | MOD2FUNC | ct v10 |
|-----|-----|-----|-----|----------------|-----------------------|
| 301 | 295 | 189 | 41 | IFNotZeroFUNC | v10 zero Offset = 218 |
| 305 | 299 | 193 | 152 | EVENT2_TWRD | 7 msec + 133 usec |
| 308 | 302 | 196 | 16 | SETPHAS90 | CH1 v1 |
| 311 | 305 | 199 | 150 | HighSpeedLINES | RXOFF |
| 314 | 308 | 202 | 151 | EVENT1_TWRD | 10.000 usec |
| 316 | 310 | 204 | 150 | HighSpeedLINES | RXOFF TXON |
| 319 | 313 | 207 | 151 | EVENT1_TWRD | 7.000 usec |
| 321 | 315 | 209 | 150 | HighSpeedLINES | RXOFF |
| 324 | 318 | 212 | 150 | HighSpeedLINES | (void) |
| 327 | 321 | 215 | 152 | EVENT2_TWRD | 7 msec + 143 usec |
| 330 | 324 | 218 | 150 | HighSpeedLINES | (void) |
| | | | | | |

The ifzero statement generates an IFNotZeroFUNC instruction in the Acode, which performs a jump to address 218 (i.e., after the endif statement) if the variable (v10 in this case) is non-zero. The endif statement itself does not generate Acode, but it is used to determine the jump address in the IFNotZeroFUNC instruction.

The argument with endif is not used in the Acode, but it serves to associate the endif with the corresponding ifzero call: ifzero constructs can also be nested (see also the last example in Section 10.4, "Real-Time Logical Decisions," on page 106), and in such cases the argument is an easy way for the software to "know" in which IFNonZeroFUNC instruction it has to set the jump address².

Mostly in the case of such conditional events, we would like to ensure that all scans are performed with exactly the same overall timing, to avoid adding up FIDs with different phase and amplitude. In this case, the elsenz function is also used between the ifzero and the endif calls:

In this example, we are performing a 45-degree pulse with every odd scan, and with the even scan numbers this pulse is replaced by a delay. Of course, the delay length must not only include the pulse length, but also any pre- and post-pulse delay included with the rgpulse statement (beware of implicit delays in statements like pulse!)³.

In the Acode, the above construct looks as follows:

| 352 | 346 | 240 | 36 | MOD2FUNC | ct | v10 | | | |
|-----|-----|-----|-----|----------------|-------|------|--------|---|-----|
| 355 | 349 | 243 | 41 | IFNotZeroFUNC | v10 | zero | Offset | = | 268 |
| 359 | 353 | 247 | 16 | SETPHAS90 | CH1 | v1 | | | |
| 362 | 356 | 250 | 150 | HighSpeedLINES | RXOFF | | | | |

² Strictly speaking, this is not true because it is always clear which ifzero and endif belongs to, as long as the ifzero ... endif constructs are hierarchically stacked, and non-hierarchical stacking (e.g., ifzero(v1)...ifzero(v2)... endif(v1)... endif(v2)) doesn't make sense at all (the same is true for real-time looping).

³ To be accurate, we cannot simply replace a pulse with a delay of the same length, since we also need to take into account the precession during the pulse. This can be done using a correction term that shortens the delay. If pw is a 90-degree pulse (pw/2.0 therefore 45 degrees), the accurate delay length is pw/2.0 - pw/3.1416 + rof1 (see also Chapter 19, "Multidimensional Experiments," on page 215).

| 365 | 359 | 253 | 151 | EVENT1_TWRD | 10.000 usec |
|-----|-----|-----|-----|----------------|-------------|
| 367 | 361 | 255 | 150 | HighSpeedLINES | RXOFF TXON |
| 370 | 364 | 258 | 151 | EVENT1_TWRD | 3.500 usec |
| 372 | 366 | 260 | 150 | HighSpeedLINES | RXOFF |
| 375 | 369 | 263 | 150 | HighSpeedLINES | (void) |
| 378 | 372 | 266 | 20 | BRANCH | Offset 270 |
| 380 | 374 | 268 | 151 | EVENT1_TWRD | 13.500 usec |
| 382 | 376 | 270 | 150 | HighSpeedLINES | (void) |

The elsenz statement again takes the same argument as the ifzero and endif calls (the real-time variable v10 in this case). It terminates the if part of the construct with a BRANCH (jump) instruction and sets the jump address in the IFNotZeroFUNC instruction to the first instruction after the if part (after the BRANCH instruction). The endif statement in this case sets the jump address in the BRANCH instruction at the end of the if part.

Real-time decisions can also be used in real-time math for the calculation of complex phase cycles (see Chapter 10, "Phase Calculations," on page 95). A table index should not be incremented with each scan, but perhaps with every fourth scan only.

Generating the Flag Variable

In most cases, the flag variable for real-time decisions is generated using some kind of (real-time) modulo statement. To generate a flag variable with a period of *n* scans, one has to take modulo *n* of a variable that is incremented with every scan like ct, or modulo n/2 if the variable is incremented every second scan (ct/2), etc. Here are some examples for the construction of ifzero flag variables (x in the comment stands for a non-zero value):

| <pre>mod2(ct,v1); sub(one,v2);</pre> | /* /* | | | | | | | | | | | | | | | | | |
|---|----------|---|---|---|---|---|---|----|---|----|---|---|---|---|---|---|---|-------------------|
| mod4(ct,v3); | /* | 0 | x | x | x | 0 | x | x | x | *, | / | | | | | | | |
| add(ct,one,v4); mod4(v4,v4); | /* /* | | | | | | | | | | | | | | | | | |
| <pre>hlv(ct,v5); add(v5,two,v5); mod4(v5,v5);</pre> | /* | 2 | 2 | 3 | 3 | 4 | 4 | 5 | 5 | 6 | 6 | 7 | 7 | 8 | 8 | 9 | 9 | * / * / * / |
| <pre>modn(ct,three,v6);</pre> | /* | 0 | x | x | 0 | x | x | *, | / | | | | | | | | | |

To invert the logic in a simple if zero construct (without elsenz branch), simply use the elsenz branch (to simulate a non-existing ifnotzero statement):

```
add(ct,two,v7); /* 2 3 4 5 6 7 8 9 */

hlv(v7,v7); /* 1 1 2 2 3 3 4 4 5 5 6 6 7 7 8 8 */

mod4(v7,v7); /* x x x x x x 0 0 x x x x x x 0 0 */

ifzero(v7);

elsenz(v7);

...

endif(v7);
```

Chapter 16. Waveform Generators

The optional waveform generator (or programmable pulse modulator as it was called in UNITY spectrometers) is primarily a board that permits changing the amplitude and the (small-angle) phase of the output of a transmitter board simultaneously and quickly. Different from the pulse programmer, it can set both the amplitude and the small-angle phase shifts using fast status lines. You can regard it as "autonomous fast bit extension of the pulse programmer."

16.1 How Does a Waveform Generator Fit Into the System?

Figure 19 is a diagram of waveform generator circuitry in the UNITY*plus*. The waveform generator sends output through the transmitter digital control board to the (small-angle and 90-degree) phase and amplitude modulation circuitry on the (direct synthesis) transmitter board. Both the amplitude and the phase modulation are also accessible from the AP bus, but—as we will see—the waveform generator can do the same thing, but much faster. The waveform generator is tightly coupled to a particular transmitter board and can only interact with the channel on which it is installed (any channel can be equipped with a waveform generator). Reconfiguring the waveform generator for an other channel involves jumper changes on the board (the AP bus addressing is channel-specific). It is not recommended as a routine operation.

The phase control circuitry on the transmitter board was also described in Section 4.2, "How Do Pulses Work?," on page 40. It does the small-angle phase shifts (0.25 degrees resolution) and the 90-degree phase shifts in two separate, consecutive steps. It takes 42 MHz as input frequency; the output is the phase-modulated 10.5 MHz intermediate frequency (I.F.). The phase control takes digital inputs from the pulse programmer (90-degree phase shifts via fast lines, small-angle phase shifts via AP bus), from the waveform generator (small-angle and 90-degree phase shifts via fast direct lines), and from the phase modulator (just the 90-degree phase shifts via fast status lines). Each channel is equipped with a phase modulator (for preprogrammed decoupling sequences like square wave, swept square wave, noise, MLEV-16, WALTZ-16, XY-32, GARP-1), and optionally with a waveform generator (freely programmable decoupling sequences, pulse shaping). The 90-degree and small-angle phase shifts are *added* to the phase setting of the pulse programmer set via AP bus.

The amplitude modulation circuitry adjusts the amplitude of the 10.5 MHz I.F. signal in 4096 *linear* steps (12 control bits), as opposed to the (63 or 79 dB) attenuators that operate in dB (i.e., in logarithmic units). The 4096 linear steps correspond to a power range of 72 dB. The amplitude modulator takes input from either the pulse programmer (via AP bus) or from the waveform generator (via fast status lines). From the 12 control bits, the waveform generator only controls the 10 most significant bits (1024 steps, 60-dB power range). The amplitude modulation circuitry takes the phase-modulated 10.5 MHz I.F. from the phase modulation circuitry. Its output is the phase- and amplitude-modulated I.F., which is then mixed with the local oscillator frequency (which is the same as the observe frequency *plus* the I.F.) that is generated directly by a PTS frequency synthesizer and is also used in the receiver (see Chapter 18, "Acquiring Data," on page 205).

Chapter 16. Waveform Generators

The output of the mixer (the mixing product) is the observe frequency, which carries the phase- and amplitude modulation from the 10.5 MHz (I.F.) input. The observe frequency that passes the transmitter gate (which can also be controlled by the waveform generator) and the 79 dB (or 63 dB) attenuator is the amplified in the rf power amplifier and finally enters the probe. The waveform amplitude takes precedence over the pulse programmer amplitude that is set via AP bus.

Spectrometers earlier than UNITY*plus* (UNITY and VXR) used a variety of different transmitter boards:

• Only the AM/PM (amplitude modulation and phase modulation) direct synthesis transmitter board was compatible with a waveform generator, and the amplitude modulation was only accessible via the waveform generator, not via the AP bus. The AM-PM board had 0.5-degree phase resolution and 1024-step (60-dB) linear amplitude range.




Figure 19. UNITY plus waveform generator circuitry

- Previous direct synthesis boards only had phase modulation capability (0.5-degree resolution, like the AM-PM board), accessible via the AP bus (and fast lines for the 90-degree phase shifts). Due to the lack of amplitude modulation capability, these boards were incompatible with waveform generators.
- Fixed frequency and broadband-type transmitter boards did not have amplitude or phase modulation capability (except for 90-degree phase shifts) and were not compatible with a waveform generator either.

Because the AM/PM transmitter board used in UNITY spectrometers did not allow setting the modulator amplitude by AP bus and because the AM/PM board combines both the digital and the analog functions on a single board, the corresponding connection scheme (shown in Figure 20) is slightly simpler.




Figure 20. UNITY programmable pulse modulator circuitry

16.2 How Does a Waveform Generator Work?

The waveform generator is similar to a pulse programmer in that it directly controls spectrometer hardware (via fast lines only) and provides accurate timing, independent of the Acode computing. However, there are some distinct differences: the pulse programmer is built as a *flow-through* (FIFO) *buffer*, because the information it handles is *dynamic* and in essence non-repetitive (subsequent scans are normally not identical, because of phase cycling etc.). A waveform generator, on the other hand, is a *static* buffer (based on static RAM), built to handle a constant set of *generic* shapes and pattern within a given pulse sequence. The term generic means that the duration (the duration scaling) of a pulse shape may vary an arbitrary number of times within a pulse sequence, but the phase and amplitude pattern remain constant. To better understand the way a waveform generator works, we need to have a more detailed look at the way it is built and how it operates.

As shown in Figure 21, the central part of the waveform generator board consists of a 256-Kbyte memory (static RAM, i.e., random access memory that does not require

Chapter 16. Waveform Generators

constant refreshing). That memory is organized in 65536 words (64 Kwords) of 32 bits each. A *variable size* portion is used for instructions; the rest is used to store pattern and shapes (see below). The other components of the waveform generator include an interface to the AP bus, an output buffer, control and timing circuitry (20 MHz, derived from the 40 MHz input frequency), loop control circuitry, and an amplitude multiplier that is only used (and present) for shaping pulsed field gradients.




Figure 21. Waveform generator board

The waveform generator is very fast in pulse shaping. All parameters under its control (amplitude, 90-degree and small-angle phase shifts, and the transmitter gate) are changed *instantaneously and simultaneously* (unlike changing these parameters via AP bus, see Chapter 12, "AP Bus Traffic," on page 137). The *minimum time slice* for waveforms and modulation pattern is 0.2 microseconds, the *timing resolution* is 50 nanoseconds on the UNITY*plus*, and 100 nanoseconds on UNITY spectrometers.

This is the instruction timing. Apart from internal propagation delay (discussed in the section "Shaped Pulses in the Acode" on page 176), there are always some minor delays until the parameter change takes effect in real hardware, especially with respect to phase changes in UNITY spectrometer (see also Section 4.2, "How Do Pulses Work?," on page 40). For shaped pulses and modulation patterns, these minor hardware delays are irrelevant because all time slices experience the same delays.

Overall, the waveform generator is fast enough to shape even hard pulses as short as a few microseconds. On the other hand, its pattern memory is large enough to hold even the most complex pulse shapes and modulation pattern known today.

Sequence of Events in a Waveform Generator

While the internal data and information paths on the waveform generator board seem complicated, the entire functionality can be explained fairly well based on the organization and the usage of the RAM. Let's first have a look at the sequence of events during the execution of a pulse sequence that uses the waveform generator:

- The pulse sequence executable reads the shape or pattern name from the statement (shape_pulse, obsprgon, etc.) in the sequence or from the modulation sequence parameter (dseq, etc.), reads and interprets the shape or pattern file (more about this below), and generates a waveform generator data file named expn.username.PID.RF¹ in /vnmr/acqqueue. This file contains both the waveforms and the instruction blocks that will be uploaded to the waveform generators. Each data segment (pattern and instruction blocks) is preceded by a short header (8 bytes) that contains the AP bus address of the waveform generator (to which the data will be sent) plus the starting RAM address and the number of (32-bit) words in the pattern or block.
- At the same time, a file /vnmr/acqqueue/ldcontrol is generated or updated. This file contains a memory map for all waveform generators of the system. It shows where in the waveform generator memory (and in which waveform generator) the current waveform(s) and the corresponding instruction blocks are located (i.e., the starting address). This starting address is also used in the pulse sequence when generating Acode instructions that call waveform generator instruction blocks.
- Before uploading the Acode for a new experiment that uses waveform generators, Acqproc uploads the waveform generator *data* (pulse shapes and decoupling pattern) and *instruction blocks* in the file expn.username.PID.RF from /vnmr/acqqueue into the HAL memory. The acquisition also uploads the data via the pulse programmer and the AP bus into the waveform generator pattern and instruction memory.
- The pulse sequence is then initiated by Acqproc, the same as any other experiment. Whenever a waveform generator is to be started within a pulse sequence (e.g., at the beginning of a shaped pulse), the execution of a specific instruction block is initiated via AP bus. Among other information, that instruction block contains the pattern start address, the pattern length and the duration of a time slice for that particular shape call, see below. To perform a particular shaped pulse, it is sufficient to send the (2-byte) address of an instruction block to a waveform generator, plus a one-byte command that starts its execution (5 AP bus words in total).

¹ The name of the file includes the experiment name, the user name, the process-ID of the pulse sequence executable, and the suffix .RF (see also Chapter 5, "Submit to Acquisition: go," on page 55). The name alone for Acqproc allows the association of such a file with a given Acode file for a specific experiment of a particular user.

• The final execution of a pulse shape or decoupling pattern is (normally) triggered by a fast status line (there is one fast line per waveform generator, see Section 9.2, "Fast Bits," on page 88).

How Are Patterns Stored in a Waveform Generator?

To better understand the functioning of the waveform generator, let's now have a look at the organization of its RAM:

Waveforms and modulation patterns are stored in the waveform generator as a *generic pattern* that describe the principal amplitude and phase modulation and the *relative duration* of each slice of a shaped or modulated pulse sequence element. The actual length of a slice duration unit is defined in each instruction block that calls that particular shape. This way, each shape only needs to be stored once. It can be called an almost infinite number of times with *individual duration scaling* for every single call.

The shapes or pattern are stored in 32-bit words organized as follows (see also the schematic drawing in Figure 21):

- 8 bits for the duration count—any element of a shape can be one up to 255 time slices long. The duration of a time slice is not defined in the pattern, but will be given through the instruction block.
- 10 bits for the amplitude—this allows setting the amplitude in 1024 steps (1 to 1024²) on the linear modulator. On the UNITY*plus* these are the 10 most significant bits of the 12-bit linear modulator; via AP bus the full amplitude range (4096 steps) can be addressed: the smallest step on the waveform generator corresponds to four amplitude steps on the linear modulator. On UNITY spectrometers, the linear modulator was 10 bits only, and these could be fully addressed via the waveform generator only.
- 11 bits (10 bits on the UNITY) for the phase of the pattern element—out of this, 2 bits are the 90-degree phase shift, the remaining 9 bits (8 bits on UNITY) are the small-angle phase shift (360 quarter-degree steps, or 180 half-degree steps on the UNITY). On the transmitter board, this phase shift is *added* to the current phase shift (as set by the pulse programmer through two fast lines and the AP bus) such that the 90-degree and small-angle phase shifting within a shaped pulse occurs *on top of* any preexisting phase shift.
- 1 bit for the transmitter gate—this can also be regarded as an additional amplitude bit that allows setting the amplitude to zero. Note that the regular amplitude bits have a gap between the positive amplitudes and "negative amplitudes" (amplitudes with phase inversion). Taking into account 180- degree phase shifts (and there are many pulse shapes with phase inversion), the "regular" amplitude settings are 1024 to –1 and 1 to 1024; the zero is accessible only through the transmitter gate bit (the linear modulator does not have a zero amplitude).
- 2 bits are spares. These bits are not used or addressed by any regular software and hardware up to now. These could possibly be used to trigger additional gates or other devices *during* a waveform or decoupling pattern (this would currently require changing psg/wg.c via psggen mechanism, because there is no function or feature that would allow setting these bits).

 $^{^2}$ In the shape or modulation files in shapelib, these amplitude values are represented by numbers ranging from 0 to 1023 (1 less than actually obtained); this has sometimes led to confusion.

Patterns or shapes can be between 1 and almost 64 Kwords long. Note that the 64 Kword pattern memory is shared among all shapes and instruction blocks for the active experiment. A variable size part of the waveform generator memory is used for instruction blocks.

The waveform generator shape definition is incomplete insofar as it only contains a relative duration in each pattern word. In essence, this allows replacing successive shape slices with identical phase and amplitude by a single slice with a larger duration count. The vast majority of the pulse shape definitions have a (default) duration count of 1 in each slice. In modulation patterns, pulse angles are translated to duration counts; therefore, decoupler modulation pattern frequently use slices with larger duration counts, and the slice duration can be defined in units between 1 and 90 degrees.

Waveform Generator Instruction Words

Patterns are stored at the high end (0xFFFF in 32-bit words is the last memory location) of the RAM address range. Instruction blocks are stored at the bottom end, from address 0x0000 up. A typical instruction block for an rf pattern (pulse shape or modulation pattern) consists of six 32-bit words. Based on the 3 most significant bits, we can distinguish 8 different types of instruction words. The structure and the contents depend on the word type: bits 29 to 31 are the word type, bit 28 is reserved in all types for (default) transmitter gating, the rest of the block is variable. Major parts of the instruction word are often not used and contain no further information.

Using Table 9, let's look at the various kinds of instruction words (unused parts of the instruction words were left white).

| Туре | 31 | 30 | 29 | 28 | 27 | 26 | 25 | 24 | 23 | 22 | 21 | 20 | 19 | 18 | 17 | 16 | 15 | 14 | 13 | 12 | 11 | 10 | 9 | 8 | 7 | 6 | 5 | 4 | 3 | 2 | 1 | 0 |
|------------|----|----|----|----|----|----|----|----|----|----|------------------------------|----|----|-----|-----|-----|------|-----|-------|------|-----|------|------|---|---|---|---|---|---|---|---|---|
| IB_START | 0 | 0 | 0 | G | L | | | | | | | | | RA | Μ | sta | rt a | ddr | ess | | | | | | | | | | | | | |
| IB_STOP | 0 | 0 | 1 | G | | | | | | | RAM stop address delay count | | | | | | | | | | | | | | | | | | | | | |
| IB_SCALE | 0 | 1 | 0 | G | | | | | | | loop count amplitude scale | | | | | | | | | | | | | | | | | | | | | |
| IB_DELAYTB | 0 | 1 | 1 | G | | | | | | | | | | del | ay | tim | e c | our | nt, 5 | 50 1 | nse | c ui | nits | | | | | | | | | |
| IB_WAITHS | 1 | 0 | 0 | G | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| IB_PATTB | 1 | 0 | 1 | G | | | | | | | | | 1 | att | ern | tin | ne c | cou | nt, | 50 | nse | ec u | init | s | | | | | | | | |
| IB_LOOPEND | 1 | 1 | 0 | G | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| IB_SEQEND | 1 | 1 | 1 | G | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Table 9. Waveform generator instruction words

Each instruction block starts with a IB_START instruction word that specifies the 16bit starting address of the associated pattern block; the next instruction word, IB_STOP, defines the (16-bit) address of the first word behind the pattern. Note that multiple pattern are stored next to each other without start or end marks.

The IB_STOP instruction can also contain a duration count (0 to 255) for an optional delay preceding the pattern in bits 0 through 7. This duration count goes together with a IB_DELAYTB instruction that defines the duration time base in a 28-bit number (0 to 268,435,455). The time base is in 50-nanosecond units (a 20 MHz clock frequency is used for the waveform generator timer); therefore, the time base can be up to 13.42 seconds. Together with a duration count, this permits specifying a delay of up to 3422

Chapter 16. Waveform Generators

seconds, or almost an hour. The minimum duration time base is 200 nanoseconds. Bit 28 determines whether the transmitter is on during that delay or not. The IB_DELAYTB instruction (together with the duration count in the IB_STOP instruction) could be used to specify additional delays preceding and following a shaped pulse or other pattern. It is currently not used for regular shaped pulses or programmed decoupling (see below).

The IB_SCALE instruction has two purposes. It allows specifying an amplitude multiplier for shaped gradients (see Section 16.6, "Using a Waveform Generator for Shaping Gradient Pulses," on page 195), and it also defines an 8-bit loop counter for cycling with a defined loop count over shapes or pattern. For pulse shapes, the loop counter is set to 1. Programmed decoupling normally uses infinite looping (with the loop count set to 1).

The IB_WAITHS instruction makes the waveform generator to wait for a trigger signal on its dedicated pulse programmer high-speed line before starting to execute a pattern.

The IB_PATTB instruction defines the pattern time base the same way as the IB_DELAYTB for the delay time base. The pattern time base can be between 0.2 microseconds and 13.42 seconds—together with a duration count of 1 to 255 in the pattern definition, *each slice* of a shaped pulse or pattern can be between 0.2 microseconds and 3422 seconds long.

The IB_LOOPEND instruction is used in connection with pattern looping, and the IB_SEQEND instruction word terminates each instruction block.

Waveform Generator Data File

Upon typing go, all pattern and instruction blocks for the current experiment are collected in a single file /vnmr/acqqueue/expn.username.PID.RF (PID being the process-ID of the pulse sequence executable), or /vnmr/acqqueue/acqi.RF for data generated through go('acqi'). This file consists of short headers consisting of 4 unsigned integers, each followed by either an instruction block or a data block. The headers consist of the AP bus address of the target waveform generator, the starting address in the RAM of the waveform generator, the number of 32-bit words (-1) in the instruction block or pattern, and a "spare" number that has no real function in the software, except that it allows the reader to distinguish between instruction blocks (0xabcd) and data blocks (0xfedc). The file ends with a 32-bit trailer containing the pattern 0xa5b6c7d8. From the headers, the acquisition CPU "knows" which waveform generator (i.e., which AP bus address) a block needs to be sent to, and at what (RAM) address in the target device it needs to store the data (or instruction) block.

Unlike Acode (or NMR data structures), the command od is of little help in inspecting waveform generator data, because all the (32-bit) words are complex bit patterns, either instruction words or pattern data. On the other hand, pattern data have a fixed structure, and there are only eight different kinds of instruction words, and also the headers are very simple.

Therefore, with the aid of the information found in /vnmr/psg/wg.c, it is relatively easy to write a little program that decodes these data. Let's look at a simple example³:

INSTRUCTION BLOCK: AP address = 0x0c18, WG start address = 0000, 6 words _____ 0x08ffcc00 IB_START: RAM start address = 0xffcc 0x20fff001 IB_STOP: RAM stop address = 0xfff0, delay count = 1 0x40010000 IB_SCALE: loop count = 1 0x80000000 IB_WAITHS: Wait for high-speed line trigger 0xa000063f IB_PATTB: Time count = 1600 (0.00008000 sec / 80.00 usec) 0xe0000000 IB_SEQEND: End of instruction block RF DATA BLOCK: AP address = 0x0c18, WG start address = 0xffcc, 36 words _____ count amplitude phase gates ----- ------ -----31024180.00410240.0021024180.00 1024 180.00 4 2 1024 0.00 3 1024 180.00

END OF FILE

This is a data file for waveform-generator-based WALTZ-16 decoupling (see also below for a more detailed interpretation of the data contents). The file begins with an instruction block, followed by a data block. Of course, there are often several instruction and data blocks—there is at least one instruction block per data block. If the same waveform or pattern is used several times with different parameters in the same pulse sequence, there may be several instruction blocks referring to the same data block. All instruction blocks are collected at the beginning of the file, followed by the data blocks in the order of their usage in the pulse sequence (the order of the data blocks is actually not relevant).

Both blocks are for the decoupler waveform generator (AP bus address $0 \times 0 \times 18^4$). The instruction block will be stored at location 0 in the waveform generator RAM, and the data block at location $0 \times \text{fcc}$. The instruction block is 6 words long, the pattern is 36. The data block address is referred to in the IB_START instruction word, the address in IB_STOP is $0 \times \text{fcc} + 0 \times 24 = 0 \times \text{ff0}$ (65484 + 36 = 65520). The delay count in the same instruction word is the default; there is no delay in this instruction block. The instruction block ends with the instruction word IB_SEQEND.

The other parts of the instruction block and the data block will be explained in Section 16.4, "Using Waveform Generators for Programmed Modulation," on page 177 below.

³ The decoding shown here was produced using a C program wgdecode from the VNMR user library (part of userlib/bin/apdecode).

⁴ The observe channel waveform generator has AP bus address $0 \times 0 \approx 10$, $0 \times 0 \approx 48$ is the address for the second decoupler channel, $0 \times 0 \approx 40$ is for the third decoupler channel.

Executing Waveform Generator Patterns

Through the AP bus and AP interface, the acquisition CPU can not only write to waveform generator RAM, but it can also write *directly into the output register* (output buffer) and set the output status of any waveform generator. For each call, it sends the *address of an instruction block* (using Acode instruction 102, WG3) and *control codes* (using Acode instruction 101, WGCMD) into the waveform generator *control register* (see also the Acode examples in the following sections). Control code 5 is used to start a shaped pulse, control code 7 is used for programmed decoupling, and control code 1 for starting shaped gradients. Termination of the execution of an instruction block is achieved with control code 0, which causes the waveform generator to halt at the end of the current pattern.

For programmed modulation, it is normally not desirable to wait for the end of the pattern: in such a case the instruction block execution is aborted with control code 0x80 (RESET). If the instruction block contains an IB_WAITHS instruction word, the execution of an instruction block is halted until the dedicated high-speed line is set to ON by the pulse programmer. This allows for an accurate coordination of waveform generator events with events that are directly controlled by the pulse programmer.

Instruction blocks (i.e., the associated pattern) can be looped for a predefined number of cycles (the IB_STOP instruction word contains an 8-bit loop count, allowing for up to 255 loop cycles), or they can be looped infinitely, with interruption either at the end of a pattern (loop cycle) or immediately, depending on whether code 0 is sent ("soft stop") or control code 0x80 (RESET), see above. It is also possible to loop over sequences of instruction blocks (in which case all but the last instruction block have bit 27 set in the IB_START instruction).

16.3 Using Waveform Generators for Shaped Pulses

There is one statement per rf channel for performing a simple shaped pulse:

shaped_pulse(name,width,phase,rx1,rx2); decshaped_pulse(name,width,phaser,rx1,rx2); dec2shaped_pulse(name,width,phase,rx1,rx2); dec3shaped_pulse(name,width,phase,rx1,rx2);

All these statements have the same set of five arguments: name is the base name (without ".RF" extension) of a pulse shape file in /vnmr/shapelib or \$vnmruser/shapelib; the other arguments are exactly the same as for an ordinary rgpulse statement (see also Section 4.2, "How Do Pulses Work?," on page 40). It turns out that the four statements are actually macros, all of them calling the same statement genshaped_pulse (the macros are defined in /vnmr/psg/macros.h):

genshaped_pulse(name,width,phase,rx1,rx2,g1,g2,device);

The last argument is the device (OBSch, DECch, DEC2ch, and DEC3ch for the four statements above, see the footnote on page 40). The two additional arguments g1 and g2 are additional delays surrounding the shaped pulse (but with the transmitter turned on). Because they can be defined in waveform generator instruction block, both are set to 0.0 in all four macro calls listed above (no statement currently uses these delays).

There are also statements for two and three simultaneous shaped pulses for the observe and decoupler, or the observe and the two first decoupler channels, respectively:

```
simshaped_pulse(n1,n2,w1,w2,ph1,ph2,rx1,rx2);
sim3shaped_pulse(n1,n2,n3,w1,w2,w3,ph1,ph2,ph3,rx1,rx2);
```

Also these are macro calls to the following functions:

With these two statements, it is possible to perform simultaneous shaped pulses on *any* combination of two or three rf channels. Note that with the presence of the rfchannel parameter, these "gen" functions should not longer be required.

Programming Shaped Pulses: An Example

To illustrate what has been explained in the previous sections, let's now take a practical (simple) example.

The Pulse Sequence

We take the sh2pul pulse sequence, which is basically a simple s2pul sequence with the two rectangular pulses replaced by shaped pulses:

```
#include <standard.h>
```

```
pulsesequence()
{
    char plpat[MAXSTR], pwpat[MAXSTR];
    getstr("plpat",plpat);
    getstr("pwpat",pwpat);
    ....
    status(A);
    hsdelay(d1);
    status(B);
        shaped_pulse(plpat,pl,zero,rof1,rof2);
        hsdelay(d2);
    status(C);
        shaped_pulse(pwpat,pw,oph,rof1,rof2);
}
```

The Shape Definition

We set both the plpat and the pwpat parameter to 'gauss', which causes the file /vnmr/shapelib/gauss.RF to be interpreted when typing go. This file (the ".RF" extension indicating that it contains a pulse shape definition) consists of several columns, each line defining one pulse slice:

- The phase angle of the slice in degrees; phases can be both positive or negative; phases will be converted to positive values between 0 and <360 degrees (better: the corresponding positive range of phase shift units) internally.
- The amplitude, in values between 0.0 and 1023.0 (the file interpreter reads fractional numbers; therefore, shape files can use fractional numbers throughout, but this is not a requirement). Note that these numbers correspond to real amplitudes of 1 up to 1024. If this column is not specified, the amplitude defaults to 1023.0 (1024). The values are rounded off internally, hence any positive fractional number up to 1023.0 is allowed (because it may be delivered from a shape calculation program).
- The relative slice duration, in values between 1.0 and 255.0; the default for this column is 1.0.
- The gate settings: 1 turns on the transmitter gate (TXON), 2 turns on the first spare line, 4 turns on the second spare line. This defaults to 1.0 for pulses (transmitter

Chapter 16. Waveform Generators

on) and is normally omitted. The values for the gates can just be added in column 4, resulting in the possibilities shown in Table 10.

| Value | TXON | Spare #1 | Spare #2 |
|-------|------|----------|----------|
| 0 | OFF | OFF | OFF |
| 1 | ON | OFF | OFF |
| 2 | OFF | ON | OFF |
| 3 | ON | ON | OFF |
| 4 | OFF | OFF | ON |
| 5 | ON | OFF | ON |
| 6 | OFF | ON | ON |
| 7 | ON | ON | ON |

Table 10. Waveform generator gate control for pulse shapes

• The gate field (value 0.0) can be used to generate a slice with zero output, to compensate for the fact that amplitude 0.0 in reality is amplitude 1. For short (hard) shaped pulses starting with a non-zero (internal) phase shift it may be desirable to precede the shape with a slice with the phase of the first real slice, but with the gate turned off.

The file shapelib/gauss.RF defines a Gaussian pulse with 256 slices:

```
#
# Gaussian Pulse: 256 points, 5-sigma
# This pulse is amplitude modulated and selectively excites
# a bandwidth (Hz) approximately equal to 2e+6/pulse_length
# (usec).
0.0000011.3651.000000
0.0000012.1891.000000
0.0000013.0661.000000
. . .
(120 lines deleted)
. . .
0.000001016.0001.000000
0.000001018.5141.000000
0.000001020.4741.000000
0.000001021.8771.000000
0.000001022.7191.000000
0.000001023.0001.000000
0.000001022.7191.000000
0.000001021.8771.000000
0.000001020.4741.000000
0.000001018.5141.000000
0.000001016.0001.000000
. . .
(120 lines deleted)
. . .
0.0000013.9991.000000
0.0000013.0661.000000
```

0.0000012.1891.000000

The Waveform Generator Data File

If both parameters plpat and pwpat in the sh2pul sequence are set to 'gauss', we get the following waveform generator data file for p1=1024 and pw=256:

```
INSTRUCTION BLOCK:
AP address = 0x0c10, WG start address = 0000,
                                                                     6 words
0x08feef00 IB START: RAM start address = 0xfeef
0x20fff001 IB_STOP: RAM stop address = 0xfff0, delay count = 1
0x40010000 IB_SCALE: loop count = 1
0x80000000 IB_WAITHS: Wait for high-speed line trigger
0xb000004f IB_PATTB: Time count = 80 (0.00000400 sec / 4.00 usec) TXON
0xe0000000 IB_SEQEND: End of instruction block
INSTRUCTION BLOCK:
AP address = 0x0c10, WG start address = 0x0006,
                                                                      6 words
_____
0x08feef00 IB_START: RAM start address = 0xfeef
0x20fff001 IB_STOP: RAM stop address = 0xfff0, delay count = 1
0x40010000 IB_SCALE: loop count = 1
0x80000000 IB_WAITHS: Wait for high-speed line trigger
0xb0000013 IB_PATTB: Time count = 20 (0.00000100 sec / 1.00 usec) TXON
0xe0000000 IB_SEQEND: End of instruction block
RF DATA BLOCK:
AP address = 0x0c10, WG start address = 0xfeef, 257 words
_____
count amplitude phase gates
_____ ____

        1
        12
        0.00
        TXON

        1
        13
        0.00
        TXON

        1
        14
        0.00
        TXON

   (120 lines deleted)
   . . . .
           1017
                         0.00 TXON
   1

        1017
        0.00
        TXON

        1019
        0.00
        TXON

        1021
        0.00
        TXON

        1022
        0.00
        TXON

        1023
        0.00
        TXON

   1
   1
   1
   1
         1024
   1
                      0.00 TXON

        1024
        0.00
        TXON

        1023
        0.00
        TXON

        1022
        0.00
        TXON

        1021
        0.00
        TXON

        1019
        0.00
        TXON

        1017
        0.00
        TXON

   1
   1
   1
   1
   1
    . . . .
   (120 lines deleted)
            14 0.00 TXON
13 0.00 TXON
1 0.00
   . . . .
   1
   1
   0
```

END OF FILE

Because we have used the same waveform twice, we get *two instruction blocks* (at locations 0000 and 0006) *referring to the same pattern block* starting at location $0 \times feef$, ending before location $0 \times fff0$. The pulse is triggered by the high-speed line from the pulse programmer, hence the instruction IB_WAITHS. There are 256 slices (actually 257, as we will see). For a pulse duration of 256 microseconds, this results in a slice duration of 1 microsecond. Since each slice has a duration count of 1, the duration time base for the first pulse is 1 microsecond (20 counts of 50 nanoseconds each).

Chapter 16. Waveform Generators

The TXON bit is not only set in the pattern itself, but also in the IB_PATTB instruction word. The second pulse was selected four times longer; therefore, the duration time base is 4 microseconds. Of course, the duration time base is rounded off to the 50 *nanoseconds timing resolution* of the waveform generator clock (100 nanoseconds on UNITY spectrometers). For short pulses and pulses with a large number of slices or total duration counts, the software may report round-off errors if they distort the pulse length substantially.

In the pattern block, we see that the software automatically adds a slice with zero amplitude (amplitude value 1, TXON is not set). During all other slices, the transmitter is gated on. This is the default for the case that only three fields are specified in the shape definition file. It is the waveform generator that turns on the transmitter during a shaped pulse; therefore, there is no need to gate the transmitter explicitly using xmtron and xmtroff.

Shaped Pulses in the Acode

The Acode, and the timing during the experiment, differs in two aspects between UNITY and UNITY *plus* spectrometers (apart from inherent timing differences with respect to the AP bus):

- On UNITY*plus*, shaped pulses are terminated when the high-speed line goes down. On UNITY, the waveform generator needs to be stopped using a "soft stop" by sending a control code 0 (stop at the end of the pattern). Therefore, there is a 2-word AP bus delay (4.3 microseconds) after shaped pulses on a UNITY; whereas on a UNITY*plus*, there is no implicit post-pulse delay due to AP bus traffic.
- There is a propagation delay in the waveform generator relative to the pulse programmer timing. On UNITY*plus*, this propagation delay is 450 nanoseconds. On UNITY, it is 1.5 microseconds (this is an approximation and varies slightly from system to system). This delay needs to be taken care of, or else the pulse programmer might abort the waveform generator or start some other event while the shaped pulse is still executing. To correct for the delay, shaped pulses are explicitly delayed by the expected waveform generator propagation delay.

Here is the Acode for the first of the above two pulses created on a UNITY spectrometer:

| 323 | 317 | 211 | 102 | WG3 | AP addr $0x0c10$, IB addr = $0x0000$ |
|-----|-----|-----|-----|----------------|---------------------------------------|
| 326 | 320 | 214 | 101 | WGCMD | AP addr $0x0c10$, WFG cmd = $0x05$ |
| 329 | 323 | 217 | 16 | SETPHAS90 | CH1 zero |
| 332 | 326 | 220 | 150 | HighSpeedLINES | (void) |
| 335 | 329 | 223 | 150 | HighSpeedLINES | RXOFF |
| 338 | 332 | 226 | 151 | EVENT1_TWRD | 10.000 usec |
| 340 | 334 | 228 | 150 | HighSpeedLINES | RXOFF WFG1 |
| 343 | 337 | 231 | 151 | EVENT1_TWRD | 1.500 usec |
| 345 | 339 | 233 | 152 | EVENT2_TWRD | 255 usec + 1.000 usec |
| 348 | 342 | 236 | 150 | HighSpeedLINES | RXOFF |
| 351 | 345 | 239 | 151 | EVENT1_TWRD | 10.000 usec |
| 353 | 347 | 241 | 150 | HighSpeedLINES | (void) |
| 356 | 350 | 244 | 101 | WGCMD | AP addr $0x0c10$, WFG cmd = 00 |

The WG3 instruction (three AP bus words) sets the instruction block address, the WGCMD instruction (two AP bus words) sets the control register on the waveform generator, see above: control word 0×05 initiates the execution of the instruction block, control word 0 terminates it.

| 409 | 403 | 297 | 102 | WG3 | AP addr 0x0c10, IB addr = 0x0000 |
|-----|-----|-----|-----|----------------|-------------------------------------|
| 412 | 406 | 300 | 101 | WGCMD | AP addr $0x0c10$, WFG cmd = $0x05$ |
| 415 | 409 | 303 | 16 | SETPHAS90 | CH1 zero |
| 418 | 412 | 306 | 150 | HighSpeedLINES | (void) |
| 421 | 415 | 309 | 150 | HighSpeedLINES | RXOFF |
| 424 | 418 | 312 | 150 | HighSpeedLINES | RXOFF |
| 427 | 421 | 315 | 151 | EVENT1_TWRD | 10.000 usec |
| 429 | 423 | 317 | 150 | HighSpeedLINES | RXOFF WFG1 |
| 432 | 426 | 320 | 151 | EVENT1_TWRD | 450 nsec |
| 434 | 428 | 322 | 152 | EVENT2_TWRD | 255 usec + 1.000 usec |
| 437 | 431 | 325 | 150 | HighSpeedLINES | RXOFF |
| 440 | 434 | 328 | 151 | EVENT1_TWRD | 10.000 usec |
| 442 | 436 | 330 | 150 | HighSpeedLINES | (void) |

On a UNITY *plus*, the propagation delay is shorter, and the terminating WGCMD instruction is missing:

From the Acode, it can be concluded that it takes 5 AP bus words to start a shaped pulse (causing hidden delays of 10.75 microseconds on UNITY spectrometers, and 5.75 microseconds on the UNITY plus), and on the UNITY there is a hidden delay of 4.3 microseconds incurring from two more AP bus words *after* the shaped pulse.

Special Cases

To demonstrate the effect of delays that are executed as part of the instruction block, a the shaped_pulse statement was replaced by

genshaped_pulse(plpat,pl,zero,rof1,rof2,0.1,0.2,TODEV);

This is, of course, a hypothetical example. It turns out that the g1 and g2 delays in genshaped_pulse (0.1 and 0.2 seconds in this case, respectively) are in fact built into the instruction block and are not coded in the Acode:

From an Acode point-of-view, the duration of this shaped pulse is p1 + g1 + g2(p1 + 0.1 + 0.2 in this example). There is currently no code or pulse sequence that uses these two additional delays (in all the macro calls to genshaped_pulse and similar statements, the delays g1 and g2 are set to 0.0).

16.4 Using Waveform Generators for Programmed Modulation

Programmed modulation—whether for broadband decoupling or for spinlock experiments—is nothing but a special kind of shaped pulse with the addition of infinite looping (the looping is the main reason for calling it via specific statement, and not with a shaped_pulse call). Apart from this, the requirements and priorities are slightly different for programmed modulation; hence, the pattern definition differs slightly from the definition of shaped pulses.

Programming Pattern Decoupling and Spinlock Experiments

Software for programmed modulation using a waveform generator can be set up in several ways.

Decoupler and Modulation Control

The simplest way to use a waveform generator is for programmed decoupling controlled through status fields in the pulse sequence (see "Implicit Gating" on page 53). Using the waveform generator this way requires no special programming effort at all—it works with any pulse sequence (as long as a waveform generator is available).

The primary control of the transmitter gates is done through the status-related parameters dm, dm2, and dm3. This has nothing to do with any kind of modulation. The modulation mode depends on the parameters dmm (decoupler modulation mode), dmm2, and dmm3 for the three decoupler channels: if these parameters are set to 'c' in the current status field, no modulation occurs. All other settings produce some kind of phase modulation on the corresponding rf channel. This modulation is *independent* of the gate (dm, dm2, dm3) setting! Note that whenever pulses should be performed on a channel, the corresponding modulation mode parameter must be set to 'c' in the current status field; otherwise, the modulation is turned un asynchronously and produces random phases on any pulses in that same status field.

There are two possibilities for modulating a decoupler (through status fields): using a "hardware modulator," or using a waveform generator (if one is available). The hardware modulator can be regarded as "software (a given modulation pattern) that is "programmed" (wired) in hardware and cannot be changed (except by changing or modifying the modulator hardware). The standard modulators only operate on the 90and 180-degree phase shift lines (some even just on 180-degree phase shifting), smallangle phase shifting and amplitude modulation is excluded. On UNITY and earlier systems, the following modulation modes are available through a hardware modulator (for literature references, see the manual *VNMR Command and Parameter Reference*):

'f' "fm-fm" or swept square-wave modulation

| 'g' | GARP-1 decoupling (UNITY plus only) |
|-----|---|
| | (100 pulses with odd pulse angles and 180 degrees phase shifts) |

- 'm' MLEV-16 decoupling (UNITY plus only) (based on 90_x - 180_y - 90_x composite inversion pulses)
- 'n' noise modulation
- 'r' square-wave modulation (UNITY*plus* only)
- 'w' WALTZ-16 modulation (based on 90_x -180_{-x}-270_x composite inversion pulses)
- 'x' XY-32 modulation (UNITY*plus* only) (based on 90_x - 180_y - 90_x composite inversion pulses)

On the UNITY*plus*, there is also a mode 'u' for selecting external (user-supplied) hardware modulation.

Most of these methods (definitely 'g', 'm', 'w', and 'x') need to be calibrated (i.e., for a given decoupler strength, the length of the 90-degree pulse must be determined). The modulator then needs to be triggered at the speed or frequency of the 90-degree pulses. In practice, the modulation frequency (4 times the decoupler field strength, in Hz) is entered directly using the dmf, dmf2, or dmf3 parameter.

Waveform Generator Control through the status Statement

For systems with waveform generator on a decoupler channel, the character 'p' in a dmm (or dmm2, or dmm3) field activates waveform generator-based modulation. Also here, the modulation rate is defined through the dmf (dmf2, dmf3) parameter, but unlike the hardware modulator (where the dmf frequency is generated as input to the modulator), with the waveform generator dmf is converted back into a 90-degree pulse width, which is then used to determine the length of a duration unit in the waveform generator instruction word.

A major advantage of the waveform generator is that it is freely programmable. Any decoupling modulation can be programmed, including methods using small-angle phase shifting, and methods with amplitude modulation (we might even think of using shaped pulses for decoupling).

Setting dmm to 'p' does not define a modulation mode, but rather selects modulation hardware (the waveform generator). In this case, the modulation mode is defined using the string parameter dseq (or dseq2, dseq3 for the other decoupler channels). The value of this parameter is the basename of a file in /vnmr/shapelib or \$vnmruser/shapelib (just the body of the name, without ".DEC" extension). Although this file is stored in the same directory as pulse shapes (shapelib/*.RF), it has a totally different format:

WALTZ-16 Broadband Decoupling Sequence 270.0 180.0 360.0 0.0 180.0 180.0 270.0 0.0 90.0 180.0 180.0 0.0 360.0 180.0 180.0 0.0 270.0 180.0 270.0 0.0 360.0 180.0 180.0 0.0 270.0 180.0 90.0 0.0 180.0 180.0 360.0 0.0 180.0 180.0 270.0 0.0 # 270.0 0.0 180.0 360.0 180.0 0.0 270.0 180.0 90.0 0.0 180.0 180.0 360.0 0.0 180.0 180.0 270.0 0.0 270.0 180.0 360.0 0.0 180.0 180.0 270.0 0.0 180.0 90.0 180.0 0.0 360.0 180.0 180.0 0.0

180.0

270.0

Chapter 16. Waveform Generators

Different from pulse shapes, which normally are described as amplitude and phase function in the time axis, modulation modes are described more often as a sequence of pulses of a certain tip angle and a certain rf phase. The amplitude remains constant throughout the modulation pattern in the vast majority of the case. For this reason, the format of the modulation files (*.DEC) in shapelib is different from the one for shaped rf pulses, even though it also has up to four columns describing pulse length, phase, amplitude and the transmitter gate setting:

- The first column describes the *tip angle in degrees* (positive number only).
- The second column defines the phase of the slice (equivalent to column 1 of a pulse shape file).
- The third field describes the amplitude (equivalent to column 2 of a pulse shape file) and defaults to 1023.0 (full amplitude on the linear modulator). As decoupler modulation is usually constant, this column is normally omitted.
- The last field controls the waveform generator gate settings. If a waveform generator is activated through status fields, there is an implicit gating function (decon and decoff) activated through the status statement, which overrides any transmitter gate settings in the modulation pattern definition, but the gating for the two spare lines (see Table 10) remains active⁵.

For translating the above pattern file into a waveform generator file, the software needs one more parameter that translates the tip angle into duration counts. We need to specify whether the tip angle in the pattern file is to be interpreted as multiples of 90 degrees, or whether the tip, angle is in multiples of a smaller angle. The dres parameter specifies the tip angle resolution. For WALTZ and MLEV-type modulation dres should be set to 90, for GARP-type decoupling it is normally set to 1 in order to get the most accurate tip angles (see also below).

Explicit Programming

Waveform generator-based modulation can also be turned on explicitly within the pulse sequence. There are two classes of statements for this: one to switch on and off the modulation, and the other that also handles the gating and runs for a given number of cycles. The following statements are used to switch on the waveform generator-based modulation:

```
obsprgon(name,pw90,pwres);
decprgon(name,pw90,pwres);
dec2prgon(name,pw90,pwres);
dec3prgon(name,pw90,pwres);
```

where name is the base name (without ".DEC" extension) of a pattern file in \$vnmruser/shapelib or /vnmr/shapelib, pw90 is the length of a 90-degree pulse (allowing the function to calculate the absolute duration of each pattern element), and pwres is the tip-angle resolution, which is needed to determine the length of a duration unit and the number of duration units per pattern element. These four statements are actually macros that call a single C function:

prg_dec_on(name,pw90,pwres,device);

⁵ The two transmitter gate signals are combined (ORed). If either of the two devices (the pulse programmer or the waveform generator) turns the transmitter on, it will be on. In other words, gating the transmitter on explicitly in the pulse sequence deactivates (overrides) control through the pattern file. Of course, this does not apply to the two spare gates on the waveform generator.

where the device argument is the rf device (OBSch, DECch, DEC2ch, and DEC3ch for the four statements above, see the footnote on page 40). For switching the modulation off again, there is another set of four statements:

```
obsprgoff();
decprgoff();
dec2prgoff();
dec3prgoff();
```

These statements are all macros that call a single C function:

prg_dec_off(2,device);

The 2 in the first argument causes a hard abort (reset) in execution of the waveform generator pattern; a 0 in lieu of the 2 would create a soft stop, allowing the waveform generator to stop at the end of the pattern⁶.

Note that these functions do not turn on the transmitter gate, and the default transmitter gating in the waveform generator pattern is *off*. In other words, unless the transmitter is gated on explicitly in the pattern file (through an odd value in column four), the transmitter has to be gated on explicitly in the pulse sequence using statement such as xmtron and xmtroff, decon, and decoff; for example:

```
obsprgon("mlev17",1.0/dmf,dres);
xmtron();
delay(getval("mix"));
obsprgoff();
xmtroff();
```

The sequence of the statements may be relevant. If the xmtron precedes the obsprgon call, the rf is turned on unmodulated during the AP bus traffic to the waveform generator (five AP bus words or 5.75 microseconds on a UNITY plus, 10.75 microseconds on UNITY spectrometers). On the other hand, with the above coding, there is a gap between any preceding time event (e.g., a pulse) and the start of the modulation. Similar considerations apply to the end of the modulation on UNITY spectrometers, where there is a 4.3 microseconds hidden delay (there is no implicit delay behind the modulation on a UNITY plus). Which solution is best depends on the pulse sequence—no general recommendation can be given here.

The advantage of the obsprgon and obsprgoff mechanism is that the modulation can be turned on at any point in a pulse sequence, and it can be left on while other events (pulses, delays—anything that does not use the modulated channel) take place, and it can be turned off at any point later in the sequence. There are cases, however, like the TOCSY sequence, where we would just like to perform a mixing time with spinlock on one channel, and there it would seem simpler if we could just activate the waveform generator with a single command. In fact, such statements exist:

```
spinlock(name,pw90,pwres,phase,nloops);
decspinlock(name,pw90,pwres,phase,nloops);
dec2spinlock(name,pw90,pwres,phase,nloops);
dec3spinlock(name,pw90,pwres,phase,nloops);
```

Once more, these are all macros (as defined in /vnmr/psg/macros.h), calling a single C function:

genspinlock(name,pw90,pwres,phase,nloops,device);

⁶ This is usually undesirable, because the modulation would then continue asynchronously up to the end of the pattern (i.e., the pulse programmer would *not* wait for the end of the modulation before continuing with the next FIFO words).

Chapter 16. Waveform Generators

The first three arguments to these statements are the same as with obsprgon and related functions. In the phase argument, the basic 90-degree phase shift for the spinlocking pulse train can be specified (often the spinlock pulses are subjected to phase cycling), and in the nloops argument to the macros, the number of pattern loop cycles is specified (nloops must be an *integer*, not a floating point number!). This is not the parameter the spectroscopist is normally interested in. In most cases, the relevant parameter is the spinlock time. This requires calculation of the possible number of loop cycles in the pulse sequence, which is no real problem for a given modulation pattern:

ncycles = (int)(getval("mix")/(64.667*pw90)); spinlock("mlev17",pw90,90.0,zero,ncycles);

It would not make sense to define the name of the decoupling pattern as a variable, because the number of 90-degree pulse lengths per pattern depends on the pattern⁷.

The spinlock statements have some other peculiarities. Compared to the obsprgon type of statement, the default gate setting in inverted so that, by default, the transmitter gate is turned on instead of off. This gating does not happen through the pulse programmer (using xmtron, decon, etc.), but it is done in the waveform generator data file (see below) to allow for gate switching *during* the modulation, as for instance required with the "clean" spinlocking sequences that have delays between the elements of each composite pulse. This seems correct and desirable, but it can create confusion.

As long as no gating occurs within the pattern, any spinlock statement can easily be replaced by the following C construct using an obsprgon type of function—this even avoids calculating the number of loop cycles:

```
xmtron();
txphase(zero);
obsprgon("mlev17",pw90,90.0);
delay("getval("mix")");
obsprgoff();
xmtroff();
```

This is more than just equivalent—it would even allow replacing the delay with a sequence of other events: performing pulses and delays *during* the spinlock period, without interrupting the modulation.

Now consider the case where somebody wanted to perform a "clean"-type of spinlock. In this case the xmtron and xmtroff statements cannot be used, because this would override any gating done by the shape (the transmitter would be continuously on). The new construct will be:

```
txphase(zero);
obsprgon("clean_mlev17", pw90, 90.0);
delay("getval("mix"));
obsprgoff();
```

Now we need a pattern definition that includes the gating column. To specify a gating value (column 4) we also need to specify the (default) amplitude. A programmer could now save some work by specifying only the gating where the transmitter must be on (because the transmitter gate is off by default):

⁷ Having the pattern name as a variable would require parsing the pattern to find out about the pattern length, which of course would exceed the scope of normal pulse sequence programming (and probably the capabilities of many pulse sequence programmers).

| # Cl | ean | MLEV-17 | | |
|------|-----|---------|--------|-----|
| 90. | 0 | 0.0 | 1023.0 | 1.0 |
| 90. | 0 | 270.0 | | |
| 180. | 0 | 270.0 | 1023.0 | 1.0 |
| 90. | 0 | 0.0 | | |
| 90. | 0 | 0.0 | 1023.0 | 1.0 |
| 90. | 0 | 0.0 | 1023.0 | 1.0 |
| 90. | 0 | 270.0 | | |
| 180. | 0 | 270.0 | 1023.0 | 1.0 |
| 90. | 0 | 0.0 | | |
| 90. | 0 | 0.0 | 1023.0 | 1.0 |
| | | | | |

There is no real problem with this definition—as long as it is used together with obsprgon! Now suppose somebody wanted to use this same pattern definition for the spinlock statement. Suddenly the experiment would not work, because with spinlock the transmitter gate in on by default, and we are not explicitly gating the transmitter off during the gaps! The conclusion is that if gating is used *within* a pattern, specify it for *all* slices. This makes the pattern usable with both types of functions:

```
# Clean MLEV-17
90.0 0.0 1023.0 1.0
90.0 270.0 1023.0 0.0
180.0 270.0 1023.0 1.0
90.0 0.0 1023.0 0.0
90.0 0.0 1023.0 1.0
90.0 270.0 1023.0 1.0
90.0 270.0 1023.0 1.0
90.0 0.0 1023.0 1.0
90.0 0.0 1023.0 0.0
90.0 0.0 1023.0 1.0
```

The file /vnmr/psg/wg.c also defines a statement gen2spinlock that allows simultaneous (asynchronous) spinlocking and decoupling on two channels; currently no macro calls this undocumented statement. In general, it is as easy to use the obsprgon type of statement to achieve the same thing; therefore, this will not be further discussed here.

How Does Pattern Modulation Work Internally?

Let's first have a look at the Acode generated by *status*-based waveform generator modulation

Acode and Pattern File for Programmed Modulation

The following Acode is generated when using the waveform generator for modulated (broadband) decoupling using the first decoupler channel on a UNITY spectrometer:

| 300 | 294 | 188 | 102 | WG3 | AP addr $0x0c18$, IB addr = $0x0000$ |
|-----|-----|-----|-----|----------------|---------------------------------------|
| 303 | 297 | 191 | 101 | WGCMD | AP addr $0x0c18$, WFG cmd = $0x07$ |
| 306 | 300 | 194 | 150 | HighSpeedLINES | DECUP WFG2 |
| 309 | 303 | 197 | 150 | HighSpeedLINES | DECUP WFG2 |
| 312 | 306 | 200 | 151 | EVENT1_TWRD | 1000 msec |
| 314 | 308 | 202 | 150 | HighSpeedLINES | WFG2 |
| 317 | 311 | 205 | 150 | HighSpeedLINES | (void) |
| 320 | 314 | 208 | 101 | WGCMD | AP addr $0x0c18$, WFG cmd = $0x80$ |

There are three principal differences to the Acode produced for a shaped pulse:

- The (decoupler) transmitter gate is turned on explicitly (overriding the waveform generator's transmitter gating).
- The waveform generator is started using control code 0x07 (for infinite looping).
- At the end of the modulation period the waveform generator is stopped using control code 0x80 (reset or abort), causing the modulation to halt immediately.

The Acode for a UNITY plus is very similar:

| 380 | 374 | 268 | 150 | HighSpeedLINES | DEC WFG2 |
|-----|-----|-----|-----|----------------|-------------------------------------|
| 383 | 377 | 271 | 102 | WG3 | AP addr 0x0c18, IB addr = 0x0000 |
| 386 | 380 | 274 | 101 | WGCMD | AP addr $0x0c18$, WFG cmd = $0x07$ |
| 389 | 383 | 277 | 150 | HighSpeedLINES | DEC WFG2 |
| 392 | 386 | 280 | 150 | HighSpeedLINES | DEC WFG2 |
| 395 | 389 | 283 | 151 | EVENT1_TWRD | 1000 msec |
| 397 | 391 | 285 | 150 | HighSpeedLINES | DEC WFG2 |
| 400 | 394 | 288 | 150 | HighSpeedLINES | WFG2 |
| 403 | 397 | 291 | 150 | HighSpeedLINES | (void) |

The main difference is that on the UNITY *plus* the waveform generator modulation is stopped by resetting the fast status line, and no AP bus traffic is required to achieve this. The Acode was created using dm='ynn', with a d1 of 1 second.

The Acode produced by the obsprgon and obsprgoff combined with the spinlock type of statements is identical, except that these statements do not turn on the pulse programmer rf gate. For the obsprgon type of statement, this is normally done through explicit gating in the pulse sequence using xmtron, decon, etc., which would make the Acode indistinguishable from the code above.

In the case of spinlock and the related statements, the system does *not* use finite looping on the waveform generator. The waveform generator has only an 8-bit loop counter, allowing for 255 loop cycles maximum. The number of loop cycles in spin locking may be much larger than that. Instead, the waveform generator is started with *infinite* looping. From the loop count, the total duration of the spinlock time is calculated. In the Acode the system performs a delay of that length and then stops the waveform generator (using control code 0x80 on a UNITY or by resetting the fast status line on a UNITY*plus*).

Let's see the waveform generator data (instruction block and pattern data) for WALTZ-16 decoupling during a status field:

INSTRUCTION BLOCK: AP address = 0x0c18, WG start address = 0000, 6 words _____ 0x08ffcc00IB_START:RAM start address = 0xffcc0x20fff001IB_STOP:RAM stop address = 0xfff0, delay count = 10x40010000IB_SCALE:loop count = 1 0x80000000 IB_WAITHS: Wait for high-speed line trigger 0xa000063f IB_PATTB: Time count = 1600 (0.00008000 sec / 80.00 usec) 0xe0000000 IB_SEQEND: End of instruction block RF DATA BLOCK: AP address = 0x0c18, WG start address = 0xffcc, 36 words _____ count amplitude phase gates ____ _____ ____ 1024 180.00 3 1024 0.00 4
 1024
 0.00

 1024
 180.00

 1024
 0.00

 1024
 180.00

 1024
 0.00

 1024
 0.00
 2 3 1 2 4 1024 180.00 1024 1024 2 0.00 3 180.00 1024 3 0.00 1024 180.00 4 2 1024 0.00
 1024
 180.00

 1024
 0.00

 1024
 180.00
 3 1 2 1024 4 0.00 2 1024 180.00 1024 0.00 3 1024 0.00 3 1024 0.00 1024 180.00 4 1024 2 0.00 3 1024 180.00
 1024
 0.00

 1024
 180.00

 1024
 0.00
 1 2 4 1024 180.00 2 3 1024 0.00 31024180.00410240.0021024180.00 3 1024 0.00 1024 180.00 1 1024 0.00 1024 180.00 1024 ^ 0.00 2 4 2 3 1024 180.00

END OF FILE

In this case, the parameter dres was set to 90, and the duration count in each slice corresponds to the tip angle divided by the resolution (90 degrees) in each slice. The parameter dmf (i.e., the inverse 90 degrees pulse length) was set to 12500; therefore, the pattern time units is set to 80 microseconds by the IB_PATTB word in the instruction block.

The (infinite) looping is not specified in the instruction block. This is done through the control code (0×07) sent to the waveform generator at run time. The instruction and data blocks generated by obsprgon (decprgon, etc.), are identical to those

generated by status-based decoupling (assuming the same parameters and pattern selection): the transmitter gate is off by default.

The waveform generator data are different for the spinlock statement and its relatives. Here the transmitter is on by default—not through the pulse programmer (because this would override transmitter gating from the waveform generator), but through the waveform generator itself. The following data are for (decoupler based) MLEV-16 spinlocking:

```
INSTRUCTION BLOCK:
AP address = 0x0c18, WG start address = 0000,
                                                                                                                               6 words
                  _____
0x08ffbf00 IB_START: RAM start address = 0xffbf
0x20fff001 IB_STOP: RAM stop address = 0xfff0, delay count = 1
0x40010000 IB_SCALE: loop count = 1
0x80000000 IB_WAITHS: Wait for high-speed line trigger
0xa000063f IB_PATTB: Time count = 1600 (0.00008000 sec / 80.00 usec)
0xe0000000 IB_SEQEND: End of instruction block
RF DATA BLOCK:
AP address = 0x0c18, WG start address = 0xffbf, 49 words
                                                                                                                             -----
 _____
count amplitude phase gates
 ----- ------
                                                               ____

        1
        1024
        0.00
        TXON

        2
        1024
        90.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        2
        1024
        90.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        180.00
        TXON

       2 1024 270.00 TXON
       1 1024 180.00 TXON
               1024 180.00 TXON
1024 270.00 TXON
       1
       2
                  1024 270.00 TXON
1024 180.00 TXON
       1
       1
                  1024 180.00 TXON
       2 1024 270.00 TXON

        1024
        180.00
        TXON

        1024
        0.00
        TXON

        1024
        90.00
        TXON

       1
       1
       2

        2
        1024
        90.00
        TAON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        2
        1024
        90.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        1000
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        180.00
        TXON

                1024 270.00 TXON
       2
       1 1024 180.00 TXON
       1 1024 180.00 TXON
2 1024 270.00 TXON
1 1024 180.00 TXON
       1
                1024 180.00 TXON
       2
                1024 270.00 TXON
                    1024 180.00 TXON
1024 0.00 TXON
       1
       1
                   1024 0.00 TXON
1024 90.00 TXON
       2
                 1024 0.00 TXON
       1

        1
        1021
        0.00
        TXON

        1
        1024
        0.00
        TXON

        2
        1024
        90.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        0.00
        TXON

        2
        1024
        90.00
        TXON

        1
        1024
        0.00
        TXON

        1
        1024
        180.00
        TXON

        2
        1024
        270.00
        TXON

        1
        1024
        180.00
        TXON
```

| 1 | 1024 | 180.00 | TXON |
|---|------|--------|------|
| 2 | 1024 | 270.00 | TXON |
| 1 | 1024 | 180.00 | TXON |
| 1 | 1024 | 0.00 | TXON |
| 2 | 1024 | 90.00 | TXON |
| 1 | 1024 | 0.00 | TXON |
| | | | |

END OF FILE

The transmitter gate is switched on throughout the pattern, but not within the instruction block (for pulses, the transmitter would also be on during eventual delays in the instruction block). The parameter dmf was again set to 12500.

The Influence of the dres Parameter

The duration count in a modulation pattern should be proportional to the tip angle that is specified in the first column of the pattern file. The number of counts per tip angle depends on the dres parameter that defines the tip angle resolution. For modulation pattern like WALTZ-16 or MLEV-16, where all the pulses are in multiples of 90 degrees, dres can be set to 90, which gives one count for a 90-degree pulse, two counts for a 180-degree pulse, etc. (as shown in the above examples). For sequences like MLEV-17 or GARP-1 dres=90 would not work, because it only allows for tip angles in multiples of 90 degrees: all pulses in a pattern should ideally be multiples of the tip angle resolution dres; otherwise, they are *rounded to multiples of* dres in the final pattern.

Two examples: The MLEV-17 modulation pattern ends with a 60 degrees pulse, all other pulses are 90 or 180 degrees. In order to get an accurate reproduction of that pattern, dres must be set to 30, 15, 7.5 or 3.75 (the last value only for the UNITY*plus*). For reasons shown below, the preferred value 30. The GARP-1 modulation pattern consists of 100 pulses with 23 different tip angles ranging from 26 to 268 degrees:

GARP-1 Broadband Decoupling Sequence

| 31.0 | 0.0 |
|-------|-------|
| 55.0 | 180.0 |
| 258.0 | 0.0 |
| 268.0 | 180.0 |
| 69.0 | 0.0 |
| 62.0 | 180.0 |
| 85.0 | 0.0 |
| 92.0 | 180.0 |
| 135.0 | 0.0 |
| 256.0 | 180.0 |
| 66.0 | 0.0 |
| 46.0 | 180.0 |
| 26.0 | 0.0 |
| 73.0 | 180.0 |
| | |
| | |

The only way to get an accurate reproduction of such a pattern is with dres=1. This produces the following instruction and data blocks for the waveform generator:

INSTRUCTION BLOCK: AP address = 0x0cl8, WG start address = 0000, 6 words 0x08ff7a00 IB_START: RAM start address = 0xff7a 0x20fff001 IB_STOP: RAM stop address = 0xfff0, delay count = 1 0x40010000 IB_SCALE: loop count = 1 0x80000000 IB_WAITHS: Wait for high-speed line trigger 0xa0000011 IB_PATTB: Time count = 18 (0.00000090 sec / 0.90 usec) 0xe0000000 IB_SEQEND: End of instruction block

| | A BLOCK: ress = 0x0c | 18, WG s | tart address = 0 | xff7a, | 118 words |
|-------|-------------------------|----------|------------------|--------|-----------|
| count | amplitude | phase | gates | | |
| 31 | 1024 | 0.00 | | | |
| 55 | 1024 | 180.00 | | | |
| 255 | 1024 | 0.00 | | | |
| 3 | 1024 | 0.00 | | | |
| 255 | 1024 | 180.00 | | | |
| 13 | 1024 | 180.00 | | | |
| 69 | 1024 | 0.00 | | | |
| 62 | 1024 | 180.00 | | | |
| 85 | 1024 | 0.00 | | | |
| 92 | 1024 | 180.00 | | | |
| 135 | 1024 | 0.00 | | | |
| 255 | 1024 | 180.00 | | | |
| 1 | 1024 | 180.00 | | | |
| 66 | 1024 | 0.00 | | | |
| 46 | 1024 | 180.00 | | | |
| 26 | 1024 | 0.00 | | | |
| 73 | 1024 | 180.00 | | | |
| | | | | | |

The maximum duration count is 255. If a tip angle results in a greater duration count, these counts are spread over several pattern words. In the above example, the parameter dmf was left at a value of 12500. This results in a duration unit of 900 nanoseconds, as can be seen from the instruction block above. The smaller the dres parameter is set, the smaller a duration unit is obtained.

The minimum duration unit is 200 nanoseconds. Even high modulation rates (large decoupling ranges) can be performed at a high tip-angle resolution (low values for dres). It seems that all modulations can be performed with the highest possible tip angle resolution. The only price we seem to pay is that for large tip angles, additional pattern words are required, but the result should be the same.

There is a possible drawback with small values in dres, however. The smaller the duration units (which are proportional to dres and inversely proportional to dmf), the larger a timing error we can get. With a duration unit of 0.5 microseconds and a timing resolution of 0.1 microseconds (UNITY), the round-off error can be up to 10%, on a UNITY *plus* with 50-nanosecond timing resolution it can still be up to 5%. The conclusion is that the smaller a dres is selected, the greater the overall timing error. If the duration unit is rounded down by 10%, all pulses in the pattern are too short by 10%. The software checks on the timing error when constructing the instruction block, and any time base (round-off) error above 2% is reported. Of course, we can adjust dmf (or the pulse length, the second argument to the spinlock and obsprgon functions), but this just avoids the error and does not remove the error in the effective tip angle in every single pulse in the pattern⁸!

The conclusion is that the tip angle should be selected *as large as the pattern permits*, and, unless a pattern like GARP-1 requires it, it is better not to work with the minimum value in dres and related parameters.

 $^{^{8}}$ In order to keep the modulation algorithm working accurately, we would have to adjust the power to the new (executable) value of dmf: a 5% error in the pulse length requires a power adjustment of about 1 dB.

How small can the tip-angle resolution be? It would not make sense to allow for tip angle resolutions as small as 0.1 degree because such values would lead to extremely small duration units, even at moderate modulation rates, and round-off errors would be obtained very often. Also, the pattern file would be blown up excessively—a 270degree tip angle would result in 11 pattern words. In addition to that, a tip angle resolution much below one degree does not make sense anyway, because a 90 degree pulse is never calibrated with this accuracy. The software, therefore, prevents setting dres and related parameters incorrectly, and the tip-angle resolution in obsprgon and spinlock type functions in various places:

- The parameter dres (as well as dres2, dres3) has a lower limit of 1 (degree). This can be changed easily using the VNMR setlimit command.
- In /vnmr/psg/cps.c, values smaller than 1 in these parameters are reset to a minimum of 1. This means, that for status-field controlled modulation dres can't be any smaller than 1.
- The obsprgon and spinlock families of statements limit the tip-angle resolution to a minimum of 0.7 degrees (in /vnmr/psg/wg.c). Altering this would require changes to the pulse sequence overhead (i.e., using psggen to change the precompiled libraries). Using this lower limit requires using a parameter different from dres (this would be reset to 1, see above) or (assuming modified parameter limits in VNMR) using getval("dres") as third argument to obsprgon/spinlock.

16.5 What If a Waveform Generator Is Not Available

For decoupling, the main disadvantage in not having a waveform generator is that the user is limited to the built-in, pre-programmed decoupling methods.

Programmed Decoupling

On a UNITY*plus*, programmed decoupling includes GARP-1, WALTZ-16, MLEV-16, XY-32 and others—a range that should cover all of today's needs for broadband decoupling. There are more limitations on UNITY and earlier spectrometers, where in essence only WALTZ-16 is available for broadband decoupling. This can be an experimental limitation, especially in cases where a large chemical shift range needs to be decoupled, where GARP-1 permits using much lower power, therefore reducing the sample heating. To some degree, it is possible to explicitly program other, synchronous decoupling methods as explicit acquisition with pulses between the sampling points, but for moire complex methods like GARP-1 this is rather complex, to say the least.

For single-broadband UNITY systems performing indirect detection, a waveform generator on the observe channel can be a real benefit, in that it allows for continuous, asynchronous broadband X decoupling. Without waveform generator the only possibility is to explicitly program and acquisition loop with interleaved decoupling pulses. On systems with an output board (63-word FIFO), this is limited to WALTZ-4 and also may suffer from power or duty-cycle problems, because for large proton spectral windows a composite $(90_x-180_{-x}-270_x)$ pulse may easily be longer than the dwell time. On systems with an acquisition control board (1024-word FIFO), XY-32 can be used, which definitely outperforms WALTZ-4.

For spinlocking experiments (modulation on the transmitter channel), only explicitly coded modulation is possible. This is no real limitation, because even complex modulation algorithms (up to several hundred pulses, depending on the pulse programmer) can be implemented using software or hardware looping.

Shaped Pulses

With a slice width of down to 200 nanoseconds, the waveform generator is unique in its capability to shape even short, hard pulses. The range of possibilities for systems without waveform generators depends on the type of instrument:

UNITY and Earlier Spectrometers

Systems without a waveform generator do not have a linear amplitude modulator, so that the only possibility for performing shaped pulses is through rf attenuators. Several hardware conditions must be fulfilled before pulse shaping can be done successfully:

- Linear amplifiers are prerequisite, otherwise the power level cannot be regulated well enough.
- The T/R switch is also required, otherwise (with crossed diodes) the pulse shapes would be distorted dramatically at voltage levels below 0.5 volts (see Section 4.2, "How Do Pulses Work?," on page 40).
- The standard 63 or 79 dB attenuator does not provide a power range that allows shaping selective pulses (with excitation bandwidths below several hundred Hz). An additional fast, switchable (PIN-diode-based) attenuator is required for the channel on which pulse shaping should occur. Early systems did not have the ports for additional (fine) attenuators. On these, the only way out would be to take the decoupler attenuator and route it into the transmitter channel, in series with the standard attenuator. This leaves the system without decoupler, or with a decoupler with fixed attenuation. Systems with newer AP interface cards have ports for fine attenuators (normally used for solids NMR experiments). These ports can also be used (improperly) for an additional 63-dB attenuator, which then gives the power range required for shaping even very selective pulses.

If these hardware conditions are fulfilled, the include file shape_pulse.c from /vnmr/psg can be used. It contains a statement for performing Gauss- and Hermite-type shaped pulses using attenuators. The limitations are that only built-in (programmed) pulse shapes can be used, small-angle phase shifting is not provided, and the minimum slice width is 10 microseconds (which precludes shaping hard pulses). For more information, see the text of the include file⁹. Shape definitions for waveform generators cannot be used for attenuators, because attenuators are calibrated in dB, whereas the amplitude modulator on the AM/PM transmitter board is linear.

UNITYplus

UNITY*plus* users are in a more favorable position. Both the amplitude and the phase modulators are included with every rf channel and can be addressed via the AP bus, even if no waveform generator is present. This permits using the shape definitions for

⁹ More powerful functions have been submitted to the user library, allowing for free shape programming through tables, including more pre-programmed pulse shapes (such as the BURP-type of pulse shapes), or even including small angle phase shifting.

waveform generators, basically allowing for any shapes, even those using small-angle phase shifting. This has been realized in the following statements:

```
apshaped_pulse(name,width,phase,tbl1,tbl2,rx1,rx2);
decapshaped_pulse(name,width,phase,tbl1,tbl2,rx1,rx2);
dec2apshaped_pulse(name,width,phase,tbl1,tbl2,rx1,rx2);
```

which are macros calling the following C function:

gen_apshaped_pulse(name,width,phase,tbl1,tbl2,rx1,rx2,device);

The arguments are the same as for the statements and macros for waveform generators, with the exception of two *unused* phase table names (tl to t60) that must be supplied with every call (for multiple calls of these functions, new names have to be supplied with every call). These tables are used to internally store the amplitude and phase vectors that are read out of the specified *standard waveform generator shape file* (shapelib/*.RF).

The C function gen_apshaped_pulse decodes the shape file and stores the phase and amplitude vectors in auto-incrementing tables. It then performs the shaped pulse in a real-time loop. The minimum slice length is given by the time it takes to set the power on the linear modulator and the small-angle phase shift via the AP bus (2.3 + 3.45 microseconds, 5.75 microseconds in total, plus a minimum delay of 0.2 microseconds, resulting in a minimum slice length of 5.95 microseconds). The central part of this function is coded as follows¹⁰:

```
/*-----| Calculate
time spent for AP bus events within each slice
+----* /
aptime = POWER_DELAY + SAPS_DELAY;
if ( (plength - aptime) < MINDELAY)
{
  text_error("apshaped_pulse: pulse too short or too many \
        elements in shape");
  abort(1);
}
/*-----+
| set 90 degrees phase, set phase step size |
+----* /
txphase(phs);
stepsize(0.25,OBSch);
/*-----+
| gate receiver off, do rx1 delay |
+----*/
rcvroff();
if (rx1 - aptime > 0.0)
  delay(rx1-aptime);
before turning on the transmitter and entering the
pulse loop preset phase and amplitude for the first
| slice, to avoid a glitch at the start of the pulse
+----*/
```

¹⁰ This is *not* the original coding. To keep the text understandable, low-level functions have been replaced by their standard equivalents that are also used in pulse sequences. The coding shown here is written specifically for the observe channel. The "real" function uses generalized, low-level functions that work for all rf channels.

```
pwrf(pwrtbl,OBSch);
xmtrphase(phstbl);
txphase(phs);
decr(v12);
/*_____
| turn on transmitter, then in a soft loop execute one |
| slice, then set power and phase for the next slice
| after the loop do the last slice, then switch rf off
      ----* /
xmtron();
loop(v12, v13);
 delay(plength-aptime);
  pwrf(pwrtbl,OBSch);
  xmtrphase(phstbl);
  txphase(phs);
endloop(v13);
delay(plength);
xmtroff();
rlpwrf(4095.0, OBSch);
rlxmtrphase(zero);
txphase(phs);
/*-----+
| perform rx2 delay, gate receiver back on again |
+----*/
if (rx2 - aptime > 0.0)
  delay(rx2-aptime);
rcvron();
```

In this code, plength is the slice length, calculated as pulse length divided by the number of slices in the shape file. Slices with duration counts greater than 1 are translated into multiple slices with a duration count of 1. The real-time variable v12 contains the number of loop cycles. pwrtbl and phstbl are the two table names that are given as arguments; they have been "filled" with the amplitude and table values from the shape file. Note that the amplitude values from the table have to be multiplied by 4 in order to get the proper amplitude range, because the waveform generator only addresses the ten most significant bits (0 to 1023), whereas through the AP bus the full range (12 bits, 0 to 4095) can be used.

As an example, if we take the 256-step gaussian pulse that was discussed previously in this chapter, we obtain a rather impressive piece of Acode from this single function:

| 399 | 393 | 287 | 150 | HighSpeedLINES | (void | 1) | | | |
|-----|-----|-----|-----|----------------|-------|--------|---------|--------|------------|
| 402 | 396 | 290 | 39 | ASSIGNFUNC | zero | v12 | | | |
| 405 | 399 | 293 | 31 | MULtFUNC | three | three | v13 | | |
| 409 | 403 | 297 | 31 | MULtFUNC | v13 | three | v13 | | |
| 413 | 407 | 301 | 31 | MULtFUNC | v13 | three | v13 | | |
| 417 | 411 | 305 | 31 | MULtFUNC | v13 | three | v13 | | |
| 421 | 415 | 309 | 29 | ADDFUNC | v12 | v13 | v12 | | |
| 425 | 419 | 313 | 32 | DIVFUNC | v13 | three | v13 | | |
| 429 | 423 | 317 | 32 | DIVFUNC | v13 | three | v13 | | |
| 433 | 427 | 321 | 32 | DIVFUNC | v13 | three | v13 | | |
| 437 | 431 | 325 | 29 | ADDFUNC | v12 | v13 | v12 | | |
| 441 | 435 | 329 | 32 | DIVFUNC | v13 | three | v13 | | |
| 445 | 439 | 333 | 29 | ADDFUNC | v12 | v13 | v12 | | |
| 449 | 443 | 337 | 32 | DIVFUNC | v13 | three | v13 | | |
| 453 | 447 | 341 | 29 | ADDFUNC | v12 | v13 | v12 | | |
| 457 | 451 | 345 | 16 | SETPHAS90 | CH1 | oph | | | |
| 460 | 454 | 348 | 68 | PHASESTEP | CH1 | 1 | units | (0.25 | degrees) |
| 463 | 457 | 351 | 150 | HighSpeedLINES | RXOFI | F | | | |
| 466 | 460 | 354 | 150 | HighSpeedLINES | RXOFI | 7 | | | |
| 469 | 463 | 357 | 151 | EVENT1_TWRD | 4.250 |) usec | | | |
| 471 | 465 | 359 | 105 | TABLE 360 size | 256, | autoin | c 1, di | ivn_re | t 1, ptr 0 |
| | | | | | | | | | |

| | | | | 45 | 49 | 52 | 56 | 60 | 64 | 69 | 73 | |
|------------|------------|------------|-----------|--|--|---|---|---|---|---|---|-------|
| | | | | 78 | 84 | 89 | 95 | 102 | 108 | 115 | 123 | |
| | | | | 131 | 139 | 148 | 157 | 166 | 176 | 187 | 198 | |
| | | | | 210 | 222 | 235 | 249 | 263 | 277 | 293 | 309 | |
| | | | | 326 | 343 | 362 | 381 | 401 | 421 | 443 | 465 | |
| | | | | 488 | 512 | 537 | 563 | 590 | 617 | 646 | 676 | |
| | | | | 706 | 738 | 770 | 804 | 838 | 874 | 910 | 948 | |
| | | | | | | | | 1150 1523 | | | | |
| | | | | | | | | 1949 | | | | |
| | | | | | | | | 2406 | | | | |
| | | | | | | | | 2869 | | | | |
| | | | | | | | | 3302 | | | | |
| | | | | | | | | 3669 | | | | |
| | | | | | | | | 3936 | | | | |
| | | | | 4024 | 4040 | 4055 | 4067 | 4077 | 4085 | 4091 | 4094 | |
| | | | | 4095 | 4094 | 4091 | 4085 | 4077 | 4067 | 4055 | 4040 | |
| | | | | 4024 | 4005 | 3984 | 3961 | 3936 | 3909 | 3880 | 3850 | |
| | | | | | | | | 3669 | | | | |
| | | | | | | | | 3302 | | | | |
| | | | | | | | | 2869 | | | | |
| | | | | | | | | 2406 | | | | |
| | | | | | | | | 1949 | | | | |
| | | | | | | | | 1523 | | | | |
| | | | | | 1283 948 | | | 1150 | | | | |
| | | | | 986 706 | 940 676 | 910 646 | 874 617 | 838 590 | 804 563 | 770 537 | 738 512 | |
| | | | | 488 | 465 | | | 401 | 381 | 362 | 343 | |
| | | | | 326 | 309 | | 277 | 263 | 249 | 235 | 222 | |
| | | | | 210 | 198 | 187 | 176 | 166 | 157 | 148 | 139 | |
| | | | | 131 | 123 | 115 | 108 | 102 | 95 | 89 | 84 | |
| | | | | 70 | 73 | 69 | 64 | 60 | 56 | 52 | 49 | |
| | | | | 78 | 15 | 0,0 | 01 | 00 | 50 | 52 | 17 | |
| 732 | 726 | 620 | 106 | | IGN ta | | | tblı | | 52 | 17 | |
| 732 735 | 726 729 | 620 623 | 106 59 | TASS | IGN ta | able 3 | 360 | tblı | rt | | | bytes |
| | | | | TASS | IGN ta | able : APado | 360 dr 11. | tblı | rt 150, | -logi | ic, 2 | - |
| | | | | TASS APCh TABLE | IGN ta ipOUT 630 | able : APado max 4 | 360 dr 11 4095, | tblı , reg offse auto | rt 150, et 0, | -logi value , divi | ic, 2 e tbl: n_ret | - |
| 735 | 729 | 623 | 59 | TASS APCh TABLE 0 | IGN ta ipOUT 630 0 | able 3 APado max 4 size 0 | 360 dr 11, 4095, 256, 0 | tblı , reg offse auto 0 | rt 150, et 0, inc 1 0 | -logi value , divi 0 | ic, 2 e tbl: n_ret 0 | rt |
| 735 | 729 | 623 | 59 | TASS APCh TABLE 0 0 | IGN ta ipOUT 630 0 0 | Able 3 APado max 4 size 0 0 | 360 dr 11, 4095, 256, 0 0 | tbln , reg offse auto 0 0 | rt 150, et 0, inc 1 0 0 | -logi value , divi 0 0 | ic, 2 e tbl: n_ret 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS APCh: TABLE 0 0 0 | IGN ta ipOUT 630 0 0 0 | Able 3 APado max 4 size 0 0 0 | 360 dr 11 4095, 256, 0 0 | tbln , reg offse auto 0 0 | rt 150, et 0, inc 1 0 0 0 | -logi value , divu 0 0 0 | ic, 2 e tbl: n_ret 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS APCh: TABLE 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 | able 3 APado max 4 size 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 | tbln reg offse auto 0 0 0 | rt 150, et 0, inc 1 0 0 0 0 | -logi value , dive 0 0 0 0 | ic, 2 e tbl: n_ret 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS APCh: TABLE 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 | Able 3 APado max 4 size 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 | tbln , reg offse auto 0 0 0 0 | rt 150, et 0, inc 1 0 0 0 0 0 | -logi value , divi 0 0 0 0 | ic, 2 e tbl: n_ret 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS APCh TABLE 0 0 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 0 0 | able 3 APado max 4 size 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 | tbln , reg offse auto 0 0 0 0 0 0 | rt 150, et 0, inc 1 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 | ic, 2 e tbln n_ret 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 0 0 0 | Able : APado max · size 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 | tbln , reg offse auto 0 0 0 0 0 0 0 | rt 150, et 0, inc 1 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 | ic, 2 e tbln n_ret 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 | tbln , reg offse auto 0 0 0 0 0 0 0 0 0 0 | et 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 | ic, 2 e tbln n_ret 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 0 0 0 | Able : APado max · size 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 | tbln , reg offse auto 0 0 0 0 0 0 0 | rt 150, et 0, inc 1 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 | ic, 2 e tbln n_ret 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln , reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 | ic, 2 e tbln n_ret 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN ta ipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APada max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | et 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 e tbl: n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | EGN ta ipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APada max - size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | et 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 e tbl: n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, st 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbl: n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbl: n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, st 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbl: n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | Able Control APado max A size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | Able Control APado max A size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able 3 APado max max 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, et 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able 3 APado max max 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | Able Control APado max A size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | EGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able 3 APado max size 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt 150, at 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | Able Control APado max A size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | EGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able 3 APado max size 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt 150, at 0, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11, 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN table 1000 table 1 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |
| 735 | 729 | 623 | 59 | TASS: APCh: TABLE 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | IGN taipOUT 630 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | able : APado max 4 size 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | 360 dr 11 4095, 256, 0 0 0 0 0 0 0 0 0 0 0 0 0 | tbln, reg offse auto 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ct 150, inc 1 0 0 0 0 0 0 0 0 0 0 0 0 0 | -logi value , divn 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | ic, 2 = tbln n_ret 0 0 0 0 0 0 0 0 0 0 0 0 0 | rt |

| | | | | 0 0 0 0 0 0 0 | |
|------|------|-----|-----|--|---|
| | | | | 0 0 0 0 0 0 0 | |
| | | | | 0 0 0 0 0 0 0 | |
| 1002 | 996 | 890 | 106 | TASSIGN table 630 tblrt | |
| 1005 | 999 | 893 | 65 | SETPHASE CH1f tblrt | |
| 1008 | 1002 | 896 | 16 | SETPHAS90 CH1 oph | |
| 1011 | 1005 | 899 | 28 | DECRFUNC v12 | |
| 1013 | 1007 | 901 | 150 | HighSpeedLINES RXOFF TXON | |
| 1016 | 1010 | 904 | 39 | ASSIGNFUNC zero v13 | |
| 1019 | 1013 | 907 | 42 | IFMInusFUNC v12 one Offset = 937 | |
| 1023 | 1017 | 911 | 151 | EVENT1_TWRD 44.250 usec | |
| 1025 | 1019 | 913 | 106 | TASSIGN table 360 tblrt | |
| 1028 | 1022 | 916 | 59 | APChipOUT APaddr 11, reg 150, -logic, 2 byte | s |
| | | | | max 4095, offset 0, value tblrt | |
| 1034 | 1028 | 922 | 106 | TASSIGN table 630 tblrt | |
| 1037 | 1031 | 925 | 65 | SETPHASE CH1f tblrt | |
| 1040 | 1034 | 928 | 16 | SETPHAS90 CH1 oph | |
| 1043 | 1037 | 931 | 27 | INCRFUNC v13 | |
| 1045 | 1039 | 933 | 42 | IFMInusFUNC v13 v12 Offset = 911 | |
| 1049 | 1043 | 937 | 151 | EVENT1_TWRD 50.000 usec | |
| 1051 | 1045 | 939 | 150 | HighSpeedLINES RXOFF | |
| 1054 | 1048 | 942 | 65 | SETPHASE CH1f zero | |
| 1057 | 1051 | 945 | 16 | SETPHAS90 CH1 oph | |
| 1060 | 1054 | 948 | 151 | EVENT1_TWRD 4.250 usec | |
| 1062 | 1056 | 950 | 150 | HighSpeedLINES (void) | |

The purpose of the real-time math section at the beginning (not shown in the C code above) is to store the number of slices in a real-time variable without using initval (see also "New Real-Time Numeric Constants" on page 97). The most prominent feature in the following section are the two tables containing all 256 amplitude and phase values of that pulse, taking up a full Kbyte (two times 256 16-bit words) of Acode space. Of course the second table could, in theory, be reduced to a single constant zero, but the program cannot "know" that this shape file contains no phase changes.

Apart from the amount of Acode that is necessary to perform such a pulse through Acode, and some limitations with respect to the modulation rate, everything seems to be fine, and apshaped_pulse in fact is a reasonable replacement for a waveform generator for many experiments, in particular, selective excitation. There are, however, a few fundamental differences and limitations to this solution. Some of them lie in the way the waveform generator works; others are a consequence of the above coding and of the way how phase shifting works.

One problem lies in small-angle phase shifting. Any internal phase shifting on a waveform generator is performed *on top of* any current quadrature and small-angle phase shifts. This cannot be emulated in software. The internal phase shifting in gen_apshaped_pulse is set via xmtrphase or equivalent statements, with work in an absolute frame. They even reset any existing quadrature (90-degree) phase shift! The latter is added in again after the xmtrphase call (gen_apshaped_pulse therefore works properly on top of any quadrature phase shift), but this cannot be done for the small-angle phase shifting. Small-angle phase shift coaddition is not provided, and hence gen_apshaped_pulse does not work properly on top of small angle phase shifting.

The other problem lies in the fact that quadrature phase shift needs to be *reestablished* after each xmtrphase or equivalent call. The statement txphase (or its equivalents for the current channel) is called once for each slice, which implies that

autoincrementing tables cannot be used as phase variables to apshaped_pulse or equivalent macros.

Over all, what are the limitations of the apshaped_pulse approach?

- The minimum slice length is 5.95 microseconds, compared to 0.2 microseconds with the waveform generator. Real short hard pulses cannot be shaped;
- Phase and amplitude changes don't occur simultaneously. The phase and the amplitude profile of the pulse are shifted against each other by 3.45 microseconds;
- Slices with the transmitter gate switched off are simply skipped;
- apshaped_pulse and related functions don't work as expected on top of small angle phase shifting;
- The phase variable cannot be an autoincrementing table;
- Simultaneous shaped pulses are not possible;
- apshaped_pulse generates lots of Acode, which may possibly limit its use for multidimensional (3D, 4D) experiments;
- Even though the Acode size may theoretically allow for a single shape of up to about 4000 slices (see also "Acode Size Limitations, Acode Buffering" on page 83), such long shapes will most likely lead to FIFO underflow. On a UNITY or UNITY *plus* system, the FIFO is only filled at a rate of below one word per 20 μsec¹¹, which may in many cases be the real rate-limiting step.

Nevertheless: apshaped_pulse permits performing most selective excitation experiments (even demanding ones like those using shifted laminar pulses¹² or SS pulses¹³) perfectly without a waveform generator, and some of the inherent limitations can be bypassed by programming measures.

16.6 Using a Waveform Generator for Shaping Gradient Pulses

The waveform generator was first designed for shaping field gradients and rf pulses in imaging experiments. Only later it was adapted for high-resolution machines. By the addition on one component (the amplitude multiplier), a UNITY*plus* waveform generator can be converted for shaping field gradients (of course, the AP address has to be changed by reconfiguring its jumper settings).

The gradient waveform generator works along the same principles as an rf waveform generator, except that it controls the amplitude of a single field gradient instead of phase and amplitude of an rf signal and rf gates. Therefore, a gradient pattern word has a totally different layout, as shown in Table 11.

Bits 0 to 7 are again the duration count (same as for the rf pattern), bits 8 to 23 form a 16-bit field gradient amplitude value (-32768 to 32767), and the bits 24 to 31 are unused (field gradients are dc by nature, so there is no phase parameter¹⁴). Pulsed field gradients may not only be scaled in the time scale, often (mostly, in imaging

¹¹ This is under optimum conditions, without lots of Acode overhead from soft looping or table access, etc. (both are used in the apshaped_pulse function). The actual transfer rate for the apshaped_pulse function has not been determined.

¹² S.L. Patt, J. Magn. Reson . 96, 94 (1992).

¹³ S.H. Smallcombe, J. Am. Chem. Soc. **115**, 4776 (1993).

¹⁴ Instead, the amplitude is signed (it can be positive or negative), but the rf amplitude is measured and regulated in magnitude only.

| Туре | 31 | 30 | 29 | 28 | 27 | 26 | 25 | 24 | 23 | 22 | 21 | 20 | 19 | 18 | 17 | 16 | 15 | 14 | 13 | 12 | 11 | 10 | 9 | 8 | 7 | 6 | 5 | 4 | 3 | 2 | 1 | 0 |
|----------|----|----|----|----|----|----|----|----|------|----|----|----|----|----|----|-----|------|----|------|------|----|----|---|---|---|---|---|------|------|---|---|---|
| rf | S | Р | Т | | | | | I | ohas | e | | | | | | | | а | ımpl | itud | e | | | | | | | dura | tior | 1 | | |
| Gradient | | | | | | | | | | | | | | | а | mpl | itud | e | | | | | | | | | | dura | tior | ı | | |

Table 11. Comparison of waveform generator pattern words

experiments) gradients are scaled in their amplitude. For an imaging experiment, a spatial dimension can be encoded by stepping a gradient through an array of (positive and negative) values.

With the rf definition of a waveform generator, we would have to define as many shapes as there are amplitude values. To avoid this, field gradient shapes are made generic in two ways: they can be scaled in duration (the same way as rf shapes and pattern) and in amplitude. The amplitude value in the pattern is not the final one. The gradient waveform generator has a built-in multiplier that multiplies the amplitude value in every pattern word with a scalar value in the IB_SCALE instruction word of the corresponding instruction block (see also "Waveform Generator Instruction Words" on page 169).

Most imaging experiments use at least one stepped gradient for spatial encoding. Often, there are 32, 64, 128, or 256 different gradient increments. This means that in the pattern file, there are that many instruction blocks pointing to the same data (pattern) block, each of them with a different scaling factor in the IB_SCALE instruction word.

In terms of basic data handling, the gradient waveform generator is treated the same way as its rf equivalent. Here is a (hypothetical) data file for a "hsine"-shaped gradient without amplitude scaling:

```
INSTRUCTION BLOCK: Z GRADIENT
AP address = 0x0c30, WG start address = 0000, 5 words
_____
0x08ff7100 IB_START: RAM start address = 0xff71
0x20fff001 IB_STOP: RAM stop address = 0xfff0, delay count = 1
0x40010001 IB_SCALE: loop count = 1, amplitude scale = 1
0xa000b88b IB_PATTB: Time count = 47244 (0.00236220 sec / 2362.20 usec)
0xe0000000 IB_SEQEND: End of instruction block
DATA BLOCK: Z GRADIENT
AP address = 0x0c30, WG start address = 0xff71, 127 words
count amplitude
      _____
_ _ _ _ _
  1
        810
  1 1620
  1
       2429
  (57 lines deleted)
   . . .
       32704
  1
  1
       32744
  1
       32764
      32764
  1
      32744
  1
  1 32704
  (57 lines deleted)
   . . .
```

| 1 | 2429 |
|---|------|
| 1 | 1620 |
| 1 | 810 |
| 1 | 0 |

END OF FILE

The gradient shape definition files are also stored in shapelib, but with the extension ".GRD". Their format is simpler than the format for rf shapes and modulation pattern. The two columns in the file contain the gradient amplitude and the duration count (usually one count per slice):

half sine for GRADIENTS # 810 1.0 1620 1.0 2429 1.0 (57 lines deleted) 32704 1.0 32744 1.0 32764 1.0 32764 1.0 32744 1.0 32704 1.0 (57 lines deleted) 2429 1.0 1620 1.0 810 1.0 0 1.0

Chapter 16. Waveform Generators

Chapter 17. Pulsed Field Gradients

The most prominent difference between pulsed field gradients (PFG) and non-PFG pulse sequences is the addition of statements that generate the pulsed field gradients. In certain PFG techniques (like coherence pathway selection and multiple-quantum filtering using pulsed field gradients), little or no phase cycling is used. Some PFG sequences use pulses with constant phase only and may therefore look simpler than their non-PFG equivalents. Where non-PFG pulse sequences eliminate artifacts and unwanted signals by phase cycling (e.g., subtraction), PFG sequences dephase unwanted coherences (i.e., make them non-observable, PFG experiments are often less subject to spectral artifacts than their non-PFG equivalents).

The PFG accessory involves a dc and audio frequency power amplifier with linear amplitude control through the AP bus, connected to a special probe with a Z gradient coil. The field gradient coil is shielded. While generating a strong Z field gradient inside (at the sample), an external compensating field (with opposite sign) is generated, such that the total field outside the gradient coil is greatly diminished. This helps reduce eddy currents in the rest of the probe body and in the metal walls of the magnet dewar (although eddy currents can't be totally avoided).

17.1 Pulse Sequence Statements for PFG Gradient Control

The statements for gradient control in typical PFG experiments are simple. They reflect the straightforward, scalar nature of a (linear) field gradient. A single parameter determines the amplitude of the gradient, which can have values between -32768.0 and 32767.0 (the gradient amplifier uses a 16-bit DAC, negative values are permitted and usually required in most PFG experiments).

What also simplifies programming pulsed field gradients is the fact that the gradient amplifier is constructed practically noise-free; therefore, there is no need for blanking it during off-intervals. The amplifier can be put into "standby" mode by setting the VNMR parameter pfgon to 'nnn', while setting pfgon to 'nny' turns on the (Z) gradient amplifier. There is also no fast line (gating) involved with gradient control. Typical gradient pulses are on the order of milliseconds; therefore, it is more than sufficient to set the gradient amplitude with the AP bus. Even gradient shaping can be done this way (see Section 17.2, "Shaping Pulsed Field Gradients," on page 202).

The most basic statement to set the gradient amplitude is rgradient¹:

rgradient(gid,amplitude);

For PFG experiments, the first argument (gid, the gradient identifier) is 'z' or 'Z', the second argument (amplitude) is a number between -32768.0 and 32767.0 (a floating point number of type double).

¹ There is also a statement vgradient that allows defining the amplitude from real-time variables. Typical PFG sequences use a few, predefined gradient levels that do not vary within the sequence or from transient to transient; therefore, using vgradient would be an unnecessary complication. Up to now, vgradient has only been used in imaging sequences (see also Chapter 21, "(Migro)]meging Experiments" on page 230)

[&]quot;(Micro)Imaging Experiments," on page 239).

Chapter 17. Pulsed Field Gradients

A typical pulse sequence construct for a gradient pulse of length gzl and amplitude gzllvl using the rgradient function could be written as follows:

```
rgradient('Z',getval("gzllvl"));
delay(getval("gzl"));
rgradient('Z',0.0);
```

The higher-level statement zgradpulse allows reducing these lines to a single call:

```
zgradpulse(amplitude,duration);
```

Using this statement, the above three lines for a gradient pulse of duration gtl and amplitude gzllvl would be written as

zgradpulse(getval("gz1lvl"),getval("gt1"));

Although this statement has the simplicity of the statements for rf pulses (that use gating through fast lines), one should still not forget the "underlying three lines" (i.e., that there is AP bus traffic before and after the actual gradient pulse, taking up a finite time). While the time spent in AP bus traffic is negligible compared to the duration of the gradient pulse itself, it is long enough to cause considerable dephasing for large chemical shift ranges. It is, therefore, strongly recommended to compensate for this time in the pulse sequence, in particular during a refocusing interval:

```
rgpulse(pw,t1,rof1,0.0);
delay(tau-rof1);
rgpulse(pw,t1,rof1,0.0);
zgradpulse(getval("gz1lvl"), gt1);
delay(tau-gt1-2.0*GRADIENT_DELAY-rof1);
rgpulse(pw,t1,rof1,0.0);
```

Up to now we have taken gradient pulses as ideal, rectangular pulses, which obviously was a simplification because any gradient change is slowed down by eddy currents (this applies to switching both on and off). In first approximation, the time constant of the decay of these eddy currents is given by the geometry and the construction of the gradient coil and its surroundings and should not depend on the gradient strength. Corrective delays have often been implemented in pulse sequences using pulsed field gradients to compensate for gradient pulse imperfections like the finite slew rate. Such corrections (in addition to eddy current compensation) are most necessary in imaging experiments, where eddy currents are a much bigger problem than for high-resolution PFG experiments (see also Chapter 21, "(Micro)Imaging Experiments," on page 239).

The first Varian PFG sequences were written with microimaging experiment concepts in mind. A typical gradient pulse was programmed as follows:

```
rgradient('Z',gzllvl);
delay(gzl+grise);
rgradient('Z',0.0);
delay(grise);
delay(gstab);
```

This generates the following timing scheme (times not shown proportionally):



grise is supposed to be the time intervals before and after the pulse that are affected by the eddy currents (delayed gradient buildup, delayed gradient decay), and gstab is the time required to reestablish full homogeneity after a gradient pulse. It turns out that after a typical gradient pulse of a few milliseconds, it takes about 50 microseconds until the system has recovered from the gradient, such that there are no observable phase errors (to regain full amplitude may take slightly longer). It is certainly not a good idea to have a gradient pulse followed immediately by an rf pulse (or data acquisition), but as long as some delay (of around 50 microseconds) follows the gradient pulse, there should be no need for an additional delay gstab.

It can be assumed that eddy currents affect both the gradient buildup and gradient turnoff times the same way (i.e., what we lose at the beginning of the gradient pulse we regain at the end of the pulse). The gradient pulse is not quite rectangular (neither is any pulse!). All that counts is the *area* (amplitude times duration) of the gradient pulse because that determines the amount of dephasing achieved. To avoid transversal relaxation and spin diffusion, it is desirable to have short gradients. On the other hand, very strong gradients require strong amplifiers that produce lots of noise and, therefore, affect the overall homogeneity.

With PFG probes, the eddy current time constants are definitely below 10 microseconds. Even if the time constants for turning on and off the gradient pulse were slightly different, that difference would be negligible, because their magnitude is only fractions of a percent of the total pulse duration (typically a few milliseconds). Under this assumption, it is certainly *wrong* to compensate for eddy current effects using a finite delay grise, as this prolongs the gradient pulse.

Most PFG pulse sequences use a gradient pulse to refocus magnetization (coherence) that was dephased by a preceding pulse. This can only work if the two gradient pulses have very accurate and well-defined gradient areas. grise in the above scheme would distort the dephasing ratio between two gradient pulses with different areas, because

$$\frac{(gt1 \times gz11v1)}{(gt2 \times gz21v1)} \neq \frac{((gt1 + grise) \times gz11v1)}{((gt2 + grise) \times gz21v1)}$$

Still, many spectroscopists are using grise. By intuition, you might think that these delays "take care of eddy currents." In fact, effects can be seen in many experiments where below certain values for grise (typically 10, sometimes more) some pulse sequences would simply not work properly. What has happened in those cases was that the gradients may have been slightly out of balance due to the addition of grise to their lengths (if they had different lengths or amplitudes) but by an amount that is not noticeable yet. What made the experiments work is that the second grise (in the above scheme) acted as gradient recovery delay.

In many pulse sequences, an rf pulse follows a gradient pulse, and for these cases we need to insert a recovery delay; otherwise, the rf pulse generates phase errors. This is most critical in multiple-quantum filtering experiments, where phase or amplitude distortions after a gradient pulse can seriously hamper the performance of the sequence. For most sequences a recovery delay of 10 to 20 microseconds is adequate and sufficient. For multiple-quantum filtering experiments (like gradient MQCOSY or E.COSY), the recovery time should be adjusted to 50 to 60 microseconds.

In conclusion, grise delays are not recommended (as shown in the scheme above), but if an rf pulse follows a gradient pulse, implementing a gradient recovery delay is strongly recommended. This delay can be made part of an ropulse call:

```
zgradpulse(getval("gz1lvl"),getval("gt1"));
rgpulse(pw,v1,gstab,0.0);
```

17.2 Shaping Pulsed Field Gradients

It has been suggested that using trapezoidal gradients (or gradients with other shapes) would allow minimizing eddy current effects. Because PFG probes use actively shielded gradient coils, eddy current effects are minimized; therefore, it was found that it is not necessary to use shaped gradients. Also, in order to minimize losses due to transverse relaxation and spin diffusion, it is desirable to use short, strong gradient pulses. A shaped gradient by definition has a lower duty cycle and is, therefore, longer than a rectangular gradient with the same "area." Thus, for most liquids applications, it is undesirable to use shaped gradients.

Still, there may be situations where shaped gradient pulses have advantages. For example, in diffusion experiments with very high gradient strengths, a gradient coil with high-inductive load may cause transition problems to the gradient amplifier when turning on the gradient. A trapezoidal gradient may alleviate such problems.

How can gradients be shaped? For imaging, a gradient control unit with waveform generator provides for an easy access to shaped gradients, but the PFG accessory does not include a waveform generator. Therefore, gradients have to be shaped through the AP bus, similar to the approach taken in the apshaped_pulse function (see Section 16.5, "What If a Waveform Generator Is Not Available," on page 189).


A function for creating a (trapezoidal) shaped gradient is shown below:

```
rampgrad(amp, length, ramp)
double amp, length, ramp;
  int i, steps;
  double iramp, initval, incr_val;
  if (length >= MINDELAY)
  {
    if (ramp > length)
    {
      fprintf(stdout,"ramp parameter larger than \
              gradient delay\n");
      abort(1);
    }
    steps = (int) (amp / 1000.0);
    if ((double) steps * GRADIENT_DELAY > ramp)
      steps = (int) (ramp / GRADIENT_DELAY);
    if (steps > 1)
      incr_val = amp / steps;
      initval = incr_val;
      iramp = ramp / (double) steps - GRADIENT_DELAY;
      for (i=0; i<steps; i++)</pre>
      {
        rgradient('Z',initval);
        delay(iramp);
```

{

```
initval += incr_val;
   }
   rgradient('Z',amp);
   if (ramp-(double)steps*(iramp + GRADIENT_DELAY) > MINDELAY)
     delay(ramp - (double) steps * (iramp + GRADIENT_DELAY));
 }
 else
  {
  rgradient('Z',amp);
  if (ramp - GRADIENT_DELAY >= 2e-7)
    delay(ramp - GRADIENT_DELAY);
 if (length - ramp >= MINDELAY) delay(length - ramp);
 if (steps > 1)
  {
   for (i=0; i<steps; i++)</pre>
   {
     rgradient('Z',initval);
     delay(iramp);
     initval -= incr_val;
   }
   rgradient('Z',0.0);
   if (ramp - (double)steps*(iramp + GRADIENT_DELAY) >= \
       GRADIENT_DELAY)
     delay(ramp - (double)steps * (iramp + GRADIENT_DELAY));
 }
 else
  {
   rgradient('Z',0.0);
   if (ramp - GRADIENT_DELAY >= MINDELAY)
     delay(ramp - GRADIENT_DELAY);
 }
}
```

This function creates a shape with a linear ramp at the beginning and at the end of the gradient pulse. The amplitude is changed in steps of 1000 ADC units or 1/32 of the maximum amplitude (or in bigger steps if the ramp duration isn't long enough). The shape is calculated such that the gradient "area" is the same as for a rectangular gradient of the same specified length:



Obviously, when the specified length is less than the duration of one ramp (specified in the third argument), this could not be maintained; therefore, an error message is produced.

}

17.3 PFG Experiments Using Homospoil Pulses

Spectroscopists who don't have access to a PFG accessory might ask whether it is possible to perform PFG experiments using the standard gradient (shim) coils. There are several problems with this approach:

- The achievable gradient strength is very limited. This would lead to very long gradient pulses (probably 10 to 100 milliseconds) causing losses due to transverse relaxations and spin diffusion.
- The shim coils are not shielded; therefore, the recovery time from a shim coil gradient pulse is several milliseconds. The gradient recovery delays would further increase the losses.
- There are no user-accessible functions that allow setting shim DACs from within a pulse sequence (although certainly such a function could be created).
- Even if a user-accessible function were created, its use would still be very unhandy. At least for one of the polarities, the achievable amplitude may be strongly limited if the standard gradient amplitude is far off from zero.

Of course, we can decide to use the homospoil pulse to avoid the above problems. The homospoil pulse can be triggered with existing statements in a pulse sequence (see "Delays With Homospoil Pulse" on page 38), but the amplitude for the homospoil pulse is not under software control. In particular, the sign of the homospoil pulse amplitude cannot be altered, which excludes a large number of PFG experiments.

Overall, it would certainly be possible to perform a limited range of (simple) PFG experiments using homospoil pulses, but there are rather serious limitations, especially for larger molecules and other samples with short transverse relaxation times. This method will therefore not be discussed any further here.

Chapter 18. Acquiring Data

A prominent feature in most VNMR pulse sequences is the absence of code dealing with data acquisition. Of course, there is a statement that does data acquisition:

acquire(points,dwell_time);

where points is the number of data points to be acquired and dwell_time is the time interval between the sampling trigger pulses, or the time that is apparently taken to acquire a complex data point¹. Because the receiver and ADC acquire data in complex points (real + imaginary), points (the number of values measured) must be a multiple of two.² At the same time, the dwell_time is equal to the inverse spectral window and determines the largest frequency that can be measured. Higher frequencies are folded in (this is also called "aliasing").

Larger spectral windows require shorter dwell times. With a 0.2-microsecond minimum time event, the maximum spectral window which the pulse programmers can possibly handle is 5 MHz^3 .

Still, most pulse sequences end with the last pulse. How then is the acquisition done?

18.1 Implicit Acquisition

One thing that acquire does (apart from coding the specified number of sampling intervals) is to increment a global variable acqtriggers that initially (before calling the function pulsesequence) is set to zero. If a pulse sequence does not contain the acquire statement (i.e., if there was no explicit acquisition, see Section 18.2, "Explicit Acquisition," on page 208), acqtriggers still is set to zero after the call to pulsesequence, and the software automatically completes the sequence with an implicit acquisition.

The following code segments are found in the function createPS from within the module psg/cps.c:

```
acqtriggers = 0;
...
pulsesequence(); /* generate Acodes from USER pulse sequence */
...
test4acquire(); /* if no acquisition done yet, do it */
write_Acodes(); /* write out generated lc, auto & instructions */
return;
```

¹ The acquire statement is coding points / 2 periods (with a duration as specified by the second argument), each starting with an ADC trigger pulse.

² If using the Output board (63-word FIFO), points must be 2 or a multiple of 64.

³ There is no need for "off" times in-between sampling intervals. Every time event can trigger a complex data point if the CTC bit (the "command to convert" or acquisition trigger) is set to high: unlike normal fast lines, the CTC is electronically reformed into a very short trigger pulse at the very beginning of each sampling interval (see also Chapter 9, "Pulse Programmers," on page 85).

Function test4acquire is found in psg/hwlooping.c; and starts with the lines:

```
if (acqtriqgers == 0)
                            /* No data acquisition yet? */
{
  if (nf > 1.0)
  {
     text_error("Number of FIDs (nf) Not Equal to One\n");
     abort(0);
  if (ap_interface < 4)
     HSgate(rcvr_hs_bit,FALSE); /* turn receiver On */
  else
     SetRFChanAttr(RF_Channel[OBSch],SET_RCVRGATE,ON,0);
  for (i = 1; i <= NUMch; i++) /* zero HS phaseshifts */</pre>
     SetRFChanAttr(RF_Channel[i],SET_RTPHASE90,zero,0);
  acqdelay = alfa + (1.0 / (beta * fb));
  /* acquire data */
  acquire(np,1.0/sw);
}
. . .
```

From these lines of code, we can extract the following information on the implicit acquisition:

- Implicit acquisition does *not* work for multi-FID experiments (nf > 1), such as multiecho imaging experiments (see Chapter 21, "(Micro)Imaging Experiments," on page 239) or sequences like COCONOESY (combined COSY—NOESY), i.e., sequences that acquire more than one FID *within* the pulse sequence.
- The receiver is gated on as part of the implicit acquisition. It should be harmless if this is forgotten within the pulse sequence.
- The quadrature phase of all transmitters is reset to zero before acquiring data. In combination with receiver phase cycling, this avoids coherent signal buildup in case of an rf leakage, which could produce a large, narrow peak ("glitch") in the center of the spectrum (with signal averaging this can only hurt in the case of the observe transmitter itself, or with decoupler transmitters operating at the frequency of the observe transmitter).
- The acquisition is preceded by a delay alfa+(1.0/(beta*fb)). Thus, the alfa and the filter group delays allow for a proper timing of the first data point (see "Considerations for the Delays Following the Last Pulse" on page 50).
- The number of points (values) acquired, and the dwell time used in the implicit acquisition are given by the two parameters np and sw (these C variables have the same name as the VNMR parameter from which they are initialized).

After the conditional branch for the implicit acquisition, the function (i.e., also for the case of an explicit acquisition) test4acquire resets any PFG or imaging gradient, adds the instruction for the housekeeping delay (see Section 18.5, "Housekeeping Delays," on page 214), and finally adds the code to jump back to the NSC (next scan) instruction in the Acode (see "The Instruction Section" on page 79).

The implicit acquisition through acquire is performed as a hardloop with a certain number of delays and with the CTC (command to convert) set to high. The number of dwell times per hardloop depends on the parameter and hardware configuration. With the pulse sequence controller board (2048 word loop FIFO), up to 1024 or 2048 dwell times can be coded per hardloop (depending on whether single- or double precision

time words are used). With a maximum of 32767 loop cycles, this allows for up to 64 or 128 million data points to be acquired, or half as much using the acquisition controller board (1024 word loop FIFO). From a parameter point of view, np is limited to multiples of 64. Residual dwell times are performed after the hardloop⁴. On systems with an output board (63-word FIFO), the acquisition is performed as hardloop with 16 or 32 dwell times (64 data points) per loop cycle, allowing for up to 1 or 2 million data points⁵.

The number of points has limitations other than those imposed by the looping capability of the pulse programmer. Primarily, the standard memory size on the HAL board limits the number of points. With the standard 2-Mbyte HAL memory, up to 1,048,576 data points can be added in single precision (16-bit) acquisitions or up to 524,288 points in double precision (32-bit) acquisitions. It is possible to expand the RAM on the HAL board: the MC-68000 CPU address space is 24 bits, or up to 16 Mbytes (some of the address space is used up by the acquisition CPU, ROM, and status registers), theoretically allowing for up to 4 or 8 million data points.

The STM board (see Chapter 7, "Digital Components," on page 65) uses a 24-bit counter to count the number of data points added up per transient. That number is compared with the number stored in the first long word (LC->np) in the LC structure of the Acode (see Section 8.2, "Looking at Acode," on page 69). If these numbers don't match, an error message "number of points acquired not equal to np" is produced. The STM counter can handle up to 16,777,216 points. This number is higher than those imposed by the memory and address space limitations—the STM board is not a limiting factor for the number of points.

The standard ADC for liquids NMR can handle data rates of up to 200,000 points per second (spectral windows of up to 100,000 Hz). The STM board can safely handle such data rates⁶. For larger spectral windows, a wideline receiver has to be used. The two wideline receiver/ADC models available can sample data at rates of up to 2 and 5 MHz (corresponding to 4 or 10 million values per second), much more than the STM board can possibly handle. Therefore, these boards are equipped with a fast on-board buffer memory, into which the data are acquired and from which the data are piped into the STM board (2-MHz digitizer) or transferred onto the HAL board (5-MHz digitizer).

The latest wideband (12-bit) ADC board allows for spectral windows of up to 5 MHz (as much as the pulse programmer can handle) and is equipped with 512 Kbyte of buffer memory, resulting in a maximum of 131,072 data points⁷. The older wideline ADC board allowed for spectral windows of up to 2 MHz and was equipped with 64 Kbyte of buffer memory, allowing a maximum of 16,384 points (8192 complex) only.

⁴ There is no technical reason why np could not be any multiple of two (even values below 64 should be permissible), but this has not been tested out. At the very least this would require changing the parameter limits for np.

⁵ With this configuration, all points must be acquired within the hardloop (no extra dwell times following the loop); therefore, the number of points acquired (np) *must* be a multiple of 64.

⁶ There is also a FIFO buffer between the ADC and the adder on the STM board.

⁷ This board has on-board STM functionality.

18.2 Explicit Acquisition

In some cases, it is necessary to code the acquisition explicitly because implicit acquisition lacks the desired functionality. These cases fall into three categories:

- Sequences where a pulse or other event should follow the acquisition, such as in the flipback experiment used sometimes in solids NMR, to recollect residual (spin-locked) proton magnetization after the acquisition time ("forced relaxation").
- Sequences that acquire multiple FIDs within the pulse sequence (nf > 1).
- Sequences that require performing pulses or other events in-between the acquisition of single data points⁸ (single-point acquisition), as used in sequences with explicit (synchronous) decoupling (see also the second example in Section 14.3, "Hardware Loops," on page 150), or for sequences with so-called multipulse line narrowing (see also the first example in the same section).

Typically, single-point acquisitions are coded using a hardloop, such as the following:

```
initval(np/2.0,v14);
starthardloop(v14);
    acquire(2.0,1.0/sw <length of other events in dwell time>);
    <additional events>
endhardloop();
```

Because two data points are acquired per loop cycle, the number of loop cycles is equal to half the number of data points. To obtain properly scaled and referenced spectra, it is essential that the time specified in the acquire statement, plus all the other events in the same dwell time (including hidden delays!), make up *exactly* the expected dwell time, 1/sw. In experiments with multipulse line narrowing, this is not so relevant because the scaling in these spectra is distorted anyway due to chemical shift scaling.

Of course, it is possible to acquire more data points per hardloop, for example, to allow changing the phase of the pulses between the data points:

```
initval(np/8.0, v14);
starthardloop(v14);
acquire(2.0,(1.0/sw-2.0*rof1-pw)/2.0);
rgpulse(pw,v1,rof1,rof1);
delay((1.0/sw-2.0*rof1-pw)/2.0);
acquire(2.0,(1.0/sw-2.0*rof1-pw)/2.0);
rgpulse(pw,v2,rof1,rof1);
delay((1.0/sw-2.0*rof1-pw)/2.0);
acquire(2.0,(1.0/sw-2.0*rof1-pw)/2.0);
rgpulse(pw,v3,rof1,rof1);
delay((1.0/sw-2.0*rof1-pw)/2.0);
acquire(2.0,(1.0/sw-2.0*rof1-pw)/2.0);
rgpulse(pw,v4,rof1,rof1);
delay((1.0/sw-2.0*rof1-pw)/2.0);
rgpulse(pw,v4,rof1,rof1);
delay((1.0/sw-2.0*rof1-pw)/2.0);
```

The number of time events in a hardloop must not exceed the loop FIFO size of the pulse programmer (otherwise FIFO errors would result at execution time). This also limits the number of points per acquire statement within a hardloop. Because nested hardloops are not possible, acquire in a hardloop codes a linear sequence of FIFO words and, in the best case, the number of points that can be collected in a single acquire statement within a hardloop is equal to the size of the loop FIFO for single

⁸ In practice, this means in-between acquiring *pairs* of data points.

precision timer words, or half that for double precision timer words. Only the second condition is checked at runtime (and aborts the *go* command with an error message if necessary). This is a rather loose test, because there are usually additional time events within the hardloop; otherwise, the entire FID could be collected in a single acquire statement without an explicit hardloop (see also Chapter 9, "Pulse Programmers," on page 85 and Section 14.3, "Hardware Loops," on page 150).

With this type of explicit acquisition (and irrespective of parameter limits), np is limited to multiples of the number of points acquired per hardloop, unless special provisions are taken in the pulse sequence to acquire any remaining points in separate acquire statements outside the hardloop. With explicit acquisition, the user must ensure that the correct (np) number of points is acquired with every transient. A mismatch in the number of points acquired leads to the error message "number of points acquired not equal to np".

An explicit acquisition can also be made part of a conditional part of the pulse sequence, as in the example of the following fragment of a proton flipback sequence:

```
rgpulse(pw,oph,rof1,0.0);
decphase(zero);
status(C);
if (dm[C] == `y')
{
   txphase(zero);
   delay(alfa+1.0/(beta*fb));
   acquire(np,1.0/getval("sw"));
   decrgpulse(pw,three,0.0,0.0);
   status(A);
}
```

In explicit acquisition, the user's responsibility is to ensure proper timing for the first data point and to avoid center glitches through rf leakage by resetting the relevant transmitter phases to zero (see also Section 18.1, "Implicit Acquisition," on page 205).

18.3 Multi-FID Sequences

}

VNMR has the ability to handle data with multiple FIDs per FID file trace (i.e., data from experiments where more than one FID was collected in a single pulse sequence transient). In such a case, all FIDs must be collected using explicit acquisition, as in the following (partial) example of a combined COSY and NOESY pulse sequence:

```
...
status(A);
hsdelay(d1);
status(B);
rgpulse(pw,v1,rof1,0.0);
delay(d2-rof1-4.0*pw/3.1416);
rgpulse(pw,v2,rof1,rof2);
status(C);
txphase(zero);
delay(alfa+1.0/(beta*fb));
acquire(np,1.0/sw);
hsdelay(mix-rof2-(alfa+1.0/(beta*fb))-(np/2.0)*1.0/sw-rof1);
rgpulse(pw,v3,rof1,rof2);
status(D);
```

```
txphase(zero);
delay(alfa+1.0/(beta*fb));
acquire(np,1.0/sw);
}
```

Of course, proper timing of the first data point must be ensured for *all* FIDs that are acquired, to avoid phasing and baseline problems throughout the experiment. Note that the restrictions in the number of points that were discussed in the previous section apply to the *sum* (np*nf) of all FIDs acquired in a pulse sequence, not just to the np parameter itself, i.e., failure to acquire nf FIDs per transient results in the (slightly misleading) error message "number of points acquired not equal to np". Other examples of multi-FID acquisitions are discussed in Chapter 21, "(Micro)Imaging Experiments," on page 239.

18.4 Receiver Phase Shifting

There are no fast lines from the pulse programmer that transmit the observe phase into some hardware, so how is receiver phase shifting done? How about small-angle receiver phase shifting? How are NMR signals detected?

Detection of NMR signals

All signals are measured relative to the reference frequency, which is constructed from the local oscillator (L.O.) and the intermediate frequency (I.F.), both of which are usually fixed in phase and frequency during an experiment.

As shown in Figure 22, the UNITY*plus* uses an I.F. of 10.5 MHz (i.e., the local oscillator is 10.5 MHz above the observe frequency). The signal from the probe is first amplified in the **preamplifier** (54 dB gain) and then passed through a first switchable




Figure 22. Detection of NMR signals

attenuator (6 and 12 dB, allowing for 0 to 18 dB attenuation, in steps of 6 dB). After that it is mixed with the local oscillator frequency in a double balanced **mixer**, resulting in a 10.5 MHz signal that is modulated with the observed signals (i.e., the difference between the transmitter frequency and the observed signals). All parts behind that first mixer operate at constant frequency, which simplifies the design of the receiver⁹.

After the mixer, the signal passes an **amplification and gain regulation** stage. Figure 22 simplifies this. The UNITY*plus* receiver contains three successive amplification stages (14 dB amplification each), each of which is preceded by switchable attenuators (14 dB each, one of them 2 + 4 + 8 dB) that allow setting the gain in steps of 2 dB between 0 and 42 dB¹⁰. In total, the UNITY*plus* preamplifier and receiver provide for a gain of 96 dB (54 + 3x14 dB) and a variable attenuation of up to 60 dB¹¹.

The amplified signal (δ + 10.5 MHz) is then split into two identical components that are fed into a pair of mixers. In these (double-balanced) **mixers**, the signals are then mixed with two 10.5 MHz I.F. components that are phase-shifted by 90 degrees against each other (they are both generated from a single, fixed 42 MHz frequency using a phase generator and frequency divider). The output from these mixers are two audio signal components that are phase-shifted against each other by 90 degrees, corresponding to the **real and imaginary components** of the audio signal.

After passing a pair of fixed filters, the two signals are then amplified to $10 V_{p-p}$ maximum¹² and finally fed into a pair of **programmable audio filters** (8-pole quasielliptical filters) that remove noise outside the spectral window (which otherwise would be folded into the observed spectral window).

These two audio signals are finally fed into the **ADC** (not shown in Figure 22), digitized and added to the current FID through the sum-to-memory (STM) board (see also Chapter 7, "Digital Components," on page 65). Because the transmitter frequency is located at the center of the observed spectral window, we need to detect both positive and negative frequencies. Sampling real and the imaginary signal components allows us to determine the sense of rotation of every signal component, and hence distinguish between positive and negative signal components. This is called *quadrature detection*.¹³

For liquids experiments (standard spectral windows of up to 100 kHz), UNITY and earlier systems used a similar scheme in their receivers, but mixed the 10.5 MHz signal (after the L.O. mixer) first with 10.0 MHz, resulting in a signal at 500 kHz. Mixing with a 500 kHz reference frequency then generated the audio signal. Reducing the number

⁹Wideband amplifiers—and rf devices in general—are much more demanding in their construction, and the receiver is certainly one of the most critical parts in the spectrometer.

¹⁰ It turns out that the "gain" in reality is controlling attenuators! We still call it gain, as the parameter runs "backwards" (to make it look like a gain parameter): 0 dB gain means maximum (18 + 14 + 14 + 8 + 4 + 2 dB) attenuation, 60 dB gain means no attenuation.

¹¹ The overall gain of the receiver chain can be varied between 36 and 96 dB.

¹² Any dc offset in the audio signals can be corrected at these amplifiers, and the gain of one of the amplifiers can be adjusted, allowing to accurately balance the two channels. Also, the relative phase shift of the two channels can be fine-adjusted down to fractions of a degree at the 42/10.5 MHz phase generator.

¹³ Errors in the relative phase of the two components, as well as any channel imbalance can lead to imperfect quadrature detection, resulting in quadrature ("mirror") images of all signals (partial folding of the signals around the center of the spectrum). If one of the two channels has zero amplitude, no quadrature detection is obtained, and the Fourier transform produces all signals both as positive and negative frequencies.


of mixing stages in the UNITY*plus* has improved performance, and the better distribution of amplification and gain attenuation at the receiver input has increased the dynamic range by avoiding overload with strong samples at intermediate levels.¹⁴

Quadrature Receiver Phase Shifts

Although dc offsets and the channel balance can be adjusted, there will always be slight differences between the two signal components. This can be observed in single-transient (or constant phase) spectra where often quadrature images and a center glitch can be observed¹⁵. Cycling the receiver through all four quadrature phases eliminates such artifacts, because both signal components have then passed through the two receiver and ADC channels (90-degree phase shifting), and dc offsets are cancelled out through phase alternation (180-degree phase shifting).

In most pulse sequences, the phase of the observed FID, as well as the phase of the receiver is therefore altered through phase cycles such as 0, 90, 180, 270, or 0, 180, 90, 270, or the like—the first example providing for fast quadrature image suppression, the second example offering faster dc offset (center glitch) cancellation¹⁶.

To better understand how receiver phase shifting works, let's look at a particular signal vector and its real and imaginary signal components at a particular point in time within the FID, during four successive transients with the phases 0, 90, 180 and 270 degrees.



If we added these four signals as shown, we would obtain exactly zero!

Proper signal addition (averaging) is easy to obtain. We simply need to properly route the two signal components before we add them to the real and imaginary parts of the stored FID. As shown in Table 12, to change the receiver phase, we should tell the STM board how to combine the stored FID with the incoming data: 90-degree receiver phase shifting is simple math and signal routing.

This is done implicitly with SETICM (set input card mode) instruction in the Acode, which tells the acquisition CPU to read the oph register from the LC structure and to transfer the observe phase information to the appropriate status register on the STM board. Such an instruction is placed ahead of the *first* acquisition instruction (be it implicit or explicit) or ahead of the first hardloop that contains an acquire statement.

¹⁴ Such overload could drive amplifiers in the receiver into saturation. At these levers they become non-linear, typically leading to intermodulation distortions in the spectrum (non-linear devices act as mixers), and typically adding mixing products of dominant NMR signals to the spectrum (ghost peaks at the distance of strong signals). This becomes most evident if more than one strong signal is present in the spectrum.

¹⁵ Under normal circumstances, these artifacts should be below 0.5% of the main signals.

¹⁶ Many sequences only do a partial receiver phase cycle (quad image suppression through 0^o and 90^o phase shifts), mainly in experiments where a residual solvent (water) signal would cover eventual center glitches (from dc offsets) anyway.

| Signal phase | Signal component | Operation |
|-----------------|------------------|------------------------------|
| 00 | 00 | add to real part |
| 0 | 90° | add to imaginary part |
| 90 ^o | 00 | subtract from imaginary part |
| | 90° | add to real part |
| 180° | 00 | subtract from real part |
| 100 | 90° | subtract from imaginary part |
| 270° | 00 | add to imaginary part |
| 270 | 90° | subtract from real part |

Table 12. Dealing with real and imaginary signal components

Note that because the SETICM instruction is placed implicitly and automatically with the first acquire statement only, it is *impossible* to change the "receiver phase" after the first acquire (oph can, of course, be changed, but this would have no effect on the way the signals are co-added. In other words, all data points within an FID are acquired with the same "receiver phase." Even in multi-FID experiments (see Section 18.3, "Multi-FID Sequences," on page 209), all FIDs are acquired with the *same receiver phase*.

Small Angle Receiver Phase Shifting

As the receiver phase can be shifted in increments of 90 degrees through mathematical operations, one could certainly imagine a more complex STM board that shifts the incoming data by *any* phase angle, using the following formula set:

 $R_i = r_i \cos \Theta + i_i \sin \Theta$ $I_i = i_i \cos \Theta - r_i \sin \Theta$

Unfortunately, the STM board does not provide such a functionality, hence proper small-angle receiver phase shifting is not possible. There are, however, a few possible workarounds for this problem:

- If *all* transients should be shifted by the *same* small angle phase (apart from quadrature phase shifts), the above transformation can be performed after the acquisition on the final FID file on the disk, either through an external program¹⁷ or from within VNMR by going through an ASCII file (writefid command), using a macro to reforming these data and creating a modified ASCII file, which can the be used to re-build a phase-shifted binary FID using the makefid command.
- On systems with older (type a or b) rf generation on the observe channel (some VXR and earlier systems only), the phaseshift statement not only shifts the transmitter phase, but it changes the phase of the entire reference frame by shifting

¹⁷ A C program with this functionality is supplied as part of the user library.

the phase of the local oscillator using the so-called "phase-pulse technique"¹⁸ (see the section "Pulse Sequence Statements: Spectrometer Control" in the manual *VNMR User Programming*). In fact, there are a few systems that can do smallangle receiver phase shifting! This technique may not work with newer rf generation types, because not all decades in the PTS frequency synthesizer switch frequencies in a phase-coherent way. In fact, the phaseshift statement is programmed to only work through the offset synthesizer used in older systems.

• Instead of shifting the receiver phase by a small angle, it is of course possible (and almost equivalent) to *shift the phase of all pulses* (at least those on the observe channel) by the same amount (using the xmtrphase statement), but *in the other direction*. Note that it is still advisable to do some quadrature receiver phase shifting to compensate for an eventual channel imbalance, quadrature phase error, or dc offset. It is a bad idea to keep the receiver phase constant for all scans while just shifting the phase of all pulses instead.

It turns out that from all pulse sequences and experiments, there are only a few rare examples that require small-angle receiver phase shifting at all, and for those, the solutions presented above are totally sufficient.

18.5 Housekeeping Delays

Before jumping back to the beginning of the pulse sequence, every pulse sequence performs an instruction 97 (HouseKEEPing), which introduces an extra delay for a variety of clean-up tasks that need to be performed before the next scan can be started:

- After every transient, the ct counter is incremented and compared with nt (to check whether the experiment is finished). Further, a check is made whether an sa (stop acquisition) signal has been received, and the number count from the STM board is checked to ensure the proper number of data points was acquired. The duration of this delay has changed over the different instrument and software generations: it currently is around 14 milliseconds. The only way to avoid this delay is to program a multi-FID pulse sequence (see Section 18.3, "Multi-FID Sequences," on page 209).
- After the first FID that acquires real data, the first 256 points of the FID are checked for ADC overflow, in cases where a fixed gain was used. For this task, the housekeeping delay is lengthened by X milliseconds. Autogain (gain='n') totally disrupts the timing after the first transient.
- Additional small delays occur at the beginning of the first "real" transient (i.e., after the steady-state transients) as well as with the start-up of every increment. The start-up delay is used to set up all the hardware (see also Section 8.2, "Looking at Acode," on page 69).
- If a diagnostics terminal is connected to the acquisition CPU and the bootup selector switch is set to a non-zero position (even if just the switch setting is the case), long housekeeping delays (around 0.2 seconds or more) occur, because the information for the diagnostics terminal needs to be prepared and sent. This was discussed in Section 13.2, "Diagnostics and Error Output," on page 146.

¹⁸ This method temporarily (for a few microseconds) changes the offset on the (observe) channel, which causes the reference frame to change its orientation.

Chapter 19. Multidimensional Experiments

It cannot be the purpose of this manual to explain the mechanistics of *n*D NMR experiments; only programming issues will be discussed here. Also, apart from the names of the parameters that are used, all indirect dimensions (t_1 up to t_{n-1} in *n*D experiments, where *n* currently is up to 4) are treated the same way; therefore, only the 2D case will really be discussed here. The following tables should allow transferring the information in this chapter to higher dimensions in *n*D experiments (all the parameters in the 2D column will be discussed in detail below):

| Parameter description | 2D | 3D | 4D |
|---|-------|--------|--------|
| Number of tn increments | ni | ni2 | ni3 |
| Spectral width in fn | swl | sw2 | sw3 |
| fn coherence selection mode parameter (optional) | phase | phase2 | phase3 |
| Flag for fn axial peak displacement (proposed/optional) | fad | fad2 | fad3 |
| Flag for inverting folded peaks in fn (proposed/optional) | f1180 | f2180 | £3180 |

Table 13. VNMR acquisition parameters used for *n*D experiments

 Table 14.
 Variables used in *n*D pulse sequences

| Variable description | 2D | 3D | 4D |
|---|--------|--------|--------|
| Evolution delay (double) | d2 | d3 | d4 |
| Spectral width in fn (double) | swl | sw2 | sw3 |
| Real-time index for evolution in tn (codeint) | id2 | id3 | id4 |
| fn coherence selection mode (integer) | phasel | phase2 | phase3 |
| Evolution time increment in tn (double) | inc2D | inc3D | inc4D |

19.1 Indirect Time Domain Incrementation

One of the prominent features of VNMR *n*D pulse sequences is the absence of any explicit coding for the looping and the evolution delay incrementation. Similarly, any arrayed experiment does not complicate the pulse sequence: the pulsesequence() function is simply called once per array element. In the case of *n*D experiments, we don't even define an array. The presence of the ni and sw1 parameters alone causes the pulse sequence software to set up an implicit array of ni elements on the parameter d2. That array is *not* shown in VNMR, but the arraydim parameter does reflect the extra dimension in that it multiplies the number of traces from explicit arrays with ni, ni2 and ni3.

Note that the implicit array is starting with the current value of d2 $(d3, d4)^{1}$, the array increment is 1/sw1 (1/sw2, 1/sw3):

```
d2<sub>implicit</sub>=d2, d2+1/sw1, ... d2+(ni-1)/sw1
d3<sub>implicit</sub>=d3, d3+1/sw2, ... d3+(ni2-1)/sw2
d4<sub>implicit</sub>=d4, d4+1/sw3, ... d4+(ni3-1)/sw3
```

For a basic *n*D experiment, all we need to do in terms of defining the sequence of pulse sequence events is, to include the evolution delay (d2 for 2D, d2 and d3 for 3D, etc.). Some sequences will use a refocusing pulse in the middle of an evolution time, which can be easily achieved with a construct like

```
delay(d2/2,0 - rof1);
rgpulse(2.0*pw90,v1,rof1,0.0);
delay(d2/2.0);
```

Other sequences might use what is sometimes called a "fixed evolution time", with a pulse that moves within that fixed delay. Also this can be realized easily:

```
pulse(pw90,v1,rof1,0.0);
delay(d2-rof1);
pulse(2.0*pw90,v2,rof1,0.0);
delay(tau-d2-rof1);
pulse(pw90,v3,rof1,0.0);
```

In order to obtain flat baselines in phase-sensitive experiments, it is furthermore important to not only compensate for any delays around the adjacent pulses but also for the precession during these pulses. The correction term is -2.0*pw90/3.14159 *per* 90 degrees pulse of length pw90 adjacent to the evolution time. Hence, for pulse sequences like COSY or NOESY the evolution period should be coded as follows:

```
pulse(pw90,v1,rof1,0.0);
if (d2 - 4.0*pw90/3.14159 - rof1 > 0.0)
    delay(d2-4.0*pw90/3.14159-rof1);
pulse(pw90,v2,rof1,0.0);
```

The *if* statement serves to suppress the error message from the fact that for the first increment the calculated delay is negative. This construct actually implies that the spacing between the first two increments is not the same as between all the other increments: all traces except for the first one (where theoretically the two pulses would have to overlap) are measured correctly, and this means that there will be a (minor) error in the first data point in the indirect dimension. Fortunately, this error can easily be corrected either by applying a dc offset correction after the Fourier transformation, or by reconstructing the first data point using linear prediction.

No correction for spin precession is required for TOCSY-type spin locking, and for ROESY spinlocks it seems easier to use an empirical correction term that is adjusted experimentally (such that no first order phase correction is required in the indirect dimension).

¹ This is a useful feature because it allows (re-) acquiring *any* part of an *n*D experiment with little effort, but it also is dangerous in that d2 (d3, d4) can inadvertently be set to some (possibly large) non-zero value, which can have serious consequences for the result of the experiment (like causing strong first-order phase shifts or the observation of noise only). Therefore, as of VNMR 5.1, it is a good idea to have the go macro issue a warning if *n*D experiments are started with non-zero evolution delays.

19.2 nD Quadrature Detection

This is not the place for an exhaustive discussion of methods for achieving quadrature detection in nD experiments; however, the most important methods will be presented here as a guideline for the implementation of new pulse sequences.

Absolute Value nD Experiments

F1 quadrature in 2D experiments is achieved by either co-adding two data sets with the "phase-relevant" pulses prior to the evolution period shifted by 90 degrees, or by independently incrementing (or decrementing²) the phase of these pulses (mostly just one pulse) in steps of 90 degrees. The first method only requires half the number of scans and is completely sufficient. This is a major cancellation step (the N+P-type spectrum is added or subtracted from the N-P type spectrum, resulting in either an N-or a P-type 2D spectrum. In other words, half the 2D spectrum—the anti-diagonal and the associated crosspeaks in homonuclear correlation spectra—is cancelled in this step) and should be one of the faster steps in the overall phase cycling. Certainly it is advisable to perform this step *before* CYCLOPS.³

One peculiar aspect of this type of absolute-value nD experiment is that when we do the f₁ coherence selection, we seem not to accumulate signal-to-noise ratio (as at that point we are in fact cancelling half the signal). This is most obvious when the phasecycling (or the cancellation efficiency) is checked in a 1D array using

ni=1 nt=1,2,4,8,16,32

where we expect a 2signal-to-noise improvement with every step. For example, if this is a simple double-quantum filtered COSY experiment where we first do a four-step double-quantum coherence selection (nt=4), followed by a two-step f₁ quadrature selection phase cycle (nt=8), followed by a four-step CYCLOPS phase cycle (nt=32), we would observe relative signal-to-noise ratios of 1, $\sqrt{2}$, 2, 2, $\sqrt{2}$ *2, 4 in the 1D trace, instead of 1, $\sqrt{2}$, 2, 2, $\sqrt{2}$ *2, 4, $\sqrt{2}$ *4 (as for standard 1D experiments). This at the same time is an easy test that indicates the minimum number of transients per increment for achieving quadrature detection in f₁⁴.

Phase-Sensitive nD Experiments: States/Haberkorn/Ruben

The above method does not allow phasing the *n*D spectrum, because the real and imaginary parts of the spectrum are not separated in the indirect frequency domain; hence, the display in absolute-value mode. Various methods have been proposed and used to separate the real and imaginary (i.e., absorption and dispersion) parts of a *n*D spectrum. The most "natural" method on Varian instruments seems to be the technique proposed by States et al.⁵, also called "hypercomplex" mode, which involves acquiring

² The sign of the phase incrementation determines whether P- or N-type spectra are obtained. Both can be processed by VNMR such as to give a "normal" presentation (using the 'ptype' argument to the ft2d command to obtain "normal" orientation for P-type spectra.

³ As the short-term fluctuations in an instrument (or its environment) are usually smaller than the long-term variations, it is advisable to perform the major cancellation steps (line multiple quantum filtering, f_1 quadrature) before performing steps like CYCLOPS for cleaning up minor artifacts.

⁴ Note that for multiple-quantum-filtered experiments, the evolution delay (d2 for 2D) must be set to a non-zero value in order to see a double-quantum-filtered signal in the 1D trace; do *not* forget to reset that delay to 0 before starting the real experiment!

⁵ D.J. States, R.A. Haberkorn & D.J. Ruben, J. Magn. Reson. 48, 286 (1982).

two separate data sets where in the second data set the phase of the (phase-relevant) pulses prior to the evolution time is incremented by 90 degrees compared to the first set. In VNMR we define an array using a parameter phase, which is set to the values 1 and 2 (phase=1, 2). This parameter is available within pulse sequences in the integer variable phase1 (note the difference in the names!)⁶. In the pulse sequence, we then can use a simple construct such as

if (phase1 == 2) incr(v1);

where v1 is the phase of the first pulse⁷. This method of achieving f_1 quadrature gives the maximum in processing flexibility. It falls in line with the simultaneous sampling in f_2 , and, compared to TPPI, it usually gives better baseline flatness. There is also a disadvantage in that any axial peaks (artifacts that often cannot be avoided, in particular with biomolecular NMR spectra) show up in the center of the spectrum, whereas with TPPI (see below) they are moved to the edge of the spectrum. This can be fixed, however, by combining the hypercomplex method with FAD (fn axial peak displacement, described below).

Note that the phase array is performed *before* incrementing d2 (i.e., the *implicit* array is array='d2, phase'), which is in line with the requirement that scans or data that are to be subtracted from each other for cancellation should be measured as close to each other in time as possible, to make the experiment less susceptible to environmental variations and give better cancellation. The same holds true for 3D and 4D experiments, where the *implicit* arrays usually are

```
array='d2,d3,phase,phase2' (for 3D), or
array='d2,d3,d4,phase,phase2,phase3' (for 4D).
```

We could argue that this method requires twice (3D) and four times (4D) as long to make the first complete plane available for viewing, but this method will definitely give better cancellation than

```
array='d3, phase2, d2, phase' (for 3D), or
array='d4, phase3, d3, phase2, d2, phase' (for 4D).
```

But of course it is entirely possible to set up such arrays explicitly from within VNMR. For example, for 3D experiments:

```
array('d2',ni,0,1/sw1)
array('d3',ni2,0,1/sw2)
ni=0 ni2=0
array='d3,phase2,d2,phase'
```

The real problem with this method (apart from giving bad cancellation efficiency) is, that such data *cannot* be processed using VNMR (unless we write extra software to rearrange the FIDs).

⁶In earlier VNMR releases this parameter had to be fetched from the parameter table using constructs such as int phase1 = (int)(getval("phase") + 0.5);

⁷ Note that early releases of VNMR also incremented the receiver phase. This results in data sets that require using a different selection of coefficient arguments with the wft2d command, like wft2d (1,0,0,0,0,0,0,-1) (i.e., typically the last two arguments need to be exchanged compared with the current "standard"). This early mode also had the disadvantage that the second data set was 90 degrees out of phase in f_2 compared to the first set.

Axial Peak Displacement (FAD)

It turns out that there is a trick to move the axial artifacts in hypercomplex nD spectra to the edge of the spectrum by inverting the phase of the (phase-relevant) pulses prior to the evolution time and the receiver phase with every even time increment, both for phase=1 and phase=2. In the vast majority of the cases this method will be used by default for hypercomplex experiments, but sometimes is may be desirable to have this extra phase inversion under flag control. For these cases we would propose the flag names fad1 (or fad), fad2, and fad3. If we disregard the flag control, the pulse sequence construct for a hypercomplex experiments will be a bit more complex than the one shown in the previous section:

```
if ((phase1 == 1) || (phase1 == 2))
{
    dbl(id2,v13);
    add(v1,v13,v1);
    add(oph,v13,oph);
    if (phase1 == 2) incr(v1);
}
```

where v1 is the phase of the (phase-relevant) pulse(s) prior to the evolution time, and id2 is a real-time variable that contains the number of evolution time increments $(0, 1, 2, ..., ni-1)^8$.

Phase-Sensitive nD Experiments: TPPI

On instruments with sequential sampling, quadrature detection is achieved by the receiver phase with every data point. An equivalent method can be applied to *n*D spectroscopy by continuously incrementing the phase of the (phase-relevant) pulses prior to the evolution time by 90 degrees with every increment. This is also called *TPPI method*, or *time proportional phase incrementation*.⁹ This separated the absorption and dispersion parts of the spectrum in a single data set, but it requires acquiring twice the spectral window and twice the number of increments in the indirect dimension to achieve the same digital and spectral resolution. The transformed spectrum then consists of two parts: one with the absorption lines and one with the dispersive contribution, separated by the axial peak (artifacts).¹⁰ Because TPPI leads to a single data set, the processing seems somewhat simpler, but it apparently is more difficult to achieve good baseline flatness.

```
int tl_counter =(int)(d2 * getval("swl") + 0.1);
(...)
if ((phase1 == 1)||(phase1 == 2))
{
    initval((double)(2*(tl_counter%2)), v13);
    add(v1,v13,v1);
    add(oph,v13,oph);
    if (phase1 == 2) incr(v1);
}
```

⁹ D. Marion & K. Wüthrich, Biochem. Biophys. Res. Commun. 113, 967 (1983).

¹⁰ The VNMR support for processing TPPI spectra is somewhat limited, in that it is not possible to discard the dispersion part of the spectrum, which means that the entire data matrix must be carried along. This can be a problem, especially with 3D processing where the transformed data matrix may require as much as four times the disk space compared to the hypercomplex method in which the imaginary parts of the spectrum can be discarded.

⁸ In earlier VNMR releases, the "increment counter" id2 did not exist and constructs like the following were used instead:

The VNMR convention is that TPPI is done with phase (phase2, phase3) set to 3, rather than 1 and 2. In the pulse sequence, we can use a simple construct like

```
if (phase1 == 3) add(v1,id2,v1);
```

where v1 is the phase of the (phase-relevant) pulse(s) prior to the evolution time, and id2 is a real-time variable that contains the number of evolution time increments $(0, 1, 2, ..., ni-1)^{11}$.

The main advantage of the TPPI method is that only four coefficients are required for the ft2d command (which accepts up to 32 coefficients). If we want to measure a genuine array of phase-sensitive 2D spectra (e.g.: an array of mixing times with NOESY or TOCSY) using the hypercomplex method, we would end up with a double array. This would lead to a long and complex set of coefficients for the ft2d command, and the number of array elements for the "genuine" array (e.g., mix) would be limited to 4. With TPPI, up to 8 array elements (e.g., 8 different mixing times) can be measured, and in addition to that there is a macro wft2dac that makes it easy to select and process individual data sets from the array (and this macro does not cope with arrayed hypercomplex data).

Phase-Sensitive nD Experiments: Arrayed TPPI

The TPPI implementation discussed above has the disadvantage of leading to frequency doubling in f_1 ; there is an alternative that doesn't have that disadvantage: *arrayed TPPI*. For this mode, the parameter convention phase=1, 4 has been used. The method involves acquiring two data sets, with the (phase-relevant) pulses, prior to the evolution time, shifted by 90 degrees in the second data set, the same as in the hypercomplex method, but in addition to that, for the second data set the evolution time is increased by half an increment:

```
if (phase1 == 4)
{
    incr(v1);
    d2 += inc2D/2.0;
}
```

As already mentioned, this method does not lead to frequency doubling in f_1 , but it requires processing the data set with real (instead of complex) FT in f_1 (procl='rft'). This again has some disadvantages in that some processing options, like linear prediction or lsfrq, are not implemented for real FT; therefore, this method can't really be recommended. In addition to that, most current sequences add FAD for the case of phase=1, 2 (hypercomplex method, see above), which conflicts with the definition of phase=1, 4 for arrayed TPPI (where FAD would unnecessarily be added for phase=1, but not for the traces with phase=4). If arrayed TPPI still is to be used, it would probably be better to change the convention to phase=4, 5 for this method, or

```
int t1_counter = (int)(d2 * getval("swl") + 0.1);
(...)
if (phasel == 3)
{
    initval((double) t1_counter),v13);
    add(v1,v13,v1);
}
```

¹¹ In earlier VNMR releases, the "increment counter" id2 did not exist and had to be created using an standard real-time variable in constructs like

alternatively, to have FAD under flag control (the former would be preferable, as the flag solution would still allow for a parameter mis-setting).

Folding in Indirect Dimensions

In hypercomplex phase-sensitive *n*D spectra, folded peaks (in an indirect domain, where no audio filter is involved) are in phase with the "normal" signals, provided the evolution time is corrected for the precession of the spins during adjacent pulses, as shown in Section 19.1, "Indirect Time Domain Incrementation," on page 215ff. In 3D and 4D experiments, the number on increments that can possibly be performed in any indirect dimension is very limited, and yet for many biomolecular 3D and 4D experiments, rather large homonuclear and heteronuclear shift ranges must be covered, which imposes severe restrictions on the achievable digital resolution.

In this situation, we would welcome using folding in one dimension, because this allows doubling the digital resolution with the same number of increments—or acquiring half as many increments only in one dimension (which also cuts the overall experiment time in half). Of course, this only works if there is no significant overlap between folded and "real" peaks, and if there is an easy and safe way of distinguishing between folded and "real" signals.

The method of choice to distinguish between the two groups of signals is to increase all t_n evolution delays by half an t_n increment, 0.5/sw1 (or 0.5/sw2, 0.5/sw3). This is the equivalent to increasing alfa by 0.5/sw for the observe dimension and leads to a 180 degrees first-order phase shift in the corresponding frequency domain, which inverts all signals that are folded in once in that domain (remember that there is no attenuation due to audio filters in the indirect domains!). If there was any first point distortion, this would lead to a half-sinoidal baseline distortion in that domain. Fortunately, there are no filter distortions in these domains, and with that addition to the evolution time we may actually be able to even sample the first data point correctly (assuming the corrections for rof1 and for the precession during adjacent pulses is less than half an evolution time increment).

There are various ways in which this could be implemented. A "minimalist approach" would be to simply set $d_{2=0.5/sw1}$ in the VNMR parameter set. We don't think this is an optimal solution, because (in order to avoid operator errors) we would rather prefer to see the evolution time delays to be set to zero for all *n*D experiments, with the exception of the case where either we want to test signal levels somewhere *within* the evolution delay or where we want to *re*acquire specific traces from an *n*D experiment (the latter case would also be made more difficult with this approach). Also, once the parameter is entered, it's value doesn't have an obvious meaning to the user, and so later it is not clear what was intended with that parameter setting.

Still we would like this feature to be under parameter control, because it certainly is not something that should be applied to all phase-sensitive *n*D experiments (in particular, not in homonuclear correlation experiments on biomolecular samples where baseline flatness is a primary requirement).

Several existing *n*D sequences use flags named f1180, f2180, and f3180 to control this feature in f_1 , f_2 , and f_3 . In the pulse sequence this can then be used as follows:

```
double d2incr;
char f1180[MAXSTR];
getstr("f1180",f1180);
if (f1180[1] == 'y')
    d2incr = inc2D/2.0;
else
    d2incr = 0.0;
(...)
pulse(pw90,v1,rof1,0.0);
if (d2 + d2incr - 4.0*pw90/3.14159 - rof1 > 0.0)
    delay(d2+d2incr-4.0*pw90/3.14159-rof1);
pulse(pw90,v2,rof1,0.0);
(...)
```

Combined Implementations

The majority of nD pulse sequences have the hypercomplex method (with FAD) and (standard) TPPI built in. This can be done with the following construct:

```
if ((phase1 == 1) || (phase1 == 2))
{
    dbl(id2,v13);
    add(v1,v13,v1);
    add(oph,v13,oph);
    if (phase1 == 2) incr(v1);
}
else if (phase1 == 3) add(v1,id2,v1);
```

where once more v1 is the phase of the (phase-relevant) pulse(s) prior to the evolution time, and id2 is a real-time variable that contains the number of evolution time increments (0, 1, 2, ..., ni-1). This should be more than sufficient for most phase-sensitive *n*D experiments.

There is a number of pulse sequences around (also in /vnmr/psglib) that combines all of the above *and* absolute-value acquisition in a single pulse sequence. This is certainly achievable, but it turns out to be somewhat non-trivial, because it requires combining two different phase cycles in one sequence. In order to do things in an optimum way, the f_1 quadrature phase incrementation isn't simply added at the end of the phase cycle, but instead it is *inserted* into the phase cycling sequence (see also "Absolute Value nD Experiments" on page 217). Where this is done with real-time calculations, it makes the phase cycling difficult to decode and understand, because for most of the calculation steps we have to consider the absolute value and phasesensitive cases separately, as shown in the following (streamlined) example from a NOESY pulse sequence:

```
/* CALCULATE PHASECYCLE */
sub(ct,ssctr,v12); /* ctss */
mod2(v12,v1); /* ct/s */
hlv(v12,v3); /* ct/2 */
dbl(v1,v1); /* 02 */
hlv(v3,v10); /* ct/4 */
hlv(v10,v10); /* ct/8 */
if (phase == 0)
hlv(v10,v10);
if (phase == 0)
{
   assign(v10,v9); /* ct/8 */
hlv(v10,v10); /* ct/16 *
mod2(v0,v10); /* t0118 *
                        /* ct/16 */
   mod2(v9,v9);
                        /* [01]8 */
}
else assign(zero,v9);
hlv(v10,v2);
                        /* ct/32 av
                            ct/16 ph */
                         /* [02]32 av
dbl(v2,v2);
                            [02]16 ph */
add(v9,v1,v1);
                        /* (02)4 (13)4 av
                             02
                                         ph */
mod2(v10,v10);
                        /* [01]16 av
                            [01]8 ph */
add(v2,v1,oph);
                        /* (02)4 (13)4 (20)4 (31)4 av
                            (02)8 (20)8 ph */
add(v3,oph,oph);
                        /* 0213203113203102 2031021331021320 av
                            0213203102132031 2031021320310213 ph*/
add(v10,oph,oph);
                       /* 0213203113203102 3102132002132031 av
                            0213203113203102 2031021331021320 ph*/
add(v10,v2,v2);
                         /* [0123]16 av
                            [0123]8 ph */
add(v10,v1,v1); /* (02)4 (13)4 (13)4 (20)4 av
                                                      ph */
                            (02)4 (13)4
add(v10,v3,v3);
                         /* [01230123 12301230]2 av
                            [01231230]2 ph */
if ((phase == 1) || (phase == 2)) /* hypercomplex + FAD */
{
  dbl(id2,v11);
  add(v1,v11,v1);
  add(oph,v11,oph);
  if (phase == 2) incr(v1);
}
else if (phase==3) add(v1,id2,v1);
                                       /* TPPI */
```

Needless to say, in the original sequence the phase cycle is explained verbally (at least!) as follows: "The first 90-degree pulse is cycled first to suppress axial peaks. This requires a two-step phase cycle consisting of $(v1, 0 \ 2)$. The third 90-degree pulse is cycled next using a four-step phase cycle (v3) designed to select both longitudinal magnetization, J-ordered states, and zero-quantum coherence (ZQC) during the mixing period. If the experiment is to collect data requiring an absolute-value display, the first pulse is next incremented by 1 to achieve f_1 quadrature. If the data are to be presented in a phase-sensitive manner, this step is not done. Next, the second 90-degree pulse is cycled to suppress axial peaks (v2). Finally, all pulse and receiver phases are incremented by 90 degrees (v10) to achieve quadrature image suppression due to receiver channel imbalance."

With phase tables, the same phase cycling can be achieved in a somewhat easier way, but still it requires altering the "division return factor" for those tables that define phase

cycles that have lower priority than the f_1 quadrature phase incrementation and are "slowed down" for the case of an absolute-value experiment.

This is the phase table definition for the same pulse sequence:

```
t1 = 0 2  /* 1st pulse */
t2 = { 0 2 }16  /* 2nd pulse (phase=0: divn_return=32) */
t3 = { 0 1 2 3 }2  /* 3rd pulse */
    /* calculate oph=t1+t2+t3 in sequence */
t4 = { 0 1 }8  /* CYCLOPS (phase=0: divn_return=16) */
t5 = { 0 1 }8  /* f1 quadrature for phase=0 */
```

In this table, definition the phase cycling elements have been separated in order to achieve simple phase tables. The division return factor for the tables ± 2 and ± 4 is altered within the pulse sequence, before the phases are combined:

```
sub(ct,ssctr,v10);
if (phase1 == 0)
{
 setdivnfactor(t2,32); /* modify tables for phase=0 */
 setdivnfactor(t4,16);
  getelem(t5,v10,v5);
}
else assign(zero,v5);
getelem(t1,v10,v1);
                        /* extract phases from tables */
getelem(t2,v10,v2);
getelem(t3,v10,v3);
                        /* f1 quadrature (phase=0) */
add(v1,v5,v1);
add(v1,v2,oph);
                        /* calculate oph */
add(v3,oph,oph);
                        /* add CYCLOPS */
getelem(t4,v10,v4);
add(oph,v4,oph);
add(v1,v4,v1);
add(v2,v4,v2);
add(v3,v4,v3);
if ((phase == 1) || (phase == 2)) /* hypercomplex + FAD */
{
 dbl(id2,v11);
 add(v1,v11,v1);
  add(oph,v11,oph);
 if (phase == 2) incr(v1);
}
else if (phase==3) add(v1,id2,v1); /* TPPI */
```

But still: combining absolute-value and phase-sensitive experiments in the same pulse sequence *complicates the phase cycling setup*—just to incorporate a rarely used option! Therefore, it is strongly recommended to instead write two different pulse sequences, where an absolute-value option really makes sense, or at the very least, we could also think of supplying two separate phase tables with one sequence:

Now a single pulse sequence could use these two files to implement both absolute value and phase-sensitive NOESY phase cycles:

```
sub(ct,ssctr,v10);
getelem(t1,v10,v1);
                                  /* extract phases from tables */
getelem(t2,v10,v2);
getelem(t3,v10,v3);
                                   /* CYCLOPS */
getelem(t4,v10,v4);
                                   /* f1 quadrature (phase=0) */
getelem(t5,v10,v5);
add(v1,v5,v1);
                                   /* f1 quadrature (phase=0) */
add(v1,v2,oph);
                                   /* calculate oph */
add(v3,oph,oph);
                                   /* add CYCLOPS */
add(oph,v4,oph);
add(v1,v4,v1);
add(v2,v4,v2);
add(v3,v4,v3);
if ((phase == 1) || (phase == 2)) /* hypercomplex + FAD */
{
  dbl(id2,v11);
  add(v1,v11,v1);
  add(oph,v11,oph);
  if (phase == 2) incr(v1);
}
else if (phase==3) add(v1,id2,v1); /* TPPI */
```

Of course, for the case of two different pulse sequences, the individual sequences could be further simplified in that the *if* statement that deals with phase-sensitive experiments could be left away for the absolute-value case, etc.

Coherence Selection through Gradients

Pulsed field gradients can be used in a number of ways in *n*D NMR. For just scrambling transverse magnetization (such as in some gradient-NOESY pulse sequences), up to multiple-quantum filtering and coherence selection in general. In the former case, there will still be a phase cycling section in the sequence, using the above mechanisms for f_1 coherence selection and (partial) artifact suppression. In other cases, the coherence selection will be done using gradients, which often dramatically simplifies the phase cycling (sometimes no phase cycling is used at all).

Using pulsed field gradients for coherence selection has the important advantage of not requiring subtraction (cancellation) through phase cycling; therefore, these experiments are much less susceptible to environmental variations (which otherwise often lead to bad cancellation). From a programming point-of-view, there is very little to add here that hasn't been discussed already.

Chapter 19. Multidimensional Experiments

Chapter 20. Solid-State NMR Experiments

Apart from the rotor synchronization feature, which consists of dedicated hardware and software for specific solid-state experiments, there is mostly a gradual difference between standard liquids and solid-state experiments. As far as software and the execution of pulse sequences are concerned, the typical spectral window in solids experiments is much larger and, because the signals themselves are often extremely wide, the pulse sequence timing becomes much more of an issue. Extra delays of even a few microseconds only can cause a severe loss of coherence, or will at the very least lead to severe phasing problems in the final spectrum.

We can not cover the entire area of solids NMR spectroscopy. In this chapter, just a few typical highlights are picked and discussed.

20.1 Cross-Polarization MAS Experiments

AP Bus Events in CP/MAS Experiments

The simplest CP/MAS experiment consists of a 90-degree pulse on protons, followed by the cross-polarization period, during which both protons and the X-nucleus are spinlocked with the same rf field (Hartmann-Hahn condition). After that, the decoupler is switched to full power to remove the dipolar line broadening during the acquisition. The following code is simplified and written for a UNITY*plus*, with AP bus control of the linear modulator:

```
status(A)
decpwrf(getval("crossp"));
delay(d1);
decrgpulse(pw,t1,rof1,0.0);
status(C);
decphase(zero);
rgpulse(getval("contact"),t2,0.0,0.0);
decpwrf(getval("dipolar"));
delay(rof2); rcvron();
```

status (C) turns the decoupler on. After the contact time, the decoupler stays on and is switched to the higher "dipolar" level using the pwrf statement, which takes 4.6 μ sec on a UNITY*plus*. It is possible and acceptable to switch the linear modulator while the rf is turned on for that channel (a 4.6 μ sec gap would be unacceptable anyway)—after all, that's what is happening all the time during a shaped pulse.¹

¹ Note that for these applications it is not desirable to switch the 63 or 79 dB attenuator with the rf turned on (even though that would take less AP bus cycles). On older (UNITY) systems, we have seen transition phenomena during the switching: there is no guarantee that all the bits of these attenuators switch at exactly the same time. On such older systems, we have also observed full power coming through for a very short period (there may also be a very small gap with less power). Even though that period lasted only a few nsec, it would be enough to make a probe arc with the 1-kW amplifiers switched on. (In liquids experiments there is much less power involved, and there should be less of a problem—after all, people have been creating soft shaped pulses using these attenuators.)

If we took an oscilloscope to see what really happens during this experiment, we would find that after the pulse, during the first three AP bus words (out of four), the power stays at the same level, because the new power level is being latched (all the bits are pre-stored and then enabled at the same time, with the last AP bus word). On a UNITY *plus*, the actual switching was found to happen 350 nsec into the last AP bus period. The first 150 nsec are the pulse programmer overhead, then it seems to take another 200 nsec for the AP bus chip to read the information and set the hardware. So, the power switching occurs 3*1.15 + 0.15 + 0.2 = 3.8 µsec after the transmitter pulse (while the entire decpwrf statement takes 4.6 µsec). For the standard CP/MAS experiment, this means that the full dipolar decoupling will occur with a slight delay that may influence the first 1 or 2 data points, but overall, it still seems acceptable.

There are experiments such as CPCOSY or CPNOESY, however, that require that the protons are spinlocked at high (dipolar decoupling) power *immediately* after the Hartmann-Hahn polarization transfer. To do this accurately, we would have to shorten the rgpulse by 3.8 µsec, then turn the transmitter on during the decpwrf statement and then turn the transmitter off *during* the last AP bus word; however, it is impossible to switch any high-speed lines *during* AP bus events—we can only switch them before or after the decpwrf call. In this case, the second option is the better approximation:

```
status(A)
decpwrf(getval("crossp"));
delay(d1);
decrgpulse(pw,t1,rof1,0.0);
status(C);
decphase(zero);
rgpulse(getval("contact") - 4.6e-6,t2,0.0,0.0);
xmtron();
decpwrf(getval("dipolar"));
xmtroff();
delay(d2);
rgpulse(pw,t3,0.0,rof2);
rcvron();
```

The only error we make with this construct is that the decoupler is switching to dipolar decoupling amplitude for a maximum of $0.8 \,\mu$ sec at the end of the contact time, but that effect is minimal and certainly acceptable.

Using a Waveform Generator in CP/MAS Experiments

A similar problem occurs when a waveform generator is to be used in such experiments (e.g., for amplitude modulated spinlocking). A construct like

```
decrgpulse(pw,t1,rof1,0.0);
decphase(zero);
decprgon(pattern,pw,5.0);
status(C);
rgpulse(getval("contact"),t2,0.0,0.0);
```

is clearly unacceptable, because the 5.75 µsec that it takes to set up the waveform generator (on a UNITY*plus*) are much too long. Immediately after the proton 90-degree pulse, the Hartmann-Hahn spinlocking must start; otherwise, the proton magnetization would be lost completely before the spinlocking starts. In this case it may not even help to play the "trick" that was discussed in the previous section (shortening the proton 90-degree pulse by WFG_START_DELAY, and then turn the decoupler back on while setting up the waveform generator via AP bus), because the

duration of WFG_START_DELAY may be longer than the proton 90-degree pulse. Also, starting the spinlocking while setting up the waveform generator is undesirable, because this would lead to $5.75 \,\mu$ sec of unmodulated cross-polarization (plus 0.45 μ sec for the waveform generator propagation delay on a UNITY*plus*). Clearly, a better method is needed here.

It would be nice if it were possible to start the waveform generator without the AP bus overhead. It turns out that this is possible! Let's first analyze what happens during the decprgon call (see also Section 16.4, "Using Waveform Generators for Programmed Modulation," on page 177):

| 300 | 294 | 188 | 102 | WG3 | AP addr $0x0c18$, IB addr = $0x0000$ |
|-----|-----|-----|-----|----------------|---------------------------------------|
| 303 | 297 | 191 | 101 | WGCMD | AP addr $0x0c18$, WFG cmd = $0x07$ |
| 306 | 300 | 194 | 150 | HighSpeedLINES | DECUP WFG2 |
| 309 | 303 | 197 | 150 | HighSpeedLINES | DECUP WFG2 |
| 312 | 306 | 200 | 151 | EVENT1_TWRD | 1000 msec |

First, the instruction block address is sent to the waveform generator (3 AP bus words, followed by the "waveform generator command" (2 AP bus words), then the high speed line for the waveform generator is set, after which the modulated time event starts (with the modulation actually being delayed by the waveform generator propagation delay of 0.45 µsec on a UNITY*plus*). So, decprgon prepares the waveform generator and sets the high speed line, but that high speed line is only going to be actuated *with the next time event*. We can use this to do what we want:

```
#define OBS_WFG 0x4
#define DEC_WFG 0x80
#define DEC2_WFG 0x1000
#define DEC3_WFG 0x20000
(...)
decprgon(pattern,pw, 5.0);
HSgate(DEC_WFG,FALSE);
decrgpulse(pw,t1,rof1,0.0);
HSgate(DEC_WFG,TRUE);
decphase(zero);
status(C);
rgpulse(getval("contact"),t2,0.0,0.0);
decprgoff();
```

This requires some explanation. After the decprgon call, the high-speed line for the waveform generator is set, as discussed above. *Immediately after that call* (before the next FIFO word is produced), we reset that fast bit to FALSE, such that it isn't actually turned on in the hardware. To do that we need the HSgate function with an address constant that we have extracted from Table 3 in Section 9.2, "Fast Bits," on page 88 for the UNITY*plus*.² During the pulse that follows, the waveform generator is ready (it has got the instruction block address and knows how to execute it), but it doesn't start yet. This only happens when we set the high-speed line *after* the decoupler pulse.

With this solution we are able to start a waveform generator without overhead. Of course, this only works for one pattern per waveform generator, and if the same should be repeated in the sequence with the same waveform generator, we would still need some time to perform the decprgon in-between. Note that on a UNITY, the statement

² The corresponding constants for UNITY systems can be taken from Table 2 in the same section of this manual. Note that the entire construct has been tested only on a UNITY*plus* and will probably not work on UNITY systems. Of course, only one value (DEC_WFG) is used in this pulse sequence fragment; the other values are only given here as reference and wouldn't have to be defined if they are not used.

degprgoff also involves AP bus traffic and hence will still introduce a delay (WFG_STOP_DELAY). On a UNITY*plus* the line

HSgate(DEC_WFG,FALSE);

is actually identical to

decprgoff();

and does not introduce an extra delay.

The only small timing error that we still have with the above construct is the waveform generator propagation delay (450 nsec on a UNITY*plus*), which will introduce 450 nsec of unmodulated spinlocking, but even that can be taken care of if we want to get a really perfect solution:

```
#define DEC_WFG 0x80
(...)
decprgon(pattern,pw,5.0);
HSgate(DEC_WFG,FALSE);
decrgpulse(pw-WFG2_OFFSET_DELAY,t1,rof1,0.0);
HSgate(DEC_WFG,TRUE);
decon(); delay(WFG2_OFFSET_DELAY); decoff();
decphase(zero);
status(C);
rgpulse(getval("contact"),t2,0.0,0.0);
decprgoff();
```

Here, we interrupt the 90 degrees decoupler pulse to enable the waveform generator high-speed line such that it *really* starts at the beginning of the spinlock period³.

20.2 Sideband Suppression in MAS Experiments

TOSS (total sideband suppression) has several variants. All form a relatively simple pulse sequence element: a series of pulses with some specific spacing (which is a function of the rotor speed). In principle, there isn't much to discuss here, but we picked TOSS as an example on how to code a series of spaced pulses and also to perhaps eliminate some common misunderstandings and help you to optimize your coding under various aspects.

An early version of TOSS coding, taken from /vnmr/psglib/xpolar.cin VNMR 4.3, is the following:

```
rcvroff();
delay((0.1226/srate)-pw);
rgpulse(2.0*pw,v3,0.0,0.0);
rcvroff();
delay((0.0773/srate)-2.0*pw);
rgpulse(2.0*pw,v4,0.0,0.0);
rcvroff();
delay((0.2236/srate)-2.0*pw);
rgpulse(2.0*pw,v3,0.0,0.0);
rcvroff();
delay((1.0433/srate)-2.0*pw);
rgpulse(2.0*pw,v4,0.0,rof2);
delay((0.7744/srate)-pw-rof2);
```

³ Instead of WFG2_OFFSET_DELAY we could also take WFG_OFFSET_DELAY as a constant, because all waveform generators in a system will generally behave the same way.

There is a misunderstanding involved in this coding. It is not necessary to switch off the receiver after every pulse, because after a first rcvroff statement the receiver is regarded to be "globally off" and is not switched back on after a pulse, but only by the next (implicit or explicit) rcvron call, as shown in Section 4.2, "How Do Pulses Work?," on page 40 (although this has no negative effect, of course, other than creating an unnecessary HighSpeedLINES call in the Acode). This has been cleaned up in a more recent coding from the user library (slightly simplified):

```
rcvroff();
(...)
if (toss[0] == `y')
{
    fprintf(stdout,", TOSS");
    delay((0.1226/srate)-pw);
    rgpulse(2.0*pw,v3,0.0,0.0);
    delay((0.0773/srate)-2.0*pw);
    rgpulse(2.0*pw,v4,0.0,0.0);
    delay((0.2236/srate)-2.0*pw);
    rgpulse(2.0*pw,v3,0.0,0.0);
    delay((1.0433/srate)-2.0*pw);
    rgpulse(2.0*pw,v4,0.0,0.0);
    delay((0.7744/srate)-pw);
}
```

Because this avoids the extra rcvroff calls, it looks easier to read, but it certainly still isn't perfect. In both versions the transmitter phase is not preset; therefore, over the first 500 nsec of each pulse (UNITY*plus*) there will be some phase transition error (there is no pre-pulse delay). However, the following construct is definitely undesirable:

```
delay((0.1226/srate)-pw-rof1);
rgpulse(2.0*pw,v3,rof1,0.0);
delay((0.0773/srate)-2.0*pw-rof1);
rgpulse(2.0*pw,v4,rof1,0.0);
delay((0.2236/srate)-2.0*pw-rof1);
rgpulse(2.0*pw,v3,rof1,0.0);
delay((1.0433/srate)-2.0*pw-rof1);
rgpulse(2.0*pw,v4,rof1,0.0);
delay((0.7744/srate)-pw);
```

This fixes the phase error but creates a new problem. A known trouble with the TOSS pulse sequence element is that it stops working above a certain rotor speed. In the expression (0.0773/srate)-2.0*pw, the value starts becoming negative as soon as 2*pw > (0.0773/srate). The above coding increases that problem, because we now also subtract rof1 from that delay.

The following, very compact coding avoids this problem *and* allows for phase-presetting:

```
rgpulse(2.0*pw,v3,(0.1226/srate)-pw,0.0);
rgpulse(2.0*pw,v4,(0.0773/srate)-2.0*pw,0.0);
rgpulse(2.0*pw,v3,(0.2236/srate)-2.0*pw,0.0);
rgpulse(2.0*pw,v4,(1.0433/srate)-2.0*pw,0.0);
delay((0.7744/srate)-pw);
```

Unfortunately, we have now traded in another problem: the dps command in VNMR hides pre- (and post-) pulse delays. Therefore, this element will be shown as four back-to-back pulses, which again is undesirable. We can, of course, explicitly pre-set the phase using the txphase function, as in the following coding⁴:

```
txphase(v3);
delay((0.1226/srate)-pw);
rgpulse(2.0*pw,v3,0.0,0.0);
txphase(v4);
delay((0.0773/srate)- 2.0*pw);
rgpulse(2.0*pw,v4,0.0,0.0);
txphase(v3);
delay((0.2236/srate)-2.0*pw);
rgpulse(2.0*pw,v3,0.0,0.0);
txphase(v4);
delay((1.0433/srate)-2.0*pw);
rgpulse(2.0*pw,v4,0.0,0.0);
delay((0.7744/srate)-pw);
```

This finally has the phases preset, and will be doing overall what we want. It also avoids unnecessary receiver gating (assuming the receiver is off globally), and it works with the dps command. Alternatively, we could use a coding that is also more efficient in terms of Acode space and execution speed (if that ever becomes an issue), because with the above coding the rgpulse function will create extra SETPHAS90 and HighSpeedLINES calls in the Acode:

```
rcvroff();
(...)
if (toss[A] == 'y')
{
   txphase(v3);
   delay((0.1226/srate)-pw);
   xmtron(); delay(2.0*pw);xmtroff();
   txphase(v4);
   delay((0.0773/srate)-2.0*pw);
   xmtron(); delay(2.0*pw);xmtroff();
   txphase(v3);
   delay((0.2236/srate)-2.0*pw);
   xmtron(); delay(2.0*pw);xmtroff();
   txphase(v4);
   delay((1.0433/srate)-2.0*pw);
   xmtron(); delay(2.0*pw);xmtroff();
   delay((0.7744/srate)-pw);
}
```

Although at first it seems odd not to use pulse or rgpulse for something that "is a pulse", this coding isn't really too complex, and it fulfills all the needs. It is readable, it does the proper thing in hardware, and it also is most efficient in terms of Acode size and execution speed. This efficiency can be relevant in cases where multiple hardloops or a large number of pulses in general is involved, or in *n*D experiments, where excessive Acode size can lead to a loss of the Acode buffering (leading to extra delays and a disruption of the steady-state between increments).

⁴Note that fprintf(stdout,...) in the previous version can always be replaced by the simpler printf(...) function. printf sends the output to stdout.

20.3 Rotor Synchronization

Rotor synchronization experiments require dedicated hardware (consisting of a counter that counts a given number of rotor periods, and a timer/counter that counts *within* the rotor period). This accessory can be used in three different ways, as described in the following sections.

Measuring the Rotor Period Duration

Measuring rotor period duration is achieved with the rotorperiod (vn); statement. In this mode, the rotor synchronization hardware *returns* the number of 100 nsec clock cycles *to the acquisition CPU*. The accessory actually counts the clock cycles for every rotor period and stores the last value in a register. rotorperiod causes the pulse programmer to *read* that register through the AP bus (which is bidirectional) and store it in one of its own registers, from where the acquisition CPU can retrieve it into the specified vn(v1 to v14) variable. This could be used to perform pulse sequence events at defined pulse angles (relative to the trigger position), using also the external trigger of the pulse programmer. In the example below, we want to perform a pulse at 0, 120, and 240 degrees rotor positioning on successive scans:

```
rotorperiod(v1);
                    /* rotor period in 100 nsec units */
/* 0 1 2 0 1 2 3 */
                    /* 50 nsec units */
dbl(v3,v3);
dbl(v3, v3);
                    /* 25 nsec units */
mult(v3,v2,v3);
                   /* 0, 120, 240 degrees, in 25 nsec */
(...)
                  /* wait for next trigger */
xgate(1.0);
ifzero(v3);
                   /* don't wait for 0 degrees */
elsenz(v3);
 sub(v3,three,v3); /* subtract 150 nsec overhead */
 sub(v3,three,v3);
                    /* perform 25 nsec time event */
 vdelay(NSEC,v3);
endif(v3);
rgpulse(...
```

After reading the rotor period duration using rotorperiod and constructing the scan-based multiplier, we calculate the desired angle resolution in 100-nsec units and then convert that to 25-nsec units. We then multiply that number by the scan-based multiplier and get the number of 25-nsec counts to perform on every scan.

We then wait for the next rotor trigger using xgate(1.0). Because we don't want to (and cannot) perform a vdelay with zero counts, we exclude that case (0 degrees angle) using a real-time ifzero statement—a pulse could directly follow. For the other cases (non-zero phase angles), we subtract the pulse programmer overhead of 6 timer units (150 nsec), and then we perform a vdelay, which is a delay with a given time base and a timer count given as real-time variable.

However, we have to be careful, in that we should consider that

- Real-time variables are 16-bit integers with a maximum value of 32767.
- Any counter on the pulse programmer is only 12 bits, allowing for a maximum value of 4096.

At 10 KHz rotor speed, the rotor period is 100 nsec, corresponding to 4000 clock cycles of 25 nsec. In this particular case, we would be dealing with a maximum of 2/3 of a full rotor cycle, which must be less than 4096*25 + 150 nsec or 102.55 µsec. This corresponds to a rotor period of 153.825 µsec, or a rotor speed of 6500 Hz. Therefore, it may therefore be safer to do the same thing using the µsec time base:

```
rotorperiod(v1);
                          /* rotor period in 100 nsec units */
modn(ct,three,v2);
divn(v1,three,v3);
add(two,three,v4);
                          /* 0 1 2 0 1 2 3 */
                          /* rotorperiod/3 */
                          /* 5 */
dbl(v4,v4);
                          /* 10 */
                          /* usec units */
divn(v3,v4,v3);
mult(v3,v2,v3);
                          /* 0, 120, 240 degrees, in usec */
(...)
xgate(1.0);
                          /* wait for next trigger */
ifzero(v3);
                          /* don't wait for 0 degrees */
  delay(0.35e-6);
elsenz(v3);
  delay(0.2e-6);
  vdelay(USEC,v3);
                          /* perform 25 nsec time event */
endif(v3);
rgpulse(....
```

In this case, we have to face round-off errors of up to 500 nsec (below 2 degrees at typical rotor speeds). To avoid extra errors from the timer word overhead of 150 nsec, we should perform a 150 nsec delay for the 0-degree case. Because this is not possible (with the minimum delay being 200 nsec), we perform a 350 nsec delay for the 0-degree case, and a 200 nsec delay where we use a vdelay. This solution gives us slightly less accuracy, but definitely no problem with rotor speeds as low as a few hundred Hz.

Waiting for Triggers

The xgate function (which was also used in the example in the previous section) loads the (12-bit) time counter on the pulse programmer with the specified number of counts and switches the counter to the external time base (the external trigger). On each trigger, the counter counts one unit down, and the next pulse sequence event starts when the count reaches zero. Often that time count will be just 1 (1.0, as the argument must be a floating point number). If in the above example the final pulse is to be performed after a longer delay, we have the options to perform a normal delay, followed by the xgate(1.0) call, or we could calculate how many rotor cycles that delay would be (this is typically done based on a VNMR parameter srate) and then perform xgate with that calculated number of rotor triggers. Note, however, that the number of rotor cycles that can be counted this way is 4096 only, because the pulse programmer uses a 12-bit counter. At typical rotor speeds of 5 to 10 KHz, this limits the "counted" delay to 0.8 to 0.4 seconds.

Rotor-Synchronized Experiments

True rotor-synchronized experiments go a bit further than the above examples. In the rotor-synchronized version of NOESY (XNOESYSYNC pulse sequence) the requirement is that the two pulses adjacent to the mixing time be performed *at the same rotor angle*, after a given number of rotor periods. This is done with the rotorsync pulse sequence element (simplified coding):

```
(...)
initval(srate*mix,v10);
rgpulse(pw,t2,0.0,0.0);
txphase(t3);
rotorsync(v10);
rgpulse(pw,t3,0.0,rof2);
(...)
```

The rotor synchronization accessory counts through every rotor period, as mentioned before. The rotorsync function stops that counter through a high-speed line. Then the number of rotor periods specified in the real-time argument variable (calculated from the srate parameter and the duration of the delay) is loaded into the accessory, which starts counting rotor (tachometer) triggers until that number is reached (it actually counts from that number down to zero). After that it takes the number of 100 nsec clock cycles, at which it was stopped initially, and counts these *down* to zero. After that, it sends a trigger to the pulse programmer, which has been set to "external timebase" mode with a trigger count of 1. Hence, this causes the next pulse sequence element (the pulse that follows) to be executed.

This complex sequence of events guarantees that when the trigger is sent back to the pulse programmer, the rotor is at exactly the same angle as when the pulse programmer stopped the 100 nsec timer through the high-speed line. Note that in this case we do not use the timer from the pulse programmer to count the rotor cycles; therefore, there is no 12-bit (4096) count limit. The limit is 32767, the largest number that can be stored in a real-time variable.

20.4 Multipulse Experiments

Multipulse line narrowing consists of a series of pulses that is performed between each pair of acquired data points. Although it very much looks like a TOSS pulse train (see Section 20.2, "Sideband Suppression in MAS Experiments," on page 230), the requirements are much different:

```
pulsesequence()
{
   double tau = getval("tau"),
         dtau = tau - pw - rof1 - rof2;
   initval(np/2.0,v9);
   . . .
   delay(d1);
   rgpulse(pw,v4,rof1,rof2);
                                     /*prep pulse */
   starthardloop(v9);
     delay(dtau);
                                     /* x */
     rgpulse(pw, v4,rof1,rof2);
     delay(dtau);
      rgpulse(pw,v3,rof1,rof2); /* -y */
      delay(tau + dtau);
```

```
rgpulse(pw,v1,rof1,rof2);
                                  /* y */
   delay(dtau);
                                  /* -x */
   rgpulse(pw, v2, rof1, rof2);
   delay(tau + dtau - 2.0e-7);
   acquire(2.0,2.0e-7);
                                  /* acquire */
   rgpulse(pw,v2,rof1,rof2);
                                  /* -x */
   delay(dtau);
   rgpulse(pw,v3,rof1,rof2);
                                  /* -v */
   delay(tau + dtau);
   rgpulse(pw,v1,rof1,rof2);
   delay(dtau);
   rgpulse(pw, v4,rof1,rof2);
                                  /* x */
   delay(tau);
endhardloop();
```

}

Unlike the TOSS pulse train, we cannot keep the receiver switched off globally; otherwise, we would not be able to see any signal! Also, the pulses here are extremely short (typically 1 to 2 μ sec), so a 500 nsec phase or amplitude glitch at the beginning of each pulse would be clearly unacceptable. Also, the amplifier blanking/unblanking needs to be timed very carefully. For good signal-to-noise, we need to minimize the pre-and post-pulse delays. On the other hand, we also need the amplifiers fully unblanked, and we need to avoid destroying the preamplifiers with the ring-down voltages from the 1 KW pulses (which rather calls for longer pre- and post-pulse delays!). Still, there will be preamplifier saturation from the residual probe ring-down; therefore, we should place the acquisition as late in a delay as possible.

Note also that in the above coding, the maximum number of points that can be acquired is 65K, because the 16-bit hardloop count limits the number of loop cycles to 32768. If more points were to be acquired, we would have to pack multiple pulse trains (and acquisition triggers) into a single loop cycle (in which case again we must set the parameter np to a number that corresponds to full loop cycles; otherwise, the error message "number of points acquired not equal to np" is obtained.

20.5 Other Line-Narrowing Techniques

A new line narrowing technique that is now gaining popularity is the so-called "frequency-switched Lee-Goldberg" (FLSG2) method⁵, which involves rf irradiation at a certain field strength with rapid frequency switching between two frequencies on either side of the proton spectrum.

The critical point about this technique is that the frequency switching must occur phase-coherently and very rapidly (actually with a 180-degree phase shift at the frequency switch). It was often thought that these requirements can be met only with special frequency synthesizers with latching (i.e., the frequency information is sent first and all the digits of the frequency change at the same time) and extremely fast (and phase coherent) frequency switching. the latter can only be obtained with overrange (i.e., with an extended range of DDS (direct digital frequency synthesis)). Such frequency synthesizers are available, but it is well known that digital frequency synthesis at higher frequencies causes increased spurious output (it is just much harder

⁵ A. Bielecki, A.C. Kolbert & M.H. Levitt, *Chem. Phys. Lett.*, **155**, 341 (1989). For the actual implementation see the article by Jun Ashida & David Rice in *Magnetic Moments*, Winter, 1996, Vol. VIII,1.

to do high frequencies/offsets digitally, rather that using analog techniques). Also, any phase error during the frequency and phase transition will lead to errors that accumulate and deteriorate the efficiency of the experiment.

It turns out that with the waveform generator, we can implement these frequency offsets as linear phase ramps (shifted laminar pulses or SLP) We used 52 steps per ramp, with 4 degrees phase shift per slice. When the frequency changes, the phase ramp changes direction (and 180 degrees are added to the phase).

This turns out to be an ideal way to implement this technique. The frequency is generated using small-angle phase shifts that are performed within 100 nsec, making frequency switching extremely rapid. And the 180-degree phase switching is done with the same hardware; therefore, there is no question that both the phase and frequency switching occur in an absolutely synchronous and accurate manner, with no accumulation of phase errors or the like.

Chapter 20. Solid-State NMR Experiments

Chapter 21. (Micro)Imaging Experiments

Although, with few exceptions, imaging experiments use the same spectrometer hardware, the pulse sequences for imaging and the associated data (parameters, macros) look totally different from any pulse sequence for liquids or solids NMR. This has a variety of reasons:

- When setting up an imaging experiment, we are dealing with the physical dimensions of the object we are looking at. On the parameter side, there is a complex relationship between the physical dimensions in the image (e.g., the field of view, the slice thickness and orientation) and the experimental parameters. Ultimately it all relates to a field gradient strength, but this gradient strength again depends on the gradient amplifiers, the gradient coil, etc. This all leads to a fairly large number of parameters involved in imaging experiments. That doesn't mean that the imaging experiment is extremely complex. We do not have to interact with all those parameters! It is just that the parameter aspect is totally different from the setup of a high-resolution NMR experiment.
- Most or all imaging experiments use multi-echo/multi-FID acquisition for an entire plane, if not for a complete 3D image.
- Some or all of the indirect dimensions are not using time incrementation but a gradient strength (sometimes a frequency offset) is varied instead. Many imaging experiments use waveform generators for field gradient pulse scaling and shaping.
- To make complex imaging sequences still readable and understandable, the pulse sequence code has been broken up into a number of C functions (such as for performing a shaped pulse on top of a gradient pulse) that are defined in the pulse sequence file itself and which are called from within the pulsesequence function. This leads to a coding style that is totally different from anything in standard liquids and solids NMR.

Another difference from the rest of the pulse sequences is that all imaging pulse sequences are now defined consistently, both in terms of the parameters involved as well as in the C coding used in the pulse sequence file itself. Overall, the *requirements* for imaging experiments are totally different from those of any other NMR experiment, which justifies the fundamental difference in the pulse sequence coding.

It is beyond the scope of this manual to explain all the special mechanisms involved in imaging experiment; however, if you look into the dedicated imaging functions, you will see that it still uses standard pulse sequence elements. The main and basic difference is the implicit coding for the nf-type multi-FID acquisition.

Chapter 21. (Micro)Imaging Experiments

Chapter 22. Role of Macros and Parameters

Most pulse sequences are more complex than the s2pul standard sequence for 1D experiments. They include more pulses, more delays, pulses on other channels, gradient pulses, etc., and all these additional pulse sequence elements will depend on parameters.

In principle, it is of course possible to "hard-code" the parameters for such new pulse sequence elements into the C program, for example:

pulse(13.5e-6,v1,10e-6,0.0);

which would avoid any extra parameter definition etc. Unfortunately, such parameters are rarely constant, because they will depend on a large number of "external" factors, such as the machine and amplifier type, the frequency, the probe type, rf and gradient amplitude settings, and even the rf routing, but also the probe tuning and the sample properties (like the susceptibility). Even though simple delays are not hardware-dependent, their optimum setting typically depends on the NMR properties of the molecules under investigation.

All this would make it necessary to change the pulse sequence and recompile it for the specific requirements of the experiment to be run, essentially for each and every acquisition. This is highly undesirable. The C coding of most pulse sequences is difficult to understand for most of the simple NMR users, and most users don't even look at the C code of a pulse sequence. Therefore, it is strongly recommended to keep anything in a pulse sequence variable that may possibly need adjustment for a given hardware setup and a given sample.

On the other hand, many of today's pulse sequences use a large number pulses on various rf channels, combined with several gradient pulses, possibly involving several waveform generator shapes and pattern. Keeping all this under variable control can lead to an enormous number of variables, which again can make a pulse sequence difficult and cumbersome to use.

Therefore, the target should be to *keep the number of variables at a minimum without giving up pulse sequence flexibility*. One obvious way to achieve this is to use *parameter dependencies* where possible (e.g., by deriving all pulse widths at a given channel and power level from a single pulse length, typically the length of a 90-degree pulse).

Assuming amplifier linearity and perfect pulse shapes (modern spectrometer hardware fulfils this criterion to a large extent), we could even extend this principle to pulses at different power levels, theoretically up to a point where all pulse widths (assuming the same pulse shape) on a given rf channel are derived from a single pulse width parameter at a reference power level. Similar considerations apply to gradient amplitudes and delay lengths.

22.1 Creating New Parameters in VNMR

As mentioned before, we want to be able to manipulate pulse sequence parameters from within VNMR. This way it is not necessary for the user to read and understand the C code, and complex experiments can be driven with a simple(r) user interface. In this scheme VNMR parameters become the primary medium for interacting with a pulse sequence.

VNMR comes with a basic collection of parameters: configuration parameters (/vnmr/conpar), global parameters (~/vnmrsys/global), in addition to the current and processed parameter trees (~/vnmrsys/expn/curpar and ~/vnmrsys/expn/procpar). For driving a pulse sequence, only the configuration parameters, plus the acquisition parameter group in both the global and the current parameter trees, are relevant. Only the parameters in use in the s2pul pulse sequence can be assumed to be present in a general case (in many cases the user will, in fact, first do a simple 1D experiment before switching to a more complex pulse sequence). Hence, for a new, complex pulse sequence, the user has to *create* the necessary new parameters, which is typically done using the create command in VNMR:

create('parameter_name'<,'type',<'tree'>>)

Typical pulse sequence parameters are created in the *current* tree (which is the default); therefore, it is normally not necessary to specify the parameter tree. From a parameter *value* point-of-view it is most important to distinguish *pulses* from any other numeric parameters, because pulses are specified in microseconds. Using a real or integer, frequency (Hz) or delay (sec) parameter for pulses by mistake could lead to pulses of several seconds¹.

From a parameter handling and security point-of-view (avoiding more subtle parameter mis-settings), it is definitely important to specify the correct parameter type, which is one of string, flag, delay, frequency, real, integer, or pulse. This leads to parameters with the properties shown in Table 15.

| Туре | Default properties | |
|-----------|--|--|
| string | character string, enumerals indicate entire string values | |
| flag | character string, enumerals are possible string elements | |
| delay | real, maximum 8190 sec, minimum 0, stepsize 100 or 25 nsec | |
| frequency | real, maximum 1e9 Hz, minimum –1e9 Hz | |
| real | real, maximum 9.99999984307e17, minimum –9.99999984307e17 | |
| integer | real, maximum 32767.0, minimum 0.0, stepsize 1.0 | |
| pulse | real, maximum 8190 µsec, minimum 0, stepsize 0.1 or 0.025 µsec | |

| Table 15. VN | MR parameter | types and | propertie |
|-----------------|--------------|-----------|-----------|
|-----------------|--------------|-----------|-----------|

¹ As of VNMR 5.1, the built-in probe protection software should prevent probes from being burnt or amplifiers from being destroyed due to such mistakes, but it is of course better not to rely on this.

The parameter limits for the pulse and delay type parameters are "indirect" limits² and are taken from the parameters parmin, parmax, and parstep in /vnmr/conpar (pulse type parameters use element number 13 in these parameters, delay type parameters use element number 14). These "global" parameter limits (in particular, the pulse and delay step size) are set by the config program and depend on the actual spectrometer hardware (pulse programmer).

The "indirect" parameter limit is a very useful feature. There are countless parameter sets around, each of them with numerous delays and pulses. By having the parameter limits for these parameters stored in one central place (/vnmr/conpar), it is not necessary to adjust the parameter limits individually when moving data or software between spectrometers with different acquisition hardware.

Note that there is no "proper" integer parameter type. All numeric parameters are stored as floating point numbers. Not all pulse sequence parameters are covered by an appropriate parameter type. Power levels (rf attenuator settings, linear modulator amplitude levels), gradient amplitudes and some others do not have their own parameter types. Many of these "missing types" are covered by integers with the appropriate parameter limits. Wherever possible, these limits are defined indirectly in /vnmr/conpar, via the parameters parmax, parmin and parstep (see Table 16).

For new parameters it is advisable to limit the parameter entry to applicable values by activating the appropriate parameter limits. Wherever such parameter limits may be hardware-dependent it is advisable to use indirect parameter limits from Table 16. This leads to the following recipe for creating new parameters:

• For *frequency offsets* (and also for parameters holding coupling constants), you can use either the default frequency parameter properties or the parameter limits³ for the tof parameter⁴:

```
create('parameter_name','frequency')
setlimit('parameter_name',7)
```

• For *delays*, and standard *rf pulses* the default parameter properties should be sufficient:

```
create('parameter_name','delay')
create('parameter_name','pulse')
```

• For *shaped or selective rf pulse lengths* and parameters that define the *duration of gradient pulses*, the default parameter maximum for pulses is too small, because selective pulses can be up to tenths of a second long. You can use either a delay parameter to define such pulse lengths, or, probably better, use a pulse parameter

² The definition is such that if protection bit 13 (value 8192, see the *VNMR Command and Parameter Reference* for the setprotect command) is set, the parameter definition does not contain the actual parameter limits, but rather indices into the three systemglobal parameters parmin, parmax, and parstep, which are numeric arrays (see also the *VNMR Command and Parameter Reference* for the setlimit command).

³ Note that as of VNMR 5.1, setlimit with two arguments automatically activates protection bit 13 (value 8192). Under earlier versions of VNMR, this would be equivalent to

setlimit('parameter_name',x,x,x)

setprotect('parameter_name','on',8192)

⁴ To be very correct, we would have to use index 8 for offsets on the decoupler channel, and indices 16 and 20 for offsets on the second and third decoupler channel, respectively. However, these limits are typically the same for all channels and not very critical. Also, with the possibility to reassign rf channels, that difference becomes irrelevant.

with modified parameter limits (*do not try* simulating millisecond parameters by multiplying pulse parameter values by 1000 or dividing delay values by 1000 in a pulse sequence, because this just complicates and obscures the issue!).

```
create('parameter_name','pulse')
setlimit('parameter_name',1e6,0,0.025)
```

• For *pulse power levels (attenuator settings)*, you should use the indirect limits used for the tpwr parameter. To be very correct, we should use the power parameter limits for the rf channel the new parameter refers to, but this is often not appropriate, because the parmax[9] (used for dpwr) is meant to be for a power level used in continuous decoupling and is usually set to 49, which is inappropriate for hard rf pulses.

```
create('parameter_name','integer')
setlimit('parameter_name',17)
```

• For *power levels used for continuous irradiation*, you can use the parameter limits also used with the dpwr parameter (where the maximum is at a safer level):

```
create('parameter_name','integer')
setlimit('parameter_name',9)
```

| Indir | rect parameter limit definition | – Parameter examples | |
|-------|---------------------------------|---------------------------|--|
| Index | Typical max / min / stepsize | | |
| 1 | 500 / 0 / 0.1 | sc | |
| 2 | 840 / 5 / 0.1 | WC | |
| 3 | 500 / 0 / 0.1 | sc2 | |
| 4 | 520 / 5 / 0.1 | wc2 | |
| 5 | 100000 / 100 / 25e-9 | sw | |
| 6 | 51200 / 200 / 200 | fb | |
| 7 | 100000 / -100000 / 0.1 | tof | |
| 8 | 99000 / -99000 / 0.1 | dof | |
| 9 | 49 / -16 / 1 | dpwr,dhp | |
| 10 | 39 / 0 / 1 | dlp | |
| 11 | 2e6 / 1 / 1 | dmf | |
| 12 | 500 / -500 / 1 | loc | |
| 13 | 8190 / 0 / 0.025 | p1,pw,pw90,rof1,rof2,alfa | |
| 14 | 8190 / 0 / 25e-9 | d1,d2,pad,dod,vtwait | |
| 15 | 1e6 / 0 / 0.025 | | |
| 16 | 100000 / -100000 / 0.1 | dof2 | |
| 17 | 63 / -16 / 1 | tpwr | |
| 18 | 63 / -16 / 1 | dpwr2 | |
| 19 | 32767 / -32768 / 1 | (shim gradient values) | |
| 20 | 100000 / -100000 / 0.1 | dof3 | |
| 21 | 49 / -16 / 1 | dpwr3 | |

Table 16. Predefined, indirect parameter limits

• For *linear modulator amplitude settings*, you must define the limits explicitly:

```
create('parameter_name','integer')
setlimit('parameter_name',4095,0,1)
```

• The same is true for *pulsed field gradient (PFG) amplitudes*⁵:

```
create('parameter_name','integer')
setlimit('parameter_name',32767,-32767,1)
```

• For *string parameters* (typically used for shape and pattern names), no parameter limits need to be defined. It may occasionally be useful to define enumerals⁶ to avoid parameter entries that don't make sense:

```
create('parameter_name','string')
setlimit('parameter_name',17)
```

• For *multi-field flags* (typically only allowing for the characters 'y' and 'n') enumerals are very useful, because they simplify the value testing within the pulse sequence (e.g., if there are only two allowed character values 'y' and 'n', we can simply test for one of the values, and if that test fails, the other value can be assumed). Also, flag enumerals are an easy way to define a large number of possible values in multifield flags⁷.

```
create('parameter_name','flag')
setenumeral('parameter_name',2,'y','n')
```

• A *single-field flag* can alternatively be defined as string parameter with enumerals, which would be more restrictive in the parameter entry:

```
create('parameter_name','string')
setenumeral('parameter_name',2,'y','n')
```

Needless to say, any new parameter should be *filled with a sensible value*. It is best to do this right after defining the parameter, because otherwise it may be forgotten! Note in particular that defining enumerals for strings and flag parameters does *not* automatically select or fill in a "legal" value: create followed by setenumeral will leave a parameter with an empty string (''), which may lead to unexpected results when testing for either 'y' or 'n' alone in the pulse sequence!

⁵ The proper limits for Performa I type gradients are 2047, -2047, 1; but still it is recommended to use the limits for Performa II gradients as shown in the text, because this makes the parameter set portable between systems with different PFG amplifiers.

⁶ See the VNMR Command and Parameter Reference for the setenumeral command. The problem with enumerals in strings is that VNMR is not very helpful in the case of an "illegal" parameter entry. An error message is issued, but no hint as to what the allowed values are is given. You have to use display('parameter_name') to see the enumerals for a parameter. For pulse shapes and WFG patterns, it is therefore better *not* to define enumerals (you will get an error message from the go command if a non-existent pattern or pulse shape is specified).

⁷ For a flag with *n* fields, allowing for the characters 'y' and 'n', the number of possible values is 2^n , or even more if we account for abbreviated versions (i.e., the fact that by convention the last character of a flag parameter is propagated to any subsequent position: 'n' stands for 'nnnnn...').

22.2 Using New Parameters in C

In order to be able to use the value of a new parameter that has been defined in VNMR, we have to use the special functions getval (for numeric parameters) and getstr (for flag and string variables) that read such parameter values off the VNMR parameter table and make them available in the pulse sequence environment.

Numeric Parameters

The getval function serves to extract numeric parameter. It returns a double, irrespective of the parameter type. Note that values from *pulse parameters* are converted and passed to the pulse sequence environment *in seconds*! Typically, the getval function is used to initialize a C variable:

```
double pwx;
pwx = getval("pwx");
```

The VNMR parameter name in the argument of the getval function and the name of the C variable don't necessarily have to be the same, but of course this is very much preferable. Because it is advisable to keep pulse sequences simple, we should always watch out for syntax simplifications. One such possibility is to combine variable declaration and initialization:

If a parameter value is used only once in a pulse sequence, it is not really necessary to define a C variable that holds the value. The function getval can also be used directly as argument to an other function:

rgpulse(2.0*getval("pwx"),zero,rof1,0.0);

And, of course, getval functions can also be used within mathematical expressions: tau = 1.0/(2.0 * getval("j"));

String Parameters

String parameters are slightly more complicated. In principle it would be possible to create a function similar to getval for strings (i.e., a function returning a (pointer to a) string). The problem with this is that the user would have to deal with pointer variables, and that was considered to be too complex for a simple pulse sequence language (apart from that, pointers are a possible source for errors that are hard to debug in C programs). Consequently, string parameters must be filled into a string C variable:

```
char xpol[MAXSTR];
getstr("xpol",xpol);
```

MAXSTR is defined as 256. Note that after the (fixed-length) string variable definition, the string is *not* automatically initialized; its contents can be assumed to be random.

Do not forget to initialize variables! This mechanism is used both for strings and flag parameters. For checking individual character fields within a status-related flag variable, see also "Checking Flag Parameters" on page 160.

22.3 Adding New Parameters to the Display

In VNMR, creating and setting the new parameters isn't just enough. You want to be able to see what the values are, and eventually print out complete parameter listings. Typing parameter_name?, as the only way to see a parameter value, isn't a good solution, because it is very easy to overlook parameters that are not shown on the screen or listed in a printout.

Editing display (dg, ap) templates is usually by paramvi ('parameter_name') or by paramedit('parameter_name'). If you have defined an environment variable vnmreditor in your ~/.login⁸ file, paramedit can be set up to use textedit instead of the vi editor, which is more convenient for editing parameters like dg and ap, which typically consist of a single, long line (in which occasional users of vi will have difficulties in setting the insertion point⁹). Entering such a long parameter value directly in VNMR (dg='.....') is difficult, if not impossible. Apart from that, you don't want to re-enter the entire template, but you selectively want to add the new parameters.

The parameter display templates are explained in Chapter 5, "Modifying Parameter Displays in VNMR," in the manual *VNMR User Programming*, Be careful with conditional displays. Referring to non-existent parameters in a condition expression causes nasty error messages!

As of VNMR 5.1. it is possible to split a display template into substrings by setting up an array of strings—typically one substring per column title. Although such split display templates are incompatible with earlier versions of VNMR, you will find the shorter substrings much more manageable when editing (particularly with vi). It even becomes feasible to substitute entire substrings, for example:

dg[2] = '1:PULSE SEQUENCE:seqfil,d1,p1(p1):1,d2(d2),pw:1;'

It is advisable to at least define a dg and an ap template (if you want, you can also use the dg template for the ap and pap commands, although that will create 4-column output and arrays will not be shown). You certainly want to make sure that all relevant acquisition parameters (in particular those that are pulse sequence-specific) are on contained in a template. For the ap template to work properly with arrays, it should be set up with two columns only.

A potential problem with the dg command and template is that the available space is very limited, and you often have to use several display templates to cover all relevant parameters. The dg command is well-behaved if there are too many parameters in a column. Extra parameters are shown on the top of the next column (except for the last one, where extra parameters will simply not be shown!). Note, however, that dg does *not* handle properly the case where the "parameter overflow" is more than one column!

 $^{^{8}}$ You can also type setenv vnmreditor textedit in a shell window before starting VNMR with the vn command.

⁹ The best way to position the insertion point in such a long line using vi is to search for a text string (/substring), rather than using the keyboard arrows.

22.4 Doing It All by Macro

Obviously, there is a fair amount of "extra setup work" related to creating a new pulse sequence. We certainly don't want to repeat that every time when we want to use such a sequence! Therefore, many years back, we started using pulse sequence-specific macros that allow us to easily switch between a simple 1D setup experiment and any more complex sequence.

Of course, we could also use the strategy of storing one complete and ready-to-use parameter set for every pulse sequence, which can simply be recalled before starting the desired sequence. The disadvantage of this approach is that sample-specific parameters that have been calibrated in the 1D setup experiment (such as transmitter offsets, spectral window, referencing parameters, probe- and sample-specific pulse widths and power levels) are lost and would have to be restored manually, which of course is a source for errors and omissions. We therefore strongly recommend creating such pulse sequence-specific macros.

Apart from the tasks described above, such a macro can carry the action a bit further and also adjust additional parameters (existing ones that need adjustment), display a manual file for the pulse sequence, etc. And of course the macro should also select the pulse sequence it is written for by setting the seqfil parameter. The ideal pulse sequence macro should be written in such a way that an inexperienced user can do a 1D setup experiment, then simply type the name of the desired pulse sequence (i.e., execute the pulse sequence-specific macro) and start the acquisition. This user should then end up with a *reasonable* spectrum using the built-in default parameter values, as set by the macro.

Macros for 1D Pulse Sequences

Possible tasks for a 1D pulse sequence setup macro can be described as follows:

- 1. Select the pulse sequence by setting the seqfil parameter.
- 2. Create new, pulse sequence-specific parameters.
- 3. Where necessary, set protection bits and limits for the new parameters.
- 4. Fill reasonable values into the new parameters.
- 5. Set other parameters to pulse sequence-specific settings, where necessary (e.g.: set nt to a multiple of the pulse sequence phase cycle length).
- 6. Set up parameter display templates (dg, ap).
- 7. Preadjust processing and display parameters, where necessary.
- 8. Display information on the new pulse sequence (typically a manual file) with further information.
- 9. Display the pulse sequence using dps (optional).

The first five tasks listed above can be done with standard MAGICAL constructs and commands listed above:

```
seqfil='dept'
create('mult','real')
setlimit('mult',2,0,0) mult=0.5,1,1,1.5
create('tau','delay')
create('pp','pulse')
```

```
create('pplvl','integer')
setlimit('pplvl','17)
....
```

There is a problem with this construct, in that create results in an error if a parameter already exists. A safer construct would be:

```
seqfil='dept'
exists('mult','parameter'):$e
if not($e) then
    create('mult','real')
    setlimit('mult',2,0,0)
endif
mult=0.5,1,1,1.5
...
```

Even though that looks complex for sequences with many parameters, it certainly is readable. Where this concept falls apart is with the parameter display templates. We can't edit a parameter from within a macro that may be called in background (we don't want to re-edit the templates upon recalling a sequence anyway), and we don't want to set the entire dg and ap parameters with constructs like

```
dg = ''
dg[1] = '1:SAMPLE:date,temp,solvent,file(file<>\'exp\');'
dg[2] = ...
```

because this would make setting up such macros quite complicated!

Instead, the following concept has been adopted by most users: We generate the new parameters once (as shown above) and fill in the appropriate default values, and we generate the modified display template parameters. Then we *save the entire parameter set* in parlib (/vnmr/parlib or ~/vnmrsys/parlib). Now, in the pulse sequence start-up macro, we simply *pick the new or modified parameters with their values* from that parameter set:

psgset('dept','mult','j','tau','presat','dg','ap')

This saves us from checking whether a parameter already exists (an existing parameter will simply be overwritten by the imported one), and it also allows setting the dg and ap parameters in a simple way (once they have been modified and stored, of course!).

The first argument to the psgset command (a macro) is not only the value that is going to be filled into the seqfil (and pslabel) parameter, it also is the (body of the) file name (without ".par" extension) of the parameter file from which the specified parameters are picked. As expected, a local parameter file takes precedence over a parameter file with the same name, but stored in /vnmr/parlib. The psgset macro can pick up to 11 parameters in one call (not counting the first argument: the maximum number of arguments is 12 with this command). To retrieve more parameters, several psgset calls are required.

Before recalling the parameters and their values (using the rtv command), psgset calls prune to remove any extra parameter in the local set that is not present in the parlib file. This avoids an excessive accumulation of parameters when calling different pulse sequence macros consecutively.

After having "imported" the new parameters plus the display templates, most pulse sequence macros are supposed to also set default values for existing parameters. This can be achieved in two different ways:

- The values can be stored in the parameter set in parlib and recalled with psgset, which will simply replace existing values with those from parlib.
- The values can be set directly, if necessary, with if constructs:
 - if d1=0 then d1=1 endif

The problem with the parlib approach is that the parameter values are well hidden away, and to change the defaults we have to recall the parameter set, change the values, and then overwrite the saved parameter set again. In the author's opinion, it is much easier (and more transparent) to change a macro than to change a parameter set, and hence it is recommended to rather use the first approach: *use* psgset *to recall and set the new and specific parameters and the display templates, but use simple MAGICAL assignments to set existing (standard) parameters.*

There is a special case. Some macros want to set nt based on the number of transients for the previous set-up experiment. In the case of the dept pulse sequence, the macro is supposed to decrease nt by a factor of 32 (taking into account the sensitivity enhancement due to polarization transfer), but ensuring that nt still is a multiple of 4 (the basic phase cycle for that sequence). Taking the current value of nt may not be appropriate, because in the preceding experiment nt may have been set to a very large number, and the actual experiment was perhaps stopped after having reached a satisfactory signal-to-noise ratio. Taking ct, on the other hand, can be wrong also, because any change in an acquisition parameter will set ct to zero! Therefore, it is a good idea to *first* capture the value of ct (note the precaution for arrayed nt!):

```
"dept - convert standard parameter set to dept"
if ct>0 then $nt=ct else $nt=nt[1] endif
psgset('dept','mult','j','tau','presat','dg','ap')
if $nt>128 then nt=$nt/128 nt=4*nt else nt=4 endif
if d1=0 then d1=1 endif
in='n' il='y' gain='y' pw=pw90 hs='nn'
dof=0-2.5*dfrq "move decoupler from 5 to 2.5 ppm"
dm='nny'
if waltz='y' then dmm='ccw' else dmm='ccf' endif
ai
wexp='deptp' array='mult'
man(seqfil) dps
```

Note that for the dept sequence, mult is an normally arrayed parameter, but using psgset and rtv to retrieve such an arrayed parameter does *not* automatically set the array parameter! Also, some parameter settings may have to be set depending on the actual hardware. Constructs like the one used for the dmm parameter (depending on the presence of WALTZ decoupler modulation capability) make a pulse sequence (macro) more portable between different hardware platforms.

Macros for 2D Pulse Sequences

An important additional task in 2D pulse sequence set-up macros is the creation of the 2D specific parameters ni and sw1. This is typically done by the set2d macro utility:

The set2d macro sets the seqfil parameter, then calls par2d(seqfil), a macro that creates the standard 2D acquisition parameters ni, swl and phase. For homonuclear 2D spectra (dn=tn), it sets

sw1=sw ni=sw1/24 rfl1=rfl rfp1=rfp

(12-Hz default digital resolution if f_1). For the heteronuclear case, it calls

psgset(seqfil,'ni','sw1')

or if no argument was specified, it sets ni=256 sw1=2*sw as "default guesses".

Back to the set2d macro. After having called par2d, the set2d macro switches to a fixed gain (i.e., that is something you will not have to do in the set-up macro!), then sets up the acquisition and processing parameters np, fn, ni, fn1 for the specified (or the default) digital resolution if f_1 and f_2 . The default is 6 Hz in f_2 , 12 Hz in f_1 (the desired digital resolutions can be specified in arguments 2 and 3).

The set2d macro then selects sinebell weighting functions in both dimensions for absolute-value spectra (dmg='av'), or a gaussian apodization for phase-sensitive 2D spectra (dmg='ph'). It is therefore important that set2d is only called when the correct display mode has been selected! It is recommended to call set2d prior to the first psgset call, and to always have it preceded by av or ph, to ensure the proper display mode selection, as in the example above.

In a last step, the set2d macro adjusts the display parameters for the 2D spectrum, and it disables the "automation parameters" wbs, wnt, wexp, and werr.

For *heteronuclear 2D* experiments, the set-up macro has the additional tasks of setting up the acquisition parameters for the "other" nucleus (in f_1). It would be nice if these parameters (dof, sw1, rfl1, rfp1) were taken from the equivalent parameters (tof, sw, rf1, rfp) in a second set-up experiment, typically a 1D experiment in an other experiment file. Starting with VNMR 5.1, this can also be an "internal" data set in the same experiment directory. This can be done by adding a lengthy construct like the following to the set-up macro (simplified, for VNMR 5.1):

```
jexp:$experiment
exists(curexp+'/subexp/H1','directory'):$internal
if $internal then
   $enumber=$experiment
else
   $enumber=0
endif
if ($#<1) then
   input('exp# containing 1H spectrum (0=none)? '):$enumber
else
   $enumber=$1
endif
if ($enumber<0) or ($enumber>9) then
```

```
write('error','%s - illegal experiment number',$0)
 return(1)
endif
if ($enumber>0) then
  "get proton parameters"
  if ($enumber=$experiment)
    if not($internal) then
      write('error','%s - no 1H data found in current exp',$0)
      return(1)
    else
      svtmp('tmp') rttmp('H1')
      $savesw=sw $savetof=tof
      $saverfl=rfl $saverfp=rfp $savewp=wp $savesp=sp
      rttmp('tmp')
    endif
  else
    jexp($enumber)
    $savesw=sw $savetof=tof
    $saverfl=rfl $saverfp=rfp $savewp=wp $savesp=sp
    jexp($experiment)
  endif
  sw1=$savesw dof=$savetof
  rfl1=$saverfl rfp1=$saverfp wp1=$savewp sp1=$savesp
else
  sw1=10*dfrg dof=0
  write('line3','a full 10 ppm 1H window will be used')
endif
fn1=sw1/3 ni=fn1/4
```

The released VNMR 5.1 hetcor macro actually uses a much more complex construct! This is, of course, a typical example for a construct that should be put into a utility macro. The only reason why this hasn't happened yet is that within the standard VNMR pulse sequences it was always only the hetcor macro that used such a construct. A number of related pulse sequences (coloc, flock, phase-sensitive hetcor, but also indirect detection sequences like hmqc, etc.) could profit from such a utility.

22.5 Switching Between Similar Sequences

The above macro scheme is acceptable for the way an NMR spectrometer is used in many organic chemistry and industrial environments, where users typically perform a series of set-up experiments (like one for ¹H, one for ¹³C, maybe also extra experiments with reduced spectral window), and then start a limited number of special experiments—mostly just one or two:

- 1 H, < 1 H (reduced sw),> COSY
- 1 H, < 1 H (reduced sw),> NOESY
- 1 H, < 1 H (reduced sw),> 13 C, DEPT, HETCOR
- ¹H, COSY, ¹³C, DEPT, HETCOR

Note that for homonuclear 2D spectra it is not necessary to acquire a spectrum with reduced spectral window. You *first* adjust the chemical shift referencing, then you save the 1D spectrum either in a subfile, using cptmp, or by moving into a different experiment (you will use the 1D spectrum for the 2D plotting). Now you can just place

the cursors around the relevant signals, type movesw, and *then* call the 2D pulse sequence set-up macro.

However, a typical spectroscopist in biological NMR will work in a totally different way. Often a sample stays in the magnet for a week or more, while a series of *similar or closely related* experiments are performed, such as HNCOCA, HNCA, etc. Under such circumstances, the above macro scheme has several disadvantages:

- The complex sequences used in this environment require a large number of parameters that have been carefully calibrated on that particular sample; psgset would simply overwrite these values with the ones from parlib, which is highly undesirable.
- Some of these macros may assume that they are called in a 1D experiment and may not work properly when called on a 2D or 3D data set (e.g., some of these macros cannot be called twice in sequence!).

Users have adopted two philosophies to cope with these deficiencies:

- Rather than calling a pulse sequence macro, they simply recall a parameter set from parlib. The referencing between different biological samples (at least those in one lab) is often almost identical, so this is not a problem. However, this will require *re-entering all calibrated parameters* for each experiment (considering the fact that many of these experiments last for several days, this seems acceptable to many users).
- Assuming all these pulse sequences use convergent parameter definitions, it may be possible to simply change the seqfil parameter to switch to a related sequence. The problem is, that it *cannot* be assumed that parameter sets are consistent between sequences in a general case (unless all sequences involved were developed locally), because there is no compulsory style guide for parameter naming and usage!

Clearly, this is an area which requires improvement. In fact, there are planned improvements:

- A new style of setup macro is being developed that can be used repeatedly and does not overwrite calibrated parameters.
- A detailed style guide covering parameter naming will be put forward in a future version of this manual.

Chapter 22. Role of Macros and Parameters

Chapter 23. Putting It All Together

In this chapter, we look at approaches for programming a new pulse sequence, testing a sequence and its associated macros and files, and submitting a sequence to the Varian user library.

23.1 Starting a New Sequence

To write your own pulse sequence, you have two principal options:

- 1. Look for a pulse sequence (in /vnmr/psglib, or in userlib/psglib, or in any other source) and modify or rework that into your desired pulse sequence. You can do the same thing with the associated parameter sets, macro, phase table, etc.
- 2. Start a new sequence "from scratch."

Programming by Modifying an Existing Pulse Sequence

Definitely, the first approach is very often much easier and is preferred simply because it is less work. It does have some disadvantages, such as:

- It is more susceptible to programming errors.
- It is more difficult to debug because the program is not developed in a "top-down" approach.
- It is much more difficult to do systematic programming and development this way.
- It tends to perpetuate bad syntax and style from existing sequences and old maybe outdated—pulse sequence programming features. The syntax tends to get worse with every transformation of a pulse sequence, unless the programmer is very disciplined.

Programming by the Top-Down Approach

We are not going to further discuss the first approach (there isn't much to say about this method anyway). On the other hand, it certainly is worthwhile giving some hints for the second approach, pulse sequence programming from the "top-down," for which (given the above) some obvious advantages exist:

- It can lead to clean, well-debugged programs.
- The resulting sequence will have a more systematic and simpler syntax.
- It probably uses the latest pulse sequence programming features and, therefore, is more up-to-date to begin with.

Although the top-down approach may be somewhat more work for the programmer, it certainly is the more valuable method from a didactic point-of-view.

The idea of the top-down approach is simple: *don't start with the details*! Start with a rough framework, then start filling in details level by level. Try keeping the syntax complete and compilable all the time, and recompile the program frequently. This way

Chapter 23. Putting It All Together

you never have too many bugs at once, and it is much easier to find the bugs than if you start compilation when the program is several pages long.

The best way to do this is to open up two UNIX windows side by side (the development of a complete, new sequence; don't block the NMR interface by calling an editor from within VNMR). In one window you edit the sequence (using your preferred editor); in the other window you *compile the sequence periodically*, whenever it is appropriate. You can leave these windows running for as long as the programming lasts (or until you want to log out, of course). Just close the windows to icons when you want to use the VNMR user interface.

Periodically save your work—don't give Murphy's law a chance to destroy several hours of working effort!

While you program, always add comments. Don't wait to insert comments until you are finished, because then you will want to work with the new sequence and the comments will never be added. Sensible comments never make a sequence more complicated; they keep it understandable.

Don't write novels into the comments. Use short comments, but add comments to every part of a complicated pulse sequence. Short sequences on the other hand (if they are well-written) should be more or less self-explanatory and usually require very little comment.

As you write the pulse sequence, maintaining a **manual file** for the new pulse sequence is recommended. Often, the manual has been taken from the first part of the comment of the pulse sequence text; it is not necessary to duplicate that text. Either write a separate manual and don't add that text also to the sequence, or just do it in the pulse sequence for simplicity (you can still move that text out of the sequence and into a separate manual file later).

23.2 Testing a Sequence and Related Files

At some point in the programming you start using new variables (delays, pulse widths, power levels, etc.). When this happens, always add the parameter definition and initialization in the pulse sequence (such that it still compiles without error message). At the same time you can start developing the **macro** and creating the **parameter set** for the pulse sequence. Reserve an experiment in VNMR for that new sequence, and successively create (and adjust) new parameters, as you add them to the pulse sequence. Also, you can keep the dg and ap templates up-to-date at the same time. In the macro, very likely you will have to add an argument to a psgset command. You also want to periodically save the macro and the parameter set, even when they aren't finished yet, just to secure your work.

Also, as you start using phase tables, you should of course maintain a **table file** for the sequence.

Continue doing this until you think the sequence is complete. When the sequence then compiles without error, you should complete the macro (of course, you will check with a similar, existing macro to ensure that you don't forget anything in the new macro).

Typing dps in VNMR will now indicate whether the parameter set is complete, and it will also give you a first "clear view" onto the features of the new sequence: are all the pulses and delays there, in the right sequence? When dps works without problems, you

should double-check the dg and ap templates. and then save the parameter set in parlib.

Now either retrieve a simple 1D spectrum or use the VNMR setup command to retrieve standard 1D parameters, and then call the pulse sequence macro and check again with dps. This will tell you whether the macro is complete and functional.

You can go even further with the testing, still without involving a spectrometer. You may want to check the pulse sequence for run-time errors and problems by typing go(`acqi'). This executes the pulse sequence and generates Acodes, at least for the first increment (in case of arrayed or *n*D experiments). It also serves as a test for the parameter checking built into the pulse sequence. If you then want, you could even use the apdecode tool (bin/apdecode from the user library) to have a look at the Acodes that were generated.

All this can be done on a data station. In order to be able to call go('acqi') you must call config and declare the system a spectrometer—but *don't* call setacq on a data station, of course!

The next steps must happen on the spectrometer. Try a 1D run (nt=1) first. Then you may want to *check the phase cycling* by performing an array with

nt=1,2,4,8,16,32,64,128

You may have to adjust the parameter defaults in the macro, as you continue with the testing.

Finally make sure the last versions of all related files are complete and saved:

- Pulse sequence (from psglib).
- Table file in tablib.
- Macro.
- Parameters in parlib.
- Manual file.
- Shape and pattern files, if required.

If the new sequence is functional and works to your satisfaction, don't you want to submit it to the Varian NMR user library and share the sequence with other users?

23.3 Submitting a Pulse Sequence to the User Library

To submit a new sequence to the user library, first fill in a submission form (in the user library the form becomes the README file, myseq.README in this example). You find a blank form in /vnmr/userlib/SUBMISSION. Then you collect all files related to your sequence, except for the compiled version (because this is specific to your software version, it would just take up lots of space in the user library—also, it is very easy and quick to recompile a sequence):

• If you want to send the submission by e-mail, you can directly create a uuencoded, compressed tar file:

Now you can send the submission form and the data to the userlib librarian:

```
cat myseq.README myseq.uu | \
    mail -s `myseq submission' rolf@nmr.varian.ch
```

• If you want to send the submission by floppy, you create a compressed tar file:

Then you store both myseq.README and myseq.tar.Z on a floppy and send it to the librarian: Rolf Kyburz, Varian International AG, Chollerstrasse 38, CH-6303 Zug / Switzerland.

Chapter 24. Syntax Guidelines

It is planned to include a complete style guide for writing pulse sequences in a future version of this manual. For the time being (as long as that style guide hasn't been defined yet), we'll just include some hints for writing pulse sequences such that they are portable to other systems and readable and easily understood by other users.

24.1 General C Syntax

Being disciplined in the use if the C syntax can be a great help in keeping a pulse sequence readable. A few points in more detail:

Comments

Adding comments is definitely recommended, but add reasonable comments: enough to explain what is not obvious, but don't pack too much theory into the comments (you can add a literature reference. Too many comments can also clutter a C program! You can have a larger comment segment at the top of the sequence, but for the actual C code it is better to use one-line comments only.

Indentation

Using proper and systematic indentation helps in two ways: it definitely makes a sequence more readable, and (even more important) it helps avoiding mistakes! The author recommends using the following indentation rules (no examples shown—there are enough examples for C coding in the rest of this manual!):

- Use indentation, but don't indent too much (such as an entire tab stop: this would lead to excessive line lengths in nested if statements etc.). For short sequences (one page or less), two spaces per indentation level should be enough. For longer sequences, it may be better to indent by three or four spaces.
- Always have corresponding opening and closing braces at the same indentation level. Do *not* add the opening brace to the end of an if line!
- Have if and the corresponding else at the same indentation level.
- Where an if or an else branch contains only one function call or statement, braces are not really necessary, and you should leave them away, to keep the coding shorter (although with beginners in C, this is a possible source for errors if additional statements are added to any of the branches.
- If the splitting into status fields is a prominent feature in a pulse sequence (most traditional pulse sequences were written this way), you can use an additional indentation after status calls.
- Additional indentation should also be used for functions between ifzero, elsenz, and endif calls, as well as for statements between loop and endloop, and between starthardloop and endhardloop.
- Place comments either behind a C line (preferably at an identical tab stop) or, for longer comments, at the same indentation level as the surrounding C statements

(starting comments at the beginning of a line would destroy the effect of the indentation).

When dealing with a program that has been coded with a "strange" indentation style, the author often uses the UNIX indent command to fix it up (see man indent for more information). indent can be used with arguments, but it is much easier if it is configured with a local ~/.indent.pro profile. A good set of indent settings would be the following line in ~/.indent.pro:

-bap -bad -bl -ncdb -nce -d0 -di2 -eei \ -nfc1 -i3 -nip -lp -npsl -sc -nsob

Using -in is the principal indentation. Use -i2, -i3, or -i4, as appropriate. The indent program usually does a very nice job, but of course it does not "know" about pulse sequence specifics (i.e., you will have to fix up the result by adding indentations after status, after ifzero and elsenz, after loop, and after starthardloop).

Variables

Try keeping things simple and short! Doing variable initialization within the variable declaration saves many lines of code, *and* it ensures that all variables are initialized! For numeric parameter values that are use only once, you don't need to define a variable at all—you can use the getval function directly as function argument.

24.2 Outdated PSG Utilities

The use of outdated PSG utilities may not affect the functionality of a pulse sequence *on your system*, and *the way you use it*; however, it makes a sequence less portable, in that it may be incompatible with different hardware or with the way an other person uses the spectrometer¹.

Device Addresses

A number of pulse sequence functions require specifying a device (rf channel) address. The device names used in earlier versions of VNMR are TODEV (observe channel), DODEV (decoupler channel), DO2DEV, and DO3DEV. As of VNMR 5.1, things have *changed*! These device addresses still exist, but they now refer to the physical, rf channels (TODEV is channel 1, DODEV is channel 2, etc.) and *don't* reflect the setting of the rfchannel parameter in VNMR (i.e., TODEV, DODEV, etc. will always point to the same physical channel).

If you ever want to use the rfchannel mechanism for re-assigning the *physical* channels to different *logical* channels (and if you don't, some other user of your pulse sequence may!), you should now use the new symbols OBSch (the logical observe channel), DECch (the logical decoupler channel), DEC2ch, DEC3ch instead. Table 17 summarizes the old and new naming conventions. A simple translation will make your pulse sequence compatible with the rfchannel mechanism.

¹ On the other hand, using new utilities may make a sequence partly incompatible with previous versions of VNMR, but as the majority of the VNMR users use the current release, "horizontal" portability (same release, but different platforms) is more important that trying to keep a sequence compatible with older versions of VNMR.

| Old symbol | New symbol |
|------------|------------|
| TODEV | OBSch |
| DODEV | DECch |
| DO2DEV | DEC2ch |
| DO3DEV | DEC3ch |

Table 17. RF channel naming convention

Functions with Device Addresses

The channel naming complication above can be avoided if functions requiring a device address are replaced by functions that addresses that logical channel directly. There are only very few exceptions (namely, some genpulse type function calls) where this isn't possible (or at least not straightforward). Table 18 gives some guidelines on what functions to use:

| Table 18. | Equivalent PSG functi | ions with and without | device address |
|-----------|-----------------------|-----------------------|----------------|
|-----------|-----------------------|-----------------------|----------------|

| Functions with device address | Equivalent functions without device address |
|-------------------------------|--|
| rlpower, power | obspower, decpower, dec2power, dec3power |
| rlpwrf, pwrf | obspwrf, decpwrf, dec2pwrf, dec3pwrf |
| offset | obsoffset, decoffset, dec2offset, dec3offset |
| stepsize | obsstepsize, decstepsize, dec2stepsize, dec3stepsize |
| genqdphase | txphase, decphase, dec2phase, dec3phase |
| gensaphase | xmtrphase, dcplrphase, dcplr2phase, dcplr3phase |
| genpulse | pulse, decpulse |
| genrgpulse | rgpulse, decrgpulse, dec2rgpulse, dec3rgpulse |
| genshaped_pulse | <pre>shaped_pulse, decshaped_pulse, dec2shaped_puls, /dec3shaped_pulse</pre> |
| gen_apshaped_pulse | apshaped_pulse, apshaped_decpulse, apshaped_dec2pulse |
| genspinlock | spinlock, decspinlock, dec2spinlock, dec3spinlock |
| prg_dec_on | obsprgon, decprgon, dec2prgon, dec3prgon |
| prg_dec_off | obsprgoff, decprgoff, dec2prgoff, dec3prgoff |

Note that, strictly speaking, power and rlpower (and also pwrf and rlpwrf) are not equivalent functions, because they require different argument types (real-time integer vs. C double). The next section covers this difference further.

The "gen" type functions are "not officially supported" (i.e., there is no manual page for these utilities). They have been used by people that wanted to implement rf channel independent pulse sequences before the rfchannel mechanism existed. This feature

Chapter 24. Syntax Guidelines

alone added considerable complexity to such sequences and made them very hard to read. Also, the "gen" type functions have the additional disadvantage of not being recognized by the dps utility. All these sequences can now be written for fixed logical rf channels, and the channel switching can be achieved much more easily using the rfchannel parameter.

The only case where "gen" type functions are still recommended is for some particular types of simultaneous pulses (shaped and rectangular). For pulses on "adjacent" logical channels (starting with OBSch), simpulse, sim3pulse, simshaped_pulse and sim3shaped_pulse are perfectly adequate. However, if simultaneous pulses are to be performed on either non-adjacent channels (e.g., OBSch, DEC2ch) or on adjacent channels not including OBSch, it is still better to use "gen" type functions. The simpler "sim" pulse functions don't handle the case the case of any of the pulse lengths being set to zero very well.² Such functions may be acceptable for testing purposes, but for a real experiment you don't want to really use these functions with any of the pulse widths set to zero (at least for the current VNMR release). As mentioned before, using "gen" type functions will unfortunately make these pulses "invisible" to dps.

Replacing power and pwrf Statements

The following construct was used in a large number of pulse sequences:

```
double pwxlvl;
(..)
pwxlvl = getval("pwxlvl");
(..)
initval(pwxlvl,v10);
(..)
power(v10,TODEV);
```

This has numerous disadvantages:

- It is complicated and lengthy;
- It requires reserving one of the few real-time variables (remember that with initval you *must not* use this real-time variable for anything else in the sequence!);
- There are often pages between the initval call and the actual use of the real-time variable containing the power value. Upon reading such sequences, it is often very hard to tell what power level is actually used, because with every power statement you have to go back in the sequence and look for the corresponding initval call.

It is much preferable to simply use

```
obspower(getval("pwxlvl"));
```

which avoids using and initializing a real-time variable and avoids an extra C variable declaration and initialization, and the VNMR parameter name occurs where it is actually used in the sequence.

² In the case of simpulse and sim3pulse, some extra (non-sense) FIFO words are produced for the pulse on the "unused" channel (i.e., it tries performing a timer word of zero length), and for simshaped_pulse and sim3shaped_pulse the waveform generator on the unused channel is still being reset (as if there was a shaped pulse of zero length).

C Constructs for Phase-Sensitive nD NMR

Many pulse sequences use outdated (rather complicated and sometimes even questionable) constructs for f_n coherence selection (States, TPPI, FAD). A preferred set of constructs is described in detail in Chapter 19, "Multidimensional Experiments," on page 215. In particular, it should be noted that the contents of the VNMR phase, phase2 and phase3 parameters are now already extracted and can be used in the *integer* variables phase1 (*not* phase!), phase2, and phase3. Also, the real-time variables id2, id3, and id4 should be used for setting up TPPI or FAD. This not only is simpler than referring the constructs used earlier, it also is safer and easier to understand.

24.3 General Considerations

This sections just provides a few additional hints of more "philosophical" nature.

Multipurpose Sequences

One of the nice features of VNMR pulse sequences is that you can have logical branchings based on VNMR parameters, allowing combining a number of related experiments into a single pulse sequence. But don't overuse this feature by trying to "pack the world" into a single sequence! You can use this feature to combine *closely related* experiments (such as HMQC, HMQC with nulling by partial inversion recovery, and HMBC). Essentially use this capability to turn on or off specific features of a single experiment, and write separate sequences where the combination would be hard to understand or would result in most of the pulse sequence being in an *if* or else branch!

Using dps

The dps command is a useful and powerful tool for quickly controlling the experiment setup before starting a real acquisition. We fully appreciate the value of this command (we use it ourselves all the time!), and we will continue to further develop dps. However, nothing is perfect, and dps has its deficiencies, mainly:

- It does not display "gen"-type pulses and simpulses. "G" type functions (not discussed in this manual) have the same problem.
- It does not display pre-pulse and post-pulse delays.
- It does not display a number of other statements that may be relevant to certain users.
- The statement zgradpulse is not displayed (this has been fixed in VNMR 5.1).
- The number of events that can be displayed is limited and may not be sufficient for very complex experiments (like certain heteronuclear 3D PFG sequences).
- Anumber of other—perhaps minor—problems have been noticed by users.

Most of these problems can be circumvented, and different approaches have been used to achieve this:

 Some users have created a modified version of a pulse sequence just for the purpose of using dps, while the "real" pulse sequence is used for the "real"

Chapter 24. Syntax Guidelines

acquisition only (with the known deficiencies of sim2shaped_pulse and sim3shaped_pulse, this may actually make real sense!).

- Others have resorted to explicit coding of pre-pulse and post-pulse delays by using explicit gating and extra delays around rf pulses, instead of relying on the features built into functions like rgpulse.
- Even extra if and else branchings were built into certain sequences (one branch for dps, one for acquisition).

In all these cases, you should ask yourself: is it worth the effort? Is it worth giving up the readability of a sequence, just to make dps work?

Some of the approaches like the explicit coding of pre-pulse and post-pulse delays may not work because they result in more pulse sequence elements than dps can display! In any case, it will be a long way until dps can display each and every detail (many users may not want that, because it can lead to excessively complex dps displays).

The question then would be: couldn't tools, such as apdecode³ that display the Acode instructions, replace certain aspects of dps (because they *really* show *every detail*)? Displaying the Acode, on the other hand, can be too complex for routine use, but it certainly can be useful for certain stages of pulse sequence debugging!

 $^{^{3}}$ apdecode is available from the user library, see also Section 8.2, "Looking at Acode," on page 69. It is a tool to display the data and instructions in the Acode, as they are executed by the acquisition CPU.

Chapter 25. Debugging a Pulse Sequence

Many of today's pulse sequences have grown into an astounding complexity, and in order to cope with the growing range of experimental demands, modern NMR instruments have become extremely complex machines. Even though Varian tries hard to keep the rf and digital scheme as simple as possible (not because this is cheaper but because in most cases simpler schemes are better and less prone to failures). On the other hand, every programmer makes mistakes¹, and in a complex instrument there is also a certain chance that something fails or does not work as expected.

Overall, there is a fair chance that a new, complex experiment does not work as expected in a first attempt. If this happens, it would be a bad idea to just call a Varian service or customer support person claiming that "my xyz sequence does not work." At the very least one should be able to tell the following:

- What is the instrument type and configuration?
- What was the exact VNMR software release?
- Is this a standard sequence? A sequence from the user library (in this case you should consider contacting the author of the sequence directly, e.g., via e-mail or telephone)? One that you modified yourself?
- What exactly are the symptoms? (Can you fax a plot illustrating the problem or send an FID by e-mail, or mail it on a floppy?)
- Was there any error message?
- What exactly were the parameters in use? (You may want to consider sending a parameter printout by fax or sending an entire parameter set by e-mail or on a floppy!)
- What did you do to further locate the problem? (The sections below should give you some hints on what you could do to track down a problem.)
- Do other (similar or simpler) experiments work without problems? (Very often it is not really clear whether one is dealing with a hardware or a software problem!)
- What was the last experiment that worked properly?
- What happened since then?

Even with all this information, you might spend considerable time on the phone or exchanging e-mail, and it is in the last two points where the user can be most helpful in diagnosing and fixing the problem, saving service time and costs, making e-mail or telephone calls more efficient, and—last but not least—speeding up the entire process of getting the experiment to work.

Given the complexity of modern experiments, this may not always be an easy task. It may require a fair amount of analytical and creative thinking to locate and fix a particular problem, and certainly knowledge about the functionality of the experiment, as well as about the internal functionality of the NMR instrument is most helpful for this task. Of course, it takes time to acquire the knowledge that is necessary to debug a pulse sequence. One motivation for creating this manual was to shorten the learning curve for advanced spectrometer operation and debugging. With the right approach, you should not only be faster in designing and debugging new pulse sequences, but you

¹ Many people go as far as claiming that there is no software without bugs!

will also be able to save service time and costs by providing a more accurate description of eventual problems.

It is difficult, if not impossible, to provide a general recipe fro debugging a pulse sequence or locating an instrumental problem. The tips below are just options. The "correct" choice and sequence of actions is depending on the nature of the problem.

25.1 Debugging the Parameters

First, you want to find out whether the problem is due to a simple parameter missetting, because this would be the easiest and fastest problem to solve.

- Check dg and dgs for obvious parameter errors.
- Check with dps whether the parameters really are what you think they are. You may have used a pulse type parameter for a delay or, much worse, a delay parameter for a pulse duration (this can burn your probe and damage your sample!).
- Important parameter may be hidden from the display. Check for the values of all parameters on which you do a getval or a getstr in the pulse sequence.

Parameter mis-settings may result from errors or omissions in the setup macro—you may want to check for problems there. Certainly, if there is a parameter problem, this should be fixed in the setup macro, to avoid the same problem next time.

25.2 Debugging the Software

Next, you want to check whether there *is* a software problem, either a faulty sequence or a malfunction in the acquisition software.

- You also may want to kill Acqproc, restart it, and retry the experiment, just to check whether the problem is reproducible. Power surges, vibrations or temperature variations can severely affect cancellation experiments.
- Analyzing the FID or the spectrum can often indicate the nature of a problem (lack of cancellation, artifacts, etc.).
- Use dps to verify that the sequence corresponds to what you intended it to do.

If you don't see any obvious source for the problem, try narrowing down the area in which the problem is located:

- Simplify the experiment either by setting selective parameters to zero or by changing flag parameters (if there are any).
- Check for phase cycling errors by using nt=1,2,4,8,16,32,64,128. In cancellation experiments, you should see whether you get cancellation, and how much. Apart from the cancellation, you should see a steady increase of the signal-to-noise ratio by a factor of 1.4 with every increment. For 2D experiments, use ni=0 (measure the first increment only)². For multiple-quantum filtered experiments (or *nD* experiments that acquire echo-type signals in general), it is also necessary to do the same test with d2=ni/(2*sw1) to check whether

 $^{^{2}}$ In non-PFG absolute-value 2D experiments, there will be one step in the nt array that shows no increase in signal-to-noise. This is the point where f₁ quadrature detection is achieved by subtracting (N+P) from (N-P) type spectra.

multiple-quantum signals are collected at all (the first increment should contain only noise if the experiment is functioning properly). For States-type phase-sensitive nD experiments, use phase=1 (and/or phase2=1).

- Use apdecode (supplied in the VNMR user library) to verify that the sequence of events in the pulse sequence is what you want, including events that are not shown by the dps command.
- Make the pulse sequence print parameter values and variables, maybe at several points in the sequence. This can reveal misconceptions such as the use of real-time variables in C calculations or any attempts to extract values from a displ.ay or processing parameter (only acquisition parameters can be accessed in a pulse sequence).
- It is easy to get mixed up in C decisions (or complex C constructs, in general). If the execution path within a complex sequence is not clear from looking at the dps or apdecode output, you may want to check points at certain places in the sequence, such as the line:

printf("starting refocusing period\n");

• If all this doesn't help, simplify the pulse sequence code by progressively removing or commenting out sections of the pulse sequence (of course, you should keep a copy of the original sequence!) to try locating the problem within the sequence.

25.3 Debugging the Hardware

- WARNING: High voltages are present inside the console!
- WARNING: Always have all cables properly connected or terminated when performing experiments, or while rf output is being generated (you could damage expensive power amplifiers!). For analyzing the output of power amplifiers (in particular for the 1 KW amplifiers used in solids NMR) you will need special power attenuators.

WARNING: Don't pull boards while the power is switched on, and never pull boards without protection against static electricity!

The one important rule for debugging hardware problems is to use *simple* tools (pulse sequences and experiments)! Complex pulse sequences may be a good measure for testing the overall performance of a spectrometer, but they are useless for locating a specific hardware problem in an NMR spectrometer.

- Is the probe tuned properly? *Can* it be tuned at all, and on all channels?
- Has the proper quarter-wavelength cable been selected?
- Is there no problem with the lock?
- Try a similar, but different experiments (a different pulse sequence), just to doublecheck whether you have a hardware problem.
- If that fails as well, try very simple experiments, like s2pul, d2pul, sh2pul, g2pul or the like, to further narrow down the possibilities.
- Double-check connections between the console, magnet leg, preamplifiers, and probe.

- Check fuses in the console.
- Periodically check the fans in the console: overheated boards can easily cause a malfunction of the spectrometer!

Further hardware troubleshooting requires hardware diagnostics tools such as an oscilloscope, BNC cables and a set of rf attenuators. You may want to seek the help of a local expert for accessing cables and connectors within the console and for analyzing power rf output.

- Check whether rf (pulses, decoupler frequency) is generated and gated properly at the output of the transmitter board.
- Then follow the signal path and check the input of the amplifiers (after the attenuators).
- Given the appropriate attenuators (**observe the warnings above!**), you may now check the power output at the amplifier outlets.
- To check the rf path in the magnet leg, you can then repeat this test at the input to the probe.
- To check for the signal path *from* the probe, you need a strong sample; otherwise, you just observe noise. First, check for the presence of the L.O. frequency at the preamplifier (without L.O. you will not even observe noise behind the mixer!). The signal at the input to the preamplifier may be difficult or impossible to see on the oscilloscope, because the signal level is very low at this point. You may have a better chance at the input to the receiver.
- Carefully analyze the FID: Do the signals look normal? Is there truncation? Do you see unnatural "spikes"? Do both channels have the same or a similar amplitude? A good and easy way to check for signal overload at the ADC is the command ddff(1) that shows the numeric contents of the FID.

Symbols

+= notation, 118 .DEC file extension, 179, 180 .RF file extension, 173 .rootmenu file, 63

Numerics

180-degree phase shift, 4163-dB attenuator, 13990-degree phase shifts, 41, 45, 163

A

aa command, 60 abort function, 155 aborting a sequence, 155 Acode, 57, 66 AP bus words, 140 contents, 71 data, 15 decoder, 59, 70 fast bits, 94 file header, 74 file structure, 69 from waveform generator modulation, 183 generation, 58 instruction section, 79 interpretation, 66, 82 interpreter, 149 loops, 147 phase calculations, 98 real-time decisions, 160 saving space, 117 segment, 60, 61 space, 111 space efficiency, 123 structure, 74 tables in instructions segment, 122 timer words, 94 acoustic ringing, 48 acq directory, 59 acqi command, 58, 59 acqi.Code file, 69 acqi.RF file, 170 acqparms.h file, 20, 156 acqpresent file, 59 acqproc account, 63 Acqproc error messages, 56 Acqproc not active, 56 Acqproc process, 15, 55, 59, 61, 63, 167, 266 acqqueue directory, 15, 15, 57, 57, 62, 69, 167 acqtriggers variable, 205 acquire statement, 151, 205, 206 acquisition explicit, 208 implicit, 205 lock, 62 multi-FID, 209 parameters, 20, 55 process, 55, 59 queue, 56

queue directory, 57 queuing, 60 status window, 62 trigger, 205 acquisition bootup selector switch, 146 acquisition control boards, 151 acquisition control parameters, 57 acquisition CPU, 15, 59, 60, 61, 66, 85, 145, 172 bus, 67 data blocks, 69 acquisition operating system, 59, 66, 69 acquisition process, 15 actively shielded gradient coils, 202 ADC (Analog-to-Digital Converter), 67, 207 ADC overflow, 214 alfa parameter, 51, 51, 52, 206 algorithm for phases, 95 algorithms for phase cycling, 124 aliasing, 205 alock parameter, 58 alternating phase cycles, 102 AM/PM transmitter board, 164 amplifier blanking, 45 amplifier stabilization, 47 amplifier timing, 41 amplitude modulation circuitry, 163 amplitude multiplier for shaped gradients, 170 analog-to-digital converter (ADC), 16 angled bracket notation, 32 angled brackets convention, 14 anti-parallel switching, 44 AP bus, 67, 67, 86, 137, 163 addressing, 163 delays, 141 indirect mode, 86 shaping gradients, 202 statements, 21 traffic, 140, 200 AP bus chip, 137 AP interface board, 88, 137, 140 AP words, 86, 137 apdecode program, 71, 75, 140, 267 apdelay.h file, 21, 142 aph command, 52 aph parameter, 50 apovrride statement, 139 apshaped_pulse statement, 191, 194, 195 aptable.h file, 21 argument type-checking, 23 arraydim parameter, 74 artifact cancellation, 95 artifacts, relative intensity, 135 as command, 25 assembly language compiler, 25 attenuator, 211 au command, 55 au('wait') command, 56 audio filter, 67 audio filter bandwidth, 139 audio frequencies, 50 audio signal, 211 AUTOD data structure, 74, 75, 76, 78 autoincrement attribute, 118, 119 autoincrementing tables, 119, 151

autolocking, 66 automated flag, 156 automatic statement, 39 automation control board, 66, 66, 78, 145 autophasing, 50, 52 autoshimming, 59, 66, 78, 80 autshm.out file, 59

В

background VNMR, 62 base counter for phase cycling, 110 baseline roll, 51 bc command, 51 BCD information, 137 Bessel filters, 51 beta constant, 51, 206 binary coded decimal (BCD) information, 137 biological NMR, 52 BIRD pulses, 122, 156, 158 blanking the amplifier, 44 blocksize transients, 96 boot PROMs, 66 bootup selector switch, 67, 214 BR-24 sequence, 150 braces notation, 117 brackets notation, 101, 117 broadband rf, 156 broadband-type transmitter boards, 165 bs parameter, 77 bsctr variable, 96 bsval constant, 77, 96 buffered acquisitio, 61 bug fixes, 29 bus decoder, 86 bus structures, 67 Butterworth filters, 51

С

C based decisions, 155 C compiler and linker, 18 C constructs, 95 C language errors, 18 C loop, 147 C preprocessor, 18, 19, 21, 22, 23 calculating complex phase cycles, 162 calfa macro, 52 cancellation experiments, 109, 266 cancellation of artifacts, 95 cancellation quality, 134 cc command, 18, 19, 24 cccc.c file, 129 CD-ROM drive, 59 celem parameter, 62 center glitch, 95, 212 change macro, 61 channel imbalance, 95 chown command, 63 CIDNP experiments, 88 class A/B linear amplifiers, 45 class C amplifiers, 45, 139, 156

COCONOESY sequence, 206 code optimizer option, 25 code section for pulse sequences, 106 code segments in Acode file, 74 codeint variable, 75 coherence pathway selection, 95, 199 coherent signal buildup, 206 coil inductance, 47, 48 coil mechanical movements, 48 combined COSY and NOESY pulse sequence, 209 combining two pulse sequences, 158 command to convert (CTC) bit, 87, 89, 205 comments for phase cycles, 101 compilation, 18, 24 conditional, 22 pulse sequences, 15, 15 complex phase cycle generation, 103 composite pulse inadequate, 129 composite pulses, 50 conditional compilation, 22 conditional processing, 62 config program, 60 configuration parameters, 15 conpar directory, 15 constants definition, 20 cp command, 31 cp parameter, 58, 76 cpp program, 18, 19, 23 cps.c file, 189, 205 CPU address space, 69 CPU boards, 66 CRAMPS experiments, 150 createPS function, 205 ct counter, 214 ct variable, 62, 75, 95, 96, 102, 109, 118, 119 curly brackets notation, 117 cycle phase flag, 76 CYCLOPS phase cycling, 130

D

d1-randomization, 135 data acquisition, 205 date inconsistencies, 31 date parameter, 63 dc offset, 51 dc offset cancellation, 212 dead times, 37, 41 debuggers, 25 debugging a pulse sequence, 265 dec2apshaped_pulse statement, 191 dec2prgoff statement, 181 dec2prgon statment, 180 dec2rgpulse statement, 48 dec2shaped_pulse statement, 172 dec2spinlock statement, 181 dec3prgoff statement, 181 dec3prgon statement, 180 dec3rgpulse statement, 48 dec3shaped_pulse statement, 172 dec3spinlock statement, 181 decapshaped_pulse statement, 191 DECch device, 40, 48, 172

decompose a complex phase cycle, 103 decoupler amplifier blanking control, 52 decoupler gating, 49 decoupler high-power level, 139 decoupler low-power attenuator, 139 decoupler modulation mode, 53, 88 decoupler modulator, 138 decoupler modulator frequency., 139 decoupler waveform generator, 171 decoupling during acquisition, 159 decprgoff statement, 181 decprgon statement, 180 decpulse statement, 21, 48 decrementing phase cycles, 102 decrgpulse statement, 48, 157 decshaped_pulse statement, 172 decspinlock statement, 181 degree of randomization, 112 delay, 37 delay constants, 21 delay for AP bus traffic, 141 delay statement, 37, 150, 235 delayer statement, 37 dephasing amount, 201 depth of the FIFO, 86 device addresses, 20 dg command, 266 dgs command, 266 diagnosing experiment problems, 265 diagnostics terminal, 66, 66, 146, 214 diffusion experiments, 202 digital components, 68 diode switches, 43 direct addressing mode, 138 direct binary information, 137 direct gating, 52 direct synthesis boards, 165 direct synthesis rf, 139, 156 displaying warning message, 156 division factor, 117, 118, 119 division return factor, 117 DLINT option, 19, 24 dm parameter, 178 dmf parameter, 178, 185 dmf, dmf2, dmf3 parameters, 179 dmm parameter, 53, 88, 140, 178 dmm, dmm2, dmm3 parameters, 179 double quotes notation, 32 double-precision flag, 76 double-precision timer words, 91 double-quantum filtered COSY, 106 dp parameter, 58, 76 dps command, 17, 18, 70, 266 dps_ps_gen program, 17 dpsdata file, 18 dps-modified file, 17 dres parameter, 180, 185, 187 dres, dres2, dres3 parameters, 189 dseq, dseq2, dseq3 parameters, 179 dumb terminals, 146 dwell time, 208 dwell_time argument, 205 dynamic binding, 26 dynamic run-time linking, 26

Ε

E.COSY experiment, 106 eddy currents, 199, 200, 202 elsenz statement, 162 enabling statements, 156 endhardloop statement, 150, 154, 236 endif statement, 161 endloop statement, 148, 149 errmsg file, 18, 24, 24 error detection, 22 error message file, 18, 24, 24 error messages, 22, 146 EVENT1_TWRD instructions, 93 EVENT2 TWRD instructions, 93 evolution phase, 110 evolution time, 143, 157 excitation point, 52 excitation pulse data, 48 execkillacqproc shell script, 63 executable file, 18 executable target file, 25 EXORCYCLE phase cycling, 131 experiment startup, 61 explicit 90-degree phase shifting, 52 explicit acquisition, 208 explicit modulation, 149 expn.username.PID file, 62 expn.username.PID.Code file, 69 expn.username.PID.RF file, 61, 167, 167, 170external declarations, 20 external definitions, 20 external phase tables, 15, 57, 116 external trigger, 88 externally defined code addresses, 20

F

fall-through time, 151 FALSE constant, 20 fast bits, 20, 76, 86, 87, 88, 160 assignment, 88, 89 speed, 137 fast status lines, 163, 168 fast switching lines, 68, 88 fb parameter, 51 FID data block, 69 FID file, 15 fid file, 62 fidpath parameter, 63, 63, 63, 63, 63 field gradient shapes, 196 FIFO (first-in-first-out) buffer, 67, 85 FIFO full message, 86 FIFO type, 59 FIFO underflow, 82 FIFO underflow message, 147 FIFO width, 87 FIFO words, 82, 86, 138, 138, 151 fifolpsize flag, 156 filter delay, 50, 52, 81 fine attenuators, 139, 190 first-order phase correction, 51 fixed frequency rf, 156

fixed frequency transmitter boards, 165 flag fields, 53 flag handling, 39 flag variable generation, 162 flags feature, 156 flipback experiment, 208 flow-through (FIFO) buffer, 165 fm-fm modulation, 178 forced relaxation, 208 foreground VNMR, 62 fuses, 268

G

G_Delay statement, 20, 37 G Power statement, 21 G_Pulse statement, 20, 21, 40, 49 G_Simpulse function, 50 ga command, 55 gain parameter, 58, 62, 214 GARP-1 decoupling, 178 GARP-1 modulation pattern, 187 gate statement, 89 gate switching, 45 gating directly, 52 gating information, 68, 76 gauss.RF file, 174 Gaussian distribution, 135 Gaussian pulse, 174 gen_apshaped_pulse statement, 191 general (object-oriented) routines, 20 general power statement, 21 generic shapes, 165 genpulse statement, 40, 48 genshaped_pulse statement, 172, 177 gensim2pulse function, 50 gensim3pulse function, 50 gensim3shaped_pulse statement, 173 genspinlock function, 181 getelem statement, 115, 118, 119, 122, 122 getval statement, 55 gf macro, 58, 58, 70, 74 Gilbert multipliers, 46 glitch peak, 206 go command, 15, 55, 170 go('acqi') command, 58, 58 go('debug') command, 71 gradient amplifier, 199 gradient amplitude, 199 gradient coil, 200 gradient control unit, 202 gradient identifier, 199 gradient MQCOSY experiment, 201 gradient pulse shape, 200 gradient recovery delay, 201 granularity, 149 grise time, 201 gstab time, 201

Η

HAL (Host-to-Acquisition Link), 15, 66, 207

hard abort, 181 hard-coded software, 69 hardware diagnostics tools, 268 hardware loops, 86, 92, 150, 151 hardware modulator, 178 hardware troubleshooting, 267 HDLOOP instruction, 154 hetcor sequence, 157 heteronuclear correlation experiment, 158 hidden AP delays, 141 hierarchically stacking, 161 high-level statement, 39 high-power amplifiers, 47 high-power decoupling, 159 HighSpeedLINES instruction, 94 HMQC pulse sequence, 156, 157 hmult parameter, 157 homospoil flag, 38 homospoil pulse, 39, 204 host computer, 59, 66 host computer peripherals, 67 host I/O bus, 145 Host-to-Acquisition Link (HAL), 15, 59, 66 hoult command, 52 Hoult delay, 51 housekeeping delay, 146, 206 HouseKEEPing instruction, 214 HS (High-Speed) bus, 68 hs parameter, 38 hsdelay statement, 38, 39, 39 HSgate statement, 39, 39 hsine-shaped gradient, 196 HSLINES instruction, 147 hst parameter, 38, 39 hwlooping.c file, 206 hypercomplex spectra, 110

I

I/O bus, 67 I/O functions, 20 IB_DELAYTB instruction, 169, 170 IB_LOOPEND instruction, 170 IB_PATTB instruction, 170, 176, 185 IB_SCALE instruction, 170, 196 IB_SEQEND instruction, 170, 171 IB_START instruction, 169, 171, 172 IB_STOP instruction, 169, 170, 171 IB_WAITHS instruction, 170, 172, 175 id parameter, 63 idelay statement, 39 ifnotzero statement, 162 IFNotZeroFUNC instruction, 161, 162 ifzero statement, 161 imaging experiments, 195, 200 imaging gradient amplitudes, 139 implicit acquisition, 81, 205 implicit decisions, 159 implicit gating statement, 53 improper refocusing, 95 inadqt.c file, 98, 124 incdelay statement, 39 include files, 21

include notation, 32 include statements, 19 incrementing phase cycles, 102 index of current FID, 76 index to phase table, 116 indirect detection, 189 inductance of rf coil. 48 initval statement, 77, 97, 98, 147 inline tables, 120 instruction block for RF pattern, 169 instrument problems, 265 interactive acquisition program, 59 interactive FID or spectrum mode, 58 interactive statements, 39 interleave parameter, 62 intermediate frequency (I.F.), 46, 210 intermodulation distortions, 212 internal periodicities, 103 internal propagation delay, 166 inverse spectral window, 205

J

J-coupling, 95

Κ

knobs information, 78

L

laser control, 88 LC data structure, 75 lc.h file, 78 ld program, 25 ldcontrol file, 57, 61, 167 lib-lc.lnfile, 23 libparam.a file, 26, 29 libpsglib.a directory, 26, 29, 30 libpsglib.so.1.0 file, 31 libpsglib.so.x.yfile, 56 linear amplifiers, 44, 138, 156, 190 stabilization, 45 linear amplitude modulator, 140, 190 linear ramp, 203 link editor, 25 link loader, 26 lint program, 17, 18, 19, 22, 22 lintfile.c file, 23, 35 llib-lpsg file, 29 llib-lpsg.lnfile, 23, 24 load parameter, 58, 78 local oscillator (L.O.), 163, 210 lock command, 61 lock file, 61, 61 lock parameters, 145 lock power, gain, and phase, 140 lock_n.primary file, 61, 61 log file, 62 look-up offsets, 47 loop count, 149

loop counter, 149 loop cycle, 208 loop FIFO, 87, 151, 156, 208 loop statement, 148 loops, 147 low-core data structure, 74, 75 low-level statement, 39 lp parameter, 51, 52

М

maclib directory, 15 macro for pulse sequence, 15 macros, 20, 21 macros.h file, 20, 21, 23, 172 macroscopic phase cycling, 130, 131 magnet leg pneumatics, 66, 145 make file, 18, 30, 33 make program, 18, 24, 30 makefid command, 213 makesuacqproc shell script, 63 makeuser command, 56 makeuserpsg file, 30, 33, 34 Masterlog file, 62 math libraries, 26, 33 math.h file, 33 mathematical algorithm for phases, 95 mbond parameter, 157 mechanical movements in the coil, 48 method parameter, 56 mirror images, 211 mixer, 211 mixing product, 164 mixing time variabtion, 112 mixvar scaling factor, 112 MLEV-16 decoupling, 178 MLEV-16 spinlocking, 186 mod2 statement, 162 mod4 statement, 162 modn statement, 162 modulation algorithms, 190 modulation files, 180 modulo function, 97 modulo statements, 162 Motorola 68000 CPU chip, 66 MREV-8 sequence, 150 multiecho imaging experiments, 146, 206 multi-FID experiments, 206, 209 multifield flags, 159 multiple hardware loops, 151 multiple-quantum filtering experiments, 199, 201, 266 multiplexer switch, 86 multipulse experiments, 150 multipulse line narrowing, 208

Ν

nesting shorthand notation, 117 newdec flag, 156 newdecamp flag, 156 newtrans flag, 156

newtransamp flag, 156 NextSCan keyword, 81 nf parameter, 206 ni parameter, 266 NOESY experiments, 110, 112 noise measurement, 80 np parameter, 75, 206, 207, 209 NSC (next scan) instruction, 206 nt parameter, 75, 266 NUMch flag, 156 numeric constants, 96

0

object archive files, 31 object dependencies, 31, 33 object file, 27 object library, 31 object library creation, 29 object module, 25 OBSch device, 40, 40, 48, 172 observe channel pulses, 40, 40 obsprgoff statement, 181 obsprgon statement, 180 obspulse statement, 40 od command, 69, 170 off-resonance effects, 95 offset statement, 138 offset synthesizers, 138, 139 one constant, 96, 97 oopc.h file, 20, 20 oph variable, 119 oscilloscope, 268 output boards, 151, 207 overflow, 67 Oxford VTC-4, 66 Oxford vttype flag, 156

Ρ

plpat parameter, 173, 175 pad parameter, 80 parameter debugging, 266 parameter tree, 55 parameters for pulse sequences, 15 parentheses notation, 101, 116 parlib directory, 15 passive diodes, 43 pattern time base, 170 pfgon parameter, 199 phase alternation, 212 phase calculations, 95, 98 phase changes, 150 phase cycles alternating, 102 creating phase math, 101 decrementing, 102 evaluation, 9 generating complex cycles, 103 incrementing, 102 length, 123 real-time construction, 106

shifted patterns, 103 shorthand syntax, 101 simple, 101 phase cycling, 95 base counter, 110 during steady-state pulses, 109 errors, 266 order, 135 real-time phase calculations, 111 refocusing periods, 110 phase generator and frequency divider, 211 phase mathematics, 95 phase modulation circuitry, 163 phase modulator, 163 phase parameter, 110, 267 phase shifting, 46 phase shifting times, 45 phase table names, 191 phase tables, 15, 115, 129 phase-pulse techniqu, 214 phase-sensitive NMR, 110 phase-settling delay, 47 phaseshift statement, 213 PIN diodes, 43, 190 pointer generation and incrementation, 118 polling rate, 61 post-pulse delay, 50, 51 power level storage, 98 power.h file, 21 preacquisition delay, 80 preamplifier, 210 preamplifier saturation, 48 precompiled files, 15 precompiled modules, 24 precompiled objects, 32 predefined variables, 96 preloop FIFO, 86, 87 preprocessor, 19 pre-pulse delay, 50 prg_dec_off function, 181 prg dec on function, 180 primary lock file, 61 probe, 164 PFG, 201 preamplifier, 210 ring-down, 48, 48 programmable audio filters, 211 programmable pulse modulator, 163 programmed decoupling, 189 programmed modulation, 177 proton flipback sequence, 209 proton frequency, 156 psg directory, 15, 15, 19, 26, 29, 31, 32, 56, 70 psggen command, 15, 29, 34 psglib directory, 15 PTS frequency synthesizer, 79, 138, 139, 163, 214 pulse break-through, 47 pulse length, 41 pulse programmer, 15, 37, 53, 66, 67, 82, 163 AP words, 138 compared to waveform generator, 165 Gemini, 151 timers, 89 pulse sequence control boards, 151, 206

pulse sequences adding acquisition parameters, 55 C errors, 18 code section, 106 compilation libraries, 29 compilation time, 32 compiling, 15 creation, 17 dead times, 37, 41 executable, 57 execution, 15, 56 generation, 15 macros, 15, 21 parameters, 15 purpose, 53 source coded, 29 troubleshooting, 265 pulse shape definition file, 173 pulse shaping speed, 166 pulse statement, 40, 49, 119, 161 pulse turn-off time, 47 pulse turn-on time, 47 pulsed field gradients (PFG), 139, 199 pulsed NMR experiments, 95 pulsesequence function, 21 pwpat parameter, 173, 175

Q

quad image suppression, 130 quadrature detection, 211 quadrature image suppression, 212 quadrature images, 95 quadrature phase shifting, 45, 103 quadrature phases, 206, 212 quality factor Q, 48 quarter-wavelength cable, 43, 43, 267 queue files, 62 queuing of acquisitions, 60 quiescent states, 89

R

rampgrad function, 202 random function, 112 random number generator, 113 random number table, 112, 135 random seed, 113 random variation of variables, 112 rc.local file, 59 rc.vnmr file, 59 rcvroff statement, 41, 41 rcvroff_flag variable, 41 read switch, 42 real-time branching, 107 real-time calculations vs. tables, 123 real-time calculations with phase tables, 129 real-time decisions, 160 real-time loops, 148 real-time math in hardware loops, 151 real-time math operands, 96 real-time math operators, 96

real-time numeric constants, 97 real-time phase calculations, 111 real-time phase math, 106 real-time pseudo-random number generator, 112 real-time variables, 96 real-time variables vs. table math, 121 receiver. 211 receiver board, 67 receiver gain, 78, 80, 140, 145 receiver gate, 45, 47, 52 receiver gating time, 40, 40 receiver phase, 115 receiver phase shifting, 212 receiver timing, 41 refocusing period, 160 refocusing pulses, 110 relay parameter, 110 relayed COSY experiment, 110 relayh pulse sequence, 147 relayh sequence, 110 relocation bits, 25 reverse synthesis, 103 revision check, 56 revision number, 26 rf channels, 156 rf channels indices, 20 rf gates, 47 rf power amplifier, 164 rf power attenuators, 139 rf ring-down, 48 rfchan_device.cfile, 119 rfchannel parameter, 40, 50, 173 rfconst.h file, 20, 89 rgpulse statement, 23, 40, 45, 49, 150, 161, 235 rgradient statement, 199 rhmon.out file, 59 rofl parameter, 40, 45 rof2 parameter, 40, 51, 52 root ownership of Acqproc, 63 round-off timing error, 91 RS-232 ports, 145 run-time linker, 56 run-time linking, 26, 29

S

s2pul pulse sequence, 21, 159 s2pul.ifile, 23 sa command, 60 sample changer, 66, 78, 145 sample heating, 159 sample macro, 61 sample-and-hold circuitry, 67 scalar math on tables, 121 scan number, 95 SCCSid string, 24 SCSI, 15 SCSI bus, 15, 59, 67 SCSI interface, 66 second FIFO, 86 secondary lock file, 61 send2Vnmr call, 62 seqfil parameter, 55

seggen command, 15, 17 seqgenmake file, 18, 24, 35 seqlib directory, 15, 18, 55, 55, 56 sequence of events, 157 serial ports, 66, 66 setacq program, 59, 60 setautoincrement statement, 119, 120 setdivnfactor statement, 118, 120 SETICM (set input card mode) instruction, 212 SETPHAS90 instruction, 94 setquadphase statement, 119 setreceiver statement, 119, 122 settable statement, 120 setuserpsg command, 29 SETVT instruction, 80 sh2pul pulse sequence, 173 shape_pulse.c file, 32, 147, 190 shaped gradients, 202 shaped pulse statements, 172 shaped_pulse statement, 172, 173, 177 shapelib directory, 15, 57, 168, 172, 179, 197 shaping field gradients, 195 shared object libraries, 29 shared objects, 26, 26 shell script, 17 shifted phase cycle patterns, 103 shim coils used for PFG, 204 shim command, 61 shim DAC values, 139 shim methods, 56 shimmethods directory, 56 shorthand notation for phase tables files, 116 shorthand syntax disadvantages, 127 shorthand syntax for phase cycles, 101 shorthand table, 117 signal measurement, 210 signal-to-noise ratio, 50 sim3pulse statement, 49 sim3shaped_pulse statement, 172 sim4pulse statement, 49, 50 simple delays, 37 simpulse statement, 49, 49, 123, 157, 158 simshaped_pulse statement, 172 simultaneous pulses, 49 simultaneous shaped pulses, 172 single-point acquisition, 208 single-precision timer words, 91 slew rate, 200 slice duration unit, 168 slice phase, 180 small-angle phase shifts, 45, 46, 139, 163 small-angle receiver phase shifting, 213 software debugging, 266 software loops, 148 solid-state NMR, 47 source code, 29 spare gates, 52 spare lines, 94 spatial dimension, 196 spatial encoding, 196 spectral artifacts, 199 spectrometer hardware, 155 spin command, 61 spin diffusion, 201

spinlock statement, 181, 183 spinlocking pulse train, 182 spinner control circuitry, 66, 145 spinner speed regulation, 66 square-wave modulation, 178 srandom function, 112 ss parameter, 58, 76, 109 ssctr real-time variable, 109 ssctr variable, 77, 96 ssval constant, 76, 96 stand-alone data stations, 56 standard include file, 32 standard.h file, 19, 32 standing wave reflection, 43 STartFIFO instruction, 82 starthardloop statement, 150, 235 static buffer, 165 static linking, 25, 26 static RAM, 165 status field constants, 20 status registers, 82 status statement, 53, 159, 160, 178, 178 status-field controlled modulation, 189 statusindex variable, 38 statusindx variable, 53 stdio.h file, 19, 20 steady-state phase cycling, 109 steady-state transients, 96 STM (Sum-to-Memory), 16, 67 STM counter, 207 StopFIFO instruction, 83 strobe command, 138 su macro, 61 submitting experiment to acquisition, 55 SUID protection bit, 63 Sum-to-Memory (STM), 16, 211 SunOS 4, 26 sw parameter, 51, 206 swept square-wave modulation, 178 switching time, 44 symbol table, 25 symbolic links, 30, 31 syntax check, 19, 22 system configuration parameters, 55 systemdir parameter, 63

T

T/R switch, 190 table addresses t1-t60, 121 table contents, 119 table index, 97, 118, 119 table math vs. real-time variables, 121 tablename, 122 tablert variable, 119 tables constants definition, 21 tables vs real-time calculations, 123 table-to-table math operation, 117 tablib directory, 15, 57, 115 target operand, 96 Televideo terminals, 146 terminal for diagnostics, 146 test4acquire function, 206

testing flag fields, 160 third rf channel, 88 three constant, 96, 97 time base, 86 time count, 86, 89 time event, 37 timers, 89 timing resolution, 87 tip angle, 180 tip-angle resolution, 189 TOCSY pulse sequence, 181 TODEV device, 40 touch command, 34 TPPI (time-proportional phase incrementation), 110 transient counter, 95, 96 transmit switch, 42 transmitter board, 68 transmitter digital control board, 163 transmitter frequency, 211 transmitter gates, 47, 52 transmitter timing, 41 transversal relaxation, 201 trapezoidal gradients, 202 traymax parameter, 145 TRUE constant, 20 tsadd statement, 121 tsdiv statement, 121 tsmult statement, 121 tssub statement, 121 T-switch, 42 ttadd statement, 117, 121 ttdiv statement, 117, 121 ttmult statement, 117, 121 ttsub statement, 117, 121 TUNE_FREQ instructions, 79 twisted ring shift register, 45 two constant, 96, 97 txphase statement, 119, 119

U

unary operators, 96 unblanking time, 45, 45 UNIX kernel, 59 user library, 71 userdir parameter, 63

types.h file, 20

V

v1, v2, . . .v14 real-time variables, 96 variable phase cycle, 110 vdelay statement, 39 vector math with tables, 121 VERSAbus, 67 vgradient statement, 199 vmunix file, 59 VT controller, 66, 80, 145 vttype flag, 156 vttype parameter, 145

W

WALTZ decoupling, 148, 149, 152, 171 WALTZ-16 decoupling, 185 WALTZ-16 modulation, 178 warning message, 156 waveform generator, 68, 88, 163 compared to pulse programmer, 165 data, 61, 170 explicit modulation, 180 gate settings, 180 initialization, 79 shapes, 57, 139 timer, 169 waveform generator board, 165 wbs parameter, 62 werr parameter, 62 wexp parameter, 62 wq.c file, 171, 183 WG3 instruction, 176 WGCMD instruction, 176, 177 wgdecode program, 171 while loop, 147 wnt parameter, 62 writefid command), 213 wshim parameter, 58, 78, 80

X

x_file prefix, 18
xmtron statement, 181
xmtrphase statement, 214
xr.conf file, 60
xr.out file, 60
xrop.out file, 59
xrxrh.out file, 59
xrxrp.out file, 59
XY-32 modulation, 178

Ζ

Z gradient coil, 199 zero constant, 96, 97 zero power RAM, 67 zero pulse length, 41 Z-filters, 112, 135 zgradpulse statement, 200